

Feel the freedom of cordless design.

TC2 Cordless Motor

SybronEndo's **EndoTouch TC2** is the next generation in cordless motors. The TC2 features a large LCD display, simple 5-key operation and a lightweight, ergonomic design to ensure comfort and ease of use during even the most delicate endodontic procedures. The auto reverse and alarm functions indicate the load is about to reach the pre-set torque level, allowing you to unload the file even before the auto reverse sets in. With the new digital operating interface, changing settings on the TC2 is much easier and more intuitive than ever before.

PART #	DESCRIPTION
815-9140	TC2 Cordless HP 110V Motor handpiece, EndoTouch push-button head 16:1, 110V power cord, battery charger, F-type spray nozzle, motorcap, quick start guide, instruction for use
815-9124	EndoTouch Push-Button Head 16:1
815-9125	EndoTouch Swing Latch Head 16:1
815-9142	M4 TC2 Handpiece Introducing the NEW M4 TC2 Handpiece for use with hand files


Small head with 6 position settings for access to entire oral cavity

Large, flush-mounted ON/OFF button for easy operation

Easy to read LCD display

Set speed, torque and other settings with ease

Easy to clean base for asepsis



TC2

Ready to use.

The TC2 supports most major brands of NiTi files along with 5 programs which can easily be set for each type of file.

Navigate with ease.

The 6 position type head allows the direction of the contra angle to be adjusted for increased comfort and easier access.

Versatile functions.

Choose to engage the "auto reverse", "auto stop" or "auto reverse off" functions based on your preference.

Plenty of working time.

Operate the TC2 continuously for approximately 1.2 hours before recharging.



Rotary NiTi files

For use with the TC2 cordless motor



Twisted not ground, for extreme flexibility. TF's efficient 3-file system allows you to quickly navigate even the most complicated canal anatomy with confidence.



SybronEndo's original NiTi file, the K3 has proven dependability featuring a unique cross section for superior strength and safety.



The newest edition to SybronEndo's NiTi family, K3XF features the flexibility of TF with the added safety and strength provided by K3's unique cross section.