






CONVENIENCE
COMFORT
HYGIENIC PROTECTION

MEDICAL CARE

Single-use infection
prevention products





Care with Confidence




TIDI focuses on you, the healthcare professional, continually evaluating your needs. So we understand that it's imperative to create an exam environment that presents a professional image and makes patients and clinicians feel comfortable. It's essential that exam products and personal protection equipment be easy, comfortable and convenient to use. And it's critical to minimize any chance for cross contamination.

In response, we engineer specific single-use infection prevention solutions to meet these needs. You can be confident in our reliable, consistent product quality, assured by our certification in ISO 9001 and EN/ISO 13485.

With these assurances, healthcare professionals confidently reach for TIDI's single-use infection prevention products to maintain a hygienic environment, put staff and patients at ease, and focus on care.

Clear Choices for Confident Care

Our simplified product organization includes three distinct levels that clearly describe our differentiated, value-added products, based on insight gleaned from healthcare professionals. This makes it easy to choose the best TIDI products for your procedures, your environment, and the comfort and convenience of you and your patients.



PRODUCT LEVELS

U Ultimate

High-performance products designed with thoughtful features for the highest-quality care, making the professional's job easier while enhancing the patient experience.

C Choice

A wide variety of convenient, high-quality products that offer comfort and protection on demand.

E Everyday

Economical single-use products that provide functional everyday protection.

TABLE OF CONTENTS

PATIENT APPAREL + DRAPING

- 4 Gowns
- 7 Capes
- 8 Shorts + Briefs + Vests + Slippers
- 9 Patient Drape Sheets
- 11 Emergency Blankets

PERSONAL PROTECTION EQUIPMENT



- Gowns
 - 14 Isolation Gowns
 - 17 Chemotherapy Gowns + Prep Mats
- 18 Facemasks
- Exam Gloves
 - 20 Nitrile Gloves
 - 21 Latex Gloves
- Personal Protective Eyewear
 - 22 Eye Shield
 - 23 Face Shield
- 24 Eyewear Packs & Dispensers
- 25 Dispensers + Dispenser Holders + TIDIShield Pillow Barriers
- 26 Cabinets + Stations
 - 28 8300 Cabinet Accessories
 - 28 Brackets
 - 29 Glove Dispensers
- 29 Blood Pressure Cuff Barriers

EQUIPMENT PROTECTION

- 30 Exam Table Barriers
- 31 Lab Countertop Barriers
- 33 Chiropractic Headrest Barriers
- 34 Equipment Drapes + Stretcher Sheets + Bed Sheets
- 35 Pillowcases
- 36 Towels + Scale Barriers + Tray Covers
- 38 Underpads
- 39 Baby Scale Barriers
- 39 Washcloths + Wipers

ACCESSORIES

- 39 Thermometer Sheaths
- 40 Wound Care
- 42 Infection Control
 - Sterilization Wraps
 - Autoclave Envelopes
 - Autoclavable Towels
 - Infectious Waste/Biohazard Bags
- 43 Patient Care
 - Specialty Patient Bibs
 - Cups
- 43 Aprons
- 44 Hospital Supply
 - Bedside/Suture Bags + Bath Mat + Bedpan Covers
 - Post-Mortem Bags + Shroud Kits
 - Veterinary Supply
- 44 Office Supply
 - X-Ray Storage Envelopes
 - Trash Can Liners
- 45 Glossary of Materials + Processes



PATIENT APPAREL + DRAPING

GOWNS

TIDI gowns are constructed to provide patient modesty and comfort while ensuring convenient access to your patients and hygienic coverage.

SHOWN

TIDI 918156

2XL 2-ply

scrim-reinforced tissue
gown with sewn shoulders
and reinforced neckline.




U Ultimate GOWNS


- The ultimate in patient modesty and comfort
- Styled and sized to fit, feel and cover like linen gowns
- Made from various high-performance materials
- Some feature liquid hold-out capability
- Multiple size options
- Optimal professional image for your practice
- Ideal for extended wear or multi-phase procedures

CODE	COLOR	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
950549	Dark Blue	70" x 45" (27" x 45")	2XL 50	Wrap-a-Gown SMS nonwoven. Generous-size, maximum coverage (a wrapped third sleeve prevents any gap without a belt). Sewn, reinforced shoulders for excellent durability. Repels liquid.
950545	Dark Blue	63" x 42" (30" x 42")	L 50	Patient Gown SMS nonwoven. Linen-like material. Sewn reinforced shoulders. Waist ties are attached.
980831	Medium Blue	54" x 42" (30" x 42")	L 50	Dexter material. Soft, pliable, leak-proof and durable. Crepe exterior absorbs and deters runoff.
980838	Light Blue	60" x 42" (30" x 42")	XL 50	2-ply scrim-reinforced tissue with sewn shoulders and sewn/reinforced neckline. Sewn, 4-ply scrim-reinforced tissue tie at waist, and nonwoven tie at neck.
918156	Light Blue	72" x 44" (36" x 44")	2XL 25	2-ply scrim-reinforced tissue with sewn shoulders and sewn/reinforced neckline. Sewn, 4-ply scrim-reinforced tissue tie at waist, and nonwoven tie at neck.
918158	Blue	80" x 44" (40" x 44")	3XL 25	2-ply scrim-reinforced tissue with sewn shoulders and sewn/reinforced neckline. Sewn, 4-ply scrim-reinforced tissue tie at waist, and nonwoven tie at neck.
920540	Light Blue	59" x 42" (30" x 42")	XL 50	3-ply scrim-reinforced tissue with glued shoulders.
920433	Light Blue	59" x 42" (30" x 42")	XL 50	2-ply scrim-reinforced tissue with glued shoulders.
918519	Sky Blue	56" x 42" (30" x 42")	XL 50	FABRICEL® poly/tissue; tissue inside fused to smooth heavyweight poly outside. Very strong, yet soft and leak-proof.
910521	Blue	63" x 42" (32" x 42")	XL 50	T/P/T, front or back opening. Leak Proof.
910540	Blue	80" x 40" (40" x 40")	3XL 25	T/P/T, PONCHO-style, sides open. Leak Proof.
910541	Blue	112" x 40" (56" x 40")	5XL 25	T/P/T, PONCHO-style, sides open. Leak Proof.
980949	Blue/transparent	62" x 48" (36" x 44")	XL 50	Poly barrier gown with ties. STAFF protection garment. Long sleeves. Head opening. Covers entire front, arms and upper back. Individually bagged.

*Girth is the material circumference; does not include any gap opening.



TIDI 950549 SMS
nonwoven Wrap-a-Gown's generous size and design offer maximum modesty. (Ultimate)




TIDI 950545 Patient Gown
SMS Nonwoven

C Choice GOWNS

Prints create a relaxed or fun environment for patient comfort and peace of mind for adults or pediatrics you provide care more easily.

Choice COLORS



Choice PRINTS



- Embossed; prevents leak through
- Adult and pediatric sizes
- Available in colors and prints
- Ample coverage, added strength and enhanced opacity
- Poly waist ties, stretch to desired length

CODE	COLOR/PRINT	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
910420	White	56" x 42" (30" x 42")	50	T/P/T, front or back opening
910520	Blue	56" x 42" (30" x 42")	50	T/P/T, front or back opening
910542	Peach	56" x 42" (30" x 42")	50	T/P/T, front or back opening
910528	Green	56" x 42" (30" x 42")	50	T/P/T, front or back opening
910536	Mauve	56" x 42" (30" x 42")	50	T/P/T, front or back opening
910510	Teal	56" x 42" (30" x 42")	50	T/P/T, front or back opening
980844	Sky Blue	56" x 42" (30" x 42")	50	Poly/tissue; tissue inside embossed to poly outside. Prevents leak through.
910416	White	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening
980830	Blue	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening

PRINTS

981436	Bugs & Things	40" x 36" (21" x 36")	50	FABRICEL® tissue/poly, pediatric , front or back opening with printed tissue outside. Very strong, yet soft and prevents leak through.
983836	Cat in the Hat [†]	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening
983936	Horton Hears a Who [†]	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening
981836	Schooltime	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening
981636	Under the Sea	40" x 36" (21" x 36")	50	T/P/T, pediatric , front or back opening
982442	Watercolors**	56" x 42" (30" x 42")	50	T/P/T, adult , front or back opening

**Matches cape, drape, table roll and towel prints

[†]Dr. Seuss Properties TM & © 2004 Dr. Seuss Enterprises, L.P.
All Rights Reserved.

E Everyday GOWNS

Everyday COLORS



- Reliable strength, coverage and opacity
- Patient modesty and hygienic security
- Poly waist ties, stretch for desired length

CODE	COLOR	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
910320	White	56" x 42" (30" x 42")	50	3-ply tissue, front or back opening
910326	Mauve	56" x 42" (30" x 42")	50	3-ply tissue, front or back opening
918507	Blue	56" x 42" (30" x 42")	50	3-ply tissue, front or back opening

*Girth is the material circumference; does not include any gap opening.

CAPE

Capes provide convenient coverage for specific needs. An ample fit offers modesty and security for patients while allowing access for examinations and procedures.

SHOWN
TIDI 981111
SMS nonwoven
Ultra Mammo Cape,
poncho-style with
convenient front snaps.



U Ultimate CAPES

- The ultimate in patient **modesty** and **comfort**
- Styled and sized to fit, feel and cover like linen
- Strength, softness and warmth
- Optimal professional image for your practice
- Ensure patients' physical comfort and peace of mind
- Ideal for extended wear or multi-phase procedures

CODE	COLOR	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
981111	Dark Blue	96" x 23" (48" x 23")	4XL 50	Ultra Mammo Cape, SMS nonwoven. PONCHO-style with convenient front snaps and open sides. Generous size for modesty.
980858	Medium Blue	56" x 21" (30" x 21")	- 100	Dexter material. Front or back opening. Soft, pliable, leak-proof and durable. Crepe exterior absorbs and deters runoff.
941113	Dark Blue	One size fits most	- 25	TIDI stretchable cape, Nonwoven. Latex free elastic band at top provides adjustable sizing. Openings on both sides of cape allows for easy access.

Capes are constructed to resist ripping and tearing. Choose poncho-style or front/back opening for coverage and convenience.




TIDI 941113 Cape, nonwoven.

C Choice CAPES

- Embossed; prevents leak through
- Available in colors and print
- Enhanced strength, better fit, softness and opacity

CODE	COLOR/PRINT	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
910415	White	56" x 21" (30" x 21")	100	T/P/T, front or back opening
910417	White	63" x 21" (32" x 21")	100	T/P/T, Plus Size XL , front or back opening
910532	Peach	56" x 21" (30" x 21")	100	T/P/T, front or back opening
910515	Blue	56" x 21" (30" x 21")	100	T/P/T, front or back opening
910517	Blue	63" x 21" (32" x 21")	100	T/P/T, Plus Size XL , front or back opening
910516	Mauve	56" x 21" (30" x 21")	100	T/P/T, front or back opening
910501	Teal	56" x 21" (30" x 21")	100	T/P/T, front or back opening
982430	Watercolors**	56" x 21" (30" x 21")	100	T/P/T, front or back opening
980854	Sky Blue	56" x 21" (30" x 21")	100	Poly/tissue, tissue <i>inside</i> embossed to poly <i>outside</i> , front or back opening
918230	White	44" x 20" (22" x 20")	200	T/P/T, PONCHO-style, open sides



Plus Size Cape is 2" wider and has a 2" wider neck opening. **Item #910517** shown in blue

Choice PRINTS



Watercolors

Choice COLORS



*Girth is the material circumference; does not include any gap opening.
**Matches gown, drape, table roll and towel prints

E Everyday

Everyday COLORS



White



Maude



Blue

CAPES

- Embossed
- Reliable coverage, strength and opacity
- Patient modesty and hygienic protection

CODE	COLOR	GIRTH* X LENGTH (INDUSTRY SIZE)	UNIT QUANTITY	DESCRIPTION
910315	White	56" x 21" (30" x 21")	100	3-ply tissue, front or back opening
910316	Maude	56" x 21" (30" x 21")	100	3-ply tissue, front or back opening
918407	Blue	56" x 21" (30" x 21")	100	3-ply tissue, front or back opening

*Girth is the material circumference;
does not include any gap opening

PATIENT APPAREL

SHORTS + BRIEFS + VESTS + SLIPPERS

SMS nonwoven material provides fabric-like thickness and softness, performance and durability for extended wear and demanding applications.



Exam/rehab vest

- Breast telemetry pocket
- Sturdy snaps
- Side vents for stress cable



Exam shorts

- Generous cut
- Elastic waist
- Heavyweight material
- Finished leg hems



Bikini brief

- All edges elastic and finished
- Elastic-finished side bands
- Comfortably fits a range of sizes
- Individually bagged



Slippers

- Crepe paper
- Economical
- Convenient

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
960350	Dark Blue	One size fits most	200 (individually bagged)	SMS nonwoven bikini
960202	Dark Blue	One size fits most	50 (5 bags of 10)	SMS nonwoven rehab vest
960400	Dark Blue	Medium (28" – 34")	100 (4 bags of 25)	SMS nonwoven exam shorts
960401	Dark Blue	Large (32" – 40")	100 (4 bags of 25)	SMS nonwoven exam shorts
960402	Dark Blue	X-Large (38" – 48")	100 (4 bags of 25)	SMS nonwoven exam shorts
960404	Dark Blue	Plus size (Up to 60")	50 (2 bags of 25)	SMS nonwoven exam shorts
950520	Brown	One size fits most	2000 (1000 pair)	Crepe paper slippers

PATIENT DRAPE SHEETS

Offer comfort, modesty and coverage with our drape sheet products, with absorbent tissue and added poly layers for additional barrier protection.

Shown:
TIDI 918316* (info on pg 10)
Mauve 2-ply tissue
patient drape sheet.

TIDI 910420 T/P/T gown.
(info on pg 6)

TIDI 919355 (info on pg 31)
FABRICEL® pillowcase.

TIDI 917406 (info on pg 33)
Mauve 3-ply tissue,
poly-backed towel.



U Ultimate PATIENT DRAPE SHEETS

- Our strongest patient drapes
- Generous coverage
- Patient modesty and comfort
- Some include poly protection that prevents leak through

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919374	Blue	40" x 48"	50	T/P/T prevents leak through
919375	Blue	40" x 60"	50	T/P/T prevents leak through
951939	Blue	40" x 72"	100	T/P/T prevents leak through
919370	Blue	40" x 90"	50	T/P/T prevents leak through
917270	White	18" x 26"	300 (6 boxes of 50)	T/P/T sterile field drape, non-fenestrated prevents leak through
917272	White	18" x 26"	300 (6 boxes of 50)	T/P/T sterile field drape, 3" fenestration prevents leak through
950334	White	40" x 48"	100	3-ply tissue, procto drape, 5" fenestration
918310	White	40" x 48"	100	3-ply tissue
918311	White	40" x 60"	100	3-ply tissue
918309	White	40" x 72"	50	3-ply tissue
918312	White	40" x 90"	50	3-ply tissue

C Choice PATIENT DRAPE SHEETS

- Extra strength and opacity
- Hygienic protection
- Range of sizes for the right coverage
- Range of colors and print for a comfortable environment

Choice PRINTS



Watercolors

Choice COLORS



White



Blue



Peach



Mauve



Teal



Lavender



Green

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
918305	White	36" x 48"	100	2-ply tissue
918302	White	40" x 48"	100	2-ply tissue
918313	Blue	40" x 48"	100	2-ply tissue
918316	Mauve	40" x 48"	100	2-ply tissue
918318	Green	40" x 48"	100	2-ply tissue
918322	Peach	40" x 48"	100	2-ply tissue
918321	Teal	40" x 48"	100	2-ply tissue
918320	Lavender	40" x 48"	100	2-ply tissue
918306	White	36" x 60"	100	2-ply tissue
918303	White	40" x 60"	100	2-ply tissue
950450	White	40" x 72"	50	2-ply tissue
918328	White	13" x 24"	500	2-ply tissue, breast drape
980823	White	24" x 40"	200	2-ply tissue, breast drape
980893	Mauve	24" x 40"	200	2-ply tissue, breast drape
982448	Watercolors**	40" x 48"	100	2-ply tissue
982460	Watercolors**	40" x 60"	100	2-ply tissue
918022	White	40" x 48"	100 (4 packs of 25)	2-ply tissue Space Saver drapes, banded condensed packs

** Matches cape, gown, table roll and towel prints



Coordinate gowns, capes, drapes, towels and exam table rolls in our Watercolors print for a professional image and a comforting exam environment. Ideal for women's health care.

E Everyday

PATIENT DRAPE SHEETS

- Ample coverage
- Dependable strength provides modesty and protection
- 100% recycled fiber

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9810824	White	40" x 48"	100	2-ply tissue
9810826	White	40" x 60"	100	2-ply tissue
9810827	White	40" x 72"	50	2-ply tissue



EMERGENCY BLANKETS

- Ideal for police, fire, emergency, and industrial safety use
- Disposable, fluid resistant, and comfortable to use
- Vibrant yellow color, maximum visual safety benefit

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
980043	Bright Yellow	58" x 90"	18	Poly-backed foam drape, insulated, individually bagged
980077	Bright Yellow	56" x 90"	24	Poly-backed 2 ply embossed tissue, individually bagged

PERSONAL PROTECTION EQUIPMENT

Personal protective equipment (PPE) is defined by the Occupational Safety and Health Administration, or OSHA, as “specialized clothing or equipment, worn by an employee for protection against infectious materials.”



Infection prevention planning is a core priority to providing the safest care for the patient and the healthcare professional. With proper care and personal protection equipment, the spread of infections by hospital personnel can be greatly reduced.


TIDI Products' P2® Personal Protection Equipment is the premier solution to providing your healthcare workers with a high level of protection from hazardous materials. P2® provides a selection of quality facemasks, gowns, gloves and storage cabinets that both assist in compliance with OSHA preparedness requirements, and provide a safe work environment for your employees and patients.



Products branded "HiRisk™" provide the highest levels of protection offered within the P2® line. These products should be used when anticipating exposure to high risk situations.




Products branded "SafetyPlus™" are great for everyday use. These products provide a great level of protection with increased comfort.



SHOWN: 9010 Facemask, 8576 Gown



SHOWN: 8573 Gown



SHOWN: HiRisk Gloves, Resposables Face Shield, 9010 Facemask

WHY PPE?

Protecting healthcare personnel and patients from infectious diseases is now more important than ever, as these groups are exposed to infections within healthcare settings daily. Recent studies reveal:

- Infections contracted in hospitals are the fourth largest killer in America. Every year in this country, two million patients contract infections¹
- The annual direct medical cost for hospital acquired infections to US hospitals ranges from \$28.4 to \$33.8 billion²
- As many as 10% of patients contract infections while visiting a hospital³
- Patients inflicted with MRSA have hospital costs that more than triple the cost of other patients⁴

¹ Source: Neu HC. The crisis in antibiotic resistance. Science 1992;257:1064- 1073

² Source: "Scott, R. Douglas. "The Direct Medical Costs of Healthcare-Associated Infections in U.S. Hospitals and the Benefits of Prevention." Mar. 2009. CDC. 27 Aug. 2009 http://www.cdc.gov/ncidod/dhqp/pdf/Scott_CostPaper.pdf.

³ Eisenstadt E, Carlton BC, Brown BJ. Gene mutation. In: Gerhardt P, Murray RGE, Wood WA, Krieg NR, eds. Methods for General and Molecular Bacteriology. Washington, DC: American Society for Microbiology; 1994:297-316

⁴ Source: Shopsis B, Mathema B, Martinez J, et al. Prevalence of methicillin-resistant and methicillin- susceptible Staphylococcus aureus in the community. J Infect Dis 2000;182:359-362

ISOLATION GOWNS

Shown:
TIDI® 8575
(info on page 15)
SafetyPlus™ Polyethylene Gown



ISOLATION GOWNS

TIDI's line of isolation gowns has expanded to include five different materials of isolation gowns. TIDI understands that not all situations call for the same isolation gown. Our gowns come in different materials to provide the right isolation gown for the task at hand.

HiRisk™ IMPERVIOUS GOWN

Our impervious gowns are made from soft embossed polyethylene which provides a level of protection against infectious materials. The design ensures comfort and protection.

Features:

- Closed loop neck for easy donning and tear off
- Open back for increased airflow
- Sized generously for more diverse shapes
- Thumbloop cuffs keep gown in place
- Latex Free
- Passes the ASTM F1671 and ASTM F1670 standards
- Longer and wider ties for a safe and secure fit
- Impervious Material prevents pass-thru of fluids

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8576	Light Blue	Universal	5 dispensers of 15 gowns; 75 gowns/case	HiRisk™ Impervious Gown. Open back, thumbloop cuff.



TIDI® 8576 HiRisk™
Impervious Gowns with
loop neck design




Features thumbloop cuffs that
provide added protection to the
vulnerable wrist areas, and make
donning gloves easy



TIDI® 8576 HiRisk™ Impervious
Gown Dispenser box is compact
and uses 75% less space than
conventional isolation gown
packaging

SafetyPlus™ POLYETHYLENE GOWN



Comfortable. Safe. Cost effective.

For those looking for the right gown at the right price, the SafetyPlus™ Polyethylene Gown is the perfect solution. This gown offers durable barrier protection from hazardous materials and infectious diseases at an economy price without sacrificing safety or comfort.

Features:

- Closed loop neck for easy donning and tear off
- Open back for increased airflow
- Sized generously for more diverse shapes
- Thumbloop cuffs keep gown in place
- Passes the ASTM F1671 and ASTM F1670 standards
- Longer and wider ties for a safe and secure fit
- Latex Free

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8575	Light Blue	Universal	5 dispensers of 15; 75 gowns / case	SafetyPlus™ Polyethylene Gown, closed neck loop, open back

TIDI® 8575 SafetyPlus™ Polyethylene

SafetyPlus™ POLY COATED GOWN



These gowns are made from a soft breathable material with an impervious coating to protect the healthcare worker. Recommended for use when the potential of moderate fluids and/or contact isolation is anticipated.

Features:

- White Polyethylene coated Polypropylene material
- Thumbloops hold gown in place
- Open back for increased airflow
- Closed loop neck for easy donning and tear off
- Latex Free



TIDI® 8578 SafetyPlus™ Poly Coated Gown Dispenser box is compact and uses 75% less space than conventional isolation gown packaging

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8578	White	X-Large	5 dispensers of 8 gowns; 40 gowns/case	SafetyPlus™ Poly Coated Gown. Open back, thumbloop cuff.

TIDI® 8578 SafetyPlus™ Poly Coated Gown

SMS OVER-THE-HEAD GOWN

NEW!




TIDI® 8573
SMS Over-the-Head Gown

Designed with medium weight SMS material and longer coverage in the back than other PPE gowns, this gown protects healthcare workers where others fall short.

Features:

- Thicker 25 gsm SMS material protects from the potential risk of fluid penetration
- Greater coverage in back, the cape comes down farther as compared to other gowns that fall just past the neck Slides on easily over the head
- Convenient thumb loops assist with donning and help gown stay in place under gloves
- Tear-off at back/neck for easy removal
- Gown naturally curves around the user



PPE storage cabinets with doors can conveniently hold up to 15 gowns per bin, while doorless PPE storage cabinets can hold up to 30!


CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8573	Yellow	X-Large	10 bags of 10 gowns; 100 gowns/case	SMS Over-the-Head Gown

SafetyPlus™ AAMI LEVEL 2 SMS GOWN

Features:

- Yellow three-layer fabric
- Ultrasonically sealed seams for increased protection
- Closed back for full protection
- Elastic cuff to prevent sleeve from moving up
- Ties at neck and waist for secure fit
- Latex Free

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
8580	Yellow	X-Large	10 bags of 10 gowns; 100 gowns/case	SafetyPlus™ AAMI Level 2 SMS Gown, closed back, elastic cuff



TIDI® 8580
SafetyPlus™ AAMI
Level 2 SMS Gown

SafetyPlus™ SMS GOWN

SafetyPlus™ SMS Gowns are made from materials that are breathable and comfortable while still providing protection for the healthcare worker. Recommended for use when the potential of light to moderate fluids and/or contact isolation is anticipated.

Features:

- Yellow three-layer fabric
- Elastic cuffs to prevent sleeve from moving up arm
- Closed back for full protection
- Ties at neck and waist for secure fit
- Latex Free



CODE	COLOR/ PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8579	Yellow	X-Large	10 bags of 10 gowns; 100 gowns/case	SafetyPlus™ SMS Gown. Closed back, elastic cuff.

TIDI® 8579 SafetyPlus™ SMS Gown


SafetyPlus™ SPUNBOND GOWN

Light and comfortable Spunbond Gown. Great for compliance and contact precautions when there is little to no fluid protection needed.

Features:

- Elastic cuff to prevent sleeve from moving up
- Closed back for full protection
- Ties at neck and waist for secure fit
- Latex Free

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
8570	Yellow	Universal	10 bags of 10 gowns; 100 gowns/case	SafetyPlus™ Spunbond Gown, closed back, elastic cuff



TIDI® 8570
SafetyPlus™
Spunbond Gown



TIDI® HiRisk™
Chemotherapy Gowns



ULTIMATE PROTECTION


- Closed back design for ultimate protection
- Back is lightweight SMS material for breathability
- Meets OSHA, ASHP and ONS recommendations for wearing gowns when handling chemotherapy drugs
- Polyethylene coated front and sleeves
- Non-linting materials
- Long sleeves with knit cuffs
- Ties at waist with Velcro closure at neck

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
8586	Blue	Universal	10 gowns/poly bag; 30 gowns/case	HiRisk™ Chemotherapy Gown
8587	Blue	X-Large	10 gowns/poly bag; 30 gowns/case	HiRisk™ Chemotherapy Gown

PREP MAT



Multiple Layers of Absorbency




SafetyPlus™ PREP MAT

- Multiple layers of absorbency
- Top — specially designed wicking layer
- Middle — three ply absorbent layer
- Bottom — poly backing prevents leak through
- Poly bagged within carton to keep mat clean from lint, fiber and dust
- Economical option
- Latex free

CODE	COLOR	SIZE	CASE QUANTITY	DESCRIPTION
983930	Light Blue	13" x 18"	250 mats/case	SafetyPlus™ Prep Mat

FACEMASKS

Shown:
Item #9010



FACEMASKS

TIDI's fluid resistant facemasks were designed from end user feedback to offer best-available protection and comfort. Manufactured in the USA using soft, breathable materials, our facemasks feature comfortable earloops for long-wear, and surpass protection standards.

Features:



- Made in the USA
- Tested using the ASTM F2100-07 test standard
- Earloop design for quick and easy donning
- Fog free strip (Item #9020 only)
- Made from lightweight and breathable materials
- Downward folding pleats to prevent pooling
- Odorless
- Latex free

	CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
HIGH BARRIER	9010	Light Blue	Universal	10 boxes of 50 500 facemasks/case	HiRisk™ Fluid Resistant Facemask. High Barrier Protection.
HIGH BARRIER	9020	Light Blue	Universal	10 boxes of 40 400 facemasks/case	HiRisk™ Fog Free Fluid Resistant Facemask. High Barrier Protection.
MEDIUM BARRIER	9030	Blue	Universal	10 boxes of 50 500 facemasks/case	SafetyPlus™ Fluid Resistant Facemask. Medium Barrier Protection.
LOW BARRIER	9040	Yellow	Universal	10 boxes of 50 500 facemasks/case	SafetyPlus™ Procedure Facemask. Low Barrier Protection.

ASTM F2100-07 PERFORMANCE CLASSES FOR FACEMASKS

TIDI facemasks are tested with the ASTM F2100-07 “Standard for Performance of Materials Used in Medical Facemasks.” This is the only specific standard available that provides performance classes for facemasks to fall into. The performance classes range from low, moderate, and high barriers.

REQUIREMENTS	LOW BARRIER STANDARDS	TIDI'S 9040 FACEMASKS	MODERATE BARRIER STANDARDS	TIDI'S 9030 FACEMASKS	HIGH BARRIER STANDARDS	TIDI'S 9010/9020 FACEMASKS
ASTM F1862 Resistance to penetration by synthetic blood (mm HG)	80	80	120	120	160	160
ASTM F2101 Bacterial Filtration Efficiency (BFE)	≥ 95	≥ 99	≥ 98	≥ 99	≥ 98	≥ 99
ASTM F2299 Particulate Filtration Efficiency (PFE)	Not required	≥ 99	≥ 98	≥ 99	≥ 98	≥ 99
Delta P Breathability	< 4.0	< 2.0	< 5.0	< 2.0	< 5.0	< 5.0
Flame Spread	Class 1	Class 1	Class 1	Class 1	Class 1	Class 1



9010



HiRisk™
Fluid Resistant Facemask
High Barrier Protection

9020



HiRisk™ Fog Free
Fluid Resistant Facemask
High Barrier Protection

9030



SafetyPlus™
Fluid Resistant Facemask
Medium Barrier Protection

9040



SafetyPlus™
Procedure Facemask
Low Barrier Protection



Perfect for use in areas such as:

- ICU
- Central Sterilization
- SPD
- Emergency Department
- Isolation



Customized packaging allows
for easy top and front access
when stored in cabinets, shelves,
and drawers.



Fog free strip prevents safety
eyewear from fogging
(Item #9020 fog free strip shown)

SPECIALTY GLOVES

TIDI's glove line includes
Latex Free Nitrile Patient
Exam Gloves



**HiRisk™
NITRILE
EXAMINATION
GLOVES**

Nitrile Examination Gloves provide a non-latex alternative to latex gloves while still providing the same strong amount of protection. Our Nitrile gloves have been tested to applicable ASTM standards as pertains to protection when dealing with chemotherapy drugs and other threats.

- Latex Free
- Powder-Free
- Tested using ASTM D6978-05, the most current and only chemotherapy test available specifically for testing gloves against chemotherapy drugs
- 8 mils at fingertips
- Ambidextrous for quick and easy donning
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection



Nitrile gloves are chemotherapy tested and are also great for the EMS setting.

PRODUCT NUMBER	SIZE	UNIT QUANTITY
940001-1	Small	10 dispenser boxes of 50 gloves; 500 gloves/case
940002-1	Medium	10 dispenser boxes of 50 gloves; 500 gloves/case
940003-1	Large	10 dispenser boxes of 50 gloves; 500 gloves/case
940004-1	X-Large	10 dispenser boxes of 50 gloves; 500 gloves/case



HiRisk™ POWDERFREE LATEX EXAMINATION GLOVES

- Tested using ASTM D6978-05, the most current and only chemotherapy test available specifically for testing gloves against chemotherapy drugs
- 18 mils thick at fingertips, eliminates need for double gloving
- Fitted for right and left hands, reduces hand fatigue
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection

PRODUCT NUMBER	SIZE	UNIT QUANTITY
932477-1	Small	10 dispenser boxes of 50 pair; 1000 gloves/case
932480-1	Medium	10 dispenser boxes of 50 pair; 1000 gloves/case
932483-1	Large	10 dispenser boxes of 50 pair; 1000 gloves/case
932486-1	X-Large	10 dispenser boxes of 50 pair; 1000 gloves/case



SafetyPlus™ POWDERFREE LATEX EXAMINATION GLOVES

- 14 mils thick at fingertips, eliminates need for double gloving
- Ambidextrous for quick and easy donning
- Textured for added grip
- Beaded cuff
- 12 inches long for additional protection
- UL Classified as meeting NFPA requirements

PRODUCT NUMBER	SIZE	UNIT QUANTITY
BS0450-1	Small	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0460-1	Medium	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0470-1	Large	20 dispenser boxes of 50 gloves; 1000 gloves/case
BS0480-1	X-Large	20 dispenser boxes of 50 gloves; 1000 gloves/case

Latex Examination Gloves deliver superior protection for **Oncology, Central Service, Decontamination, Environmental or EMS Services.** Made with natural rubber latex, all gloves have been tested to applicable ASTM standards.

**NOW
Chemotherapy Tested**



EXTRA LONG PROTECTION

Measures 12 inches from tip to cuff for full wrist protection

PERSONAL
PROTECTIVE
EYEWEAR






**FLIPEASE®
EYE SHIELD**

- No assembly required
- FLAT for compact storage, FLIP to wear
- Extremely lightweight and comfortable
- Wrap-around design for front and side eye protection
- Optical-grade, anti-static, lens

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9300FE-50	5 assorted colored frames/clear lens	One size fits most	50 (2 boxes of 25)	One piece disposable eyewear; stores flat and assembles with one simple step

FlipEase® Single-use Eyewear



You dispose of your gloves and gowns every use, why not eyewear?

RESPOSABLES® EYE SHIELD®

- Weighs only 3/10 of an ounce!
- "Wrap-around" design provides both front and side eye protection
- Available in pre-assembled, ready-to-wear option or as two separate components (lens and frame)
- Optical-grade replacement lenses; available in clear, UV-tinted and amber-tinted (blue-light absorbent)
- Separately available replacement frames and lenses for maximized cost-effective replenishment



CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9210-50	Clear	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210-100	Clear	One Size	Dispenser box of 100 lenses	Disposable lens for use with 9211-100
9210-250	Clear	One Size	250 (10 envelopes of 25)	Disposable lens for use with 9211-100
9210LT-50	Bronze	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210LT-250	Bronze	One Size	250 (10 envelopes of 25)	Disposable lens for use with 9211-100
9210CL-50	Amber	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210CL-200	Amber	One Size	200 (20 envelopes of 10)	Disposable lens for use with 9211-100
9211-50	Assorted frames	One Size	50 (10 bags of 5 frames)	Reusable frames for use with 9210's
9211-100	Assorted frames	One Size	100 (10 bags of 10 frames)	Reusable frames for use with 9210's
9452	OFFICE PACK	One Size	10 frames and 25 lenses	Disposable eye shields
9210A-50	Assorted color frames/clear lenses	One Size	50 assembled eye shields	Assembled eye shields

RESPOSABLES® FACE SHIELD®

- Weighs less than 1.5 ounces!
- Designed to maximize airflow and prevent fogging
- Provides over-the-top, side and front face protection
- Disposable, optical grade, anti-static replacement shields
- Includes a clear, snap-in nosepiece for support
- Separately available replacement visors and lenses for maximized cost-effective replenishment



Resposables Face Shield-
Patent #6,564,804

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
2210-50	Clear	13" x 7"	50	Disposable optical grade shield; for use with 2211-100
2210-100	Clear	13" x 7"	100	Disposable optical grade shield; for use with 2211-100
2211-50	Assorted	One size fits most	50 (5 bags of 10 visors)	Reusable visor; for use with 2210-100
2211-100	Assorted	One size fits most	100 (10 bags of 10 visors)	Reusable visor; for use with 2210-100
2452	Assorted Visors/Clear Lenses	One size fits most	3 visors; 15 lenses	Reusable visor and optical grade face shield

DISPOSABLE FACE SHIELD • Preassembled • Adjustable strap for comfortable fit • Latex free

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
FS850-30	Clear	Standard 11.5" wide x 9" long	30 (1 box of 30)	One piece face shield with foam headband, elastic strap and optical-grade lens
FS850-120	Clear	Standard 11.5" wide x 9" long	120 (4 boxes of 30)	One piece face shield with foam headband, elastic strap and optical-grade lens
FS9605-100	Clear	Full Size 14" wide x 11.5" long	100	One piece face shield with foam headband, elastic strap and optical-grade lens

ACCESSORIES

CODE	COLOR	SIZE	UNIT QUANTITY
9505-50	Neck Strap 5 Assorted Colors	One size	50 Neck Straps
9505-100	Neck Strap Assorted Colors	One size	100 Neck Straps
9810-50	Side Shield Clear	One size	50 Side Shields (25 pairs)
9902-100	Nose Clips Assorted Colors	One size	100 (50 per bag/2 bags per box) Disposable Nose Clip
9902-1000	Nose Clips Assorted Colors	One size	1000 (50 per bag/20 bags per box) Disposable Nose Clip

*Trademarks of TIDI Products, LLC.

EYE SHIELD DISPENSERS



Responsables Eye Shields® feature reusable frames and disposable lenses*



The New TIDI Tower™ dispenses assembled TIDISHIELD™ Disposable Eye Shields

PERSONAL PROTECTIVE EYEWEAR IN CONVENIENT, COMPACT DISPENSERS

Clean eye protection is not always readily available to healthcare providers. To solve this problem, TIDI® created convenient, portable Eye Shield Dispensers that can be easily located at point-of-use.

READILY AVAILABLE AT POINT-OF-USE FOR QUICK RESPONSE

- Eye Shield Dispensers are easily portable so eye protection is always available
- Eyewear is ready to use—no need to search for lenses and frames
- Position the dispenser at point-of-use using a plastic dispenser holder (two options available)

PACKS IDEAL FOR USE IN:

- Emergency Rooms
- Operating Rooms
- Labor & Delivery
- Burn Units



The NEW Eyewear Packs from TIDI® provide protective eyewear when and where you need it.

EYE SHIELDS



TIDI TOWER™ (patent pending)
Containing **TIDISHIELD™**
DISPOSABLE EYE SHIELDS
9210A-100

25 assembled eye shields per dispenser, 4 dispensers per case
Fits in Dispenser Holder 7003
6.75" W x 17" H x 5.75" D



TIDISHIELD™ DISPOSABLE EYE SHIELD LENSES & FRAMES
9209

100 frames and 100 lenses per dispenser, 4 dispensers per case
Fits in Dispenser Holder 7001
9.0" W x 8.25" H x 5.75" D

DISPENSERS


RESPOSABLES* EYE SHIELD®
500 LENS DISPENSER
9210-500DSP

500 lenses per dispenser
 Fits in Dispenser Holder 7001
 9.0" W x 2" H x 2.75" D


RESPOSABLES* EYE SHIELD®
250 LENS DISPENSER
9210-250DSP

250 lenses per dispenser
 Fits in Dispenser Holder 7001
 9.0" W x 2" H x 5.625" D


FLIPEASE® EYE SHIELD DISPENSER
9300FE-50

50 individually bagged eye shields
 (2 boxes of 25)
 Fits in Dispenser Holder 7001
 9.25" W x 7.625" H x 4.625" D


RESPOSABLES* EYE SHIELD®
FRAMES DISPENSER
9211-400

100 frames per dispenser,
 4 dispensers per case
 Fits in Dispenser Holder 7001
 9.0" W x 8.25" H x 5.75" D

DISPENSER HOLDERS


TIDI® EYE SHIELD DISPENSER HOLDER
7001

2 dispenser holders per case
 Adheres with 3M Command Strips®
 9.75" W x 10.25" H x 6" D


TIDI® EYE SHIELD
TOWER DISPENSER HOLDER
7003


2 dispenser holders per case
 Adheres with 3M Command Strips®
 6" W x 12" H x 7.5" D

TIDISHIELD™ DISPOSABLE PILLOW BARRIERS

NEW!



Uniquely designed, the premium layered poly-backed material helps keep patients protected. In addition, the durable material means your pillows are cleaner and last longer so you spend less on laundering costs.

**REDUCE RISK. SAVE TIME & MONEY**

- Non-corrugated packaging
- Easy to grab, easy to apply
- Extends life of the pillow and saves on laundering costs

TIDISHield™ Disposable Pillow Barriers

Product Code: 6355

Color: White

Size: 21" x 30"

Unit Quantity: 2 boxes of 50 / 100 per case

IMPROVE PATIENT SAFETY & COMFORT

- Effective barrier against bacteria²
- Fluid-resistant material
- Soft, linen-life feel
- Single-use



CABINETS + STATIONS



Shown:
TIDI® P2® 8468
ICU Isolation Station



TIDI® P2® 8569
ICU Isolation Station



STORAGE UNITS

- Keep personal protection and patient care equipment in a clean and organized space at point of use
- Constructed with ABS plastic which is durable, mar-resistant, and easy to clean
- Storage units are all compatible with Glove dispensers. Dispensers provide added space and easy accessibility to gloves

P2® ICU ISOLATION STATION

for glass, door, or wall mounting (for door mounting, add item # 8566)



CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8468	20.5" H x 6.75" D x 17.75" W	1 station/case	ICU Isolation Station containing: 1 glass-mount bracket with adhesive strips 4 adjustable shelves 2 pockets in door 2 storage bins 1 ambu bag shelf 1 glove dispenser

P2® ISOLATION STATION

for glass, wall, or over-the-door mounting



CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8569	20.5" H x 6.75" D x 17.75" W	1 station/case	Isolation Station containing: Over-the-door brackets 4 adjustable shelves 2 pockets in door 2 storage bins 1 ambu bag shelf 1 glove dispenser

P2® STATION

for glass, wall, or over-the-door mounting (for glass mounting, add item #8461)

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8565	20.5" H x 6.75" D x 17.75" W	1 station/case	Station containing: 4 adjustable shelves 2 pockets in door 1 storage bin 1 ambu bag shelf

**P2® CABINET**

for wall mounting

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8560	12.75" H x 5.5" D x 13.5" W	1 cabinet/case	Cabinet containing: 3 adjustable shelves 2 pockets in door

**P2® CABINET**

for glass, door or wall mounting (see accessories for mounting options)

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8300	20.5" H x 4.75" D x 19.75" W	1 per case	Made of quartz high impact plastic Meets fire codes Made in USA Door, wall, and glass mount options available

**P2® PROTECTION ORGANIZER ACCESSORIES**

TIDI P2® glove box dispenser with shelf attaches to protection organizer item #8300

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8301	5.5" H x 4" D x 22.5" W	3 per case	Holds one box of gloves with additional shelf to hold medical supplies Attaches to bottom of Protection Organizer Made of quartz ABS plastic



8300 CABINET ACCESSORIES

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8330	17.75" H x 34" W	2 per case	Fire Door Hanger Holds 8300 cabinet. Attaches to top of fire door. Made of quartz powder coated metal.
8331	4.5" H x 22" W	2 per case	J Hook Glass Hanger Holds 8300 cabinets on a wall or on glass. Back of hanger contains double-sided, semi-permanent tape. Made of quartz powder coated metal.
8332	2.5" H x 5.5" W	2 per case	J Hook Wall Hanger B-Hanger accessory for walls or doors. Mounting holes. Made of quartz powder coated metal.
8333	12" H x 5.5" W	2 per case	Door Hanger 12" Door hanger accessory for use with 8300 cabinets Hooks over the top of a door. Made of quartz powder coated steel.
8334	11.25" H x 8.5" W	2 per case	Sign Holder Holds signage. Clips on to front of protective organizer. Made of PETG plastic.
8301	5.5" H x 4" D x 22.5" W	3 per case	Glove Box Dispenser with Shelf Holds one box of gloves with additional shelf to hold medical supplies. Attaches to bottom of Protection Organizer. Made of quartz ABS plastic.

8330



8332



8333



8334



8301



8331



Printed material available
on tidproducts.com

BRACKETS

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8563	17.75"	1 set of 2 vertical extension brackets/case	Vertical Extension Bracket for Isolation Station Extension hardware accommodates taller hospital doors so PPE is more easily reached by staff. For use with door mount isolation station #8569.
8564	36"	1 horizontal extension bracket/case	Horizontal Extension Bracket for Isolation Station Extension hardware accommodates most hospital doors with built-in door damper slots. For use with door mount isolation station #8569.
8566		1 door hang bracket	Door Hang Bracket for Isolation Station Bracket converts #8468 and #8565 stations to allow for door mounting.
8461	3" H x 12" W	2/case	Glass-Mount Bracket 2 glass-mount brackets with adhesive strips. For use with wall mount station #8565.

8563



8564



8566



8461



GLOVE DISPENSERS

CODE	SIZE	UNIT QUANTITY	DESCRIPTION
8555-SA	7" H x 4.5" D x 12" W	10 glove dispensers/case	P2® Side Loaded Glove Dispenser For use with any station or cabinet. Can be stacked for multiple glove dispensing.
8562	7" H x 4.5" D 12" W	10 glove dispensers/case	P2® Glove Dispenser For use with any station or cabinet or as stand-alone dispenser with double glove bracket #8554-B.
8554-B	3.5" H x 4.5" x 12" W	5 glove brackets/case	P2® Double Glove Bracket Stand-alone bracket for double glove dispensing with dispenser #8562.

8555-SA



8562



8554-B



BLOOD PRESSURE CUFF BARRIERS



U Ultimate BLOOD PRESSURE CUFF BARRIERS

- An effective and convenient way to minimize the risk of transmitting infections via blood pressure cuffs
- 99.9% protection against bacteria pass-through (independent tests)
- Easy to apply
- Will not alter BP readings
- Made with soft FABRICEL® tissue fused to leak-proof polyethylene film
- Available in three sizes

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
916113	White	13" x 5-3/8," Pediatric	300 (6 boxes of 50)	FABRICEL® tissue/poly
916114	White	18-1/8" x 6-3/4," Adult Regular	300 (6 boxes of 50)	FABRICEL® tissue/poly
916115	White	22-1/2" x 8," Adult Large	150 (3 boxes of 50)	FABRICEL® tissue/poly

EQUIPMENT PROTECTION

EXAM TABLE BARRIERS

Provide comfortable crepe surfaces, in widths that cover exam table sizes from pediatrics to bariatrics.

Shown:
TIDI 981918
Schooltime printed
crepe exam table rolls.

TIDI 981836
Schooltime printed
T/P/T pediatric gown.



U Ultimate EXAM TABLE BARRIERS

- The ultimate in both **strength** and **softness**
- Designed to provide maximum patient comfort and protection
- Embossed or creped material provides quiet softness for patients
- **Poly-backing** provides extra strength, leak-proof protection, and reduced clean-up time for your staff
- **Perforated** options add convenience and neatness
- 99.9% protection against bacteria pass-through (independent tests)

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
916218	White	18" x 100'	12	Embossed FABRICEL® tissue fused with a layer of polyethylene film on the back. PERFORATED every 30" for neat and convenient tear-off and dispensing.
916221	White	21" x 100'	12	Embossed FABRICEL® tissue fused with a layer of polyethylene film on the back. PERFORATED every 30" for neat and convenient tear-off and dispensing.
916118	White	18" x 125'	9	Crepe tissue poly-coated on the back. PERFORATED every 30" for neat and convenient tear-off and dispensing.
916121	White	21" x 125'	9	Crepe tissue poly-coated on the back. PERFORATED every 30" for neat and convenient tear-off and dispensing.
916018	White	18" x 125'	12	Crepe tissue poly-coated on the back. Non-perforated.
916021	White	21" x 125'	12	Crepe tissue poly-coated on the back. Non-perforated.



Leak-proof FABRICEL® and poly-coated rolls reduce cleanup time.



Perforated rolls for neat and convenient tear-off.

Thick, absorbent, Double Re-Crepe (DRC) rolls and sheets. Pebble-embossed for added cushioning. Poly-backed for fluid holdout.

LAB COUNTERTOP BARRIERS

980980	White	20" x 150'	4	DRC rolls
980982	White SHEETS	12" x 20"	600	DRC flat SHEETS
980984	White SHEETS	18" x 20"	400	DRC flat SHEETS

C Choice EXAM TABLE BARRIERS

- Assortment of durable options
- Available in a variety of widths and roll lengths
- Smooth or crepe (quiet, softness) materials
- Adult and pediatric prints; aesthetics and patient peace of mind
- Rolls and sheets
- 99.5% protection against bacteria pass-through (independent tests)

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
913142	White	14" x 225'	12	Smooth rolls
913182	White	18" x 225'	12	Smooth rolls
913212	White	21" x 225'	12	Smooth rolls
913183	White	18" x 260'	12	Smooth rolls
913213	White	21" x 260'	12	Smooth rolls
911243	White	24" x 260'	8	Smooth rolls
950906	White	18" x 550'	6	Smooth rolls
912214	White	21" x 550'	6	Smooth rolls

C Choice

Prints create a relaxed or fun environment for patient comfort and peace of mind for adults or pediatrics—you provide care more easily.

Choice PRINTS



Watercolors



Bugs & Things



Cat in the Hat[†]



Horton Hears a Who[†]



Schooltime



Under the Sea



Pastel Vines



Raining Cats & Dogs



Pediatric Combo

EXAM TABLE BARRIERS *Continued*

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
916143	White	14" x 125'	12	Crepe rolls
916183	White	18" x 125'	12	Crepe rolls
916213	White	21" x 125'	12	Crepe rolls
914243	White	24" x 125'	12	Crepe rolls
918703	White SHEETS	15" x 20"	2000	Crepe SHEETS, pre-cut, flat
918701	White SHEETS	18" x 24"	1000	Crepe SHEETS, pre-cut, flat
918702	White SHEETS	20" x 30"	1000	Crepe SHEETS, pre-cut, flat

PRINTS

982718	Pastel Vines	18" x 225'	6	Smooth rolls
982721	Pastel Vines	21" x 225'	6	Smooth rolls
982618	Pastel Vines	18" x 125'	6	Crepe rolls
982621	Pastel Vines	21" x 125'	6	Crepe rolls
982518	Watercolors**	18" x 225'	12	Smooth rolls
982521	Watercolors**	21" x 225'	12	Smooth rolls
982418	Watercolors**	18" x 125'	12	Crepe rolls
982421	Watercolors**	21" x 125'	12	Crepe rolls
983818	Cat in the Hat [†]	18" x 225'	6	Smooth rolls
983918	Horton Hears a Who [†]	18" x 225'	6	Smooth rolls
981514	Bugs & Things	14" x 225'	6	Smooth rolls (contains printed measure in inches)
981518	Bugs & Things	18" x 225'	6	Smooth rolls (contains printed measure in inches)
981521	Bugs & Things	21" x 225'	6	Smooth rolls (contains printed measure in inches)
981414	Bugs & Things	14" x 125'	6	Crepe rolls (contains printed measure in inches)
981418	Bugs & Things	18" x 125'	6	Crepe rolls (contains printed measure in inches)
981421	Bugs & Things	21" x 125'	6	Crepe rolls (contains printed measure in inches)
981718	Under the Sea	18" x 225'	6	Smooth rolls
981618	Under the Sea	18" x 125'	6	Crepe rolls
982018	Schooltime	18" x 225'	6	Smooth rolls
981918	Schooltime	18" x 125'	6	Crepe rolls
981814	Raining Cats & Dogs	14" x 125'	6	Crepe rolls
981818	Raining Cats & Dogs	18" x 125'	6	Crepe rolls
982114	Pediatric Combo	14" x 125'	12	Crepe rolls
982118	Pediatric Combo	18" x 125'	12	Crepe rolls

**Matches gown, cape, drape, and towel prints

[†] Dr. Seuss Properties TM & © 2004 Dr. Seuss Enterprises, L.P.
All Rights Reserved.

E Everyday EXAM TABLE BARRIERS

- Economical strength, comfort and convenience
- Reliable infection barrier protection
- 99% protection against bacteria pass-through (independent tests)

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
980912	White	18" x 225'	12	Smooth rolls
980914	White	21" x 225'	12	Smooth rolls
9810891	White	18" x 200'	12	Smooth rolls
9810892	White	21" x 200'	12	Smooth rolls
981002	White	18" x 125'	12	Crepe rolls
981004	White	21" x 125'	12	Crepe rolls

CHIROPRACTIC HEADREST BARRIERS

Durable protective coverage in soft crepe or smooth finishes that offer liquid holdout. Pre-cut sheets are available with or without facial slits.



C Choice

CHIROPRACTIC HEADREST BARRIERS Strong, tear-resistant barrier protection.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
ROLLS				
980899	White	8-1/2" x 225'	25	Smooth rolls
980879	White	12" x 225'	12	Smooth rolls
980898	White	8-1/2" x 125'	25	Crepe rolls
SHEETS				
980881	White	12" x 12"	1000	Smooth sheets, pre-cut, flat with facial slit
980883	White	12" x 24"	1000	Smooth sheets, pre-cut, flat with facial slit
980880	White	12" x 12"	1000	Smooth sheets, pre-cut, flat <i>without</i> facial slit
980882	White	12" x 24"	1000	Smooth sheets, pre-cut, flat <i>without</i> facial slit

E Everyday

CHIROPRACTIC HEADREST BARRIERS Barrier protection in a lighter-weight material.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
980900	White	8-1/2" x 225'	25	Smooth rolls

EQUIPMENT DRAPE SHEETS

Our equipment drape sheets, including stretcher sheets and bed sheets, provide strength and fluid holdout for effective infection barrier protection and reduced cleanup time.



U **Ultimate** EQUIPMENT DRAPE SHEETS STRETCHER SHEETS

- High tear-strength nonwoven materials
- Soft; strong
- Fluid holdout for dependable protection
- Available in fitted sheets

Fitted sheets
950556
950562



CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
950552	Light Blue	25" x 72"	50	Fitted, polypropylene nonwoven with 6" sewn gusset pocket for better fit on thicker pads; resists fluids
950551	Light Blue	25" x 72"	50	Fitted, polypropylene nonwoven with 6" sewn pocket; resists fluids
950453	Dark Blue	40" x 90"	50	SMS nonwoven, repels fluid
980939	Blue	40" x 90"	25	Double Re-Crepe (DRC) with poly back, soft, thick, strong, highly absorbent and prevents leak through
950556	Light Blue	30" x 84"	50	Fitted, Polypropylene nonwoven, resists fluid
950562	Light Blue	40" x 89"	25	Fitted, Polypropylene nonwoven, resists fluid

C **Choice** EQUIPMENT DRAPE + STRETCHER + BED SHEETS

- FABRICEL® tissue fused with a layer of polyethylene film
- Super strong, tear-resistant, prevents leak through
- Quiet, with optimal softness

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
919372	White	40" x 72"	50	FABRICEL stretcher sheet
919380	White	40" x 84"	50	FABRICEL stretcher sheet
919396	White	60" x 96"	25	FABRICEL bed sheet

E **Everyday** EQUIPMENT DRAPE + STRETCHER + BED SHEETS

- Strength and fluid holdout
- Reliable coverage and infection control
- Embossed for softness

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
918211	White	30" x 48"	100	Tissue/poly, poly-backed
918213	Blue	30" x 48"	100	Tissue/poly, poly-backed
980924	Blue	40" x 48"	100	Tissue/poly, poly-backed, heavier, stronger tissue

E Everyday

EQUIPMENT DRAPE + STRETCHER + BED SHEETS *Continued*

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
980926	Blue	40" x 60"	100	Tissue/poly, poly-backed, heavier, stronger tissue
918372	Blue	40" x 72"	50	Tissue/poly, poly-backed
980940	White	60" x 96"	25	Tissue/poly, poly-backed, heavier, stronger poly
918248	White	40" x 48"	100	2-ply tissue/poly, poly-backed, extra absorbency and strength
918272	White	40" x 72"	50	2-ply tissue/poly, poly-backed, extra absorbency and strength
9810455	White	40" x 90"	50	Tissue/poly, poly-backed
980925	Blue	40" x 72"	50	Tissue/poly, poly-backed
980928	Blue	40" x 90"	50	Tissue/poly, poly-backed
980927	Blue	40" x 72"	50	Tissue/poly, poly-backed, heavier tissue and film for extra strength
980929	Blue	40" x 90"	50	Tissue/poly, poly-backed, heavier tissue and film for extra strength

PILLOWCASES

Create a comforting environment with reliable barrier protection.



Pillowcases combine barrier protection and patient comfort.

Choice COLORS



White



Blue



Mauve

C Choice PILLOWCASES

- Very soft; patient comfort
- Strong, tear-resistant
- Prevents leak through; hygienic barrier protection

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
919353	Blue	21" x 30"	100	FABRICEL®, tissue/poly, prevents leak through
919355	White	21" x 30"	100	FABRICEL®, tissue/poly, prevents leak through
919356	Mauve	21" x 30"	100	FABRICEL®, tissue/poly, prevents leak through
980915	White	21" x 30"	100	Spunbond polypropylene, nonwoven, resists fluids

E Everyday PILLOWCASES

- Comfortable, reliable hygienic barrier

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
919365	White	21" x 30"	100	Tissue/poly
919363	Blue	21" x 30"	100	Tissue/poly
919350	White	21" x 30"	100	Embossed, poly, strong, prevents leak through

TOWELS + SCALE BARRIERS + TRAY COVERS

Our towels offer convenience and absorbency to perform while providing a soft, gentle touch to patients' skin. Multi-purpose, poly-backed options effectively protect surfaces, such as scales and trays.

Wash, wipe and protect with our multi-purpose towels.



Ultimate TOWELS

- Our softest, strongest, most absorbent 3-ply tissue towels
- Available in colors to add a warm touch to any décor
- Poly-back feature prevents leak-through
- Available in Diamond and Waffle Emboss
- **Diamond-emboss** provides greater bulk and thickness

Ultimate TEXTURE



Waffle



Diamond

Choice COLORS



White Green Blue



Mauve Teal Yellow



Grey Beige Peach



Lavender

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
919101	White	13" x 18"	500	3-ply tissue, diamond-embossed
919102	Green	13" x 18"	500	3-ply tissue, diamond-embossed
919103	Blue	13" x 18"	500	3-ply tissue, diamond-embossed
919106	Mauve	13" x 18"	500	3-ply tissue, diamond-embossed
919110	Teal	13" x 18"	500	3-ply tissue, diamond-embossed
918101	White	13" x 18"	500	3-ply tissue, waffle-embossed
918102	Green	13" x 18"	500	3-ply tissue, waffle-embossed
918103	Blue	13" x 18"	500	3-ply tissue, waffle-embossed
918106	Mauve	13" x 18"	500	3-ply tissue, waffle-embossed
918110	Teal	13" x 18"	500	3-ply tissue, waffle-embossed
918111	White	17" x 18"	500	3-ply tissue, waffle-embossed
919401	White	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919402	Green	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919403	Blue	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919404	Yellow	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919405	Gray	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919406	Mauve	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919407	Beige	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919408	Peach	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919410	Teal	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
919400	Lavender	13" x 18"	500	3-ply tissue, poly-backed, diamond-embossed
917401	White	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917402	Green	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917403	Blue	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917404	Yellow	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917405	Gray	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917406	Mauve	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917407	Beige	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917408	Peach	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917410	Teal	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917400	Lavender	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917411	White	17" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
917413	Blue	17" x 18"	500	3-ply tissue, poly-backed, waffle-embossed

C Choice TOWELS

- Enhanced 2-ply and 3-ply absorbencies and strength for better performance
- A range of pastel colors and colorful prints
- Poly-back feature prevents leak-through
- Available in Diamond and Waffle Emboss
- **Diamond-emboss** provides greater bulk and thickness

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
918161	White	13" x 18"	500	2-ply tissue, waffle-embossed
919461	White	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919462	Green	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919463	Blue	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919464	Yellow	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919465	Gray	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919466	Mauve	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919467	Beige	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919468	Peach	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919470	Teal	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
919459	Lavender	13" x 18"	500	2-ply tissue, poly-backed, diamond-embossed
917461	White	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917462	Green	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917463	Blue	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917464	Yellow	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917465	Gray	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917466	Mauve	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917467	Beige	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917468	Peach	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917470	Teal	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917459	Lavender	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917471	White	17" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
917473	Blue	17" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
982413	Watercolors**	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
982433	Watercolors**	13" x 18"	500	3-ply tissue, poly-backed, waffle-embossed
982434	Watercolors**	13" x 18"	500	3-ply tissue, waffle-embossed
983813	Cat in the Hat [†]	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
983910	Horton Hears a Who [†]	10" x 13"	250	2-ply tissue, poly-backed, waffle-embossed
983913	Horton Hears a Who [†]	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
981410	Bugs & Things	10" x 13"	250	2-ply tissue, poly-backed, waffle-embossed
981413	Bugs & Things	13" x 18"	250	2-ply tissue, poly-backed, waffle-embossed
981610	Under the Sea	10" x 13"	250	2-ply tissue, poly-backed, waffle-embossed
981613	Under the Sea	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
981913	Schooltime	13" x 18"	500	2-ply tissue, poly-backed, waffle-embossed
918189	TIDI Toes	13" x 18"	500	3-ply tissue, PODIATRY towel, waffle-embossed
917489	TIDI Toes	13" x 18"	500	2-ply tissue, poly-backed PODIATRY towel, waffle-embossed

**Matches gown, cape, drape, and table roll prints

[†]Dr. Seuss Properties TM & © 2004 Dr. Seuss Enterprises, L.P. All Rights Reserved.

Prints create a relaxed or fun environment for patient comfort and peace of mind for adults or pediatrics—you provide care more easily.

Choice TEXTURE



Waffle



Diamond

Choice COLORS



White



Green



Blue



Mauve



Teal



Yellow



Grey



Beige



Peach



Lavender

Choice PRINTS



Watercolors



Bugs & Things



Cat in the Hat[†]



Horton Hears a Who[†]



Schooltime



Under the Sea



TIDI Toes

E Everyday TOWELS

- Reliable 2-ply tissue
- Poly-back strength and fluid holdout for effective hygiene
- 100% recycled fiber

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9810865	White	13" x 18"	500	2-ply tissue, poly-backed, diamond embossed

TOWEL DISPENSER

7400	Clear	7" x 5" x 6"	1	Acrylic countertop dispenser for all 13" x 18" towels
------	-------	--------------	---	---



Convenient acrylic dispenser fits all TIDI 13" x 18" wipers and towels.

KAY-PEES® PROFESSIONAL TOWELS + DISPENSER BOX

- For sensitive skin
- Lint-free
- Low friction

KAY-PEES towels in handy individual dispenser boxes for portability and product protection.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
911898	White	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed
911900	Blue	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed



DURAWICK™ COUNTER TOWEL

Specially created towel to provide superior protection for surfaces and counter tops. Multi-layer design; top layer wicks away moisture, middle 3-ply tissue layer absorbs fluids, and the poly back prevents leak through. Durawick offers a durable infection prevention barrier.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
983914	Blue	13" x 18"	100	Multi-layer wicking towel



ABSORBENT UNDERPADS

- Soft nonwoven top ply provides non-irritating skin protection from wetness
- Highly absorbent tissue-filled or fluff-and-tissue-filled center
- Waterproof poly backing

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
16653	Blue	12" x 17"	500 (10 bags of 50)	3-ply tissue-filled
16650	Blue	17" x 24"	300 (6 bags of 50)	3-ply tissue-filled
16652	Blue	23" x 24"	200 (4 bags of 50)	3-ply tissue-filled
16651	Blue	24" x 36"	150 (3 bags of 50)	3-ply tissue-filled
16660	Blue	17" x 22"	300 (6 bags of 50)	Fluff and 2-ply tissue filled
16661	Blue	22" x 24"	200 (4 bags of 50)	Fluff and 2-ply tissue filled
16662	Blue	23" x 36"	150 (6 bags of 25)	Fluff and 2-ply tissue filled

BABY SCALE BARRIERS

C Choice BABY SCALE BARRIERS

- Protect your smallest patients from the transfer of infections
- Two customized sizes provide optimal coverage on most baby scales
- 2-ply tissue offers softness for patient comfort and reliable absorbency
- Poly liner prevents leak-through
- Printed with an inviting, baby-friendly pattern

CODE	COLOR/ PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
981225	Just for Baby	13" x 25"	250	2-ply tissue, poly liner
981222	Just for Baby	13" x 22"	250	2-ply tissue, poly liner

Choice PRINT




Just for Baby

WASHCLOTHS + WIPERS

C Choice WASHCLOTHS + WIPERS

- Double Re-Crepe (DRC) material
- Superior absorbency, strength and softness

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
950750	White	10" x 13"	500 (10 packs of 50)	DRC hygienic washcloths/wipers
950753	White	10" x 13"	1000 (20 packs of 50)	DRC hygienic washcloths/wipers
950755	White	13" x 13"	500 (bulk)	DRC hygienic washcloths/wipers
950754	White	13" x 13"	800 (16 packs of 50)	DRC hygienic washcloths/wipers



Durable, absorbent materials make TIDI washcloths a reliable choice.

ACCESSORIES

With a variety of products to meet your needs, we help you cover all the details, keeping your focus where you want it: on your patients.

CODE	UNIT QUANTITY	DESCRIPTION
------	------------------	-------------

THERMOMETER SHEATHS

20524	5000 (10 boxes of 500)	Mercury oral sheath
20521	5000 (50 boxes of 100)	Mercury oral sheath
20538	11000 (Bulk)	Mercury oral sheath
20710	5000 (50 boxes of 100)	Pre-lubricated mercury rectal sheath
20631	11000 (Bulk)	Digital oral sheath
20633	5000 (50 boxes of 100)	Digital oral sheath
20634	5000 (10 boxes of 500)	Digital oral sheath



Thermometer sheaths provide reliable infection control and simple convenience for care professionals.



COTTON + GAUZE + NONWOVEN WOUND CARE

GAUZE SPONGES + PADS, STERILE + NON-STERILE

Sterile 1's are packaged with a straight seal, while Sterile 2's use chevron packaging, which allows the user to open the package without touching its contents.

- Ideal for general cleaning, dressing, prepping and packing wounds.

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
VENTURE™	908282	2" x 2"	5000 (25 sleeves of 200)	8-ply non-sterile gauze sponges
VENTURE	908294	2" x 2"	8000 (40 sleeves of 200)	12-ply non-sterile gauze sponges
VENTURE	908206	3" x 3"	4000 (20 sleeves of 200)	8-ply non-sterile gauze sponges
VENTURE	908293	3" x 3"	4000 (20 sleeves of 200)	12-ply non-sterile gauze sponges
VENTURE	908292	4" x 4"	4000 (20 sleeves of 200)	8-ply non-sterile gauze sponges
VENTURE	908295	4" x 4"	2000 (10 sleeves of 200)	12-ply non-sterile gauze sponges
VENTURE	908291	4" x 4"	2000 (10 sleeves of 200)	16-ply non-sterile gauze sponges
VENTURE	908200	4" x 8"	2000 (20 sleeves of 100)	12-ply non-sterile gauze sponges
VENTURE	908112	2" x 2" Sterile	3600 (36 boxes of 100)	12-ply Sterile 1's gauze pads
VENTURE	908110	3" x 3" Sterile	2400 (24 boxes of 100)	12-ply Sterile 1's gauze pads
VENTURE	908115	4" x 4" Sterile	1200 (12 boxes of 100)	12-ply Sterile 1's gauze pads
VENTURE	908271	3" x 3" Sterile	2400 (30 boxes of 40)	12-ply Sterile 2's gauze pads
VENTURE	908276	4" x 4" Sterile	1200 (24 boxes of 25)	8-ply Sterile 2's gauze pads
VENTURE	908272	4" x 4" Sterile	1200 (24 boxes of 25)	12-ply Sterile 2's gauze pads

COTTON-FILLED GAUZE SPONGES, STERILE + NON-STERILE

- Cotton filled sponges feature a woven gauze facing with an absorbent cotton fill
- Available in sterile and non-sterile packaging

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
VENTURE	919000	2" x 2"	5000 (25 bags of 200)	8-ply, Cotton-filled non-sterile gauze sponges
VENTURE	908224	4" x 4"	2000 (20 bags of 100)	8-ply, Cotton-filled non-sterile gauze sponges
VENTURE	908223	3" x 3"	4000 (40 bags of 100)	8-ply, Cotton-filled non-sterile gauze sponges
TIDI	919010	2" x 2" Sterile	5000 (25 bags of 200)	8-ply, Cotton-filled sterile gauze sponges. Made in the USA.

NONWOVEN SPONGES, NON-STERILE

- Nonwoven sponges are highly absorbent
- Low-linting
- General cleaning, prepping and applying ointments

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
OMNISORB®	942022	2" x 2"	5000 (25 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
OMNISORB	942033	3" x 3"	4000 (20 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
OMNISORB	942044	4" x 4"	2000 (10 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
VENTURE™	912022	2" x 2"	5000 (25 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
VENTURE	912033	3" x 3"	4000 (20 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
VENTURE	912044	4" x 4"	2000 (10 bags of 200)	4-ply, closed weave, non-sterile, nonwoven sponges
VENTURE	260002	2" x 2"	8000 (40 bags of 200)	4-ply, open weave, non-sterile, nonwoven sponges
VENTURE	260003	3" x 3"	4000 (20 bags of 200)	4-ply, open weave, non-sterile, nonwoven sponges
VENTURE	260004	4" x 4"	2000 (10 bags of 200)	4-ply, open weave, non-sterile, nonwoven sponges

DRESSING SPONGES, STERILE + NON-STERILE

- Multi-layered post-op sponges
- Excellent absorption and dispersal

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
VENTURE	908242	3" x 3"	4000 (40 boxes of 100)	4-ply, Post-op non-sterile dressing sponges
VENTURE	908244	4" x 4"	2000 (20 boxes of 100)	4-ply, Post-op non-sterile dressing sponges

ABDOMINAL PADS + ROLLS

- Made with a combination of highly absorbent materials and a moisture resistant barrier
- Open edges offer greater absorbency
- Single-sheet design on pads reduces linting and stray particles
- Non-absorbent rolls for secondary wound dressing.

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
TIDI	979553	8" x 10"	432	Non-absorbent, non-sterile abdominal pads. Made in the USA.
TIDI	979555	10" x 12"	288	Non-absorbent, non-sterile abdominal pads. Made in the USA.
TIDI	979556	12" x 16"	144	Non-absorbent, non-sterile abdominal pads. Made in the USA.
TIDI	979562	12" x 24"	144	Non-absorbent, non-sterile abdominal pads. Made in the USA.
TIDI	979560	8" x 20 yards	12	Non-absorbent, non-sterile abdominal rolls. Made in the USA.
TIDI	979540	8" x 20 yards	6	Double thick, non-absorbent, non-sterile abdominal rolls. Made in the USA.
TIDI	979580	12" x 30 yards	1	Non-absorbent, non-sterile abdominal rolls. Made in the USA.

COTTON + RAYON BALLS

- Manufactured for maximum absorbency and uniform consistency, shape and size.

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
TIDI	969153	Medium	4000 (2 bags of 2000)	Cotton balls. Made in the USA.
TIDI	969152	Large	2000 (2 bags of 1000)	Cotton balls. Made in the USA.
TIDI	969164	X-Large	2000 (4 bags of 500)	Cotton balls. Made in the USA.
TIDI	969163	Medium	4000 (2 bags of 2000)	Rayon balls. Made in the USA.
TIDI	969162	Large	2000 (2 bags of 1000)	Rayon balls. Made in the USA.



INFECTION CONTROL

STERILIZATION WRAP

- Durable, water-repellant material allows penetration for various sterilization methods.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
960172	Blue	12" x 12"	1000	Sterilization wrap
960173	Blue	15" x 15"	1000	Sterilization wrap
960174	Blue	18" x 18"	500	Sterilization wrap
960175	Blue	20" x 20"	500	Sterilization wrap
960176	Blue	24" x 24"	500	Sterilization wrap
960177	Blue	30" x 30"	400	Sterilization wrap
960178	Blue	36" x 36"	200	Sterilization wrap

STERILIZATION POUCHES

- Self-sealing and have external color indicators for peace of mind.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
5549	Clear one side	2-3/4" x 10"	300	Sterilization pouches, dual indicators, latex free
5550	Clear one side	3-1/4" x 6-1/2"	200	Sterilization pouches, dual indicators, latex free
5551	Clear one side	3-1/2" x 10"	200	Sterilization pouches, dual indicators, latex free
5552	Clear one side	5-1/4" x 11"	200	Sterilization pouches, dual indicators, latex free
5553	Clear one side	7-1/2" x 12"	200	Sterilization pouches, dual indicators, latex free
5543	Blue transparent	2-1/2" x 8-1/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5544	Blue transparent	3-1/2" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5545	Blue transparent	6" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
30203	White	2-1/2" x 1-1/2" x 10-1/2"	1000	Autoclave envelopes, contain latex in adhesive strip

AUTOCLAVE TOWELS

- A poly layer, sandwiched between two layers of tissue
- Allows for use in kits and in autoclaves while maintaining the integrity of the towels

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
8253	White	13" x 18"	500	T/P/T autoclavable towels
8252	White	19" x 22"	300	T/P/T autoclavable towels
8251	White	19" x 30"	300	T/P/T autoclavable towels

INFECTIOUS WASTE BAGS

- Generously sized, durable waste bags make disposal of infectious waste easier.



Assure proper disposal with durable, generously sized infectious waste bags.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
8653	Red	24" x 24" (10 gal)	250	High density, 12 microns, infectious waste (Biohazard) bag
8654	Red	24" x 33" (16 gal)	250	High density, 12 microns, infectious waste (Biohazard) bag
8656	Red	33" x 40" (33 gal)	250	High density, 17 microns, infectious waste (Biohazard) bag
8657	Red	40" x 48" (40 gal)	250	High density, 17 microns, infectious waste (Biohazard) bag

PATIENT CARE

SPECIALTY PATIENT BIBS

Convenience for long-term care environments.

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
919398	White	20" x 30"	250	FABRICEL® tissue/poly slipover bib
980966	Blue	20" x 29"	500	XL tissue/poly bib, die-cut neck
917673	Blue	18" x 25"	250	Chainless bib with ties
917676	Mauve	18" x 25"	250	Chainless bib with ties
980964	White	16" x 32"	300	Regular tissue/poly bib with ties
920461	White	18" x 30"	150	Tissue/poly bib with ties
920462	White	18" x 30"	150	Economy tissue/poly slipover bib
920463	White	18" x 30"	150	Standard tissue/poly slipover bib
920861	White	18" x 32"	300	Thrifty tissue/poly bib with ties
920862	White	18" x 32"	300	Thrifty tissue/poly slipover bib
920962	White	18" x 30"	500	Economy tissue/poly slipover bib
950443	White	16" x 24"	200	Adult size poly bib with ties
980962	White	16" x 24"	500	Embossed poly bib, front full-width pocket, with ties
305	White	16" x 24"	500	FOODCARE™ adult bib with ties

PATIENT CUPS

Choose resilient poly or basic waxed paper.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9211	White	5 oz.	1000	Multi-use ribbed plastic drinking cups
9212	Green	5 oz.	1000	Multi-use ribbed plastic drinking cups
9213	Blue	5 oz.	1000	Multi-use ribbed plastic drinking cups
9214	Yellow	5 oz.	1000	Multi-use ribbed plastic drinking cups
9215	Gray	5 oz.	1000	Multi-use ribbed plastic drinking cups
9216	Mauve	5 oz.	1000	Multi-use ribbed plastic drinking cups
9217	Beige	5 oz.	1000	Multi-use ribbed plastic drinking cups
9218	Peach	5 oz.	1000	Multi-use ribbed plastic drinking cups
9210	Lavender	5 oz.	1000	Multi-use ribbed plastic drinking cups
9241	White	3.5 oz.	1000	Multi-use ribbed plastic drinking cups
9243	Blue	3.5 oz.	1000	Multi-use ribbed plastic drinking cups
9226	Jazz	5 oz.	1000	Waxed paper cups
9240	Jazz	3 oz.	5000	Waxed paper cups

CUP COLORS



CUP PRINTS



Jazz

CUP DISPENSERS

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9246	Chrome	3 or 5 oz.	1	Adjustable cup dispenser

PREMIUM POLY APRONS

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
10412	White	24" x 42"	1000 (10 boxes of 100)	Individual poly-bagged ultra-lightweight embossed poly aprons
10401	White	28" x 46"	1000 (10 boxes of 100)	Individual poly-bagged ultra-lightweight embossed poly aprons
10403	White	28" x 46"	500 (5 boxes of 100)	Individual poly-bagged regular weight embossed poly aprons
10416	White	28" x 46"	500 (5 boxes of 100)	Lightweight taffeta-embossed poly aprons
10430	Clear	28" x 46"	500 (5 boxes of 100)	Lightweight taffeta-embossed poly aprons
10418	White	28" x 46"	500 (5 boxes of 100)	Regular weight taffeta-embossed poly aprons
10419	White	28" x 46"	500 (5 boxes of 100)	Heavy weight taffeta-embossed poly aprons



TIDI 10419 heavyweight embossed poly apron. All TIDI aprons are designed as shown.

HOSPITAL SUPPLIES

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
------	-------------	------	---------------	-------------

BEDSIDE/CHAIRSIDE/SUTURE BAGS

950232	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, bedside/chairside
950231	White, Blue Floral Print	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, bedside/chairside
950251	White, Blue Soft Shell Print	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, bedside/chairside
950252	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, bedside/chairside
950234	White, Blue Floral Print	8" x 3-3/4" x 12"	2000 (20 bags of 100)	Non-flame-retardant poly bag, bedside/chairside

BATH MAT

950255	Soft Shell	14-1/2" x 20"	500	Blue print on white paper
--------	------------	---------------	-----	---------------------------

BEDPAN COVERS

950260	White	11-1/2" x 3-1/2" x 19"	1500 (15 packs of 100)	Bed pan cover, closed-back, flushable
950261	White	11-1/2" x 3-1/2" x 19"	1250 (5 packs of 250)	Bed pan cover, open-back, flushable

POST MORTEM BAGS + SHROUD KITS

950259	White	36" x 90"	10	Vinyl post-mortem bag with curved zipper
950359	White	36" x 90"	10	Vinyl post-mortem bag with straight zipper
950358	White	36" x 90"	5	Contamination kit. Heavy vinyl cadaver bag with straight zipper, chin strap, ID tags, cavity pad and personal effects bag.
950361	White	54" x 108"	10	Shroud kit, adult size (ID tags, chin strap, cavity pad and ties)
950362	White	54" x 72"	10	Shroud kit, child size (ID tags, chin strap, cavity pad and ties)
950365	White	72" x 108"	10	Shroud kit, XL, adult size (ID tags, chin strap, cavity pad and ties)
950363	White	54" x 108"	25	Shroud sheet, adult size, embossed poly
950371	White	72" x 108"	10	Shroud sheet, XL, embossed poly

VETERINARY SUPPLY

919595	Blue	100 yds x 38"	1 roll	Disposable wrap, Dexter material. Soft, pliable, leak proof and durable
--------	------	---------------	--------	---

OFFICE SUPPLY

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
------	-------	------	---------------	-------------

STORAGE ENVELOPES

950222	Kraft	8-1/2" x 10-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950221	Kraft	10-1/2" x 12-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950218	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, heavyweight
950220	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, lightweight
950219	Kraft	15" x 18"	50/case, bulk packed	X-ray mailer envelopes with insert

TRASH CAN LINERS

950211	Brown	8-3/4" x 7-3/4" x 17"	600 (12 packs of 50)	Waxed can liners, 16-20 quarts
8613	Mint Green	10 gal	500	Trash can liners, heavy-duty poly
8614	Mint Green	16 gal	500	Trash can liners, heavy-duty poly
8625	Mint Green	30 gal	250	Trash can liners, heavy-duty poly
8627	Mint Green	40 gal	200	Trash can liners, heavy-duty poly

GLOSSARY

Materials + Processes

SMS: Spunbond/Meltblown/Spunbond

This multi-layer nonwoven material is made up of fine synthetic fibers for superior softness, strength and liquid repellant.

T/P/T: Tissue/Poly/Tissue

A layer of polyethylene film, sandwiched between two sheets of tissue, adds strength and resistance against fluid strike-through to the soft, absorbent tissue.

Tissue (2-ply, 3-ply and 4-ply)

TIDI offers tissue products in configurations to suit a range of needs. Multiple tissue plies add absorbency and strength, while two-ply tissue provides economical softness and absorbency that is well suited for shorter procedures.

Smooth/Crepe Exam Table Barriers

Choose smooth paper with added liquid holdout, or the softer, more comfortable surface feel of crepe.

Poly

A layer of premium octane linear polyethylene provides high tensile strength and reliable barrier properties.

DRC

Double Re-Crepe material provides a soft, cloth-like feel for durable absorbency.

FABRICEL®

A proprietary material that fuses a sheet of soft, supple tissue and a layer of polyethylene film. This heat lamination process creates a very soft yet super-strong, tear-resistant material that prevents soak-through.

Scrim-reinforced Tissue

We fortify our high-quality tissue by fusing nylon woven netting directly into the fiber matrix, for a soft, yet highly durable and tear-resistant fabric that is ideal for extended procedures.

Nonwoven

Synthetic materials are bonded into a web of fibers. Nonwoven garments are ideal for longer and more involved procedures as they provide patient comfort, strength and durability. TIDI also offers nonwovens in our wound care selections.

Spunlace

Also called hydroentanglement, this process uses water to create an interlocked web of fibers. The resulting fabric is absorbent with a textile-like drape and softness for patient comfort.

Sustainable Solutions

TIDI Products is focused on reducing resource consumption, minimizing environmental impact and acting in an environmentally responsible manner, in order to provide effective solutions without compromising the well-being of future generations.

Online Information

Explore our broad range of product areas and learn more about how our solutions can enhance the care experience. Visit **TIDIproducts.com**.

Product Samples

Experience the softness of a gown, the absorbency of a towel, the coverage of a drape, or the performance of any TIDI product with a no-charge sample pack.

Call **1.800.521.1314**, visit **TIDIproducts.com**, or e-mail us at **Excellence@TIDIproducts.com**.

Value-Added Packaging

Newly designed packaging provides clear, concise, instant communication. Quickly and simply verify product identity and description for ease in ordering and selecting the right TIDI products for your needs.



Certified in ISO 9001 and EN/ISO 13485






TIDI Products, LLC

570 Enterprise Drive
Neenah, WI 54956 USA

Customer Service

800.521.1314
+1 920.751.4138
800.837.7770 FAX

www.tidiproducs.com



MEDCAT1113