

PORTER

PRODUCTS AND PARTS CATALOG JULY 2006



POLICIES

- PRICES** ALL PRICES ARE F.O.B. HATFIELD, PA, AND SUBJECT TO CHANGE WITHOUT NOTICE. ALL MERCHANDISE WILL BE INVOICED AT PRICES IN EFFECT AT TIME OF SHIPMENT. PRICES DO NOT INCLUDE INSURANCE, FREIGHT, TAXES OR SPECIAL HANDLING. THESE CHARGES, IF APPLICABLE, WILL BE SHOWN SEPARATELY ON INVOICE. **MINIMUM ORDER \$60.00**
- PAYMENT TERMS** NET 30 DAYS AFTER INVOICE DATE. ALL INVOICES PAST DUE ARE SUBJECT TO A FINANCE CHARGE OF 1 1/2 % PER MONTH. (18% ANNUAL RATE).
- SHIPMENTS** SHIPMENT OF MERCHANDISE SHALL AT ALL TIMES BE SUBJECT TO CREDIT APPROVAL AND WILL BE CONTINGENT UPON FIRES, ACCIDENTS, EMERGENCIES, ACTS OF GOD OR ANY OTHER CAUSES WHICH ARE BEYOND PORTER INSTRUMENT DIVISION'S CONTROL.
- CANCELLATIONS** NO CANCELLATIONS WILL BE ACCEPTED ON NON-STANDARD OR SPECIAL MERCHANDISE, EXCEPT BY PAYMENT OF FULL PURCHASE PRICE. IF BUYER REQUESTS CANCELLATION OF ANY ORDER OR PART THEREOF, AND IS AGREED TO BY PORTER INSTRUMENT DIVISION IN WRITING, BUYER WILL BE SUBJECT TO CANCELLATION CHARGES TO COVER THE COST OF MATERIAL AND/OR FABRICATION INCURRED BY PORTER INSTRUMENT DIVISION TO DATE OF CANCELLATION.
- CHANGES OF ORDER** A MINIMUM OF 90 DAYS NOTICE IS REQUIRED ON ALL CHANGES TO ORDERS AND WILL BE SUBJECT TO RE-SCHEDULING AS A NEW ORDER AT PORTER INSTRUMENT DIVISION'S DISCRETION
- RETURNS** NO RETURNS WILL BE ACCEPTED UNLESS AUTHORIZED IN WRITING BY PORTER INSTRUMENT DIVISION, AND ACCOMPANIED BY A PROPERLY COMPLETED RETURN GOODS AUTHORIZATION. ALL RETURNS ARE SUBJECT TO A RE-STOCKING FEE AND POSSIBLE REWORK CHARGES TO BE DETERMINED BY PORTER INSTRUMENT DIVISION.

CERTIFICATE OF WARRANTY

THIS WARRANTY IS GIVEN IN PLACE OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OTHERWISE.

Under no circumstances shall Parker Hannifin Corporation be liable for incidental or consequential damages as those terms are defined in the uniform commercial code.

Parker Hannifin Corporation, Porter Instrument Division warrants that each product or part shall be free from defects in workmanship and materials, under normal use and with appropriate maintenance, for one (1) year from the date of delivery to customer unless otherwise specified in writing. All rubber and plastic parts and accessories are warranted under the same conditions for a period of ninety (90) days from date of purchase.

No statement or claim about the product by any employee, agent, representative, or dealer of Parker Hannifin Corporation shall constitute a warranty by Parker Hannifin Corporation or give to rise to any liability or obligation of Parker Hannifin Corporation.

Parker Hannifin Corporation shall not be liable for any damage, injury or loss arising out of the use of the product, whether as a result of a defect in the product or otherwise, if, prior to such damage, injury or loss, the product was (1) damaged or misused; (2) repaired, altered or modified by persons other than Parker Hannifin Corporation; (3) not installed in strict compliance with applicable codes and ordinances; or (4) not installed by an authorized Parker Hannifin Corporation dealer. Parker Hannifin Corporation's obligation for breach of this warranty, or for negligence or otherwise, shall be strictly and exclusively limited to the repair or replacement of the product or part. This warranty shall be void on any product on which the serial number has been altered, defected or removed.

ORDERS All orders are to be made through authorized Parker Hannifin Corporation distributors. All billing will be done through said distributors. Direct orders will be handled through the authorized local dealer as determined by Parker Hannifin Corporation.

PRODUCTS

Automatic Vacuum Systems	8
Cabinet Mount Flowmeters	4
Emergency Oxygen/Resuscitation Systems	16
Flowmeters (only)	20
Flowmeters Systems	2,3,4,7
Flowmeter Mounts and Stands	14,15
Gas Supply Tubing	13
Manifold/Tankroom Systems	
Sentinel	10
Vanguard.....	11
Manifold System Installation Kit (Nitrogen Test Fixture).....	10
Nitrogen Equipment	17
Office Communication Systems.....	18,19
Oral Surgery Flowmeter Systems	6
Outlet Stations	12
Porter Disposable Liners.....	7
Porter Navicom Communication System	18
Porter Scavenger Retrofit Kits	
for other manufacturer's Flowmeters	9
Porter Scavenger Systems	8
Porter Scavenger Systems Parts Breakdown	22
Porter Vacuum Control Kits	8
Zone Valves	10,11


PARTS

Cabinet Mount Accessories	5,21
Flowmeter Parts	20
Gas Supply Fittings Only.....	21
Oxygen Resuscitation Parts	29
Oxygen Therapy	29
Meter Mounts and Stand Parts	26,27
Manifold System Parts	
Sentinel.....	25
Vanguard.....	24
Manifold/Tankroom Parts	24
Oral Surgery Accessories	29
Outlet Station Replacement Parts	28
Porter Scavenger Retrofit Adapters	23
Porter Flowmeter & Bagtee Mounts	27
Sterilizer System Parts (See Website)	30
Wall Cut-Outs.....	31



3400AV	MXR-Flowmeter	\$3,245.00
	with Porter Scavenger Breathing Circuit	
	with New Automatic Vacuum Switch	
3000	MXR-Flowmeter	
4100-3NL	Three Liter Bag (No-Latex)	
5155-3	Porter Scavenger Breathing Circuit	
AVS-5000	Automatic Vacuum Switch	

3400 substitutes 5155-1 Breathing Circuit \$2,908.00
See page 8 (no AVS)




3420AV	MXR Package (3400 AV) plus	\$3,570.00
2020	Solid Wall Mount	
	Telescoping	

3420 substitutes 5155-1 Breathing Circuit \$3,232.00
See page 8 (no AVS)


3435-1AV	MXR Package (3400 AV) plus	\$3,602.00
	Cabinet Mount	
2035-1	Slide Away	

3435-1 substitutes 5155-1 Breathing Circuit \$3,265.00
See page 8 (no AVS)



3436AV	MXR Package (3400 AV) plus	\$3,637.00
2036-3	Cabinet Mount	
	-Slide	


3436 substitutes 5155-1 Breathing Circuit \$3,303.00
See Page 8 (no AVS)



3440AV	MXR Package (3400 AV) plus	\$3,655.00
2040	Compact Mobile Stand	

3440 substitutes 5155-1 Breathing Circuit \$3,318.00
See page 8 (no AVS)

NOTE: When Scavenger is connected to a 6256-3 outlet with vacuum, order a 5602 found on page 13.
Alternative vacuum packages can be found on page 8. Necessary Gas Supply Tubing can be found on page 13.




3445-2AV **MXR Package (3400AV) plus** \$5,468.00
 2045-2 Tall Mobile Stand
 for 2+2 E Tanks
 (Top of meter telescopes from 43" to 50")

3445-2 substitutes 5155-1 Breathing Circuit \$5,133.00
 See page 8 (no AVS)
 (Top of meter telescopes from 43" to 50")


3445E-2AV **Same As Above plus** \$6,054.00
 1506-1 Positive Pressure/Demand Valve

3445E-2 substitutes 5155-1 Breathing Circuit \$5,717.00
 See page 8 (no AVS)



3445-SAV **MXR Package with under cabinet** \$5,509.00
 mobile stand for 2+2 "E" Tanks
 31" to top of stand in storage position
 (Top of meter telescopes from 31" to 40")

3445-S substitutes 5155-1 Breathing Circuit \$5,174.00
 See page 8 (no AVS)
 31" to top of stand in storage positions
 (Top of meter telescopes from 31" to 40")



3445E-SAV **Package with under cabinet** \$6,095.00
 Mobile stand for 2+2 "E" Tanks
 31" to top of stand in storage position
 1506-1 positive pressure/demand valve

3445E-S substitutes 5155-1 Breathing Circuit \$5,759.00
 See page 8 (no AVS)
 31" to top of stand in storage position

3465WF-AV MXR-1 Flushmount Flowmeter \$4,528.00

MXR-1 Flushmount Flowmeter Package
(Universal connections) One model with
connections that are easily switched in the
field from back to bottom

*** Bezel for wall installs not included. Order separately
Bezel part number (B-2625) see page 21**



3465WF (NO AVS) \$4,192.00

Same package as above with 5155-1
substituted for the AVS package



Package with AVS

Package without AVS

Package includes:

Package includes:



A-1679 included
with meter package



2037-1 Under cabinet
mount




Bag tee Mounting
Options for
Flushmount
Flowmeters

2036-4





3 1/2'' wide

2036-2
5'' wide




For optional mount pricing see page 14






- 2400-OS Oral Surgery MXR-Flowmeter** \$2,990.00
with Porter Scavenger Breathing Circuit
 includes:
 2000-OS MXR Flowmeter (includes P2407A oral surgery threaded bag tee w/adapter)
 4100 Three Liter Bag
 5155-1 Porter Adult Inhaler
 Scavenger Assembly, complete
 Special Features:
 Higher Maximum flow rates
 Oral Surgery bag tee for 2210K Y valve




- 2444-OS Oral Surgery MXR-Flowmeter with Mobile Stand System** \$4,309.00
 includes:
 2000-OS MXR flowmeter (includes P2407A oral surgery threaded bag tee w/adapter)
 4100 Three Liter Bag
 5155-1 Porter Adult Inhaler
 Scavenger Assembly, complete
 2210K Directional Y Valve
 1597 Non-rebreathing valve
 4200 Corrugate tube
 2044 Oral surgery Star-Base stand top of meter telescopes form 46" to 55"



- 2200-OS Oral Surgery Accessory Package** \$847.00
 includes:
 2210K Directional Y Valve
 1597 Non-rebreathing valve
 1575-1 Universal size full face, clear mask
 4200 Corrugated tube


Note: Must have threaded bag tee to attach 2210K Part # P2407A




- B-5129-000 Flowmeter Rubbergoods Hook** \$70.00

STANDARD FLOWMETER SYSTEMS

PORTER



1400AV	Standard Flowmeter	\$3,124.00
	with Porter Scavenger Breathing Circuit	
	(with Vacuum Control Block Kit)	
1100	Standard Flowmeter	
4100	Three Liter Bag	
5155-3AV	Porter Adult Inhaler & Automatic Vacuum Switch	






1420AV 2020	Same as above with Wall Mount	\$3,449.00
1435-1AV 2035-1	Same as above with Under Counter Slide Mount	\$3,483.00
1436AV 2036-3	Same as above with Slide Mount	\$3,519.00
1440AV 2040	Same as above with Mobile Stand	\$3,536.00
1445-2AV 2045-2	Same as above with Tall E Stand	\$5,349.00









DISPOSABLE LINERS FOR PORTER SCAVENGER MASKS



THE MARKET LEADER IN SCAVENGING TECHNOLOGY		
5053 AD12	Porter Adult Disposable Liners	\$34.00
5053 AD12C	with Citrus Scent	\$34.00
One bag of (12) Adult Liners. Not Returnable		
5053 AD144	Porter Adult Disposable Liners	\$395.00
5053 AD144C	with Citrus Scent	\$395.00
One box of (12) bags. Each bag contains (12) Adult Liners Contents returnable when original box is unopened only.		
5053 PD12	Porter Pado Disposable Liners	\$34.00
5053 PD12C	with Citrus Scent	\$34.00
One bag of Pado Liners. Not Returnable		
5053 PD144	Porter Pado Disposable Liners	\$395.00
5053 PD144C	with Citrus Scent	\$395.00
One box of (12) bags. Each bag contains (12) Pado Liners Contents returnable when original box is unopened only.		

	<div> AVS-5000 Automatic Vacuum Switch Sytem \$414.00 </div> <div> Retrofits to ALL PORTER Manufactured Flowmeters (except 2000-0S) </div>
	<div> 5155-3AV Porter Adult Scavenger \$899.00 </div> <div> with Automatic Vacuum Switch </div> <div> 5155-4AV Porter Pedit Scavenger \$899.00 </div> <div> with Automatic Vacuum Switch </div>
	<div> 5155-1* Porter Adult Scavenger \$613.00 </div> <div> Breathing Circuit with In-Line Vacuum Control Kit </div> <div> 5155-2* Porter Pedit Scavenger \$613.00 </div> <div> Breathing Circuit with In-Line Vacuum Control Kit </div> <div> 5155-3* Porter Adult Scavenger \$516.00 </div> <div> Breathing Circuit (no Vacuum Control Block) </div> <div> 5155-4* Porter Pedit Scavenger \$516.00 </div> <div> Breathing Circuit (no Vacuum Control Block) </div> <div> * Note: Packages do NOT include 5602 Vac Quick Conn or 4100 Bag. Please order separately. </div>
	<div> 5054A Porter Adult Hood Kit - w/3 liners \$242.00 </div> <div> 5054B Porter Pedit Hood Kit - w/3 liners \$242.00 </div> <div> 5054-1 Porter Adult Hood Kit - pkg of 3 \$49.00 </div> <div> 5054-2 Porter Pedit Hood Kit - pkg of 3 \$49.00 </div>

OPTIONS FOR VACUUM HOOK-UP

 <p>Part No. AVS-5000 Automatic Vacuum Switch</p>	<div>  <p>High Volume Evacuation (HVE) Attachment</p> </div> <div>  <p>Tee Connection in the Vacuum Line</p> </div> <div>  <p>PVC Vacuum Line Hook-Up</p> </div>
 <p>Part No. 5501 RK Vacuum Control Block</p>	<div>  <p>High Volume Evacuation (HVE) Attachment</p> </div> <div>  <p>Tee Connection in the Vacuum Line</p> </div> <div>  <p>PVC Vacuum Line Hook-Up</p> </div>



All kits include the Porter Scavenger with In-Line Vacuum (specify "A" for Adult Hood or "P" for Peds Hood), Three Liter Breathing Bag, Flowmeter Adapters and Bag Tee where necessary.

5800-1A (Adult)

\$651.00

5800-1P (Peds)

\$651.00

This one kit retrofits all these flowmeters:

Matrix - Quantiflex, MDM

Ohmeda

Fraser Harlake

Fraser Sweatman

MDT-McKesson - *Analor Four/Compalor*

Parkell-Sedatron (*not pictured*)

Belmed

Accutron (*not pictured*)

*Also Replaces: Porter Clean-Air™ Rubber Goods
Non-Scavenging System*



Matrix (Quantiflex, MDM)

Ohmeda

Fraser Harlake

Fraser Sweatman



MDT
McKesson
Analor Four
Compalor



Belmed
PC7



Coastal

NCG Chemetron



Hampton
Simplaire Micro



McKesson
Analor
I & II



Starflite
Denetron



McKesson
Analor 3C



5450 Vac Block
&
5300 Canister

5800-2A Coastal

NCG Chemetron

\$738.00

5800-4A Hampton

Simplaire Micro

\$750.00

5800-5A McKesson

Analor I&II

Starflite Denetron

\$700.00

5800-6A McKesson

Analor 3C

\$750.00

5800-9A Porter Clean-Air™

Breathing Circuit & Canister
Scavenging System

\$597.00

Note: *It is not recommended to retrofit Porter Scavenger Breathing Circuit to Adec, Veriflo or Spica machines.*

**For 2 Tanks Oxygen and 2 Tanks Nitrous Oxide**

3222CX **Manifold System, 2-O₂, N₂O,** \$3,687.00
with wall model sentinel, regulators
hoses, and tank restraints

3222DX **Manifold System, 2-O₂, N₂O,** \$3,687.00
with desk model sentinel, regulators
hoses, and tank restraints

For 2 Tanks Oxygen and 1 Tanks Nitrous Oxide

3221CX **Manifold System, 2-O₂, 1-N₂O,** \$3,396.00
with wall model sentinel, regulators, hoses, and
tank restraints

3221DX **Manifold System, 2-O₂, 1-N₂O,** \$3,396.00
with desk model sentinel, regulators, hoses, and
tank restraints

3222CX-N/R&H **Manifold System** \$2,530.00
with wall sentinel;
WITHOUT regulators and hoses

3222DX-N/R&H **Manifold System** \$2,530.00
with desk sentinel;
WITHOUT regulators and hoses

3222MX **Sentinel Manifold ONLY** \$1,879.00
5251 **Wall Sentinel** \$725.00
(includes 75' cable & power supply)
5251 **Desk Sentinel** \$725.00
(includes 75' cable & power supply)

Manifolds without Sentinel DO NOT meet the requirements of the NFPA 99. Porter Instrument Division accepts no responsibility for the installation and/or use of manifolds without Sentinel alarm.

ZONE VALVES WITH CROSS-PROTECTION

9000-1 **Concealed Zone Valve** \$1,148.00
For use with 3222MX Sentinel Manifold

Note: NFPA 99-GAS SHUT-OFF VALVES Where the central supply is remote from the medical gas system use points, the main supply line shall be provided with a shut-off valve (Zone Valve) is located in the single treatment facility as to be accessible from use-point locations in an emergency.

MANIFOLD INSTALLATION KIT

B-2065-000 **Nitrogen Test Fixture** \$700.00
(includes adapter hoses)



with Automatic Changeover

4222CXA **Vanguard Manifold Systems** \$4,710.00
with automatic changeover for
2-O₂, 2-N₂O, wall model alarm, regulators,
hoses and tank restraints

4222DXA **Vanguard Manifold Systems** \$4,710.00
with 6201A Desk Model

4222CXALPI **Vanguard Manifold Systems** \$4,137.00
for 2+2
wall alarm WITHOUT Pre-Install Kit

4222DXALPI **Vanguard Manifold Systems** \$4,137.00
for 2+2
Desk alarm WITHOUT Pre-Install Kit

4222MXVA **Vanguard Manifold Only** \$3,016.00

6251A **Vanguard Wall** \$1,136.00
Alarm Only

6201A **Vanguard Desk** \$1,136.00
Alarm Only

ZONE VALVES WITH CROSS + PROTECTION



9000-2 **Concealed Zone Valve** \$1,148.00
For use with 4222 Vanguard Manifold


Note: NFPA 99-GAS SHUT-OFF VALVES Where the central supply is remote from the medical gas system use points, the main supply line shall be provided with a shut-off valve (Zone Valve) is located in the single treatment facility as to be accessible from use-point locations in an emergency.



6100-1	Single Oxygen Outlet Station Requires Hose, O ₂ , diss/Porter listed on page 13	\$268.00
6200-1	Dual Outlet Station, O₂/N₂O	\$343.00
6250K	Exposed Cover for 6100-1/6200-1	\$41.00
6300-1	Floor Mounted Outlet Station	\$578.00




6256-3	Triple Outlet O ₂ /N ₂ O/Vacuum	\$538.00
6257-3	Exposed Triple Outlet O ₂ /N ₂ O/Vacuum	\$566.00
6255-1	Compact Wall Triple Outlet O ₂ /N ₂ O/Vacuum	\$538.00










6258-1	Compact Triple Outlet for Cabinets (no hoses or shutoffs)	\$538.00
6258-3	Compact Triple Outlet for Cabinets with 3 foot hoses and shutoffs	\$932.00
6258-5	Compact Triple Outlet for Cabinets with 5 foot hoses and shutoffs	\$932.00



6258-1B	Compact triple outlet w/back connection for Cabinets (no hoses or shutoffs)	\$538.00
6258-3B	Compact triple outlet w/back connection with 3 foot hoses and shutoffs	\$932.00
6258-5B	Compact triple outlet w/back connection with 5 foot hoses and shutoffs	\$932.00
















5600-3	Vacuum Outlet Station	\$308.00
---------------	------------------------------	----------



	<table><tr><td>8065</td><td>Hose, O₂/N₂O, diss/Porter, 5ft.</td><td>\$215.00</td></tr><tr><td>8055*</td><td>Hose, O₂, diss/Porter, 5ft.</td><td>\$173.00</td></tr><tr><td>8005*</td><td>Hose, O₂, diss/diss, 5ft.</td><td>\$115.00</td></tr><tr><td>8505*</td><td>Hose, N₂O, diss/diss, 5ft.</td><td>\$115.00</td></tr></table>	8065	Hose, O ₂ /N ₂ O, diss/Porter, 5ft.	\$215.00	8055*	Hose, O ₂ , diss/Porter, 5ft.	\$173.00	8005*	Hose, O ₂ , diss/diss, 5ft.	\$115.00	8505*	Hose, N ₂ O, diss/diss, 5ft.	\$115.00
8065	Hose, O ₂ /N ₂ O, diss/Porter, 5ft.	\$215.00											
8055*	Hose, O ₂ , diss/Porter, 5ft.	\$173.00											
8005*	Hose, O ₂ , diss/diss, 5ft.	\$115.00											
8505*	Hose, N ₂ O, diss/diss, 5ft.	\$115.00											
	<table><tr><td>8015*</td><td>Hose, O₂, diss/NCG (Chemetron), 5ft.</td><td>\$186.00</td></tr><tr><td>8515*</td><td>Hose, N₂O, diss/NCG (Chemetron), 5ft.</td><td>\$186.00</td></tr></table>	8015*	Hose, O ₂ , diss/NCG (Chemetron), 5ft.	\$186.00	8515*	Hose, N ₂ O, diss/NCG (Chemetron), 5ft.	\$186.00						
8015*	Hose, O ₂ , diss/NCG (Chemetron), 5ft.	\$186.00											
8515*	Hose, N ₂ O, diss/NCG (Chemetron), 5ft.	\$186.00											
	<table><tr><td>8025*</td><td>Hose, O₂, diss/Puritan, 5ft.</td><td>\$186.00</td></tr><tr><td>8525*</td><td>Hose, N₂O, diss/Puritan, 5ft. for use with McKesson & Fraser Harlake outlet stations</td><td>\$186.00</td></tr></table>	8025*	Hose, O ₂ , diss/Puritan, 5ft.	\$186.00	8525*	Hose, N ₂ O, diss/Puritan, 5ft. for use with McKesson & Fraser Harlake outlet stations	\$186.00						
8025*	Hose, O ₂ , diss/Puritan, 5ft.	\$186.00											
8525*	Hose, N ₂ O, diss/Puritan, 5ft. for use with McKesson & Fraser Harlake outlet stations	\$186.00											
	<table><tr><td>8035*</td><td>Hose, O₂, diss/Schrader, 5ft.</td><td>\$207.00</td></tr><tr><td>8535*</td><td>Hose, N₂O, diss/Schrader, 5ft.</td><td>\$207.00</td></tr></table>	8035*	Hose, O ₂ , diss/Schrader, 5ft.	\$207.00	8535*	Hose, N ₂ O, diss/Schrader, 5ft.	\$207.00						
8035*	Hose, O ₂ , diss/Schrader, 5ft.	\$207.00											
8535*	Hose, N ₂ O, diss/Schrader, 5ft.	\$207.00											
	<table><tr><td>8045*</td><td>Hose, O₂, diss/Oxequip, 5ft.</td><td>\$186.00</td></tr><tr><td>8545*</td><td>Hose, N₂O, diss/Oxequip, 5ft.</td><td>\$186.00</td></tr></table>	8045*	Hose, O ₂ , diss/Oxequip, 5ft.	\$186.00	8545*	Hose, N ₂ O, diss/Oxequip, 5ft.	\$186.00						
8045*	Hose, O ₂ , diss/Oxequip, 5ft.	\$186.00											
8545*	Hose, N ₂ O, diss/Oxequip, 5ft.	\$186.00											
	<table><tr><td>8075*</td><td>Hose, O₂, diss/Ohio Medical/ Matrx, 5 ft.</td><td>\$207.00</td></tr><tr><td>8575*</td><td>Hose, N₂O, diss/Ohio Medical/ Matrx, 5 ft.</td><td>\$207.00</td></tr></table>	8075*	Hose, O ₂ , diss/Ohio Medical/ Matrx, 5 ft.	\$207.00	8575*	Hose, N ₂ O, diss/Ohio Medical/ Matrx, 5 ft.	\$207.00						
8075*	Hose, O ₂ , diss/Ohio Medical/ Matrx, 5 ft.	\$207.00											
8575*	Hose, N ₂ O, diss/Ohio Medical/ Matrx, 5 ft.	\$207.00											
<div>Note Hose Sizes</div>	<div>*ALL STOCK NUMBERS LISTED ABOVE ARE FOR A 5 FOOT HOSE ASSEMBLY. Last digit indicates length in feet. Example: 8053 Hose, O₂, diss/Porter - 3 feet long.</div> <div>Note: Hoses longer than 5 feet are custom made. Add \$6.00 per foot ALL SALES ARE FINAL</div>												
	<table><tr><td>5602</td><td>Vacuum Quick Connect</td><td>\$102.00</td></tr></table>	5602	Vacuum Quick Connect	\$102.00									
5602	Vacuum Quick Connect	\$102.00											

FLOWMETER MOUNTS

PORTER

	<div> <div>2020</div> <div> Solid Wall Mount (Telescoping) Extends to 18" -OR- </div> <div>\$350.00</div> </div> <div> <div>2022</div> <div> Wall Mount (Telescoping) with Cut-Out for Porter Triple Outlet Station Model 6255-1 (only) </div> <div>\$379.00</div> </div>
	<div> <div>2023</div> <div>Non-Telescoping Wall Bracket</div> <div>\$173.00</div> </div>
	<div> <div>2030</div> <div>Cabinet Mount-Swing Away</div> <div>\$485.00</div> </div>
   	<div> <div>2035-1</div> <div>Cabinet Mount, Slide Away</div> <div>\$388.00</div> </div> <div> <div>2036-2</div> <div> Cabinet Slide Mount For Bag Tee and Porter AVS-5000 (5" wide) </div> <div>\$410.00</div> </div> <div> <div>2036-3</div> <div> Cabinet Slide Mount For Flowmeter </div> <div>\$420.00</div> </div> <div> <div>2036-4</div> <div> Cabinet Slide Mount For Bag Tee and Porter AVS-5000 (3" wide) </div> <div>\$410.00</div> </div> <div> <div>2037-1</div> <div> Under Cabinet Slide Mount For Bag Tee and Porter AVS-5000 </div> <div>\$410.00</div> </div>
	<div> <div>5335</div> <div> Chair/Light Post Mount (BagTee or AVS) </div> <div>\$235.00</div> </div>

	<p>2040 Compact Mobile Stand \$420.00 (Top of Meter Telescopes from 24" to 36")</p> <p>2042 Tall Mobile Stand \$430.00 (Top of Meter Telescopes from 43" to 50")</p>
	<p>2045-2 Tall Mobile Stand \$2395.00 for 2 "E" Tanks of O2 and 2 "E" Tanks of N2O (Top of meter telescopes from 43" to 50") <i>Gas Supply Tubing included</i></p>
	<p>2045-Short Under Cabinet Stand \$2506.00 for 2 "E" Tanks of O2 and 2 "E" Tanks fo N2O (Top of meter telescopes from 31" to 41") <i>Gas Supply Tubing included</i></p>
	<p>C-349-2 E-Block Assembly \$2000.00</p>
	<p>A-3492-000 Tank Restraints \$89.00</p>



1536-1 Positive Pressure/Demand Valve \$1019.00 with Universal Mask; 10 foot Hose; Mobile Cart; Regulator and Yoke; and Porter Connection

The Porter/O-TWO Demand Valve Resuscitator provides a fast, simple and effective means to ventilate non-breathing patients and provides 100% oxygen on Demand.

Designed to meet the varying requirements of the Emergency Medical Services/Rescue Personnel, it is the only Demand Valve Resuscitator that comes standard with an Audible Pressure Relief System and sensitive Manual Ventilation Button, conveniently located on the front of the resuscitator. The Manual Ventilation Button helps to ensure a patient airway by enabling both hands to remain on the mask to facilitate a superior mask to face seal (without having to deal with cumbersome trigger levers or inconveniently located top buttons).

The Manual Ventilation Flow rate is fixed at 40 lpm in line with the 1992 American Heart Association Recommendations and Guidelines for Resuscitation.

During spontaneous respirations, the Demand Valve Resuscitator is designed to provide 100% oxygen to a breathing patient with minimal respiratory effort. The oxygen demand flow is controlled by the patients inspiratory effort, with flow rates up to 160 lpm, to accommodate their own respiratory distress situation.

1538 Mobile Cart Only \$166.00 (Tank not included)

1506-1 Positive Pressure/Demand Valve \$593.00 with Universal Mask, 10 foot hose and Porter Connection


1540-2 Oxygen Demand Valve \$729.00

The Porter/O-TWO Oxygen Demand Valve provides trained individuals with a safe and effective means of delivering 100% oxygen, by inhalation, in line with the flow rate demanded by the patient.

Lightweight, portable and extremely durable, the Oxygen Demand Valve is designed for the demands of the pre-hospital and in-hospital environment and can be used anywhere a regulated oxygen supply is available.

The oxygen Demand Valve is supplied with a B-V-M Oxygen Supply Valve Adapter. This allows the Oxygen Supply Valve to be connected to the intake valve on a Bag-Valve-Mask resuscitator and refills the bag with 100% oxygen when the bag squeeze is released, thereby replacing the standard reservoir system.


Components on Page 29



8201 Nitrogen Manifold System \$1,825.00
Includes one regulator & one 3ft. hose


8202 Nitrogen Manifold System \$2,335.00
Includes two regulators & two 3ft. hoses

8203 Nitrogen Manifold System \$2,843.00
Includes three regulators & three 3ft. hose



6400 Nitrogen Outlet Station \$507.00
with Schrader Quick Connect


6410 Nitrogen Outlet Station \$507.00
with Imperial Quick Connect



PARTS ONLY

1.7400-1 Nitrogen Regulator Only \$398.00

2. 8803 Nitrogen Hose, diss/diss, 3ft. \$115.00
Note: Stock number listed is for 3 ft. Last digit indicates length.
Example: 8805 hose N₂, diss/diss=5ft. long




3. 8200 Nitrogen Manifold Only \$1,380.00

4. PA-1140-000 Nitrogen Check Valve Assy. \$80.00

5. PA-1143-000 Nitrogen Gauge \$62.00

PORTER NAVICOM® Communication System




PNC-30W	Navicom Wall Station	\$807.00
PNC-30D	Navicom Desk Station	\$883.00
PNC-PCWJ	Computer Wall Jack Kit w/cd software	\$528.00
PNC-PS1	Power Supply	\$743.00
PV-100	Round Cable w/PVC Jacket (100 feet)	\$204.00
PV-250	Round Cable w/PVC Jacket (250 feet)	\$500.00
TF-100	Round Cable w/Teflon Jacket (100 feet)	\$416.00
TF-250	Round Cable w/Teflon Jacket (250 feet)	\$1003.00



24 BUTTON

OCS-24W 24 Button Wall Unit \$613.00


OCS-24D 24 Button Desk Unit \$758.00
Includes: wall outlet (OCS-WJ)
and 6 ft. cable



POWER SUPPLY

OCS-PS1 Power Supply Unit \$743.00
contains ON/OFF switch and 110 volt power
plug
UL approved

*One power supply is needed for every 20 stations and/or
1000 feet of cable.*



CABLE

PV-100 Round Cable \$204.00
w/PVC Jacket (100 feet)
Four conductor, 16 gauge, with ground wire
shield

PV-250 Round Cable \$500.00
w/PVC Jacket (250 feet)
Four conductor, 16 gauge, with ground wire
shield

TF-100 Round Cable \$416.00
w/Teflon Jacket (100 feet)
Four conductor, 16 gauge, with ground wire
shield


TF-250 Round Cable \$1003.00
w/Teflon Jacket (250 feet)
Four conductor, 16 gauge, with
ground wire shield



*The average installation requires approximately 40 feet of
cable per station. This amount will vary according to your
specific installation.*

OCS-MR Pre-Installation Bezel \$21.00
for wall unit location in dry wall.
Provided at no charge when requested with
complete OCS order.

Note: See Porter Instrument drawing FM-634 for complete
installation instructions.




MXR and Standard Flowmeters

STOCK NO .	DESCRIPTION	
1. PA-236-PL	Plastic Window Covers both O₂/N₂O tubes	\$34.00
PA-236-001	Plastic Window for MXR 3000 (white), & 1100 (chrome)	\$34.00
2. P2400	Sterile Knob Kit	\$77.00
3. A-2134	Emergency Air Intake Valve Assembly	\$96.00
4. A-2146	Non-rebreathing Check Valve Assembly	\$49.00
5. 4650	Flapper Valve Only	\$11.00
6. P1407A	Bag Tee Assembly	\$220.00
7. P2407A	Oral Surgery Bag Tee Assembly (threaded)	\$297.00
A-1890-TF4	Teflon Washer for Oral Surgery Bag Tee Only	\$9.00
8. A-3379-000	Cannula/Bag Tee Adapter	\$80.00
PA-296-001	Viton Sealing Washer for Bag Tees	\$4.00
9. A-1762	On/Off Switch Assembly	\$80.00
9A.PA-128-PL3	Green Knob for On/Off Switch	\$12.00
10.A-4059-000	On/Off Switch Rebuilt Kit	\$22.00

FLOWMETERS (ONLY)

3000	MXR-Flowmeter Head Only (includes Bag Tee)	\$2,502.00
2000-OS	Oral Surgery Meter Only (includes Bag Tee)	\$2,582.00
1100	Standard Flowmeter Head Only (includes Bag Tee)	\$2,381.00
3065W	MXR-1 with Universal Connection	\$3,017.00

Cabinet Flowmeter Flush Mount (Bag Tee Not Included)



FLEXIBLE INSTALLATION

STOCK NO.	DESCRIPTION	
PA-450-1-001	1/4" ID Green O₂ Hose, 20ft.	\$58.00
A-2633-000	Through Wall Adapter Assy. for Bag Tee	\$151.00
-OR-		
A-1679	Remote Surface Mount for Bag Tee Assembly	\$151.00
B-2460-001	O₂ Shut-Off Valve, diss fittings with check valves	\$122.00
B-2461-001	N₂O Shut-Off Valve, diss fittings with check valves	\$122.00
A-2508-000	Rubber Goods Hook	\$31.00
P1407A	Bag Tee Assembly	\$220.00
8005	O₂ diss/diss Hose, 5ft.	\$115.00
8505	N₂O diss/diss Hose, 5ft.	\$115.00

ACCESSORIES

STOCK NO.	DESCRIPTION	
B-2625	Bezel Mounting Bracket	\$88.00
	Can be used with either flexible or hard copper installation	
PA-128-PL3	Green O₂ Flush Knob	\$12.00
	Complete with set screw	


GAS SUPPLY FITTINGS ONLY

A-1675	Porter O₂ diss/Quick Disconnect Fitting	\$93.00
PA-866-000	O₂ diss/NCG (Chemetron) Fitting	\$117.00
PA-864-000	N₂O diss/NCG (Chemetron) Fitting	\$117.00
PA-867-000	O₂ diss/Puritan (McKesson & Fraser Harlake) Fitting	\$117.00
PA-865-000	N₂O diss/Puritan (McKesson & Fraser Harlake) Fitting	\$117.00
PA-927-000	O₂ diss/Schrader Twist Fitting	\$126.00
PA-928-000	N₂O diss/Schrader Twist Fitting	\$126.00
PA-871-000	O₂ diss/Oxequip Fitting	\$126.00
PA-872-000	N₂O diss/Oxequip Fitting	\$126.00
PA-868-000	O₂ diss/Ohio Medical (Matrx) Fitting	\$126.00
PA-869-000	N₂O diss/Ohio Medical (Matrx) Fitting	\$126.00

GAS SUPPLY TUBING

B-1329-5	Duplex Hose, O₂/N₂O, Male diss/Porter, 5ft.	\$254.00
	used on NCG Flowmeter.	







SPECIAL ORDER ALL SALES FINAL





ITEMS BELOW IN GRAY BOXES ARE STEAM AUTOCLAVEABLE!

STOCK NO.	DESCRIPTION	
1. 5054A	Porter Adult Hood Complete with 3 Liners	\$242.00
1a. 5054-1	Adult Inner Hood Liners Package of 3	\$49.00
1b. 5054C	Porter Adult Hood with 1 Liner	\$215.00
2. 5054B	Porter Pedo Hood Complete with 3 Liners	\$242.00
2a. 5054-2	Pedo Inner Hood Liners Package of 3	\$49.00
2b. 5054D	Porter Pedo Hood with 1 Liner	\$215.00
3. 5056	Coaxial Inhaler Tubing 30"	\$182.00
4. 5057	Tension Slide Adjuster	\$18.00
5. 5058	Fresh Gas Y Connector	\$21.00
6. 5059	Vacuum Hose (8 feet)	\$41.00
PA-679-000	Vacuum Hose by the foot	\$7.00
7. 5060-3	Fresh gas Corrugated Tubing (3 feet)	\$38.00
5060-6	Fresh gas Corrugated Tubing (6 feet)	\$60.00
8. 1571-22	22mm Right Angle Adapter	\$23.00
1570-24	24mm Right Angle Adapter	\$34.00
9. 5061	Mask to Tubing Plastic Connectors (set of 2)	\$25.00
10. 5062	Vacuum Line Y Connector	\$23.00
11. 5501-RK	In-Line Vacuum Kit (can be utilized as pictured above)	\$96.00
12. 4100	3 Liter Bag-Latex	\$36.00
4100-3NL	3 Liter Bag Non-Latex	\$36.00
4100-2NL	2 Liter Bag Non-Latex	\$36.00
13. 5063	T Adapter for In-Line Vacuum Block 1/2"	\$11.00
5068	T Adapter for In-Line Vacuum Block 5/8"	\$11.00
14. 5064	Straight Adapter for In-Line Vacuum Block	\$11.00
15. 5065	Vacuum Tuber Holder	\$21.00
16. 5066	Porter Liner Valve Only (package of 3)	\$21.00
17. A-3679-000	HVE Adapter	\$27.00
5602	Vacuum Quick Connect	\$102.00

Disposable Liners found on page 7.


	<p>B-1396-000 7/8" Delrin Adapter \$84.00 Manufacturer and Model: FRASER SWEATMAN FRASER HARLAKE OHMEDA MATRIX QUANTIFLEX MDM ACCUTRON BELMED MDT-McKESSON ANALOR FOUR PARKELL-SEDATRON TRIDENT</p>
	<p>P1407A Bag Tee Standard \$220.00 A-3379-000 Canula Adapter \$80.00 with PB-102-7-BN Buna O-ring Manufacturer and Model: COASTAL NCG CHEMETRON</p>
	<p>P1407A Bag Tee Standard \$220.00 A-3378-000 Canula Adapter \$124.00 with PB-102-116BN Buna O-ring Manufacturer and Model: HAMPTON SIMPLAIRE MIRCO</p>
	<p>P1407B Special Bag Tee \$235.00 Manufacturer and Model: McKESSON ANALOR I&II STARFLIGHT DENETRON</p>
	<p>P1407A Bag Tee Standard \$220.00 A-3380-000 Canula Adapter \$124.00 with PB-102-44-BN Buna O-ring Manufacturer and Model: McKESSON ANALOR 3C</p>
	<p>5651 Right Angle Adapter \$49.00 Manufacturer and Model: Porter Instrument Kit Replaces Porter Clean-Air™ Rubber Goods (mask with plastic sweeps and canister)</p> <p>Note: It is not recommended to retrofit Porter Scavenger Breathing Circuit to Adec Veriflo or Spica machines.</p>

Above Adapters in addition to the PORTER Scavenger Rubber Goods with In-Line Vacuum and 3L Breathing Bag make up a complete PORTER Retrofit Kit found on page 9 in this catalog.




STOCK NO.	DESCRIPTION	
A. A-2521-001	Electric Pneumatic Solenoid	\$117.00
B. A-2554-000	O ₂ Gauge 0-150-psi	\$62.00
C. A-2554-001	N ₂ O Gauge 0-150-psi	\$62.00
D. B-1322	O ₂ diss Check Valve Assembly	\$68.00
E. B-1323	N ₂ O diss Check Valve Assembly	\$68.00
F. B-1324	Pressure Relief Valve Assembly	\$86.00
G. C-394-001	Pressure Switch Assembly	\$230.00
H. A-2144	Diss/Copper Pigtails (straight) Oxygen - 1/2 O.D. copper x diss O ₂ female Nitrous Oxide - 3/8 O.D. copper x diss N ₂ O female	\$113.00
1. B-1780-000	Molded Plastic Manifold Cover	\$128.00
2. B-1805-002	Power Supply	\$258.00
3. KA18028-020	75ft PVC Cable, #22 AWG, 4 Conductor for Vanguard	\$115.00
KA18028-021	75ft Teflon Cable #22 AWG, 4 Conductor for Vanguard	\$2.87/ft.
KA18028-000	PVC Cable, #22 AWG, 4 Conductor for Vanguard-PER FOOT	\$3.00/ft.
KA18028-001	Teflon Cable, #22 AWG, 4 Conductor for Vanguard-PER FOOT	\$4.00/ft.
4. B-1811-000	Concealed Interconnect Assembly for Wall Alarm	\$144.00
5. 5700-10	10ft Modular Cable (10ft)	\$30.00
7. A-2522-000	Signal Interconnect Assembly for Desk Alarm	\$144.00
8. A-3381-000	Vanguard Manifold PC Board Assembly Only	\$969.00
8a. A-3381-001	Vanguard Manifold PC Board Assembly Only (model "A" only)	
3653-1	Vanguard Wall Pre-Install	\$573.00
3653-2	Vanguard Desk Pre-Install	\$573.00
7000-2	Oxygen Regulator, Hi pressure gauge only, 50 PSI	\$369.00
7500-2	Nitrous Oxide Regulator, Hi pressure gauge only, 50 PSI	\$369.00
PA-514-000	Oxygen Pipeline Label, vinyl	\$2.00
PA-515-000	Nitrous Oxide Pipeline Label, vinyl	\$2.00
3525	Set of two (2) - O ₂ -N ₂ O Restraints	\$134.00

**NOTE: MODULAR CABLE SHOULD NOT BE CUT OR SPLICED.
THE VOLTAGE DROP IS TOO GREAT AND IT WILL AFFECT THE SYSTEM OPERATION.**




STOCK NO.	DESCRIPTION	
A. A-2554-000	O₂ Gauge, 0-150 psi	\$62.00
B. A-2554-001	N₂O Gauge, 0-150 psi	\$62.00
C. B-1322	O₂ diss Check Valve Assembly	\$68.00
D. B-1323	N₂O diss Check Valve Assembly	\$68.00
E. B-1324	Pressure Relief Valve Assembly	\$86.00
F. C-394-000	Pressure Switch Assembly	\$230.00
G. A-2144	diss/Copper Pigtails (straight)	\$113.00
	Oxygen - 1/2 O.D. copper x diss O₂ female	
	Nitrous Oxide - 3/8 O.D. copper x diss N₂O female	
H. A-2041	Sentinel Interconnect Board Assembly	\$68.00
I. D-160-PL	Molded Plastic 3222MX Cover	\$60.00
1. 5775	75' Sentinel Modular Connector Cable (standard)	\$146.00
5700-150	150' Sentinel Modular Connector Cable (special order)	\$198.00
5700-200	200' Sentinel Modular Connector Cable (special order)	\$223.00
2. 5700	50' Sentinel Modular Connector Cable	\$128.00
3. 5700-1T	100' Teflon Sentinel Modular Connector Cable	\$287.00
4. KC13161-001	Power Supply for Wall Sentinel	\$119.00
6. C-637-000	Power Supply for Desk Sentinel	\$119.00
7. A-2083	Modular Cable Connector Box	\$139.00
PA-636	Sentinel Light Bulb All Glass	\$7.00
PA-571	Sentinel Light Bulb Old Bayonet Style	\$9.00



MOBILE STAND


20. A-3369-000	Mounting stud w/Washer for all Porter Mobile Stands	\$41.00
21. A-3372-000	Inner Post Assy. with mounting studs	\$119.00
22. A-3399-000	Yoke Washers	\$2.00
23. B-1883-000	Compact Stand Post/Tube Assy. Only (15" to 27" tall)	\$205.00
24. B-1883-001	Tall Stand Post/Tube Assy. Only (33" to 41" tall)	\$249.00
25. PA-1185-000	Tank Wrench with chain	\$10.00
26. A-3492-000	E-stand Tank Restraint Kit	\$89.00
27. PA-232-1	19" Base	\$157.00
28. PA-234-002	Caster	\$33.00
29. PA-4060-000	"T" Handle	\$26.00
30. A-3279K	Check Valve Kit	\$64.00



UNDER CABINET MOBILE STAND

31. B-4669-000	S Bracket Assy	\$86.00
32. A-4052-000	Inner Post Assembly	\$123.00
33. A-3399-000	Yoke Washers	\$2.00
34. A-2508-002	Rubber Goods Hook	\$31.00
35. B-4668-000	Post/Tube Assy.	\$249.00
36. PA-1185-000	Tank Wrench w/chain	\$10.00
37. A-3492-000	E-Stand Tank Restraint Kit	\$89.00
38. PA-232-1	19" Base	\$157.00
39. PA-234-002	Caster	\$33.00
40. PA-4060-000	"T" Handle	\$26.00
41. 8082	O₂ Hose 90° elbow*	\$147.00
42. 8582	N₂O Hose 90° elbow*	\$147.00

*Not Shown


Wall/Post Mounts

1. A-1613-AL1	Mounting Stud Assembly for 2020, 2022, 2033	\$64.00
2. A-1614-000	Telescoping Arm Assy. for 2020 & 2022	\$213.00
3. B-810-000	Solid Back Plate for 2020	\$126.00
4. B-805-000	Cut Out Back Plate for 2022	\$154.00
5. 2023	Non-Telescoping Wall Bracket	\$173.00
6. 2030	Cabinet Mount-Swing Away	\$485.00
7. 2033	Remote Bag Tee Mount Kit	\$166.00
8. 5335	Chair/Light Post Mount (Bag Tee or AVS)	\$235.00
9. A-2633-000	Through Wall/Cabinet Mount (Bag Tee or AVS)	\$151.00
10. A-2800-000	Knob for 2030	\$27.00
11. A-1037-AL2	Meter Mounting Stud for 2030, 2036-3	\$64.00
12. B-2075-000	Mounting Block for 3000 MXR	\$75.00
13. A-1679	Remote Bag Tee Mount	\$151.00
14. PA-450-1-BR	1/4" ID Green O ₂ Tubing (specify length)	\$5.00/foot


OUTLET STATION REPLACEMENT PARTS

PORTER

For Models 6100-1, 6200-1, 6300-1, 6256-1, 6256-3, 6258-1, 6255-1

PORTER Instrument advises replacing both N₂O and O₂ check valves when servicing Outlet Stations.

Numbers 3 through 6 represent the Secondary Check Valve Assembly. This valve is only in use when the Primary Check Valve is removed for servicing and will NOT remedy a leaking outlet station.



STOCK NO.	DESCRIPTION	
6200K	Outlet Station Repair Kit Includes N₂O, O₂ primary check valves and O-rings for 6200, 6256, 6300	\$144.00
6200-1K	Outlet Station Repair Kit Includes N₂O, O₂ primary check valves and O-rings for 6200-1, 6256-1, 6256-3, 6300-1, 6258-1, 6255-1	\$109.00
1. A-2689-000	O₂ Primary Check Valve Assembly (slotted back)	\$58.00
A-2690-000	N₂O Primary Check Valve Assembly (slotted back)	\$58.00
*A-1655-A	O₂ Primary Check Valve Assembly for 6100, 6200, 6256, 6300	\$77.00
*A-1655-B	N₂O Primary Check Valve Assembly for 6200, 6256, 6300	\$77.00
2. PB-102-104NE	014 O-ring for O ₂ or N ₂ O Primary Check Valve, outside	\$4.00
3. B-2237-000	O₂ Cartridge Body	\$29.00
B-2275-000	N₂O Cartridge Body	\$29.00
4. PB-102-215	016 O-ring	\$4.00
5. A-2699-000	Secondary Check Valve Assembly with O-ring	\$21.00
6. PB-115-63SS	Compression Spring	\$4.00
7. B-1231-003	Vacuum Check Valve Assembly	\$64.00
A-3356-000	5600-3 Vacuum, Face Plate Kit	\$51.00
8. B-2568-000	6200-1 Dual, Face Plate Kit	\$45.00
9. LB-119-000	Vacuum Pipeline Label	\$2.00
10. B-2568-000	6256-3 Triple, Face Plate Kit	\$93.00
B-2737-004	6255-1 Triple Face Plate Kit	\$93.00
11. PA-514-000	O₂ Pipeline Label Vinyl	\$2.00
12. PA-515-000	N₂O Pipeline Label Vinyl	\$2.00







***Note: Primary and secondary check valve assemblies marked (*) above are for outlet stations manufactured prior to February 1986 only.**



5602-RK	In-Line Vacuum to 6256-1 Retrofit Kit includes:	\$170.00
1. 5602	Vacuum Quick Connect	\$102.00
2. B-1231-003	New Knotched Vacuum Check Valve Assembly	\$64.00
3. LB-169	O₂/N₂O Label	\$4.00
4. LB-343-000	Vacuum Label (not pictured)	\$4.00

ORAL SURGERY

PORTER

	2210K	Directional Y Valve for instant selection of nasal inhaler or full face mask.	\$568.00
	A-1890-TF4	Teflon Sealing Washer for 2210K	\$9.00
	1597	Non-Rebreathing Valve	\$220.00
	1550	Pedo Size Full Face, Clear Mask	\$77.00
	1575-1	Universal Size Full Face, Clear Mask	\$77.00
	4200	Corrugated Tube (Latex)	\$53.00
	1585	Plastic Elbow (15mm - 22mm)	\$11.00
	2044	Oral Surgery Star-Base Stand	\$474.00


OXYGEN THERAPY METERS AND OXYGEN RESUSCITATION PARTS



A-3919-000	O₂ Flowmeter 0-8 Liter with Porter Connector	\$167.00
A-3919-001	O₂ Flowmeter with diss connection	\$111.00
A-1675	Porter O₂ diss/Quick Connect	\$93.00
1560-1	Regulator and Yoke w/diss Male Connection	\$434.00
1566-1	Regulator and Yoke w/Porter Imperial Connection	\$434.00

PARTS LIST FOR SERIAL NUMBERS 1000 TO 3400

Stock No.	Description	
KA13250-RK	15 Amp Circuit Breaker Kit (replaces fuseholder)	\$38.00
PA-926-000	Air Vent Assembly	\$52.00
KA18144-000	Thermostat Assembly (SES-2000 ONLY)	\$115.00
B-2921-000	Fusible Link	\$51.00
A-3354-000	Chamber Fittings Seal Package	\$38.00
KA13206-000	SES-2000 2.0 Amp Fuse for PC Board	\$6.00
*A-3184-000	Door Sealing Gaskets -Package includes 3 different seals	\$38.00
KC13311-001	Heater	\$215.00
A-3897-000	Water level Sensor	\$31.00
A-3338-000	Safety Relief Valve Assembly , 38 psig.	\$122.00
BX-267	Shipping Box with foam inserts	\$37.00
A-3224-RK	100% drainable Reservoir	\$218.00






**All Sales Final for
All Sterilizer Parts**

Additional Parts List Available on Website


STOCK NO.	DESCRIPTION	
SES2000E		
1. PB-138-025	C Clip	\$4.00
2. PB-123-010	Bowed Washer	\$5.00
3. A-2959-000	Anti-Friction Washer	\$5.00
*4. A-3184-000	Door Sealing Gaskets Package includes 3 different seals	\$38.00
5. B-2915-000	Door Clamp Disc	\$27.00
6. PB-892-000	Silicone Washer	\$5.00
7. A-2961-000	Stainless Steel Washer	\$5.00
8. PB-891-000	Airtight Nut	\$7.00
SES2012E		
*4. A-3314-000	Door Sealing Gaskets -Package includes 2 different seals	\$37.00
6. PA-1037-001	Silicone Washer	\$5.00
7. PB-103-050	Stainless Steel Washer	\$5.00
8. PA-1042-000	Airtight Nut	\$7.00

***We recommend that all door gaskets be changed at the same time; therefore, a package of 1 each of the different seals is supplied.**



OUTLET STATIONS



WALL JACK CUT-OUT



NAVICOM WALL STATION CUT-OUT



Use same dimensions for the following:

Vanguard	6251-A
Sentinel	5251
OCS	OCS-24W

WALL STATION CUT-OUT



PORTER

**Website: www.porterinstrument.com
Fax#: for orders 215-723-5106**

**Parker Hannifin Corporation
Porter Instrument Division
245 Township Line Road
P.O. Box 907
Hatfield, PA 19440-0907 USA
Phone (215) 723-4000
Fax (215) 723-5106**