



Better Health. Better Future.

Pandemic Preparedness

Products that provide protection for healthcare professionals working in infectious conditions.



**Enhancing family &
community health
is the way to a
better future**

Gowns

Isolation Gowns	4
CPE Gowns	6
AAMI Level 1 Gowns.....	8
AAMI Level 2 Gowns.....	10

Disposable Scrubs and Coveralls

Coversoft™ Disposable Apparel.....	12
Disposable Scrub Wear	13
Disposable Coveralls	14

Caps and Shoe Covers

Bouffant Caps.....	15
Shoe Covers.....	16

Lab Coats

Lab Coats.....	17
----------------	----

Patient Exam Wear

Patient Exam Gowns.....	18
Patient Exam Shorts.....	18

Face Masks

General Use Face Mask	19
Surgical Masks	20
Face Shield	22

Diagnostics and Instruments

Digital Thermometer	23
Disposable Stethoscope	24
Bowles Stethoscope.....	25
Single Head Stethoscope	26
Dual Head Stethoscope.....	27
Liquid Dispenser.....	28
Non-Latex Tourniquet	29
Quick Release Tourniquet	30
Finger Cots	31
Nitrile Finger Cots.....	32
IV Stands.....	33

Post Mortem

Body Bag.....	34
---------------	----



Isolation Gowns

Made of high quality spunbonded polypropylene

Style Guide



Wrist Style



Neck Style

- Each gown features elastic cuffs with waist and neck tie closures
- Available in multiple colors to allow easy detection between departments or functions



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
301	Isolation Gown, Yellow	10/bg, 5 bg/cs
301BL	Isolation Gown, Blue	10/bg, 5 bg/cs



Poly Coated Isolation Gowns

Made of high quality spunbond polypropylene with a polyethylene coating

Style Guide



Wrist Style



Neck Style

- Each gown features elastic cuffs with waist and neck tie closures
- Available in multiple colors to allow easy detection between departments or functions



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
303	Isolation Gown, coated, Yellow	10/bg, 5 bg/cs
303BL	Isolation Gown, coated, Blue	10/bg, 5 bg/cs



CPE Isolation Gowns

Polyethylene film provides enhanced fluid protection

- Over-the-head protective cover gown with thumb loops
- Tear away design enables gowns to be removed quickly
- Thumb loop holes provide easy glove donning
- Open back style provides comfortable fit
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
307R	CPE Gown, Regular	15/bg, 5 bg/cs
307XL	CPE Gown, XL	15/bg, 5 bg/cs

Style Guide



Wrist Style



Neck Style



AAMI LEVEL 1

Isolation Gowns

Made of multi-ply SMS material

- Meets AAMI Level 1 requirements for fluid protection
- Provides a slight barrier against small amounts of fluid penetration
- Elastic wrist with thumb loops to provide a more secure fit under gloves
- Tape tab neck with waist ties
- Used in minimal risk situations
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
308	Level 1 Isolation Gown, Universal, Yellow	10/bg, 10 bg/cs
308XL	Level 1 Isolation Gown, XL, Yellow	10/bg, 10 bg/cs

Style Guide



Wrist Style



Neck Style



AAMI LEVEL 2

Isolation Gowns

Made of multi-ply SMS material

- Meets AAMI Level 2 requirements for fluid protection
- Proprietary over-the-head design with easy tear away removal
- Elastic wrist with thumb loops to provide a more secure fit under gloves
- Used in low risk situations
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
309	Level 2 Isolation Gown, Universal, Yellow	10/bg, 10 bg/cs
309XL	Level 2 Isolation Gown, XL, Yellow	10/bg, 10 bg/cs

Style Guide



Wrist Style




Neck Style



CoverSoft™ Disposable Apparel

Made of durable and extra soft SMS material that feels like real fabric

Fabric Detail



- Extra soft fabric provides optimal comfort
- Scrub pants have a drawstring tie waist, one front pocket, and one back pocket
- Shirts have two pockets (one top left, one bottom right)
- Dark blue for maximum coverage
- Not made with natural rubber latex



Item #	Description	Packaging
Non-Sterile		
376S	Shirt, Small	10/bg, 5bg/cs
376M	Shirt, Medium	10/bg, 5bg/cs
376L	Shirt, Large	10/bg, 5bg/cs
376XL	Shirt, X-Large	10/bg, 5bg/cs
376XXL	Shirt, XX-Large	10/bg, 5bg/cs
381S	Pants, Small	10/bg, 5bg/cs
381M	Pants, Medium	10/bg, 5bg/cs
381L	Pants, Large	10/bg, 5bg/cs
381XL	Pants, X-Large	10/bg, 5bg/cs
381XXL	Pants, XX-Large	10/bg, 5bg/cs

Disposable Scrub Wear

Made of durable spunbonded polypropylene

- Scrub pants have a drawstring tie waist, one front pocket, and one back pocket
- Shirts have two pockets (one top left, one bottom right)
- Dark blue
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
375S	Shirt, Small	10/bg, 5bg/cs
375M	Shirt, Medium	10/bg, 5bg/cs
375L	Shirt, Large	10/bg, 5bg/cs
375XL	Shirt, X-Large	10/bg, 5bg/cs
375XXL	Shirt, XX-Large	10/bg, 5bg/cs
380S	Pants, Small	10/bg, 5bg/cs
380M	Pants, Medium	10/bg, 5bg/cs
380L	Pants, Large	10/bg, 5bg/cs
380XL	Pants, X-Large	10/bg, 5bg/cs
380XXL	Pants, XX-Large	10/bg, 5bg/cs



Disposable Coveralls

Made of durable spunbonded polypropylene

Style Guide



Wrist Style

- Elastic wrist and ankle cuffs
- Front zip closure
- Available in two different colors



- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
382XL	White, XL	5/bg, 5bg/cs
382XXL	White, XXL	5/bg, 5bg/cs
382XXXL	White, XXXL	5/bg, 5bg/cs
392XL	Blue, XL	5/bg, 5bg/cs



Bouffant Caps

Made of a comfortable spunbonded polypropylene



- Cool and can be worn for hours
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
321B	21", Blue	100/bx, 5 bx/cs
321-10	21", Accordion Style, Blue	100/bg, 10 bg/cs
324B	24", Blue	100/bx, 5 bx/cs
324-10	24", Accordion Style, Blue	100/bg, 10 bg/cs

Shoe Covers

Made of spunbonded polypropylene



- Non-skid to provide better traction
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
350	Shoe Covers, Non-Skid, Blue	100/bx, 3 bx/cs
350-10	Shoe Covers, Non-Skid, Blue	100/bg, 10 bg/cs

Lab Coats

Multi-layer SMS material
to provide quality and comfort

- Full length with knit collar and cuff
- Convenient hang tab and size tag sewn into back of coat
- White
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile - Without Pockets		
340	Small	10/bg, 5 bg/cs
341	Medium	10/bg, 5 bg/cs
342	Large	10/bg, 5 bg/cs
343	XL	10/bg, 5 bg/cs
344	XXL	10/bg, 5 bg/cs
345	XXXL	10/bg, 5 bg/cs
Non-Sterile - With Pockets		
340P	Small	10/bg, 5 bg/cs
341P	Medium	10/bg, 5 bg/cs
342P	Large	10/bg, 5 bg/cs
343P	XL	10/bg, 5 bg/cs
344P	XXL	10/bg, 5 bg/cs
345P	XXXL	10/bg, 5 bg/cs





Patient Exam Wear

Made of spunbonded polypropylene

- Made of spunbonded polypropylene
- Gown has open back with belt
- Gown also available in a deluxe version with neck tie
- Shorts have elastic waist and no pockets
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
304	Exam Gown, Dark Blue	10/bg, 5 bg/cs
305	Deluxe Exam Gown, Dark Blue	10/bg, 5 bg/cs
360	Shorts, Unisize, Dark Blue	10/bg, 5 bg/cs

General Use Face Mask

Lightweight mask that provides comfortable face protection



- 3-ply, pleated with earloop
- Not for medical use
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1525	General use, ear loop (not for medical use)	50/bg, 40 bg/cs

ASTM LEVEL 1

Surgical Masks

3-ply pleated lightweight,
breathable material



- BFE of $\geq 99\%$
- Meets ASTM Level 1 standards
- Available in earloop or tie on with flexible nosepiece
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1540	Tie, Blue	50/bx, 6 bx/cs
1541	Ear Loop, Blue	50/bx, 6 bx/cs

Performance Characteristics	Test Standard Methods	ASTM Level 1 Requirements
Bacterial Filtration Efficiency Performance (BFE) (%)	ASTM F2101	$\geq 95\%$
Particulate Filtration Efficiency Performance (PFE) (%)	ASTM F2299	$\geq 95\%$
Differential Pressure (Delta-P) (mm H ₂ O/cm ²)	MIL-M-36945C EN 14683	< 5.0
Fluid Resistance Performance (mmHg)	ASTM 1862	80 mmHg
Flammability	16 CFR part 1610	Class I No Flame Spread

ASTM LEVEL 2

Surgical Masks

3-ply pleated lightweight,
breathable material



- BFE of $\geq 99\%$
- Meets ASTM Level 2 standards
- Available in earloop or tie on with flexible nosepiece
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
CM-101	Ear Loop, Blue	100/bg, 6 bg/cs
CM-102	Tie, Blue	50/bx, 6 bx/cs

Performance Characteristics	Test Standard Methods	ASTM Level 2 Requirements
Bacterial Filtration Efficiency Performance (BFE) (%)	ASTM F2101	$\geq 98\%$
Particulate Filtration Efficiency Performance (PFE) (%)	ASTM F2299	$\geq 98\%$
Differential Pressure (Delta-P) (mm H ₂ O/cm ²)	MIL-M-36945C EN 14683	< 5.0
Fluid Resistance Performance (mmHg)	ASTM 1862	120 mmHg
Flammability	16 CFR part 1610	Class I No Flame Spread

Protective Plastic Face Shield

Clear polyethylene for optimal visibility

- Double-sided anti-fog
- Protection from splashing or spraying
- Foam headband for gentle comfort
- Single Use
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
CM-105	Face Shield	5/bg, 40 bg/cs



Digital Thermometer

Quick 10 Second Read with easy to read large LCD display and low battery indication

- Reads both Centigrade and Fahrenheit; Range: 32° – 43.9°C, 90° – 109.9° F
- Alarm/Beeper when reading is complete
- Battery life: 200 hours of continuous operation
- Auto shut off after 9 minutes of non use
- Built in memory enabling display of the previously measured value



Item #	Description	Packaging
3210D-360	Dual Scale, 10 second read	12 ea/bx, 30 bx/cs

Disposable Stethoscopes

Single Head Stethoscope with
22" Y Piece Tubing

- Overall Length: 32"
- Individually poly bagged
- Single Patient Use
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1115	Yellow	10/bx, 10 bx/cs



Bowles Stethoscope

Chrome plated chest piece with 22" black PVC "Y" tubing and adjustable chrome plated binaural

- Classic blood pressure reading stethoscope
- Lifetime Warranty
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1211	Black	1/bx



Single Head Stethoscopes

Features 22" PVC tubing, anodized aluminum chest piece, and chrome plated brass binaural

- Basic single head stethoscope for blood pressure readings
- Color coded non chill ring
- One year warranty
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1100	Black	1/bx
1100BL	Blue	1/bx
1100GR	Green	1/bx
1100GY	Grey	1/bx
1100R	Red	1/bx



Dual Head Stethoscopes

Features 22" PVC tubing, anodized aluminum chest piece, and chrome plated brass binaural

- Basic dual head stethoscope
- One year warranty
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
1200	Black	1/bx
1200BL	Blue	1/bx
1200GR	Green	1/bx
1200GY	Grey	1/bx
1200R	Red	1/bx



Liquid Dispenser

Dispensing design produces steady stream of liquid controlled by finger pressure

- Size: 4-3/4" x 2-1/2"

Item #	Description	Packaging
Non-Sterile		
4026	9 oz. Lotion Style Pump	12 ea/bx, 12 bx/cs




Non-Latex Tourniquets

Used for constricting or compression

- 18"x1" textured nitrile
- Hypoallergenic
- Superior strength
- Convenient packaging designed for kit assembly and resale
- Not made with natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4408	1" x 18", Powder Free	10/bg, 50 bg/bx



Quick Release Tourniquet

Single handed use tourniquet with push button to release

- Universal in size and adjustments will work with child, adult, and thigh size
- Reusable
- Not made of natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4418	Red	12/bx



Latex Finger Cots

Pre-rolled and non-powdered

- Provide user protection and to help prevent contamination and infection
- Approved for medical use and intended to protect health care professionals and patients

Item #	Description	Packaging
Non-Sterile		
4423S	Small	144/bx
4423M	Medium	144/bx
4423L	Large	144/bx
4423XL	X-Large	144/bx



Nitrile Finger Cots

Pre-rolled natural white, powder free, and anti-static

- For medical exam use
- Minimum of .08mm thickness
- Not made of natural rubber latex

Item #	Description	Packaging
Non-Sterile		
4404S	Small	144/bx
4404M	Medium	144/bx
4404L	Large	144/bx
4404XL	X-Large	144/bx



IV Stands

Available in chrome plated steel
or stainless steel

- Adjustable from 47" to 84"

Item #	Description	Weight Capacity	Base Diameter	Packaging
4353	Four Caster, 2 Hook	15 pounds	23"	1/ea
4355	Five Caster, 2 Hook	22 pounds	20"	1/ea
4356	Five Caster, 4 Hook	22 pounds	20"	1/ea
4356-SS	Stainless Steel, 4 Hook	22 pounds	20"	1/ea



4353



4355



4356



4356-SS

Body Bag

Made of strong polyvinyl chloride (PVC) fabric

- Bottom of bag is leak proof
- 350lb weight capacity
- Three toe tags



Toe Tags



Item #	Description	Packaging
Non-Sterile		
400B	Body Bag, White	10/cs



Dukal Corporation

2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800

Toll Free: 1 800 243 0741

Email: customerservice@dukald.com

To Place Orders: orders@dukald.com

Find our full line of products at www.Dukal.com