

3M Infection Prevention

Surgical Masks and Respirators



N95 Healthcare Respirators




Surgical & Specialty Masks



Procedure Masks



Product Name	Breathability (ΔP : mm H ₂ O/cm ²)	Filtration BFE ¹ PFE ²	Fluid and Splash-Resistant	Additional Features
3M™ N95 Particulate Respirator and Surgical Mask — cone 1860 1860S (smaller size) 	Meets or exceeds NIOSH ⁴ specifications for ΔP	— —	Yes	Standard size surgical / laser mask (1860) Double elastic band
3M™ N95 Particulate Respirator and Surgical Mask — three-panel, flat-fold 1870 	Meets or exceeds NIOSH ⁴ specifications for ΔP	— —	(160 mm Hg) ASTM ³	Soft inner liner, Individually packaged surgical / laser mask Single dispensing, Double elastic band
3M™ Molded Surgical Mask 1800+NL 	<0.5	>98% n/a	—	Collapse-resistant Design
3M™ Standard Tie-on Surgical Mask 1816 	<2.5	>99% >98%	—	—
3M™ Tie-on Surgical Mask 1818 	<2.1	>99% >99%	—	Manufactured with 3M Filtration Technology and soft, comfortable rayon
3M™ Tie-on Surgical Mask 1818FSG 	<2.1	>99% >99%	—	Premium anti-fog, anti-reflective face shield
3M™ Ultra-soft Fog-Free Surgical Mask for Sensitive Skin 1833 	<2.0	>99% >99%	—	Soft inner liner for sensitive skin, white, all dyes or colorants removed medical grade non-woven anti-fog strip
3M™ Anti-fog Surgical Mask with Film 1834 	<2.3	>99% >99%	(80 mm Hg) ASTM ³	High medical grade anti-fog film strip
3M™ High Fluid-Resistant Surgical Mask 1835 	<2.5	>99% >99%	(160 mm Hg) ASTM ³	—
3M™ High Fluid-Resistant Surgical Mask 1835FS 	<2.5	>99% >99%	(160 mm Hg) ASTM ³	Anti-fog face shield, black anti-glare strip
3M™ High-Performance Tie-On Surgical Mask 1838R 	<3.6	>99% >99%	(120 mm Hg) ASTM ³	Fog-free design, off-the-face comfort
3M™ High-Performance Tie-On Surgical Mask 1838FSG 	<3.6	>99% >99%	(120 mm Hg) ASTM ³	Premium anti-fog, anti-reflective face shield
3M™ Filtron™ Duckbill Tie-On Surgical Mask* 1900 	<1.6	>99% >99%	—	Fog-free design, off-the-face comfort, Duck bill style, tie-on (vertical) Manufactured with 3M Filtration Technology and soft, comfortable rayon.

* Filtron 1900 surgical mask is manufactured in the United States of 89 percent U.S. components, 11 percent Austrian components and is latex-free.

¹BFE = Bacterial Filtration Efficiency; ²PFE = Particulate Filtration Efficiency meets the new classification system for mask standard ASTM F2100-04. PFE test results are provided for comparing performance of surgical masks under similar laboratory conditions. This does not reflect expected levels of filtration efficiency in actual use conditions nor replace the need to use NIOSH certified respirators if exposed to hazardous airborne contaminants. ³ASTM = American Society of Testing Materials F1862. ⁴NIOSH 42 CFR 84 Certified.



3M Procedure Masks

Procedure Masks

3M™ Procedure Mask

1820
1820FS



	1820	1820FS
Breathability	<2.1	<2.1
Filtration	BFE >99%	BFE >99%
	PFE >98%	PFE >99%
		Black anti-glare strip

3M™ Standard Procedure Mask

1826



	1826
Breathability	<2.5
Filtration	BFE >99%
	PFE >98%

3M™ High Fluid-Resistant Procedure Mask

1840



1840FS



	1840	1840FS
Breathability	<2.5	<2.5
Filtration	BFE >99%	BFE >99%
	PFE >99%	PFE >99%
		Anti-fog face shield black anti-glare strip



3M Health Care
3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957
www.3M.com/healthcare

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

1 800 364-3577

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
70-2009-4534-6