

Generator+LED

AirLight®

M800 Surgical 45

It provides excellent accessibility, and visibility to make stress-free treatment a reality. M800 Surgical 45 features 22W of cutting power for faster preparation and rear exhaust to help increase patient safety. With its uniquely angled head, enhanced access is simple. Experience for yourself, the difference between ordinary and extraordinary.




Accessibility

45° Head design for perfect access to the molars.



Visibility

Stream water (not mist) for perfect bur cooling and visibility.



Hygiene

Rear air exhaust to prevent contamination into cutting area.



Power

Strong chuck holding force to prevent slippage and improve power with 22W torque.

Technical Data

Colour temperature(K)	5,500
Brightness(LUX)	25,500
Head Size(Width x Height)	10.8mm x 13.0mm
Drive Air Pressure(PSI)	40-60 PSI
Power(WATT)	22
Speed(RPM)	410,000
Rotary Instruments	FG burs according to EN ISO 1797-1:1995
Ø / max. length(MM)	1.6 / 25
Chuck Type	Push Button
Water Port	Triple Water Jet

Generator+LED



Order Info.

Midwest 4 Hole Connection



Code **HP3048**
Model **M800-45/M4**

Borden 2 Hole Connection



Code **HP3052**
Model **M800-45/B2**

Beyes® PD Connection



Code **HP6038**
Model **M800-45/PD**

Kavo® Connection



Code **HP3037**
Model **M800-45/K**

NSK® Connection



Code **HP3044**
Model **M800-45/N**

GLED

Beyes®

The Next Standard®

Beyes Dental Canada Inc.

30 Casebridge Court
Toronto, ON., M1B 3M5
Canada

Tel: 1-416-281-1888

Email: info@beyes.ca

Website: www.beyes.ca

Your Dealer for Beyes Products



Printed in Canada.
B111EN.REV0X1707

© Beyes Canada Inc 2017. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.