

# EZ-Swivel® Couplers

Customer satisfaction is the main priority under the Beyes quality policy. EZ-Swivel Couplers have been developed, manufactured and subjected to final inspection according to legal regulations, quality norms and industry standards.

## KEY SPECIFICATION

<b>Hose Type</b>	ISO 9168 type 1, ISO 9168 type 3
<b>Recommended Voltage Range</b>	3.2 ± 0.1 V AC/DC
<b>Temperature for Storage</b>	-40 °C to +70 °C (-40 °F to +158 °F)
<b>Air humidity for storage</b>	8 % to 80 % (relative), non-condensing
<b>Ambient temperature during operation</b>	+10 °C to +35 °C (+50 °F to +86 °F)
<b>Air humidity during operation</b>	15 % to 80 % (relative), non-condensing

## ORDER INFORMATION

### QD Couplers



#### 4 Holes

**Code:** HP3018

**Model:** EZ-Swivel QD M4

#### 2 Holes

**Code:** HP3019

**Model:** EZ-Swivel QD B2

### PD Couplers



#### 4 Holes

**Code:** HP6008

**Model:** EZ-Swivel PD M4

#### 2 Holes

**Code:** HP6002

**Model:** EZ-Swivel PD B2

### K Couplers



#### 4 Holes

**Code:** HP2063

**Model:** EZ-Swivel K M4

#### 2 Holes

**Code:** HP2064

**Model:** EZ-Swivel K E6

**Beyes**®  
The Next Standard®

**Beyes Dental Canada Inc.**  
23-595 Middlefield Road  
Toronto, Ontario, M1B 3M5  
Canada  
**Tel:** 1-855-603-1888  
**Email:** info@beyes.ca  
**Web:** www.beyes.ca



**Lotus NL B.V.**  
Koningin Julianaplein 10, le Verd,  
2595AA, The Hague,  
Netherlands  
**Tel:** +31645171879 (English)  
+31626669008 (Dutch)  
**Web:** www.beyes.ca

## Your Dealer for Beyes Products

Scan the QR code for more information



Printed in Canada  
File #: EN8032  
Rev.1 / 06.11.20

© Beyes Canada Inc 2020. All Rights Reserved.  
Specifications are subject to change without notice.  
All brands are holders of their respective trademarks.  
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.  
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDCAS).