

CareMaster®

Automatic Handpiece Maintenance System

- Maximize the life span of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces for safety
- Care for any brand of handpieces
- Disposable mist filtration system
- Detachable door for easy cleaning
- Easy replenishment of solution tank



Handpieces are one of the significant investments in the dental office. It is crucial to maintain their effectiveness, so the practitioner will have the cutting efficiency and durability that they need. CareMaster is a solution for the dynamic clinical needs.



Individual controlled

S: Short Cycle for HS* (Approx. 35s)
L: Long Cycle for LS* (Approx. 45s)
E: Extra Long Cycle for SR* (Approx. 50s)



Air Regulator

Built-in air regulator to ensure proper maintenance.



No mess

Spray mist is generally absorbed within the disposable filter sponges that can be changed.

*HS: High Speed Handpiece, LS: Low Speed Handpiece, SR: Speed Reduction Handpiece

Key Specifications

Rated Voltage	AC100V-240V, 50/60Hz
Input Power	35W
Working Pressure	0.35-0.40Mpa (3.5-4.0KGf/cm2)(50-60psi)
Capacity	350ml
Dimensions	282x300x370mm
Kit Content	CareMaster Unit, 1000ml CareMaster Solution Plus, Mist Filter Set (pack of 12 pieces), Oil Absorber Sheets (Pack of 20), Oil Replenish Funnel with Filter, Power Cord, Air Tube

Order Information



2 Turbine Handpieces
1 E-Type Attachment
CODE: HPL2010



1 Turbine handpiece
2 E-Type Attachments
CODE: HPL2011



BeyesLube®
Aerosole Lubricant
CODE: LC2010

Beyes®
The Next Standard

Beyes Dental Canada Inc.
23-595 Middlefield Road
Toronto, Ontario, M1V 3S2
Canada
Tel: 1-416-281-1888
Email: info@beyes.ca
Website: www.beyes.ca

Your Dealer for Beyes Products

Scan the QR code for more information



Printed in Canada
Document #: ENB022
Rev.1 /01.10.20

© Beyes Canada Inc. 2020. All Rights Reserved.
Specifications are subject to change without notice.
All brands are holders of their respective trademarks.
All Beyes products are manufactured in strict compliance with Health Canada regulations ensuring quality.
Manufacturing facility is ISO 13485 certified under Canadian Medical Devices Conformity Assessment System (CMDACS).