

Use PowerScrub™
with your NXT Hg5
and we'll extend
your warranty to
THREE years!

Solmetex™ **POWERSCRUB™** VACUUM LINE CLEANER

Vacuum Line Cleaning Solutions pH and Hg5® Performance

The EPA “Dental Rule” on the proper handling of amalgam waste along with the adoption of the ADA’s Best Management Practices (BMPs) went into effect nationwide on July 14, 2017. The rule states that only non-chlorine or non-oxidizing disinfectants and neutral vacuum line cleaners with a **pH >6 and <8** are to be used.

Solmetex exclusively recommends PowerScrub Vacuum Line Cleaning Solution for use with the Hg5 Series of Amalgam Separators



Solmetex™ **POWERSCRUB™** VACUUM LINE CLEANER

- SmartScrub Technology — Continuous cleaning action
- Microbial formulation optimizes collection container capacity
- Cleans vacuum line debris on contact
- Safe for the clinician, the patient and the environment — No harsh chemicals
- Releases microbial cleaning enzymes at any water temperature
- Meets all EPA “Dental Rule” requirements for pH >6 and <8
- Compatible with both dry and wet ring vacuums
- Compatible with all amalgam separators