



## Nexus™ RMGI

*Resin-Modified Glass Ionomer Cement*



### One-Peel Cleanup

Enjoy quicker and easier cleanups with Nexus RMGI. Advanced Nexus technology gives this cement an optimal gel state and tack-cure capability for a more predictable clinical experience.



### Self-Adhesive with Exceptional Bond Strength

Powered by Kerr's innovative dual-adhesive technology, Nexus RMGI delivers outstanding bond strength to tooth structures and common dental substrates, including zirconia-based ceramics. Nexus RMGI's bond strength surpasses that of other RMGI products on the market.



### Smart Response Ion Technology Helps Prevent Secondary Caries

Nexus RMGI offers long-term sustained fluoride release, as well as calcium ion release and zinc ion release in response to acidic conditions, to help protect teeth from secondary caries.




### Exceptional Esthetics

No other RMGI product comes close to delivering the esthetic results that Nexus RMGI does. A patented redox initiator system enables Nexus RMGI to deliver unrivaled translucency and color stability for the most esthetically pleasing outcomes.

# Kerr™

# Formulated for Stronger Bonds


## Common Substrates




\*Data on file.

■ Nexus™ RMGI ■ RelyX™ Luting Plus ■ GC FujiCem® Evolve

## Dentin/Enamel




## Zn<sup>2+</sup> Cumulative Release



Nexus RMGI offers significantly higher zinc ion release than 3M RelyX Luting Plus and GC FujiCem 2.

\*Data on file.

## Calcium Release of Nexus RMGI



Nexus RMGI releases higher levels of calcium during episodes of elevated acidity in the oral environment.

**Easy  
One-Peel™  
Cleanup**



**Scan with  
Your Camera**  
to Learn More

## Ordering Information

35640 Nexus RMGI Kit Contains: 3 dual syringes (5g): translucent white, 24 1:1 automix tips (regular), Directions for Use, Technique Guide

33361 Automix Regular Tips 1:1 (pkg of 50)

35650 Automix Wide Tips 1:1 (pkg of 8)

33880 Curved Intraoral Tips (pkg of 8)

33881 Curved Root-Canal Tips (pkg of 8)

Nexus is a trademark of Kerr Corporation. All other trademarks are properties of their respective owners.