



**Works Better.
Lasts Longer.
Costs Less.**



EVACUATION



“Works Better”

Made from a biodegradable, non-chlorine, non-corrosive formulation that contains fast acting protease enzymes. The solution will greatly assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents.

“Lasts Longer”

These easy-to-use enzymatic tablets dissolve in a gallon of water in under a minute to create a powerful and effective ultrasonic cleaning solution. Can also be used as a highly effective evacuation system cleaner, these tablets have a neutral working PH making them safe for use with amalgam separators.

“Costs Less”

The powerful formula has many uses and is less expensive than the competitors in this category.

Ultrasonic Enzymatic Tablets

These easy-to-use enzymatic tablets dissolve in a gallon of water in under a minute to create a powerful and effective ultrasonic cleaning solution. Made from a biodegradable, non-chlorine, non-corrosive formulation that contains fast acting protease enzymes. The solution will greatly assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. Virtually eliminates handling and storage of heavy liquid concentrates. Contains a pleasant fragrance and is a neutral pH. May also be used as an instrument pre-soak or as an evacuation system cleaner.

Packaging:	Item No:	Compares to*:
64 per box <i>makes 32 gallons</i>	#UT-1000	Efferzyme® Effervescent Cleaning Tablets, Cleanlets™ Ultrasonic Cleaning Tablets & Ultrasonic Enzymatic Tablets – Henry Schein®

*Legal disclaimers: Efferzyme® Brush is a registered trademark of Crosstex. Cleanlets™ is a registered trademark of Sultan Healthcare. Use of these trademarks does not imply any affiliation with or endorsement by these companies. Mydent International is not responsible for any inaccuracy or typographical errors in this document.

Request FREE Samples at: www.defend.com/freesamples



80 Suffolk Court, Hauppauge, NY 11788
800-275-0020 • FAX 631-434-7750
www.defend.com • EMAIL: sales@defend.com



Follow us on Facebook.com/MydentDefend
Follow us on Twitter.com/Mydent_Defend