

# Qnala Dental Products

## plastics

### Qnala Set-Up Tray

Flat, B-size, (13 $\frac{3}{8}$ " x 9 $\frac{5}{8}$ " x  $\frac{7}{8}$ ") trays are available in four colors: Blue, White, Beige and Green.

#### FEATURES

Wide edges  
Non-Slip feet  
High gloss surface  
Clear Plastic Tray Cover  
(sold separately)

#### BENEFITS

Easy to grip  
Provides stability  
For stain resistance and easy cleaning  
Helps prevent cross-contamination



### Qnala Denture & Retainer Boxes

Denture boxes are available in a 12-pack of assorted colors: Blue, Pale Yellow, Mauve, Beige, Teal, Purple (two each).

Retainer boxes are 1" deep and can be used for retainers, mouth guards, whitening trays, etc. Available in a 12-pack of assorted colors: Teal, Neon Pink, Neon Purple, Neon Blue, Neon Yellow, Neon Green (two each).

#### FEATURES

One-Piece plastic construction  
Contoured corners  
Hinged lid with latch closure  
Imprintable

#### BENEFITS

Strong; durable  
For easy cleaning  
Lid stays securely closed  
Customized for the dental office



### Qnala Disposable Screen Traps

Available in seven different configurations. 25/bag.

#### FEATURES

For use with dental evacuation systems  
Color-coded  
Disposable plastic

#### BENEFITS

Improved performance  
Easy identification  
Reduces possibilities of infection and contamination



### Qnala Disposable Cuspidor Screens

Available in seven different configurations. 25/bag.

### Qnala Saliva Ejector Screens

Universal size. 25/bag.



# Qnala Dental Products

## plastics

CODE	DESCRIPTION
	<b>B Size Flat Setup Tray:</b> Tray size: 13 $\frac{3}{8}$ " x 9 $\frac{5}{8}$ " x $\frac{7}{8}$ "
Q9007	Blue
Q9008	White
Q9009	Beige
Q9010	Green
Q9011	<b>Flat Tray Covers:</b> Non-Locking, Clear, Cover size: 14" x 10 $\frac{5}{8}$ " x $\frac{5}{16}$ "
Q9015	<b>Retainer Box:</b> 1" deep, Assorted Colors <i>*Imprinting Available</i>
Q9020	<b>Denture Box:</b> Assorted Colors <i>*Imprinting Available</i>
	<b>Disposable Screen Traps:</b>
Q9030	Blue, 2 $\frac{1}{8}$ " diameter <i>A-dec w/plastic canisters; Belmont, DCI, Forest, Marus &amp; others</i>
Q9031	Mauve, 1 $\frac{1}{4}$ " diameter <i>A-dec Cascade and Performer</i>
Q9032	Green, 1 $\frac{1}{8}$ " diameter <i>Belmont, Chayes, Den-Tal-Ez, Dentech, Forest, Pelton &amp; Crane, Proma and others</i>
Q9033	Yellow, 2 $\frac{1}{4}$ " diameter <i>A-dec model w/metal canisters</i>
Q9034	Coral, 2 $\frac{5}{8}$ " diameter <i>Pelton &amp; Crane Spirit w/plastic canisters</i>
Q9035	Beige, 2 $\frac{7}{8}$ " diameter <i>Den-tal-ez with black canisters</i>
Q9036	Gray, 2 $\frac{3}{4}$ " diameter <i>Midwest 210</i>
Q9037	<b>Saliva Ejector Screens:</b> Universal
	<b>Cuspidor Screens:</b>
Q9038	Small, White, 1 $\frac{5}{8}$ " <i>Fits A-dec Cascade and Performer</i>
Q9039	Large Universal, Gray, 2 $\frac{1}{2}$ " <i>Fits A-dec, Ampco, Belmont, Dabi Atlante, Dansereau, DCI, Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter</i>



**Qnala**  
DENTAL PRODUCTS  
www.qnala.com

#QSS003R1

**Distributed by:**