

INSTRUCTIONS FOR USE


Treatment Accessories



HVE Tips

Instructions for use:

1. Attach tip to High Volume Evacuator (HVE) hand piece, ensuring that the preferred working end is not inserted into the valve.
2. Dispose properly when procedure is finished.



Ref.#'s ZET, ZETNV, ZES50



Do not reuse.

Not made with natural rubber latex.



Crosstex International, Inc., a Cantel Medical Company
10 Ranick Road, Hauppauge, NY 11788-4209 USA
www.crosstex.com 800.276.7783 631.582.6777



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

