

INSTRUCTIONS FOR USE

Treatment Accessories



MaxVac Plus® Combo Tip®

Instructions for use:

1. Use when needed to remove large matter from oral cavity.
2. Attach tube to suction hand piece.
3. Dispose properly when procedure is finished.



Ref.#'s ZETMVBL, ZETMVGR,
ZETMVPK, ZETMVLV, ZETMVEY, ZETMVP

MaxVax Plus® and Combo Tip® are registered trademarks of Crosstex International, Inc.



Do not reuse.



Not made with natural rubber latex.



CROSSTEX[®]
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in USA

Manufactured by: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

ZETMV-IFU
Rev. 0 – 11/16