

Quala Disposable Screen Traps
Cross Reference Guide

Quala Item #	Pinnacle Item #	Product Name	Color	Diameter	Fits the following:	UOM
Q9030	5501	Disposable Screen Traps	Blue	2-1/8"	A-dec w/plastic canisters; Belmont, DCI, Forest, Marus & others	25/bag
Q9031	5500	Disposable Screen Traps	Mauve	1-1/4"	A-dec Cascade and Performer.	25/bag
Q9032	5503	Disposable Screen Traps	Green	1-7/8"	Belmont, Chayes, Den-Tal-Ez, Dentech, Forest, Knight, Pelton & Crane, Proma and others	25/bag
Q9033	5502	Disposable Screen Traps	Yellow	2-1/4"	A-dec model w/metal canisters.	25/bag
Q9034	5505	Disposable Screen Traps	Coral	2-5/8"	Pelton & Crane Spirit w/plastic canisters.	25/bag
Q9035	5506	Disposable Screen Traps	Beige	2-7/8"	Den-tal-ez with black canisters.	25/bag
Q9036	5507	Disposable Screen Traps	Gray	2-3/4"	Midwest 210.	25/bag
Q9038	6200-C	Cuspidor Screens-Small	White	1-5/8"	A-dec, Cascade and Performer Cuspidor bowls.	25/bx
Q9039	6100-C	Cuspidor Screens-Large Universal	Gray	2-1/2"	A-dec, Ampco, Belmont, Dabi Atlante, Dansereau, DCI, Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter.	25/bx
Q9037	5508-S	Saliva Ejector Screens-Universal	Clear		All Saliva Ejectors	25/bag