

INSTRUCTIONS FOR USE


Treatment Accessories



Comfort Plus® Premium Saliva Ejector

Instructions for use:

1. Attach Ejector head onto the end of the low volume evacuation tube.
2. Place in mouth and bend to desired position.
3. After use, dispose of the ejector.



Ref.#'s ZWWCP, ZCWCP,
ZWPCPB, ZWGCPM

Comfort Plus® is a registered trademark of Crosstex International, Inc.



Do not reuse.



Not made with natural rubber latex.



CROSSTEX
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in USA

Manufactured by: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

ZWCP-IFU
Rev. 0 – 11/16