


INSTRUCTIONS FOR USE

Treatment Accessories



Hygoformic® Saliva Ejectors



Ref.#'s ZFF, ZFFAD

Instructions for use:

1. Bend Saliva Ejector/Tongue Retractor to the desired position.
2. Shape for posterior or anterior use.
3. Unlimited flexibility.
4. Adjust for desired chin length.
5. Ideal for pedodontic procedures.
6. Adjust coil over evacuation slots to prevent suctioning of soft tissue.

Hygoformic® is a registered trademark of J.H. Orsing AB.



Do not reuse.



Not made with natural rubber latex.



CROSSTEX
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in Sweden

Manufactured for: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

ZFF-IFU
Rev. 0 – 11/16