

INSTRUCTIONS FOR USE

Treatment Accessories



Surgical Aspirator Tips


Surgical Aspirator Tips: with small, standard or large openings and color-coded according to size. They can fit 11 mm suction holders.

Indication:

Blood and fluid suction during surgeries.

Instructions for use:

1. Connect the Surgical Aspirator Tip to the high volume evacuation hose.
2. Discard aspirator after use.



Ref.#'s ZSABL, ZSAWH, ZSAGR

Blue–Small Tip (.0625" / .16 cm)

White–Standard Tip (.125" / .32 cm)

Green–Large Tip (.25" / .64 cm)



Do not reuse.



Not made with natural rubber latex.



CROSSTEX[®]
A CANTEL MEDICAL COMPANY



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Made in China

Manufactured for: Crosstex International, Inc.
10 Ranick Road | Hauppauge, NY 11788 | 631.582.6777 | 888.276.7783 | Crosstex@crosstex.com

ZSA-IFU
Rev. 0 – 11/16