

Get the perfect mix of precision, speed and accuracy.

With the touch of a button, you can mix and dispense impression material. Doctors get more usable working time while assistants have less stress and hand fatigue. In the end, your material will be homogeneous and void-free which means the finals will fit better.

The Pentamix 3 mixing unit mixes and dispenses all your 3M impression materials—even putty—faster than any automatic mixing unit on the market. With the working time you gain, you'll even be able to take advantage of quick-setting materials.



Even easier to operate.

The 3M™ Pentamix™ 3 Automatic Mixing Unit operates with push-button from either side. This makes it easy to use whether you are right- or left-handed. Effortless cartridge exchange and easy-to-attach mixing tips also increase procedure efficiency.



Easy to attach mixing tip without opening device cover.



Easy cartridge exchange.



Push-buttons on left and right side allow for easy operation.




Product name and working/setting times are visible, even when the device cover is closed.

Speed isn't its only advantage.

The unit has the smallest footprint on the market.

That helps you make the most of valuable office space. Its sleek, high-tech design keeps your office looking good, and the wall-mount design option lets you free up your counter completely.



Ordering information

Item #	Product Information
77875	3M [™] Pentamix [™] 3 Automatic Mixing Unit Universal Cartridge sold separately
71545	Universal Cartridge for 3M™ Pentamix™ 3
77601	3M™ Pentamix™ 3 Wall Mount Kit
77949	3M [™] Penta [™] Mixing Tips — Red, Refill Pack
P3792	3M [™] Impregum [™] Penta [™]
P3788	3M™ Impregum™ Penta™ Soft MB
P3794	3M™ Impregum™ Penta™ Soft HB
P3796	3M™ Impregum™ Penta™ Soft LB
P3776	3M [™] Impregum [™] Penta [™] Soft QS MB
P3774	3M [™] Impregum [™] Penta [™] Soft QS HB
77784	3M [™] Imprint [™] 4 Penta [™] HB
77785	3M™ Imprint™ 4 Penta™ SQ HB
77786	3M [™] Imprint [™] 4 Penta [™] Putty
71527	3M [™] Imprint [™] 4 Preliminary Penta [™]
71528	3M [™] Imprint [™] 4 Preliminary Penta [™] SQ
71615	3M [™] Directed Flow Impression Tray Refills, S Upper (10)
71616	3M [™] Directed Flow Impression Tray Refills, M Upper (10)
71617	3M™ Directed Flow Impression Tray Refills, L Upper (10)
71618	3M [™] Directed Flow Impression Tray Refills, S Lower (10)
71619	3M [™] Directed Flow Impression Tray Refills, M Lower (10)
71620	3M™ Directed Flow Impression Tray Refills, L Lower (10)
77780	Plunger Discs
71225	12 Syringe Tips
71200	10 3M [™] Elastomer Syringe Barrels
77720	Syringe Knurled Hub Cap
71210	3M [™] Elastomer Syringe Set Kit includes: 1 syringe; 4 tips; 4 barrels; 1 brush

3M

3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental 3M Canada Post Office Box 5757 London, Ontario N6A 4T1 Canada

Phone 1-888-363-3685

Customer Care Center: 1-800-634-2249 3M.com/greatimpressions

3M, Imprint, Impregum, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved. 70-2013-0069-9 (Rev. 8/2016)