


Paradigm™

VPS Impression Materials

NOW
AVAILABLE
IN PENTA™
DELIVERY



Benefits:

- Easy extrusion and easy to use
- Hydrophilic
- Great flow, yet does not drip or slump
- High tear strength

The exceptional value in dentistry.

Paradigm™ VPS Impression Materials are designed to capture detail and margins well, so you can count on consistently accurate impressions. They flow smoothly into the tray or mouth, yet do not drip or slump, making them easy to use.

They're hydrophilic, which helps produce void-free impressions — even with moisture present, like saliva or blood.

High tear strength and recovery from deformation help protect the impression from damage and distortion as it's removed from the mouth.

Their versatility and price make it great for all indications, including:

- Crowns and bridges
- Bleaching trays
- Study models
- Denture impressions
- Provisional impressions


Great value
from a company you trust.



Bring the control and precision of automixing to your practice.


3M ESPE's NEW affordable and portable mixer — Pentamix™ Lite Automatic Mixing Unit.

- Economical and exact dispensing of the amount of material needed
- Fast, time-saving procedure and push-button activation
- Homogeneous and void-free mixing results of reproducible quality
- Easily carry from operator to operator with the lightest device on the market



COMPACT DESIGN

Saves valuable space in the upright standby position.



PORTABLE DESIGN

Easy to carry from operator to operator thanks to light weight and integrated handgrip.

Always use 3M™ ESPE™ Penta™ Mixing Tips – Red in your Pentamix device to ensure proper mixing and impression results.

Ordering Information

PARADIGM™ VPS IMPRESSION MATERIALS

Item No. Product Information

Refills — Tray Material

36716	Paradigm™ Penta™ VPS Heavy Body 2 – 300ml base pastes; 2 – 60ml catalysts; Step-by-step card; Instructions for use
36717	Paradigm™ Penta™ VPS Heavy Body Fast Set 2 – 300ml base pastes; 2 – 60ml catalysts; Step-by-step card; Instructions for use
5313	Paradigm™ VPS Heavy Body Refill 2 – 50ml cartridges
5513	Paradigm™ VPS Heavy Body Fast Set Refill 2 – 50ml cartridges

Refills — Wash Material

Each Wash Material refill includes: 2 – 50ml cartridges

5312	Paradigm™ VPS Monophase Refill
5315	Paradigm™ VPS Regular Body Refill
5515	Paradigm™ VPS Regular Body Fast Set Refill
5311	Paradigm™ VPS Light Body Refill
5511	Paradigm™ VPS Light Body Fast Set Refill

PARADIGM™ VPS IMPRESSION MATERIALS

Item No. Product Information

Jumbo Pack Refills — Tray Material

Each Jumbo Pack Tray Material refill includes: 30 – 50ml cartridges

5313J	Paradigm™ VPS Heavy Body Jumbo Refill
5513J	Paradigm™ VPS Heavy Body Fast Set Jumbo Refill

PENTAMIX™ LITE AUTOMATIC MIXING UNIT AND ACCESSORIES

Item No. Product Information

77907	Pentamix™ Lite Automatic Mixing Unit – 120 V/60 Hz Includes: 1 Penta™ Cartridge for Pentamix™ Lite
77944	1 Penta™ Cartridge for Pentamix™ Lite
77949	Penta™ Mixing Tips – Red Includes: 50 Mixing Tips

Customer Care Center: 1-800-634-2249 www.3MESPE.com/GreatImpressions



3M ESPE Dental

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

3M, ESPE, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2015. All rights reserved.
Please recycle. Printed in USA.
70-2013-0469-1 (Rev. B)