Paradigm™ VPS Impression Materials Overview Sheet

		Viscosity			Indications	
Product	Color	low high	Working Time (Including Mixing)	Setting Time (From Start of Mix)	Crown and Bridge Impressions	Alginate Replacement
Paradigm VPS Heavy Body			2:00	4:00	•	
Paradigm VPS Heavy Body Fast Set			1:15	2:30	•	
Paradigm VPS Monophase	•		2:00	4:00	•	•
Paradigm VPS Regular Body			1:00	4:00		
Paradigm VPS Regular Body Fast Set			0:40	2:30	•	
Paradigm VPS Light Body			1:00	4:00		
Paradigm VPS Light Body Fast Set			0:40	2:30	•	



Paradigm™

VPS Impression Material

The NEW value in dentistry.

- Easy extrusion and **easy-to-use**
- **Great flow,** yet does not drip or slump
- **■** Hydrophilic
- Good tear strength, and recovery from deformation
- Great for all indications

Great value from a company you trust.



