

Quala Dental Products

impression materials



Vinyl Polysiloxane (VPS)

Super hydrophilic formulation offers excellent dimensional stability and tear strength for crisply detailed impressions every time. No odor and no irritation provide greater patient comfort. Tips sold separately.

FEATURES

State-of-the-art hydrophilic technology

Uniform consistency

High dimensional stability

Pleasantly neutral taste and odor

New cartridge delivery system

BENEFITS

Provides great "wetability"; it won't "run away" from prepared impression area

Assures crisply detailed impressions

Maintains the accuracy of the impression over time, even after disinfection by immersion

Increases patient acceptance, helping to reduce gag reflexes during impressions

Easier dispensing of material

Contains: 4 x 50 ml cartridges. Tips sold separately.

	WORK TIME	INTRAORAL SET TIME
Regular Set	2:30	1:00
Fast Set	1:00	1:00

CODE	DESCRIPTION
Light Body (Yellow)	
32821	Vanilla Crème Fast Set
32822	Unflavored Fast Set
32823	Unflavored Regular Set
Medium Body/Monophase (Purple)	
32824	Vanilla Crème Fast Set
32825	Unflavored Fast Set
32826	Unflavored Regular Set
Heavy Body (Rose)	
32827	Vanilla Crème Fast Set
32828	Unflavored Fast Set
32829	Unflavored Regular Set