

**Medicom® SafeWear™ Form-Fit Isolation Gown™**
**Intended Use**

Non-Sterile Procedure Isolation Gown with knitted cuffs and collar and neck tie, designed to protect health professionals and patients from cross infection. Single use and disposable isolation gowns made of polypropylene. Latex-Free.

**Country of Origin**

China

**Packaging**

12 gowns / bag. 60 gowns per case

**Storage**



Store under ambient temperature and dry conditions

**Composition**

100% Spunbond Polypropylene (PP). Latex Free

**Basis Weight**

25 g/m<sup>2</sup>

Description & Color	Color Reference	Size & Item Numbers	
SafeWear™ Form-Fit Isolation Gown™ Sunny Yellow		Regular	XL
		8113	8113-D
SafeWear™ Form-Fit Isolation Gown™ Bright Blue		Regular	XL
		8114	8114-D

Dimensions		
Description	Regular (cm)	Extra Large Size (cm)
Underarm Width (A)	71	78
Gown Length (B)	110	116
Sleeve Length (C)	58	61
Arm Height (D)	28	32
Open Gown Width (E)	142	156
Cuff (F)	5	5
Neckline (G)	11	11
Belt Length	183	183

