



TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116

Intended Use

Procedure SMS Lab Jacket designed to protect health professionals and patients from cross contamination. Single use and disposable gowns made of polypropylene.

Characteristic	Specification		
Description	Jacket made of non-woven material with front closure with snap buttons. Features knitted collar and cuffs, three (3) pockets and a hanging loop inside the jacket at the back of the collar.		
Material Composition	Spunbond Meltblown Spunbond (SMS) Polypropylene. Latex-Free		
Packaging	Jackets/Bag	Bags/Case	Jackets/Case
	10	5	50

Color and Size Tables

Color	Size & Item Code					
	Small	Medium	Large	X-Large	2X-Large	3X-Large
Deep Blue	8101-A	8101-B	8101-C	8101-D	8101-E	8101-F
Pretty Pink	8102-A	8102-B	8102-C	8102-D	N/A	N/A
Purple	8103-A	8103-B	8103-C	8103-D	N/A	N/A
Soft Blue	8104-A	8104-B	8104-C	8104-D	N/A	N/A
White	8105-A	8105-B	8105-C	8105-D	N/A	N/A
Teal	8115-A	8115-B	8115-C	8115-D	N/A	N/A
Poppy Pink	8116-A	8116-B	8116-C	8116-D	N/A	N/A

TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116


Dimension Table

Dimensions/Size (cm)						
Description	Small	Medium	Large	X-Large	2X-Large	3X-Large
Gown Length (A)	71	75	82	86	89	91
Underarm Width (B)	55	60	65	70	78	86
Sleeve Length (C)	54	54	57	58	59	62
Arm Height (D)	28	29	30	33	33	34
Shoulder Width (E)	17	17	21	22	26	27
Knitted Cuff Length (F)	6	6	6	6	6	6
Sleeve Opening Width (Stretched) (L)	15	15	15	15	15	15
Knitted Cuff Width (M)	6	6	6	6	6	6
Button Spacing (G)	11,5	11,5	11,5	11,5	11,5	11,5
Lower Front Pocket	Height (I)	18	18	18	18	18
	Width (H)	18	18	18	18	18
Upper Pocket	Height (K)	14	14	14	14	14
	Width (J)	14	14	14	14	14
Upper Pocket Distance From Opening Edge (N)	8	9	10	11	13	15
Lower Pocket Distance From Opening Edge (P)	5	5	7	9	10	12
Distance between Upper and Lower Pocket (S)	12	14	16	20	23	23
Lower Pocket Distance From Bottom Edge (O)	5	7	13	13	13	13
Number of Buttons	5	5	5	5	5	5

TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116

Schematic representation of gown measurements



Item Code	Description	Color	Specific Attributes
8101-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8101-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8101-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8101-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8101-E	Medicom® SafeWear™ Hipster Jacket™, 2X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8101-F	Medicom® SafeWear™ Hipster Jacket™, 3X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Deep Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile



TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116

Item Code	Description	Color	Specific Attributes
8102-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8102-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8102-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8102-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8103-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8103-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8103-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8103-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8104-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8104-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8104-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8104-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile



TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116

Item Code	Description	Color	Specific Attributes
8105-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	White	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8105-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	White	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8105-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	White	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8105-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	White	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8115-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8115-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8115-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8115-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8116-A	Medicom® SafeWear™ Hipster Jacket™, Small, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Poppy Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8116-B	Medicom® SafeWear™ Hipster Jacket™, Medium, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Poppy Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8116-C	Medicom® SafeWear™ Hipster Jacket™, Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Poppy Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile
8116-D	Medicom® SafeWear™ Hipster Jacket™, X-Large, 12 Jackets/Bag, 5 Bags/Case, 60 Jackets/Case	Poppy Pink	Knitted Cuff, Knitted Collar, Three Pockets, Non-Sterile



TECHNICAL DATA SHEET

Medicom® SafeWear™ Hipster Jacket™ Items # 8101 to 8105 and, 8115 to 8116

Regulatory

FDA Class 1, FDA Product Code OEA
Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

China