



TECHNICAL DATA SHEET

Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Intended Use

Procedure SMS Lab Coat designed to protect health professionals and patients from cross contamination. Single use and disposable gowns made of polypropylene.

Characteristic	Specification		
Description	Coat made of non-woven material with front closure and snap buttons. Features knitted collar and cuffs, three pockets, 2 slits in the back and a hanging loop inside the coat at the back of the collar.		
Material Composition	Spunbond Meltblown Spunbond (SMS) Polypropylene. Latex-Free		
Packaging	Coats/Bag	Bags/Case	Coats/Case
	10	5	50

Color and Size Tables

Color	Size & Item Code					
	Small	Medium	Large	X-Large	2X-Large	3X-Large
Deep Blue	8107-A	8107-B	8107-C	8107-D	8107-E	8107-F
Pretty Pink	8108-A	8108-B	8108-C	8108-D	N/A	N/A
Plum Purple	8109-A	8109-B	8109-C	8109-D	N/A	N/A
White Frost	8110-A	8110-B	8110-C	8110-D	N/A	N/A
Soft Blue	8112-A	8112-B	8112-C	8112-D	N/A	N/A
Tropical Teal	8117-A	8117-B	8117-C	8117-D	N/A	N/A

TECHNICAL DATA SHEET

Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Dimension Table


Dimensions/Size (cm)						
Description	Small	Medium	Large	X-Large	2X-Large	3X-Large
Gown Length (A)	100	112	114	120	116	120
Underarm Width (B)	55	60	67	70	76	90
Sleeve Length (C)	52	54	57	60	60	64
Arm Height (D)	28	29	30	31	33	36
Shoulder Width (E)	18	18	20	23	26	29
Knitted Cuff Length (F)	6	6	6	6	6	6
Sleeve Opening Width (Stretched) (L)	15	15	15	15	15	15
Knitted Cuff Width (M)	6	6	6	6	6	6
Button Spacing (G)	15.5	15.5	15.5	15.5	15.5	15.5
Lower Front Pocket	Height (I)	18	18	18	18	18
	Width (H)	18	18	18	18	18
Upper Pocket	Height (K)	14	14	14	14	14
	Width (J)	14	14	14	14	14
Upper Pocket Distance From Opening Edge (N)	6	9	9	9	10	11
Lower Pocket Distance From Opening Edge (P)	5	6	6	8	8	10
Distance between Upper and Lower Pocket (S)	27	30	30	31	31	31
Lower Pocket Distance From Bottom Edge (O)	18	26	27	32	29	33
Back Slit Height (R)	30	30	30	30	30	30
Back Panel Width (Q)	43	48	55	58	64	78
Number of Buttons	5	5	5	5	5	5

TECHNICAL DATA SHEET


Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Schematic representation of gown measurements

Front View



Back View



TECHNICAL DATA SHEET

Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Item Code	Description	Color	Specific Attributes
8107-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8107-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8107-C	Medicom® SafeWear™ High Performance Lab Coat™, Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8107-D	Medicom® SafeWear™ High Performance Lab Coat™, X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8107-E	Medicom® SafeWear™ High Performance Lab Coat™, 2X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8107-F	Medicom® SafeWear™ High Performance Lab Coat™, 3X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Deep Bleu	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8108-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8108-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8108-C	Medicom® SafeWear™ High Performance Lab Coat™, Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8108-D	Medicom® SafeWear™ High Performance Lab Coat™, X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Pretty Pink	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile

TECHNICAL DATA SHEET

Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Item Code	Description	Color	Specific Attributes
8109-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8109-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8109-C	Medicom® SafeWear™ High Performance Lab Coat™, Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8109-D	Medicom® SafeWear™ High Performance Lab Coat™, X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Plum Purple	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8110-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	White Frost	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8110-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	White Frost	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8110-C	Medicom® SafeWear™ High Performance Lab Coat™, Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	White Frost	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8110-D	Medicom® SafeWear™ High Performance Lab Coat™, X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	White Frost	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8112-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8112-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile



TECHNICAL DATA SHEET

Medicom® SafeWear™ High Performance Lab Coat™ Items # 8107 to 8110, 8112 and 8117

Item Code	Description	Color	Specific Attributes
8112-C	Medicom® SafeWear™ High Performance Lab Coat™, Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8112-D	Medicom® SafeWear™ High Performance Lab Coat™, X-Large, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Soft Blue	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8117-A	Medicom® SafeWear™ High Performance Lab Coat™, Small, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile
8117-B	Medicom® SafeWear™ High Performance Lab Coat™, Medium, 12 Coats/Bag, 5 Bags/Case, 60 Coats/Case	Tropical Teal	Knitted Cuff, Knitted Collar, Three Pockets, Two Back Slits, Non-Sterile

Regulatory

FDA Class 1, FDA Product Code OEA
Health Canada Class 1

Storage Conditions

Store under ambient temperature and dry conditions

Country of Origin

China