

3M Science.
Applied to Life.™

Choose the Right Protection For You

3M™ Face Masks



Dental professionals should wear appropriate personal protective equipment—including masks, which adhere to Standard Precautions for protection of the mouth and nose—when splashes or sprays of blood and body fluids are likely.^{1,2}

3M dental masks designed for procedures where there is potential exposure to blood and body fluids:

Triple Layer Molded Face Mask

Catalog Number:
2042F and
2042FP (Petite)



- Fluid resistant³
- Helps reduce exposure to blood and body fluids³
- BFE⁴ ≥ 98%
- EN 14683: Type IIR⁵
- Middle layer is designed to attract and trap airborne particles
- Malleable nose piece

Tie-on Face Mask

Catalog Number: 1818



- Flat, pleated
- Helps reduce exposure to blood and body fluids⁶
- BFE⁴ ≥ 98%
- EN 14683: Type II⁵
- Soft, lightweight, comfortable to wear
- Malleable nose piece

Ear-loop Face Mask

Catalog Number: 2000F



- Flat, pleated
- Helps reduce exposure to blood and body fluids⁶
- BFE⁴ ≥ 98%
- EN 14683: Type II⁵
- Lightweight, soft, easy to put on and remove
- Malleable nose piece

Utility masks:

Molded Face Mask

Catalog Number: 1942FB (Blue) and 1942FR (Rose)



- Designed for procedures that do not produce fluid, spray or aerosols
- Malleable nose piece

Ordering Information



Item # Product Information

Flat Face Masks

1818 **3M™ Tie-On Face Masks**
50 Masks per Box

2000F **3M™ Ear Loop Face Mask**
50 Masks per Box



Molded Face Masks

1942FB **3M™ Molded Face Masks, Blue**
50 Masks per Box

1942FR **3M™ Molded Face Masks, Rose**
50 Masks per Box

2042F **3M™ Triple Layer Molded Face Masks**
50 Masks per Box

2042FP **3M™ Triple Layer Molded Face Masks, Petite**
50 Masks per Box

3M™ Face Masks do not contain components made from natural rubber latex.

Also available from 3M

3M is a leader in providing quality infection prevention solutions that help protect you and your patients against the spread of infectious diseases. We share your commitment to doing what's right for your patients—and for you. We offer a family of reliable infection prevention products to meet the needs of dental professionals around the world.



3M™ SecureFit™ Protective Eyewear



3M™ Comply™
Lead Free Autoclave Indicator Tape



3M™ Avagard™ D
Instant Hand
Antiseptic With
Moisturizers



3M™ Attest™
Biological Monitoring
System

¹Centers for Disease Control and Prevention. Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care. Atlanta, GA: US Department of Health and Human Services, Centers for Disease Control and Prevention, National Center for Chronic Disease Prevention and Health Promotion, Division of Oral Health; March 2016.

²Centers for Disease Control and Prevention. Guidelines for infection control in dental health-care settings—2003. MMWR Recomm Rep 2003;52(RR-17):1–61. Available at: www.cdc.gov/mmwr/PDF/rr/rr5217.pdf.

³Fluid resistance determined using test method ASTM F1862:2013

⁴BFE = Bacterial Filtration Efficiency using test method EN 14683:2005

⁵Meets EN 14683:2005 requirements

⁶Hydrostatic pressure test method

www.3M.com



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

3M, Attest, Avagard, Comply, ESPE and SecureFit are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2017. All rights reserved.
70-2013-0889-0