



## *Product Guide*



---

*Product Solutions You Trust*

***Quality goes beyond consistently good products.  
True quality permeates an organization. It breeds excellence.***

Welcome. Innovative Healthcare Corporation is committed to taking care of you. Our people and our products are the best in the world, and we prove it to our customers day in and day out, month after month, and year after year.

We have one of the highest customer retention rates in our industry. Why? Because we genuinely believe in, and practice our **One Corporate Policy**:

***“Do what the customer wants done”***

We've learned over the years that the partnerships we bond are far more valuable than the dollar value of the invoices. In fact, our customers return value to us every day by recommending us to their associates. A fact which naturally makes us very proud.

We don't have a "Customer Care" Department, because Customer Care isn't a department. Customer Care is what comes naturally from a company whose employees are committed to building successful partnerships.

***At Innovative Healthcare, we're committed to the kind  
of service that creates successful, long-term partnerships.***



We'll provide the best we have to offer in every contact you have with us. Our people are professional and knowledgeable. Our glove line is one of the broadest and most comprehensive in the industry. All of our disposable products are top quality and priced to be competitive.

We're here to help you find the most appropriate product for the job.



**Intertek**  
ISO 9001:2008 Certified



**Intertek**  
ISO 13485:2003 Certified

Innovative Healthcare has been a leading supplier of high quality gloves to the medical, dental and industrial markets in North America since 1991. We've developed a wide range of products designed to maximize user performance, efficiency and safety.

Our products include some of the most trusted and widely distributed gloves and disposables in healthcare. Names like DermAssist®, NitriDerm®, DermaGold®, ProDerm™ and Pulse®. Products that, through their consistent growth, continue to prove that quality does make a difference.

And we continue to develop new products that satisfy the ever-evolving requirements of our industry. Products like our new Prestige® Polyisoprene surgeon's gloves utilizing the latest in synthetic polymer compounding technology for incredibly soft, ultra low modulus comfort.

No matter what segment of the healthcare industry you serve, we have a glove product with characteristics well suited to it's specific tasks or requirements. Acute care, primary care, ambulatory care, extended care, dentistry, laboratory, corrections, EMS and many others. Each has a unique need, and you'll find many of those product solutions on the pages that follow.




But while we're proud of our products, we're even more proud of the confidence our customers place in us every day. We've built our company around our customers and their needs. Nothing is more gratifying than the knowledge that we're exceeding their expectations.



#### PRODUCT SOLUTIONS YOU TRUST

Visit us: [www.ihcsolutions.com](http://www.ihcsolutions.com) Contact us: [info@ihcsolutions.com](mailto:info@ihcsolutions.com)


Tel 800.272.1533 Fax 800.650.1014



■ Latex Exam Gloves - Powder-Free.....	4
■ Latex Exam Gloves - Powdered.....	5
Vinyl Synthetic Exam Gloves - Powdered.....	5
Vinyl Synthetic Exam Gloves - Powder-Free.....	5
Chloroprene Synthetic Exam Gloves - Powder-Free.....	6
Colloidal Oatmeal Active Therapeutic System Exam Gloves - Powder-Free.....	7
Nitrile Synthetic Exam Gloves - Powder-Free.....	8-9
Latex Exam Gloves - <i>Extra Protection</i> - Powder-Free.....	10-11
Nitrile Synthetic Exam Gloves - <i>Extra Protection</i> - Powder-Free.....	11
Latex Sterile Surgical Gloves - Powdered/Powder-Free.....	12
Latex Specialty Sterile Surgical Gloves - Powder-Free.....	12-13
Synthetic Sterile Surgical Gloves - Powder-Free.....	13-14
Latex - Sterile Ambi Exam Gloves - Singles and Pairs.....	14
Nitrile - Sterile Ambi Exam Gloves - Singles & Pairs, 9 1/2" & 12" Length.....	14-15
Wipes Disposable - Adult - Thermalbond and Spunlace.....	15
Wipes Disposable - Baby - Spunlace.....	16
Wipes Disposable - Incontinence.....	16
Patient Cleansing Bath Systems.....	17
Protective Apparel - Isolation Gowns.....	18
Protective Apparel - Bouffant Caps.....	18
Protective Apparel - Shoe Covers.....	18
Protective Apparel - Disposable Patient Footwear.....	18
Masks - Procedure.....	18
Product By Series Index.....	19
Chemotherapy Drug Permeation Guide.....	19



## LATEX EXAM GLOVES - POWDER-FREE



### DermAssist® COATS™

- Textured Finish
- Low-Modulus
- Low Protein (<50µg)

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	124050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	124100	Finger Thickness	avg. 4.5 mil (0.11mm)
M	(7 ½ - 8)	124200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	124300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	124350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:** 100 Gloves/Box, 10 Boxes/Case

**Characteristics:** COATS™ barrier technology is a totally powder free, all-natural physical barrier between the latex surface and the wearer's skin, utilizing the anti-irritant and anti-inflammatory properties of colloidal oatmeal. In our COATS™ gloves, the colloidal oatmeal barrier not only shields the skin surface, it also creates a moisturizing and pH balancing environment that helps dry, damaged skin to repair itself.



### Pulse®

- Textured Finish
- Powder-Free
- Online-Chlorination
- Low Protein (<50µg)

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	151050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	151100	Finger Thickness	avg. 4.5 mil (0.11mm)
M	(7 ½ - 8)	151200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	151300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	151350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Engineered to be slightly thinner for superior tactile sensitivity. Online chlorination provides superior gripping ability and results in very soft gloves with a clean white color. Excellent donning properties and comfort.



### DermaGold® PF

- Smooth Finish
- Powder-Free
- Polymer Bonded
- Low Protein (<50µg)

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	153050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	153100	Finger Thickness	avg. 5.5 mil (0.13mm)
M	(7 ½ - 8)	153200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	153300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	153350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Polymer bonded interior provides excellent donning without the use of chlorination. Smooth surface provides good gripping characteristics without tackiness. A very soft glove with excellent elasticity.



### ProDerm™


- Textured Finish
- Powder-Free
- Polymer-Bonded
- Low Protein (<50µg)

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	155050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	155100	Finger Thickness	avg. 4.5 mil (0.11mm)
M	(7 ½ - 8)	155200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	155300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	155350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Polymer bonded interior and lightly textured surface, engineered to be slightly thinner for incredible tactile sensitivity and comfort. Unique on-line polymer process creates a glove that's very easy to don without the use of chlorination and very comfortable for extended wear.

## LATEX EXAM GLOVES - POWDERED



### DermAssist®


- Smooth Finish
- Powdered

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	115050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	115100	Finger Thickness	avg. 5.0 mil (0.12mm)
M	(7 ½ - 8)	115200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	115300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	115350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Excellent donning properties, low modulus and exceptional strength. Consistently low powder levels.

## VINYL SYNTHETIC EXAM GLOVES - POWDERED



### DermAssist®

- Smooth Finish
- Powdered
- Non-Latex

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	119050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	119100	Finger Thickness	avg. 5.5 mil (0.13mm)
M	(7 ½ - 8)	119200	Tensile Strength Before Aging	avg. 18 MPa
L	(8 ½ - 9)	119300	Tensile Strength After Aging	avg. 17 MPa
XL	(9 ½ - 10)	119350	Elongation Before Aging (%)	avg. 410
			Elongation After Aging (%)	avg. 400

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Synthetic material contains no natural rubber latex proteins. Pure PVC polymer is soft and semi-elastic. Light powder for easy donning. Low cost choice for non-latex requirements.

## VINYL SYNTHETIC EXAM GLOVES - POWDER-FREE



### DermAssist®

- Smooth Finish
- Powder-Free
- Polymer Coated
- Non-Latex

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	161050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	161100	Finger Thickness (mm)	avg. 5.5 mil (0.13mm)
M	(7 ½ - 8)	161200	Tensile Strength Before Aging	avg. 18 MPa
L	(8 ½ - 9)	161300	Tensile Strength After Aging	avg. 17 MPa
XL	(9 ½ - 10)	161350	Elongation Before Aging (%)	avg. 410
			Elongation After Aging (%)	avg. 400

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Powder-Free synthetic contains no natural rubber latex proteins. Soft, flexible, PVC polymer is bonded with a polyurethane interior, creating a wonderfully comfortable non-latex glove.



### DermAssist® Stretch Vinyl

- Smooth Finish
- Powder-Free
- Polymer Coated
- Non-Latex

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	162050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	162100	Finger Thickness (mm)	avg. 5.5 mil (0.13mm)
M	(7 ½ - 8)	162200	Tensile Strength Before Aging	avg. 18 MPa
L	(8 ½ - 9)	162300	Tensile Strength After Aging	avg. 17 MPa
XL	(9 ½ - 10)	162350	Elongation Before Aging (%)	avg. 410
			Elongation After Aging (%)	avg. 400


**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** The next generation of vinyl synthetic has arrived. This unique PVC polymer blend has a soft, elastic feel with a more snug fit at the cuff than traditional vinyl gloves. Pure polymer interior for easy donning. Tan latex-like color.



## - Thinner, Better Chloroprene Exam Gloves

Made from non-latex, polychloroprene synthetic polymer, these new DermAssist® CR exam gloves offer all of the benefits of natural rubber latex, with none of the allergen problems associated with latex. Thicker than most nitrile gloves, yet incredibly elastic with very high tensile strength. Textured finish provides excellent wet or dry gripping ability.



- Textured Finish
- Powder-Free
- Ultra Low Modulus



### **DermAssist® CR Pink** Powder-Free Chloroprene Exam Gloves

#### SIZES

XS	(5 ½ - 6)
S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)

#### ORDER NO

191050
191100
191200
191300
191350

#### PRODUCT SPECIFICATIONS

Length	avg. 9.5" (240mm)
Finger Thickness (mm)	avg. 5.1 mil (0.13mm)
Tensile Strength Before Aging	avg. 15.9 MPa
Tensile Strength After Aging	avg. 20.9 MPa
Elongation Before Aging (%)	avg. 817
Elongation After Aging (%)	avg. 732



- Textured Finish
- Powder-Free
- Ultra Low Modulus



### **DermAssist® CR** Powder-Free Chloroprene Exam Gloves

#### SIZES

XS	(5 ½ - 6)
S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)

#### ORDER NO

192050
192100
192200
192300
192350

#### PRODUCT SPECIFICATIONS

Length	avg. 9.5" (240mm)
Finger Thickness	avg. 5.1 mil (0.13mm)
Tensile Strength Before Aging	avg. 15.9 MPa
Tensile Strength After Aging	avg. 20.9 MPa
Elongation Before Aging (%)	avg. 817
Elongation After Aging (%)	avg. 732

### Features

- Non-Latex
- Textured Finish
- Look and Feel of Latex

### Benefits

- No Risk of Latex Allergens
- Improved Wet/Dry Grip
- Fit and Comfort of Latex without Protein Allergens




**Tested for use with Chemotherapy Drugs**



## - Gloves that protect you and your skin.

Our COATS™ barrier technology is a totally powder free, all-natural physical barrier between the glove surface and the wearer's skin. The colloidal oatmeal barrier not only shields the skin surface from the glove, it also creates a moisturizing, and pH balancing environment that helps dry, damaged skin to repair itself, and helps to keep healthy skin healthy.




- Fingertip Textured
- Powder-Free
- Ultra Low Modulus



## NITRIDERM® COATS™ POWDER-FREE NITRILE EXAM GLOVES

SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	125052	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	125102	Finger Thickness (mm)	avg. 3.5 mil (0.09mm)
M	(7 ½ - 8)	125202	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	125302	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	125352	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480



- Textured Finish
- Powder-Free
- Low Protein (<50µg\*)

## DERMASSIST® COATS™ POWDER-FREE LATEX EXAM GLOVES

SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	124050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	124100	Finger Thickness	avg. 4.5 mil (0.11mm)
M	(7 ½ - 8)	124200	Tensile Strength Before Aging	avg. 28 MPa
L	(8 ½ - 9)	124300	Tensile Strength After Aging	avg. 24 MPa
XL	(9 ½ - 10)	124350	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

## COATS™ - Colloidal Oatmeal Active Therapeutic System

- COATS™ is an innovative new barrier coating comprised of all natural colloidal oatmeal USP which is an FDA recognized skin protectant drug.
- Applied to the don side of the glove through a proprietary process, COATS™ is a powder-free process.
- COATS™ technology reduces irritation and hydration dermatitis by creating a physical barrier between glove film and skin; absorbing perspiration, urea, and salt excretions; providing anti-irritant and anti-inflammatory benefits; effectively pH buffering the skin; and naturally moisturizing the surface of the skin.

**And healthy skin is the body's best defense against infection.**

***What a naturally great idea!***

## NITRILE SYNTHETIC EXAM GLOVES - POWDER-FREE



### NitriDerm® COATS™

- Low-Modulus
- Fingertip Textured
- Powder-Free
- Non-Latex

SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	125052	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	125102	Finger Thickness (mm)	avg. 3.5 mil (0.09mm)
M	(7 ½ - 8)	125202	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	125302	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	125352	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

**All Sizes:** 200 Gloves/Box, 10 Boxes/Case

**Characteristics:** COATS™ barrier technology is a totally powder free, all-natural physical barrier between the nitrile surface and the wearer's skin, utilizing the anti-irritant and anti-inflammatory properties of colloidal oatmeal. In our COATS™ gloves, the colloidal oatmeal barrier not only shields the skin surface, it also creates a moisturizing and pH balancing environment that helps dry, damaged skin to repair itself.



### NitriDerm® Ultra Blue

- Fingertip Textured
- Powder-Free
- Non-Latex




**Tested for use with  
Chemotherapy Drugs**

SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	157050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	157100	Finger Thickness (mm)	avg. 4 mil (0.10mm)
M	(7 ½ - 8)	157200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	157300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	157350	Elongation Before Aging (%)	avg. 650
2XL	(10 ½ - 11)	157400	Elongation After Aging (%)	avg. 480

**Sizes XS-XL:**  
100 Gloves/Box, 10 Boxes/Case

**Size 2XL:**  
80 Gloves/Box, 10 Boxes/Case

**Characteristics:** Exceptionally soft, low modulus nitrile synthetic compound in a traditional weight. Tough yet very comfortable for extended periods of wear. Excellent donning properties. Textured and slightly tacky exterior finish provides unequalled gripping ability. Light blue color.




### ProDerm™ XR Nitrile

- Thinner and Lighter
- Fingertip Textured
- Powder-Free
- Non-Latex
- 250 Glove per Box

SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	158050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	158100	Finger Thickness (mm)	avg. 3.0 mil (0.08mm)
M	(7 ½ - 8)	158200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	158300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	158350	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

**All Sizes:**  
250 Gloves/Box, 10 Boxes/Case

**Characteristics:** New ProDerm XR Nitrile Exam Gloves have been re-designed from the ground up to meet the constantly evolving needs of today's health care practitioners. This exciting new nitrile compounding technology allows for a thinner, lighter glove with extraordinary tactile sensitivity, yet equally extraordinary strength. Easy to don, soft, comfortable and providing all of the required barrier properties for a wide range of health care applications. Sky blue color and packaged 250 Gloves/Box, 10 Boxes/Case



### NitriDerm® 100

- Fingertip Textured
- Powder-Free
- Non-Latex


SIZES		ORDER NO	PRODUCT SPECIFICATIONS	
XS	(5 ½ - 6)	159050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	159100	Finger Thickness (mm)	avg. 3.5 mil (0.09mm)
M	(7 ½ - 8)	159200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	159300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	159350	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

**All Sizes:**  
100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Thin film compounding technology allows for a glove with excellent tactile sensitivity while maintaining strength and durability. This very low modulus nitrile glove has good donning characteristics for all-day comfort. The traditional 100 count box makes this the preferred option for those that prefer a smaller package size. Blue color.



## NITRILE SYNTHETIC EXAM GLOVES - POWDER-FREE



### NitriDerm® Ultra White

- Fingertip Textured
- Powder-Free
- Non-Latex
- Ultra Low Modulus

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	167050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	167100	Finger Thickness	avg. 3.5 mil (0.09mm)
M	(7 ½ - 8)	167200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	167300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	167350	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

#### All Sizes:

100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Our exceptionally soft and comfortable thin nitrile in a white version which is particularly popular with dental end users. White color eliminates the problem of colored gloves reflecting light during tooth shade matching procedures. Textured for improved grip in both wet and dry environments. White color.



### Pulse® Nitrile

- 200 Gloves Per Box\*
- Finger-Tip Textured
- Powder-Free
- Non-Latex
- Ultra Low Modulus

✓ **Tested for use with  
Chemotherapy Drugs**

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	177052	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	177102	Finger Thickness	avg. 3.5 mil (0.09mm)
M	(7 ½ - 8)	177202	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	177302	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	177352	Elongation Before Aging (%)	avg. 650
2XL	(10 ½ - 11)	177402	Elongation After Aging (%)	avg. 480

#### XS-XL:

200 Gloves/Box, 10 Boxes/Case

#### Size 2XL:

180 Gloves/Box\*, 10 Boxes/Case

**Characteristics:** This revolutionary nitrile compounding technology allows for a thinner glove with incredible elasticity. The finely textured finger tips create an unequalled tactile sensitivity and excellent gripping properties. Excellent donning properties and one of the most comfortable nitrile gloves available. Dark lavender blue color.



### NitriDerm® Ultra Black

- Fingertip Textured
- Powder-Free
- Non-Latex

✓ **Tested for use with  
Chemotherapy Drugs**



	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	187050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	187100	Finger Thickness	avg. 4.3 mil (0.11mm)
M	(7 ½ - 8)	187200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	187300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	187350	Elongation Before Aging (%)	avg. 650
2XL	(10½ -11)	187400	Elongation After Aging (%)	avg. 480
3XL	(11½ -12)	187450		


#### Sizes XS-XL:

100 Gloves/Box, 10 Boxes/Case

#### Size 2XL-3XL:

90 Gloves/Box, 10 Boxes/Case

**Characteristics:** This popular black color masks a variety of stains and is a true medical grade exam glove. Engineered with a medium modulus and a high tensile strength, this glove is first choice for a wide range of end users including law enforcement and corrections. Textured for excellent grip. Black color.



### NitriDerm® Ultra Pro

- Non-Latex
- Textured Finish
- 5.1 milFinger Thickness

✓ **Tested for use with  
Chemotherapy Drugs**

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	188050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	188100	Finger Thickness	avg. 5.1 mil (0.13mm)
M	(7 ½ - 8)	188200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	188300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	188350	Elongation Before Aging (%)	avg. 650
2XL	(10½ -11)	188400	Elongation After Aging (%)	avg. 480

#### Sizes XS-XL:

100 Gloves/Box, 10 Boxes/Case


#### Size 2XL:

90 Gloves/Box, 10 Boxes/Case

**Characteristics:** This powder-free nitrile synthetic glove has a 5.1 mil Finger Thickness to provide extra protection in high risk situations. This glove is rated for use with chemotherapy drugs. Excellent donning properties in a soft, low modulus formula. Textured, slightly tacky surface for a better grip. Cobalt blue color.



## NITRILE SYNTHETIC EXAM GLOVES - POWDER-FREE



**NitriDerm®**  
**Scents**

**Great Grape™**


- Fingertip Textured
- Low Modulus
- Non-Latex
- Grape Color and Scent

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	196050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	196100	Finger Thickness	avg. 4 mil (0.10mm)
M	(7 ½ - 8)	196200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	196300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	196350	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

**All Sizes:**

100 Gloves/Box, 10 Boxes/Case

**Characteristics:** Designed for dentistry, our Great Grape™ gloves are colored to match their light and enjoyable scent. The exceptionally soft, low modulus nitrile synthetic compound is very comfortable for extended periods, perfect for dentistry! Not only are they patient fun and friendly, they are excellent nitrile gloves with very good donning properties and a slightly tacky exterior finish which provides excellent gripping properties.



**NitriDerm®**  
**AquaBLUE™**

- Fingertip Textured
- Powder-Free
- Non-Latex

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	197050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	197100	Finger Thickness	avg. 4 mil (0.10mm)
M	(7 ½ - 8)	197200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	197300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	197350	Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

**Sizes XS-XL:**


100 Gloves/Box, 10 Boxes/Case

**Size 2XL:**

90 Gloves/Box, 10 Boxes/Case

**Characteristics:** Tested for use with 12 common chemotherapy drugs, with all drug permeation data printed on the box. This medium weight nitrile glove is extremely elastic and comfortable while providing the tensile strength required for today's demanding procedures. Polymer inner surface for effortless donning. Textured, slightly tacky finish for excellent gripping ability even in the presence of moisture. Aqua blue color.

✓ **Tested for use with  
Chemotherapy Drugs**



**NitriDerm®**  
**Ultra Orange™**

- Low Modulus
- Non-Latex
- Textured Finish
- Orange Color

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
XS	(5 ½ - 6)	199050	Length	avg. 9.5" (240mm)
S	(6 ½ - 7)	199100	Finger Thickness	avg. 4.3 mil (0.11mm)
M	(7 ½ - 8)	199200	Tensile Strength Before Aging	avg. 32 MPa
L	(8 ½ - 9)	199300	Tensile Strength After Aging	avg. 27 MPa
XL	(9 ½ - 10)	199350	Elongation Before Aging (%)	avg. 650
2XL	(10½ - 11)	199400	Elongation After Aging (%)	avg. 480
3XL	(11½ - 12)	199450		

**Sizes XS-XL:**

100 Gloves/Box, 10 Boxes/Case

**Size 2XL-3XL:**

90 Gloves/Box, 10 Boxes/Case

**Characteristics:** Orange - The universal symbol for Caution - and ideal for EMS, police and fire, haz mat functions, spills, etc. This medium weight nitrile glove is extremely elastic and comfortable while providing the tensile strength required for these demanding procedures. Polymer interior for rapid donning. Textured finish for an excellent grip even in the presence of moisture. Color: CAUTION Orange.



PROTECTIVE SINGLE-USE EXAMINATION GLOVES FOR EMERGENCY MEDICAL OPERATIONS IN ACCORDANCE WITH THE NATIONAL FIRE PROTECTION ASSOCIATION STANDARD ON PROTECTIVE CLOTHING FOR EMERGENCY MEDICAL OPERATIONS, NFPA 1999-2008

## LATEX EXAM GLOVES - EXTRA PROTECTION - POWDER-FREE



**DermaGold® EP**

- Textured Finish
- Powder-Free
- Low Protein (<50µg\*)
- 10 Mil

	<u>SIZES</u>	<u>ORDER NO</u>	<u>PRODUCT SPECIFICATIONS</u>	
S	(6 ½ - 7)	180100	Length	avg. 11.5" (290mm)
M	(7 ½ - 8)	180200	Finger Thickness	avg. 10 mil (0.25mm)
L	(8 ½ - 9)	180300	Tensile Strength Before Aging	avg. 28 MPa
XL	(9 ½ - 10)	180350	Tensile Strength After Aging	avg. 24 MPa
2XL	(10½ - 11)	180400	Elongation Before Aging (%)	avg. 850
			Elongation After Aging (%)	avg. 820

**All Sizes:**

50 Gloves/Box, 10 Boxes/Case


**Characteristics:** Thicker and longer than ordinary exam gloves for extra protection in high risk situations, yet providing exceptionally good tactile sensitivity. Ultra low modulus compound is very comfortable for extended periods.



PROTECTIVE SINGLE-USE EXAMINATION GLOVES FOR EMERGENCY MEDICAL OPERATIONS IN ACCORDANCE WITH THE NATIONAL FIRE PROTECTION ASSOCIATION STANDARD ON PROTECTIVE CLOTHING FOR EMERGENCY MEDICAL OPERATIONS, NFPA 1999-2008

✓ **Tested for use with  
Chemotherapy Drugs**

## LATEX EXAM GLOVES - EXTRA PROTECTION - POWDER-FREE



### DermAssist® EP BLUE®

- Textured Finish
- Powder-Free
- Low Protein (<50µg\*)
- 15 Mil

✓ **Tested for use with  
Chemotherapy Drugs**



PROTECTIVE SINGLE-USE EXAMINATION  
GLOVES FOR EMERGENCY MEDICAL  
OPERATIONS IN ACCORDANCE WITH  
THE NATIONAL FIRE PROTECTION  
ASSOCIATION STANDARD ON PROTECTIVE  
CLOTHING FOR EMERGENCY MEDICAL  
OPERATIONS, NFPA 1999-2008

#### All Sizes:

50 Gloves/Box, 10 Boxes/Case

**Characteristics:** For maximum protection in high risk environments. Chlorinated and coated for rapid donning in critical situations. Textured finish with a medium surface tack for positive grip. High tensile strength for added protection. Blue color.

#### SIZES

S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)
2XL	(10½ - 11)


#### ORDER NO

181100
181200
181300
181350
181400

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 15 mil (0.38mm)
Tensile Strength Before Aging	avg. 28 MPa
Tensile Strength After Aging	avg. 24 MPa
Elongation Before Aging (%)	avg. 850
Elongation After Aging (%)	avg. 820

## NITRILE EXAM GLOVES - EXTRA PROTECTION - POWDER-FREE



### NitriDerm® EP

- Textured Finish
- Powder-Free
- Non-Latex
- 5 Mil

✓ **Tested for use with  
Chemotherapy Drugs**



PROTECTIVE SINGLE-USE EXAMINATION  
GLOVES FOR EMERGENCY MEDICAL  
OPERATIONS IN ACCORDANCE WITH  
THE NATIONAL FIRE PROTECTION  
ASSOCIATION STANDARD ON PROTECTIVE  
CLOTHING FOR EMERGENCY MEDICAL  
OPERATIONS, NFPA 1999-2008

#### Sizes XS-XL:

100 Gloves/Box, 10 Boxes/Case

#### Size 2XL:

80 Gloves/Box, 10 Boxes/Case

**Characteristics:** This powder-free nitrile synthetic glove is a full 12 inches in length to provide extra protection in high risk situations. Excellent donning properties in a soft, low modulus formula. Textured, slightly tacky surface for a better grip.

#### SIZES


S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)
2XL	(10½ - 11)

#### ORDER NO

182100
182200
182300
182350
182400

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 5 mil (0.13mm)
Tensile Strength Before Aging	avg. 32 MPa
Tensile Strength After Aging	avg. 27 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 480



### NitriDerm® EP Pro™

- Textured Finish
- Non-Latex
- 5 Mil Finger Thickness

✓ **Tested for use with  
Chemotherapy Drugs**

#### SIZES

S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)
2XL	(10½ - 11)

#### ORDER NO

185100
185200
185300
185350
185400

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 8.0 mil (0.20mm)
Tensile Strength Before Aging	avg. 32 MPa
Tensile Strength After Aging	avg. 27 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 480

#### Sizes XS-XL:

100 Gloves/Box, 10 Boxes/Case

#### Size 2XL:

80 Gloves/Box, 10 Boxes/Case

**Characteristics:** This powder-free nitrile synthetic glove has an 8 mil Finger Thickness to provide extra protection in high risk situations. This glove is rated for use with chemotherapy drugs. Excellent donning properties in a soft, low modulus formula. Textured, slightly tacky surface for a better grip. Light blue color.



### NitriDerm® EP Orange®

- NFPA & UL Approved
- Non-Latex
- Textured Finish



PROTECTIVE SINGLE-USE EXAMINATION  
GLOVES FOR EMERGENCY MEDICAL  
OPERATIONS IN ACCORDANCE WITH  
THE NATIONAL FIRE PROTECTION  
ASSOCIATION STANDARD ON PROTECTIVE  
CLOTHING FOR EMERGENCY MEDICAL  
OPERATIONS, NFPA 1999-2008

#### Sizes XS-XL:

100 Gloves/Box, 10 Boxes/Case

#### Size 2XL-3XL:

90 Gloves/Box, 10 Boxes/Case

**Characteristics:** Orange - The universal symbol for Caution - and ideal for EMS, police and fire, haz mat functions, spills, etc. This powder-free nitrile synthetic glove is a full 12 inches in length to provide extra protection in high risk situations. Excellent donning properties in a soft, low modulus formula. Textured, slightly tacky surface for a better grip.

#### SIZES

S	(6 ½ - 7)
M	(7 ½ - 8)
L	(8 ½ - 9)
XL	(9 ½ - 10)
2XL	(10½ - 11)
3XL	(11½ - 12)


#### ORDER NO

189100
189200
189300
189350
189400
189450

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 6.5 mil (0.17mm)
Tensile Strength Before Aging	avg. 32 MPa
Tensile Strength After Aging	avg. 27 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 480

## STERILE SURGICAL - LATEX POWDERED & POWDER-FREE



### DermAssist® Powdered

50 Pairs/Box, 4 Boxes/Case

- Bisque Finish
- Powdered

SIZES	ORDER NO
5.5	132550
6.0	132600
6.5	132650
7.0	132700
7.5	132750
8.0	132800
8.5	132850
9.0	132900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 8 mil (0.21mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** Manufactured at an incredible 0.65 AQL level, these gloves exceed the standards for all general surgical applications. Bisque surface provides a sure grip whether wet or dry.



### DermAssist® Powder-Free

50 Pairs/Box, 4 Boxes/Case

- Bisque Finish
- Powder-Free
- Low Protein (<50µg\*)

SIZES	ORDER NO
5.5	133550
6.0	133600
6.5	133650
7.0	133700
7.5	133750
8.0	133800
8.5	133850
9.0	133900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 9.5 mil (0.24mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** These low-protein, powder-free gloves eliminate air-borne particles of latex protein associated with powdered surgical gloves. Ultra-low modulus is designed for greater comfort, less hand fatigue during long surgical procedures.

## STERILE SURGICAL - LATEX SPECIALTY POWDER-FREE



- Bisque Finish
- Powder-Free
- Low Protein (<50µg\*)
- Greater Puncture Resistance

### DermAssist® PRESTIGE® Orthopaedic Powder-Free

25 Pairs/Box, 4 Boxes/Case

SIZES	ORDER NO
6.0	131600
6.5	131650
7.0	131700
7.5	131750
8.0	131800
8.5	131850
9.0	131900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 13.8 mil (0.35mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** 30% thicker at the fingertip for improved protection during rigorous orthopaedic surgery. Powder-Free to reduce risk of powder induced allergies, adhesion and granulomas. Bisque finish improves grip.



- Bisque Finish
- Powder-Free
- Low Protein (<50µg\*)
- Extra Sensitive

### DermAssist® PRESTIGE® Microsurgical Powder-Free

25 Pairs/Box, 4 Boxes/Case

SIZES	ORDER NO
5.5	137550
6.0	137600
6.5	137650
7.0	137700
7.5	137750
8.0	137800
8.5	137850
9.0	137900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 7 mil (0.18mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** Thinner than normal surgical gloves for specialists who demand extra sensitivity in the fingertips for micro and ophthalmic procedures. Strong, yet highly sensitive latex barrier is roughened for improved control.

## STERILE SURGICAL - LATEX SPECIALTY POWDER-FREE



### DermAssist® PRESTIGE® DHD Powder-Free

50 Pairs/Box, 4 Boxes/Case

- Bisque Finish
- Powder-Free
- Damp Hand Don
- Low Protein (<50µg\*)


SIZES	ORDER NO
5.5	139550
6.0	139600
6.5	139650
7.0	139700
7.5	139750
8.0	139800
8.5	139850
9.0	139900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 10 mil (0.26mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** Damp don gel interior provides effortless donning. Curved finger design increases comfort. Powder-Free to reduce risk of powder induced allergies, adhesion and granulomas.

\*These gloves contain no more than 50 micrograms of protein per gram of glove.



### DermAssist® 18½" Elbow Length

25 Pairs/Box, 4 Boxes/Case

- Micro Textured
- Powder-Free
- Low Protein (<50µg\*)
- 480mm Length

SIZES	ORDER NO
6.5	141650
7.5	141750
8.5	141850

#### PRODUCT SPECIFICATIONS

Length	avg. 18.5" (490mm)
Finger Thickness	avg. 9.5 mil (0.24mm)
Tensile Strength Before Aging	avg. 18 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 650
Elongation After Aging (%)	avg. 500

**Characteristics:** This elbow length sterile surgical glove is 18 1/2 inches making it essential for deep cavity procedures. The DermAssist Elbow Length is low-protein, and powder-free eliminating air-borne particles of latex protein associated with powdered surgical gloves. Ultra-low modulus is designed for greater comfort, and less hand fatigue during long surgical procedures.

\*These gloves contain no more than 50 micrograms of protein per gram of glove.

## STERILE SURGICAL - SYNTHETIC POWDER-FREE



### Nitriderm® Nitrile Synthetic Powder-Free

25 Pairs/Box, 4 Boxes/Case

- Textured Finish
- Powder-Free
- Non-Latex

SIZES	ORDER NO
6.0	135600
6.5	135650
7.0	135700
7.5	135750
8.0	135800
8.5	135850
9.0	135900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 7 mil (0.18mm)
Tensile Strength Before Aging	avg. 14 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 500
Elongation After Aging (%)	avg. 400

**Characteristics:** These nitrile synthetic surgical gloves contain no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort, fit and protection, while offering exceptional tactile sensitivity.



### DermAssist® PRESTIGE® Chloroprene Powder-Free 25 Pairs/Box, 4 Boxes/Case

- Polymer Coated
- Bisque Finish
- Unique Polychloroprene Formulation
- Latex-Free



SIZES	ORDER NO
5.5	134550
6.0	134600
6.5	134650
7.0	134700
7.5	134750
8.0	134800
8.5	134850
9.0	134900

#### PRODUCT SPECIFICATIONS

Length	avg. 11.5" (290mm)
Finger Thickness	avg. 7.5 mil (0.19mm)
Tensile Strength Before Aging	avg. 14 MPa
Tensile Strength After Aging	avg. 14 MPa
Elongation Before Aging (%)	avg. 500
Elongation After Aging (%)	avg. 400

**Characteristics:** Made from non-latex, polychloroprene synthetic polymer, these new surgical gloves offer all of the benefits of natural rubber latex, with none of the allergen problems associated with latex. The unique formulation allows them to be worn for extended periods without hand fatigue. Slightly roughened finish provides the ideal combination of gripping ability and tactility.



## STERILE SURGICAL - LATEX SPECIALTY POWDER-FREE




### DermAssist® PRESTIGE® Polyisoprene Powder-Free 25 Pairs/Box, 4 Boxes/Case

SIZES	ORDER NO	PRODUCT SPECIFICATIONS	
5.5	140550	Length	avg. 11.5" (290mm)
6.0	140600	Finger Thickness	avg. 9.8 mil (0.25mm)
6.5	140650	Tensile Strength Before Aging	avg. 18 MPa
7.0	140700	Tensile Strength After Aging	avg. 14 MPa
7.5	140750	Elongation Before Aging (%)	avg. 650
8.0	140800	Elongation After Aging (%)	avg. 500
8.5	140850		
9.0	140900		



- Textured Finish
- Powder-Free
- Non-Latex
- Damp Hand Don

**Characteristics:** This revolutionary new synthetic material has all of the preferred physical characteristics of natural rubber latex, but contains no NR latex proteins. Aqueous gel interior provides superior donning even when hands are damp.



### DermAssist® PRESTIGE® Green Powder-Free 25 Pairs/Box, 4 Boxes/Case

SIZES	ORDER NO	PRODUCT SPECIFICATIONS	
5.5	142550	Length	avg. 11.5" (290mm)
6.0	142600	Finger Thickness	avg. 9.8 mil (0.25mm)
6.5	142650	Tensile Strength Before Aging	avg. 18 MPa
7.0	142700	Tensile Strength After Aging	avg. 14 MPa
7.5	142750	Elongation Before Aging (%)	avg. 650
8.0	142800	Elongation After Aging (%)	avg. 500
8.5	142850		
9.0	142900		



- Textured Finish
- Powder-Free
- Non-Latex
- Damp Hand Don

**Characteristics:** This revolutionary new synthetic material has all of the preferred physical characteristics of natural rubber latex, but contains no NR latex proteins. Aqueous gel interior provides superior donning even when hands are damp.

## LATEX EXAM GLOVES - STERILE

### DermAssist® Powdered & Powder-Free Sterile Singles & Pairs Series 101/102/103/104



SIZES	PWD Singles	PWD PAIRS	PF Singles	PF PAIRS	PRODUCT SPECIFICATIONS	
S (6 1/2 - 7)	101100	102100	103100	104100	Length	avg. 9.5" (240mm)
M (7 1/2 - 8)	101200	102200	103200	104200	Finger Thickness	avg. 4.5 mil (0.11mm)
L (8 1/2 - 9)	101300	102300	103300	104300	Tensile Strength Before Aging	avg. 28 MPa
XL (9 1/2 - 10)	101350	102350	103350	104350	Tensile Strength After Aging	avg. 24 MPa
					Elongation Before Aging (%)	avg. 850
					Elongation After Aging (%)	avg. 820

- Sterile
- 9 1/2" Length
- Ambidextrous

**All Sizes:** Singles 100 Gloves/Box, 4 Boxes/Case, Pairs 50 Pairs/Box, 4 Boxes/Case

**Characteristics:** Exceptional strength with excellent donning properties. Inner wrapped singles or pairs in ozone resistant multi-vac packaging.

## NITRILE EXAM GLOVES - STERILE

### NitriDerm® Powder-Free Sterile Singles & Pairs, 9 1/2" Length Series 109/106



SIZES	SINGLES	PAIRS	PRODUCT SPECIFICATIONS	
S (6 1/2 - 7)	109100	106100	Length	avg. 9.5" (240mm)
M (7 1/2 - 8)	109200	106200	Finger Thickness	avg. 4.0 mil (0.10mm)
L (8 1/2 - 9)	109300	106300	Tensile Strength Before Aging	avg. 32 MPa
XL (9 1/2 - 10)	109350	106350	Tensile Strength After Aging	avg. 27 MPa
			Elongation Before Aging (%)	avg. 650
			Elongation After Aging (%)	avg. 480

- Powder-Free
- Ambidextrous
- Non-Latex
- 9 1/2" Length

**All Sizes:** Singles 100 Gloves/Box, 4 Boxes/Case, Pairs 50 Pairs/Box, 4 Boxes/Case

**Characteristics:** Nitrile synthetic contain no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Inner wrapped singles or pairs in ozone resistant multi-vac packaging.

## NITRILE EXAM GLOVES - STERILE

### NitriDerm® Powder-Free Sterile Singles & Pairs, 12" Length

Series 111/112



- Powder-Free
- Ambidextrous
- Non-Latex
- 12" Length

SIZES		SINGLES	PAIRS	PRODUCT SPECIFICATIONS	
S	(6 ½ - 7)	111100	112100	Length	avg. 11.5" (290mm)
M	(7 ½ - 8)	111200	112200	Finger Thickness	avg. 5.1 mil (0.13mm)
L	(8 ½ - 9)	111300	112300	Tensile Strength Before Aging	avg. 32 MPa
XL	(9 ½ - 10)	111350	112350	Tensile Strength After Aging	avg. 27 MPa
				Elongation Before Aging (%)	avg. 650
				Elongation After Aging (%)	avg. 480

All Sizes: Singles 100 Gloves/Box, 4 Boxes/Case, Pairs 50 Pairs/Box, 4 Boxes/Case

**Characteristics:** 12" Nitrile synthetic gloves contain no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Inner wrapped singles or pairs in ozone resistant multi-vac packaging.

### NitriDerm® Powder-Free Sterile Pairs CHEMO, 12" Length

Series 110/114



- Powder-Free
- Ambidextrous
- Non-Latex
- CHEMO Rated
- 797 Compliant

SIZES		PEEL	MULTI-VAC	PRODUCT SPECIFICATIONS	
S	(6 ½ - 7)	110100	114100	Length	avg. 11.5" (290mm)
M	(7 ½ - 8)	110200	114200	Finger Thickness	avg. 7.9 mil (0.20mm)
L	(8 ½ - 9)	110300	114300	Tensile Strength Before Aging	avg. 32 MPa
XL	(9 ½ - 10)	110350	114350	Tensile Strength After Aging	avg. 27 MPa
				Elongation Before Aging (%)	avg. 650
				Elongation After Aging (%)	avg. 480

All Sizes: Pairs 50 Pairs/Box, 4 Boxes/Case

**Characteristics:** Sterile nitrile pairs, tested for use with chemotherapy drugs. Meets all requirements of USP 797: Pharmaceutical Compounding – Sterile Preparations. Designed to allow easy double gloving. Pair packaged, inner wrapper in either ozone resistant multi-vac packaging or traditional peel-open pouch.

✓ Tested for use with Chemotherapy Drugs\*\*

## WIPES - DISPOSABLE - ADULT INCONTINENCE & BABY

### DermAssist® ADULT WIPES - Premium Spunlace



Tub and Softpack with "stay-open" lid

STYLE	ORDER NO	SIZE	PACKAGING
Adult Wipes Tub	80-300	9 x 13"	50 Wipes/Tub, 12 Tubs/Case
Adult Wipes Softpack	80-301	9 x 13"	50 Wipes/Pack, 12 Packs/Case

**Characteristics:** These ultra-soft spunlace washcloths are pre-moistened with our hypoallergenic, alcohol-free formula. Our premium spunlace material is less bulky and more cloth-like than other materials allowing for more thorough cleaning, particularly between skin folds.

### DermAssist® ADULT WIPES - Economical Thermal Bond



STYLE	ORDER NO	SIZE	PACKAGING
Adult Wipes Softpack	80-201	9 x 13"	50 Wipes/Pack, 12Packs/Case

**Characteristics:** These soft, strong wipes are pre-moistened with a gentle cleansing lotion, providing effective patient perineal care, yet they are gentle enough for sensitive skin.

### DermAssist® ADULT WIPES - Economical Spunlace - Fragrance Free



STYLE	ORDER NO	SIZE	PACKAGING
Adult Wipes Softpack	80-504	8 x 12"	64 Wipes/Pack, 12Packs/Case

**Characteristics:** Our most economical spunlace wipes are soft, yet strong and pre-moistened with a gentle cleansing lotion. They deliver effective patient perineal care and are gentle enough for sensitive skin.



## WIPES - DISPOSABLE - ADULT INCONTINENCE & BABY

**DermAssist®**

**ADULT WIPES - Premium Spunlace - Low Profile & Unique Pop-Up Dispenser**



Tub and Softpack  
with "stay-open" lid

<u>STYLE</u>	<u>ORDER NO</u>	<u>SIZE</u>	<u>PACKAGING</u>
Adult Wipes Softpack	80-501	9 x 13"	64 Wipes/Pack, 12Packs/Carton
Adult Wipes <b>Pop-Up Tub</b>	80-502	9 x 13"	64 Wipes/Tub, 12Tubs/Carton

**Characteristics:** Our ultra-soft spunlace washcloths in a low-profile, easy-to-store tub or soft-pack with flip-top lid. Pre-moistened with our hypoallergenic, alcohol-free formula. Our premium spunlace material is more cloth-like than other materials allowing for more thorough cleaning, particularly between skin folds.

**DermAssist®**

**BABY WIPES - Premium Spunlace**

<u>STYLE</u>	<u>ORDER NO</u>	<u>SIZE</u>	<u>PACKAGING</u>
Baby Wipes Tub	80-400	7.1 x 7.9"	100 Wipes/Tub, 12Tubs/Carton
Baby Wipes Softpack	80-401	6.7 x 7.7"	100 Wipes/Pack, 12Packs/Carton

**Characteristics:** Premium cloth-like spunlace wipes are pre-moistened with a gentle cleansing lotion containing aloe and lanolin. Hypoallergenic and pH balanced for baby's sensitive skin. Tubs have a gentle fragrance, softpacks are fragrance free. Convenient, stay-open lid.



Tub with "stay-open" lid

**DermAssist®**

**ADULT WIPES - Spunlace**

<u>STYLE</u>	<u>ORDER NO</u>	<u>SIZE</u>	<u>PACKAGING</u>
Adult Dry Wipes	80-600	9 x 13"	50 Wipes/Pack, 10Packs/Carton
Adult Dry Wipes	80-601	12 x 13"	50 Wipes/Pack, 20Packs/Carton

**Characteristics:** Only DermAssist® offers cloth-like, premium spunlace material in a dry wipe. Stronger when wet than common dry wipes made of air laid material, these wipes are perfect for bathing and incontinent cleanup, as well as general housekeeping applications.



Outer Carton

## WIPES - DISPOSABLE - INCONTINENCE & PERINEAL CARE

**TotalBath®** - Perineal Care - Fragrance Free

<u>STYLE</u>	<u>ORDER NO</u>	<u>SIZE</u>	<u>PACKAGING</u>
Perineal Care Cloths	85-202	8 x 8"	8 Cloths/Pack, 24 Packs/Carton

**Characteristics:** These extra-thick cleansing barrier cloths contain 3% dimethicone for proper perineal care. These cloths are a one step solution for cleaning, moisturizing, and protecting the skin, and fragrance free. Gentle enough for the most sensitive skin. Extra soft with a touch of Aloe Vera.

**DermAssist®** - Adult Wipes with 3% Dimethicone - Perineal Care - Fragrance Free

<u>STYLE</u>	<u>ORDER NO</u>	<u>SIZE</u>	<u>PACKAGING</u>
Adult Wipes <b>Pop-Up Tub</b>	80-503	9 x 13"	64 Wipes/Tub, 12 Tubs/Carton
Adult Perineal Care Wipes	80-505	9 x 13"	24 Wipes/Pack, 24 Packs/Carton
Adult Perineal Care Wipes	80-506	9 x 13"	48 Wipes/Pack, 12 Packs/Carton

**Characteristics:** Premium disposable cleansing wipes for perineal care. Contain 3% dimethicone and are hypoallergenic and fragrance free. Gentle enough for the most sensitive skin. Extra soft with a touch of Aloe Vera.



## TotalBath® - Standard Weight



<u>ORDER NO</u>	<u>SIZE</u>	<u>DESCRIPTION</u>	<u>PACKAGING</u>
85-100	8 x 8"	Antibacterial, Deodorant, Clean Scent, Moisturizing	8 Cloths/Pack, 30 Packs/Carton
85-101	8 x 8"	Antibacterial, Deodorant, Fragrance Free, Moisturizing	8 Cloths/Pack, 30 Packs/Carton

**Characteristics:** These soft, premoistened cleansing cloths provide an effective rinse-free bathing experience that leaves skin feeling moisturized. These cloths are the lowest cost alternative to basin bathing. Available in fresh scent and fragrance free

## TotalBath® - Medium Weight



<u>ORDER NO</u>	<u>SIZE</u>	<u>DESCRIPTION</u>	<u>PACKAGING</u>
85-200	8 x 8"	Antibacterial, Deodorant, Clean Scent, Moisturizing	8 Cloths/Pack, 30 Packs/Carton
85-201	8 x 8"	Antibacterial, Deodorant, Fragrance Free, Moisturizing	8 Cloths/Pack, 30 Packs/Carton

**Characteristics:** Our medium weight cleansing cloths provide more cleansers and moisturizers and a better absorbency than the standard cloths. These cloths are not only cost-effective, but also high in quality. Available in fresh scent and fragrance free.

## TotalBath® - Heavy Weight



<u>ORDER NO</u>	<u>SIZE</u>	<u>DESCRIPTION</u>	<u>PACKAGING</u>
85-300	8 x 8"	Antibacterial, Deodorant, Clean Scent, Moisturizing	8 Cloths/Pack, 30 Packs/Carton
85-301	8 x 8"	Antibacterial, Deodorant, Fragrance Free, Moisturizing	8 Cloths/Pack, 30 Packs/Carton

**Characteristics:** Our heavy weight cleansing cloths provide the ideal bathing experience. These cloths resemble traditional wash cloths in thickness, but are much softer. The moisturizing properties of these cloths help promote healthy skin. Available in fresh scent and fragrance free.

## TotalBath® - Shampoo Cap




<u>ORDER NO</u>	<u>SIZE</u>	<u>DESCRIPTION</u>	<u>PACKAGING</u>
85-600	One Size	Standard Weight with Conditioner, Clean Scent	1 ct pack; 30 Packs/Carton

**Characteristics:** These rinse-free shampoo and conditioner caps are an easy alternative to conventional hair washing methods. These caps clean and condition the hair and scalp in one easy step.

## PROTECTIVE APPAREL - GOWNS, CAPS, SHOE COVERS

### DermAssist® Isolation Gowns - Disposable, Spun Bond - Yellow

<u>SIZES</u>	<u>PIECES/CASE</u>	<u>ORDER NO</u>	<u>GSM</u>
Regular Isolation Gowns	50	75-400	19
Fluid Resistant Isolation Gowns	50	75-500	28
Poly Coated Barrier Gowns	50	75-600	35
XL Poly Coated Barrier Gowns	50	75-601	35
2XL Poly Coated Barrier Gowns	50	75-602	35



**Characteristics:** Full back style in a durable, but comfortable spun bond non-woven material. Also available in a polyethylene coated version for greater fluid protection. Fully sewn.

### DermAssist® Bouffant Caps - Disposable Traditional Style, Lightweight - Blue

<u>SIZES</u>	<u>BOXES/CASE</u>	<u>ORDER NO</u>
21" (Fully Stretched)	100 caps/box, 5 boxes/case	76-210
24" (Fully Stretched)	100 caps/box, 5 boxes/case	76-240

**Characteristics:** Traditional style in a comfortable and lightweight 14 gsm weight.

### DermAssist® Shoe Covers - Disposable Durable, Non-Woven - Blue


<u>SIZES</u>	<u>PACKING</u>	<u>ORDER NO</u>
Universal	50 pairs/box, 3 boxes/case	77-101
Universal Non-Skid	50 pairs/box, 3 boxes/case	77-102
X-Large	50 pairs/box, 3 boxes/case	77-201
X-Large Non-Skid	50 pairs/box, 3 boxes/case	77-202

**Characteristics:** Durable 40 gsm non-woven blue material stands up to wear.

### DermAssist® SOX Disposable Footwear

**Characteristics:** Sox provide lightweight, comfortable, yet durable patient footwear, with gripper soles. Multidirectional stretchability, graduated sizing and elastic cuffs offer easy donning and comfort, even in cases of swelling. Easy to distinguish color-coded sizes.

- Treads
- 91% Polyester
- Graduated Sizing
- Elastic Cuffs
- Color Coded Sizes




<u>ORDER NO</u>	<u>SIZES</u>	<u>COLOR</u>
72-100	Small	Yellow (not shown)
72-200	Medium	Blue
72-500	Unsize	Red
72-600	Unsize	Green

All Sizes: 48 Pairs/Case

## MASKS - PROCEDURE AND SURGICAL


### DermAssist® PREMIER



<u>STYLE</u>	<u>COLOR</u>	<u>ORDER NO</u>	<u>QTY</u>
Pleated, Ear Loop	Blue	78-100	50 Masks/Disp, 10 Disp/Case
Pleated, Ear Loop	Yellow	78-200	50 Masks/Disp, 10 Disp/Case
Pleated, Ear Loop	Pink	78-300	50 Masks/Disp, 10 Disp/Case

Comfortable and highly breathable. Appropriate for all general patient care situations.

- Latex Free, 3-Ply Construction
- High Filtration Meltblown Filter
- BFE > 99% @ 3.0µm
- Easy breathing; Delta-P < 2.2
- Blue, Yellow or Pink



# INDEX

Series	Description	Page	Series	Description	Page
72	DermAssist <sup>®</sup> SOX Disposable Footwear	18	139	DermAssist <sup>®</sup> Prestige <sup>®</sup> DHD, Pwdr-free Latex Surgical Gloves	13
75	DermAssist <sup>®</sup> Disposable Isolation Gowns	18	140	DermAssist <sup>®</sup> Prestige <sup>®</sup> Polyisoprene, Pwdr-free Surgical Gloves	14
76	DermAssist <sup>®</sup> Bouffant Caps	18	141	DermAssist <sup>®</sup> Powder-Free Latex Elbow Length Surgical Gloves	13
77	DermAssist <sup>®</sup> Disposable Shoe covers	18	142	DermAssist <sup>®</sup> Prestige <sup>®</sup> Green Powder-Free Surgical Gloves	14
78	DermAssist <sup>®</sup> Premier Procedure Masks	18	151	Pulse <sup>®</sup> Powder-Free Latex Exam Gloves	4
80	DermAssist <sup>®</sup> Disposable Wipes - Adult, Baby	15	153	DermaGold <sup>®</sup> PF Powder-Free latex Exam Gloves	4
85	TotalBath <sup>®</sup> Patient Cleansing Bath Systems	16	155	ProDerm <sup>™</sup> Powder-Free Latex Exam Gloves	4
101	DermAssist <sup>®</sup> Powdered Latex Sterile Exam Gloves -Singles	14	157	NitriDerm <sup>®</sup> Ultra Blue Powder-Free Nitrile Synthetic Exam Gloves	8
102	DermAssist <sup>®</sup> Powdered Latex Sterile Exam Gloves -Pairs	14	158	ProDerm <sup>™</sup> XR Powder-Free Nitrile Synthetic Exam Gloves	8
103	DermAssist <sup>®</sup> Powder-Free Latex Sterile Exam Gloves- Singles	14	159	NitriDerm <sup>®</sup> 100 Powder-Free Nitrile Synthetic Exam Gloves	8
104	DermAssist <sup>®</sup> Powder-Free Latex Sterile Exam Gloves- Pairs	14	161	DermAssist <sup>®</sup> Powder-Free Vinyl Synthetic Exam Gloves	5
106	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves-Pairs, 9 1/2"	14	162	DermAssist <sup>®</sup> Stretch Powder-Free Vinyl Synthetic Exam Gloves	5
109	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves-Singles, 9 1/2"	14	167	NitriDerm <sup>®</sup> Ultra White Powder-Free Nitrile Synthetic Exam Gloves	9
110	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves- Chemo 12"	15	177	Pulse <sup>®</sup> Powder-Free Nitrile Synthetic Exam Gloves	9
111	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves- 12" Singles	15	180	DermaGold <sup>®</sup> EP Pwdr-free Latex Exam Gloves - Extra Protection	10
112	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves- 12" Pairs	15	181	DermAssist <sup>®</sup> EP Blue <sup>™</sup> Pwdr-Free Latex Exam Gloves - Extra Protection	11
114	NitriDerm <sup>®</sup> Powder-Free Nitrile Sterile Exam Gloves- Chemo 12"	15	182	NitriDerm <sup>®</sup> EP Pwdr-free Nitrile Synthetic Exam Gloves - Extra Protection	11
115	DermAssist <sup>®</sup> Powdered Latex Exam Gloves	5	185	NitriDerm <sup>®</sup> EP Pro <sup>™</sup> Nitrile Synthetic Exam Gloves - Extra Protection	11
119	DermAssist <sup>®</sup> Powdered Vinyl Synthetic Exam Gloves	5	187	NitriDerm <sup>®</sup> Ultra Black Pwdr-free Nitrile Synthetic Exam Gloves	9
124	DermAssist <sup>®</sup> COATS <sup>™</sup> Powder-Free Latex Exam Gloves	6	188	NitriDerm <sup>®</sup> Ultra Pro <sup>™</sup> Pwdr-free Nitrile Synthetic Exam Gloves	9
125	NitriDerm <sup>®</sup> COATS <sup>™</sup> Powder-Free Nitrile Exam Gloves	8	189	NitriDerm <sup>®</sup> EP ORANGE <sup>™</sup> Pwdr-free Nitrile Synthetic Exam Gloves	11
131	DermAssist <sup>®</sup> Prestige <sup>®</sup> Orthopaedic Pwdr-free Latex Surgical Gloves	12	191	DermAssist <sup>®</sup> P CR Pink Powder-Free Synthetic Exam Gloves	6
132	DermAssist <sup>®</sup> Powdered Latex Sterile Surgical Gloves	12	192	DermAssist <sup>®</sup> P CR Powder-Free Synthetic Exam Gloves	6
133	DermAssist <sup>®</sup> Powder-Free Latex Sterile Surgical Gloves	12	196	NitriDerm <sup>®</sup> Scents Great Grape <sup>™</sup> Pwdr-free Nitrile Exam Gloves	8
134	DermAssist <sup>®</sup> Prestige <sup>®</sup> Chloroprene, Powder-Free Surgical Gloves	13	197	NitriDerm <sup>®</sup> AquaBLUE <sup>™</sup> Pwdr-free Nitrile Synthetic Exam Gloves	8
135	NitriDerm <sup>®</sup> Powder-Free Nitrile Synthetic Surgical Gloves	13	199	NitriDerm <sup>®</sup> Ultra Orange <sup>™</sup> Pwdr-free Nitrile Synthetic Exam Gloves	8
137	DermAssist <sup>®</sup> Prestige <sup>®</sup> Microsurgical Pwdr-free Latex Surgical Gloves	12			

## Chemotherapy Drug Permeation

(Minimum breakthrough detection time in minutes, per ASTM D6978-05)

Chemotherapy Drugs	Breakthrough Detection Time					
	Series 110	Series 157	Series 177	Series 187	Series 195	Series 197
Carmustine (3.3 mg/ml)	Do Not Use	35.8*	Do Not Use	Do Not Use	Do Not Use	Do Not Use
Cisplatin	>240	>240	>240	>240	>240	>240
Cyclophosphamide (20 mg/ml)	>240	>240	>240	>240	>240	>240
Dacarbazine (10.0 mg/ml)	>240	>240	>240	>240	>240	>240
Doxorubicin Hydrochloride (2.0 mg/ml)	>240	>240	>240	>240	>240	>240
Etoposide (20 mg/ml)	>240	>240	>240	>240	>240	>240
5-Fluorouracil (50 mg/ml)	>240	>240	>240	>240	>240	>240
Ifosfamide (50 mg/ml)	Not Tested	Not Tested	Not Tested	Not Tested	Not Tested	>240
Methotrexate (25 mg/ml)	Not Tested	>240	>240	Not Tested	Not Tested	>240
Mitomycin C (0.5 mg/ml)	Not Tested	>240	>240	Not Tested	Not Tested	>240
Mitoxantrone (2.0 mg/ml)	Not Tested	Not Tested	Not Tested	Not Tested	Not Tested	>240
Paclitaxel (Taxol) (6.0 mg/ml)	>240	>240	>240	>240	>240	>240
Thio-Tepa (10 mg/ml)	Do Not Use	85.48	Do Not Use	Do Not Use	Do Not Use	Do Not Use
Vincristine Sulfate (1.0 mg/ml)	>240	>240	>240	Not Tested	>240	>240

\* Testing showed an average breakthrough time of 35.8 minutes with Carmustine. Double gloving is recommended when handling this drug.



**INNOVATIVE HEALTHCARE CORPORATION**

8220 Charles Page Blvd ▼ Sand Springs, OK 74063

[www.ihcsolutions.com](http://www.ihcsolutions.com)

**Tel** 800.272.1533 **Fax** 800.650.1014

productguide05052015

