

BAGETTE USAGE CHART

PART CODE	SIZE	POUCHES/BOX	SUITABLE CONTENTS
IMS-1346	2-1/4" x 4" (57 mm x 102 mm)	200	Rubber dam clamps or burs
IMS-1347	2-3/4" x 9" (70 mm x 229 mm)	200	1 handpiece or 1-3 hand instruments
IMS-1236	3-1/2" x 9" (89 mm x 229 mm)	200	2-5 hand instruments
IMS-1348	5-1/4" x 10" (133 mm x 254 mm)	200	Handpiece or multiple instruments
IMS-1237	7-1/2" x 13" (190 mm x 330 mm)	100	• Exam Cassette – 5-Instrument; IM605X • Double-Decker® Cassette – 8 instrument; IM908X • Large Handpiece Cassette; IM60HX • Small Ultrasonic Insert; IM30HUX • Symmetry IQ® Handpiece Cassette; IM60IQ8, IM30IQ8 • 1/4 DIN Cassette – 5 instruments; IM14DIN5X
IMS-1238	10" x 14" (254 mm x 356 mm)	100	• Small Cassette – 8 instruments; IM508X • Small Cassette – 10 instruments; IM510X • Double-Decker Cassette–14 instrument; IM914X • Radiography Cassette; IM7000 • Space Saver Seven Cassette; IMSS7X • Small Thin Ortho Cassette; IMORTHSM8 • Small Resin Cassette – 7-16 instruments; IMS-3107, IMS-3112, IMS-3116 • 1/2 DIN Cassette – 8-10 instruments; IM12DINXXX
IMS-1239	12" x 17" (305 mm x 432 mm)	100	• Large Cassette – 16 instruments; IM416X • Large Cassette – 20 instruments; IM420X • Double-Decker Cassette – 20-instrument; IM920X • Signa-Stat Cassette – 12-instrument; IM812X • Large Thin Ortho Cassette; IMORTHLG8 • Large Resin Cassettes – 12-24 instruments; IMS-1118, IMS-1124, IMS-1112, IMS-1116 • DIN Cassette – 13-20 instruments; IMDINXXX
IMS-1345	13" x 20" (330 mm x 508 mm)	100	• Oral Surgery Cassette; IMOSX • Large Custom Cassette; I18

How the best perform

