



Needle & Syringe Solutions
Safety, Standard, Dental & Sharps Disposal Range



Contents

• Safety Needle & Syringe Products	5
Magellan Safety Products	6
Monoject Safety Range	8
Oral Medication Syringes	9
Angel Wing Blood Collection & Transfer Sets	10
Needleless Vial & IV Access Range	12
• Standard Needles & Syringes	13
Standard Hypodermic Needles	14
Filter & Medication Needles	15
Standard Syringes	16
Blunt Cannula	16
• Dental Products	17
Dental Needles	18
Endodontic Irrigation	19
Dental Injectors	20
Oral Evacuators	21
Aspiration Tip	21
Saliva Ejector	22
Curved Tip Syringes	22
• Sharps Disposal	23
Phlebotomy Containers	24
SharpStar Containers	25
Multi Purpose Containers	26
Large Volume Containers	27

Numerous scientific studies reveal that administering injections and drawing blood with single-use syringes is not only a risk for medical personnel but also a high risk for ancillary staff, such as cleaners and laundry workers, laboratory and other technical support employees

The risks of a health care worker who has been pricked by a needle contaminated by an infected patient, are:

- HBV - 1 out of 3
- HCV - 1 out of 30
- HIV - 1 out of 300

80% of needlestick injuries can be prevented by:

- Eliminating unnecessary use of needles where safe and effective alternatives are available.
- Increasing the use of devices with safety features.
- Promoting education and safe work practices for handling needles and related systems within health care facilities.

In today's healthcare environment, nothing is more important than safety. Protecting clinicians and patients has always been a priority whereas now, it's rapidly becoming a necessity. Healthcare organisations must assess risk and prevent workers' exposure to biological threat and are responsible of adapting their operations to technical progress by using new technologies that can reduce or eliminate risk.

Covidien now offers a range of needle and syringe products to hospitals and healthcare organisations that includes a focus on safety products, as well as the standard and dental ranges of needles & syringes and a range of sharps disposal containers.



Safety Needles & Syringes



Magellan™ Safety Needle

ALL MAGELLAN™ SAFETY PRODUCTS INCLUDE MANY INNOVATIVE FEATURES:

- A choice of easy, one-handed activation methods that minimises training
- An integral, robust safety shield that completely encases the needle
- Bevel orientation so that if the safety shield is up, the needle bevel is up
- Color coded hub to correspond with ISO standards
- A full range of gauges and lengths

MAGELLAN™ HYPODERMIC SAFETY NEEDLES - **Sterile**

- Needle-based safety technology
- Intuitive, one-handed activation
- Ultra-sharp, tri-beveled, anti-coring needle
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Self-leveling needle sheath

	18 Pink
	19 Cream
	20 Yellow
	21 Green
	22 Black
	23 Blue
	25 Orange

All Magellan Safety Hypodermic Needles are ISO colour-coded.

ITEM CODE	NEEDLE SIZE (G X INCH)	NEEDLE SIZE (Ø X LENGTH MM)	COLOUR CODE	BOX/CASE
5551-850558	25g x 5/8"	0.5 x 16mm	Orange	50/500
5551-850510	25g x 1"	0.5 x 25mm	Orange	50/500
5551-850358	23g x 5/8"	0.6 x 16mm	Blue	50/500
5551-850310	23g x 1"	0.6 x 25mm	Blue	50/500
5551-850210	22g x 1"	0.7 x 25mm	Black	50/500
5551-850215	22g x 1.5"	0.7 x 40mm	Black	50/500
5551-850158	21g x 5/8"	0.8 x 16mm	Green	50/500
5551-850110	21g x 1"	0.8 x 25mm	Green	50/500
5551-850115	21g x 1.5"	0.8 x 40mm	Green	50/500
5551-850010	20g x 1"	0.9 x 25mm	Yellow	50/500
5551-850015	20g x 1.5"	0.9 x 40mm	Yellow	50/500
5551-850910	19g x 1"	1.1 x 25mm	Cream	50/500
5551-850915	19g x 1.5"	1.1 x 40mm	Cream	50/500
5551-850810	18g x 1"	1.2 x 25mm	Pink	50/500
5551-850815	18g x 1.5"	1.2 x 40mm	Pink	50/500



MAGELLAN™ INSULIN AND TUBERCULIN SAFETY SYRINGES

- Sterile

- Permanently attached needle
- Needle-based safety technology
- Intuitive, one-handed activation
- Ultra-sharp, tri-beveled, anti-coring needle
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Low dead space

CODE	DESCRIPTION	BOX/CASE
INSULIN		
8881-892930	3/10mL 29g x 13mm	50/500
8881-892950	1/2mL 29g x 13mm	50/500
8881-892910	1mL 29g x 13mm	50/500
8881-893030	3/10mL 30g x 8mm	50/500
8881-893050	1/2mL 30g x 8mm	50/500
8881-893010	1mL 30g x 8mm	50/500
TUBERCULIN		
8881-882712	1mL 27g x 13mm	50/500
8881-882812	1mL 28g x 13mm	50/500



Magellan™ Insulin Safety Syringes



Magellan™ Tuberculin Safety Syringes

MAGELLAN™ SAFETY BLOOD COLLECTION DEVICE - Sterile

- Ready to use 3-in-1 safety design
- One piece safety needle and blood tube holder combination
- Intuitive, one-handed activation
- Integral safety shield
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Single use to help comply with OSHA regulations

CODE	DESCRIPTION	CASE
8881-227121	21g x 1 1/4"	300
8881-227122	22g x 1 1/4"	300



Magellan™ Safety Blood Collection Device

MONOJECT™ INSULIN SAFETY SYRINGES - Sterile

- Soft pack packaging
- Low dead space
- Ultra-Comfort, anti-coring, permanently attached needle
- Accu-Tip flat plunger tip
- For U-100 Insulin

CODE	DESCRIPTION	BOX/CASE
8881-511136	1/2mL Insulin Safety Syringe 29g x 13mm	100/500
8881-511110	1mL Insulin Safety Syringe 29g x 13mm	100/500



Insulin Safety Syringes

MONOJECT™ TUBERCULIN SAFETY SYRINGES - Sterile

- Soft pack packaging
- Low dead space
- Ultra-Comfort, anti-coring, permanently attached needle
- Accu-Tip flat plunger tip
- Accu-Scale graduations: 0.01mL

CODE	DESCRIPTION	BOX/CASE
8881-511235	1mL Tuberculin Safety Syringe 25g x 16mm	100/500
8881-511201	1mL Tuberculin Safety Syringe 28g x 13mm	100/500



Tuberculin Safety Syringes

ALL MONOJECT™ SAFETY SYRINGES OFFER THE FOLLOWING FEATURES AND BENEFITS:

- Barrel-based safety technology
- Robust safety shield that covers the entire needle during transport and after usage
- Transport position that allows clinicians to safely transport medication and eliminates the need to recap the needle
- Rigid pack packaging that may be autoclaved or gas sterilized (Insulin and TB are available only in soft pack packaging)



Monoject™ Safety Syringes

MONOJECT™ SAFETY SYRINGES - Sterile

- Rigid pack packaging
- Bold graduations: 0.1mL / 0.2mL
- Ultra-sharp, tri-beveled, anti-coring needle
- Luer lock syringe tip

CODE	DESCRIPTION	BOX/CASE
8881-566007	6mL Safety Syringe Luer Lock Tip	50/300
8881-522000	12mL Safety Syringe Luer Lock Tip	50/300

MONOJECT™ ORAL MEDICATION SYRINGES - Non-Sterile

- Bulk packaged
- Permanent, heat-etched dual graduations
- Scale markings in milliliters and teaspoons
- Tip will not accept hypodermic needle
- Ribbed tip caps packaged separately
- Amber color available for light-sensitive drugs

CODE	DESCRIPTION	GRADUATIONS	BOX/CASE
8881-901014	Clear 1mL	0.01ml/0.5 minims	100/500
8881-903002	Clear 3mL	0.1ml/0.25 tsp	100/500
8881-906104	Clear 6mL	0.2ml/0.25 tsp	100/500
8881-907102	Clear 10mL	0.2ml/0.25 tsp	100/500
8881-901006	Amber 1mL	0.01ml/0.5 minims	100/500
8881-903010	Amber 3mL	0.1ml/0.25 tsp	100/500
8881-906005	Amber 6mL	0.2ml/0.25 tsp	100/500
8881-907003	Amber 10mL	0.2ml/0.25 tsp	100/500



Monoject™ Oral Medication Syringes

ANGEL WING BLOOD COLLECTION SETS WITH FEMALE LUER - Sterile

- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Female luer connector compatible with all blood drawing equipment
- Latex free, DEHP free
- Colour coded wings to coincide with ISO standards

CODE	NEEDLE GAUGE & LENGTH	TUBING LENGTH	CASE
8881-225174	19g x 19mm	12"	50
8881-225182	21g x 19mm	12"	50
8881-225390	23g x 19mm	3"	50
8881-225190	23g x 19mm	12"	50
8881-225208	25g x 19mm	12"	50



ANGEL WING BLOOD COLLECTION SETS WITH MULTI-SAMPLE LUER ADAPTER - Sterile

- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Luer adapter for multi-sample collection
- Latex free, DEHP free
- Colour coded wings to coincide with ISO standards

CODE	NEEDLE GAUGE & LENGTH	TUBING LENGTH	CASE
8881-225281	19g x 19mm	12"	50
8881-225339	21g x 19mm	3"	50
8881-225299	21g x 19mm	12"	50
8881-225707	23g x 19mm	7"	50
8881-225307	23g x 19mm	12"	50
8881-225315	25g x 19mm	12"	50



ANGEL WING BLOOD COLLECTION SETS WITH HOLDER ATTACHED - Sterile

- Unique self-sheathing needle system
- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Blood tube holder attached
- Colour coded wings to coincide with ISO standards

CODE	NEEDLE GAUGE & LENGTH	TUBING LENGTH	CASE
8881-225221	21g x 19mm	12"	50
8881-225265	23g x 19mm	12"	50
8881-225273	25g x 19mm	12"	50



ANGEL WING LUER LOCK COLLECTION SET - Sterile

- Eliminates the need for a needle when collecting blood from valve based IV ports
- Fast, effective and safe
- Latex free

CODE	DESCRIPTION	CASE
8881-225119	Luer Lock Collection Set	50

**Angel Wing Luer Lock Collection Set****Angel Wing Multi-Sample Transfer Set****ANGEL WING MULTI-SAMPLE TRANSFER SET (DOME) - FEMALE - Sterile**

- Needle-shielded holder for transferring blood specimens from syringes into blood culture bottles
- Protective cap

CODE	DESCRIPTION	CASE
8881-225240	Holder and Safety Cap (Female)	50

**Angel Wing Bact/Alert™ Safety Collection Device****ANGEL WING BACT/ALERT™ SAFETY COLLECTION DEVICE MALE - Sterile**

- Single use for safe transfer of blood samples
- Adapter insert for evacuated tubes
- Designed specifically to fit Bact/Alert blood culture bottles

CODE	DESCRIPTION	CASE
8881-225245	Male Safety Collection Device	50

ANGEL WING TUBE HOLDER WITH FEMALE ADAPTER - Sterile

- Needle-shielded holders for transferring blood specimens from syringes in blood culture bottles
- Protective cap

CODE	DESCRIPTION	CASE
8881-225232	Tube Holder (Female)	50

**Angel Wing Bact/Alert™ Tube Holder with Female Adapter**

MONOJECT™ SMARTIPT™ NEEDLELESS MED PREP CANNULA

- Sterile

- Plastic cannula specifically designed to access rubber stopper medication vials
- 16g flow rate
- Distinctive colour and appearance
- Anti-coring/anti-particulate design

CODE	DESCRIPTION	BOX/CASE
8881-540111	Needleless Vial Access Cannula	100/1000



SmartTip™ Needleless Med Prep Cannula



BlunTip™ Safety I.V. Access Cannula

MONOJECT™ BLUNTIP™ SAFETY I.V. ACCESS PRODUCTS LATEX FREE SAFETY I.V. ACCESS CANNULA - Sterile

- Compatible with a variety of needleless I.V. systems
- Distinctive colour
- Rigid, autoclavable hardpack
- Product identification printed on rigid pack
- Easy-grip design
- 15 gauge x 1/2 cannula

CODE	DESCRIPTION	BOX/CASE
8881-540101	I.V. Access Cannula	100/1000

MONOJECT™ BLUNTIP™ SAFETY I.V. ACCESS CANNULA WITH VIAL ACCESS PIN - Sterile

- Designed to provide needleless access to drug vials with rubber closures
- SoftPack packaging
- Includes a vial-piercing pin and I.V. access cannula

CODE	DESCRIPTION	BOX/CASE
8881-540120	BlunTip with Vial Access Pin	100/1000



BlunTip™ Safety I.V. Access Cannula with Vial Access Pin

Standard Needles & Syringes



MONOJECT™ HYPODERMIC NEEDLES - Sterile

- Polypropylene/Aluminium luer lock hub
- Ultra-sharp, tri-beveled, anti-coring
- May be autoclaved or gas sterilised for use in O.R./procedure packs

CODE	DESCRIPTION	BOX/CASE
WITH POLYPROPYLENE HUB		
8881-250040	Hypodermic Needle 18g x 1" A	100/1000
8881-250339	Hypodermic Needle 20g x 1/2"	100/1000
8881-250123	Hypodermic Needle 20g x 1" A	100/1000
8881-250248	Hypodermic Needle 22g x 3/4" A	100/1000
8881-250222	Hypodermic Needle 22g x 1" A	100/1000
8881-250230	Hypodermic Needle 22g x 1" B	100/1000
8881-250206	Hypodermic Needle 22g x 1 1/2" A	100/1000
8881-250305	Hypodermic Needle 25g x 1" A	100/1000
8881-250313	Hypodermic Needle 25g x 5/8" A	100/1000
8881-250321	Hypodermic Needle 26g x 1/2" IN	100/1000
NZ ONLY		
8881-250016	Hypodermic Needle 18g x 1 1/2" A	100/1000
8881-250107	Hypodermic Needle 20g x 1 1/2" A	100/1000
8881-250289	Hypodermic Needle 23g x 3/4" A	100/1000
WITH ALUMINIUM HUB		
8881-200029	Hypodermic Needle 15g x 1 1/2" B	100/1000
8881-200805	Hypodermic Needle 14g x 1" IN	100/1000
8881-200011	Hypodermic Needle 14g x 1 1/2" IN	100/1000
8881-200573	Hypodermic Needle 14g x 2" IN	100/1000
8881-200755	Hypodermic Needle 16g x 5/8" B	100/1000
8881-200060	Hypodermic Needle 16g x 3/4" IN	100/1000
8881-200052	Hypodermic Needle 16g x 1" A	100/1000
8881-200037	Hypodermic Needle 16g x 1 1/2" A	100/1000
8881-200045	Hypodermic Needle 16g x 1 1/2" B	100/1000
8881-200672	Hypodermic Needle 18g x 5/8" B	100/1000
8881-200714	Hypodermic Needle 18g x 1" A	100/1000
8881-200078	Hypodermic Needle 18g x 1 1/2" A	100/1000
8881-200169	Hypodermic Needle 19g TW*x 1" A	100/1000
8881-200136	Hypodermic Needle 19g TW*x 1 1/2" A	100/1000
8881-200219	Hypodermic Needle 20g x 1" A	100/1000
8881-200185	Hypodermic Needle 20g x 1 1/2" A	100/1000
8881-200235	Hypodermic Needle 20g x 3" A	100/1000
8881-200383	Hypodermic Needle 23g x 1" A	100/1000
8881-200441	Hypodermic Needle 25g x 2" A	100/1000
8881-200508	Hypodermic Needle 27g x 1 1/4" A	100/1000



Monoject™ Hypodermic Needles with Polypropylene Hub



Monoject™ Hypodermic Needles with Aluminium Hub

"A" Bevel- Regular bevel for intramuscular and subcutaneous injections

"B" Bevel- Short bevel for intravenous or intra-arterial applications and nerve blocks

"IN" Bevel- Intermediate bevel for intradermal injections/skin testing

*Thin wall needle

MONOJECT™ FILTER NEEDLES - Sterile

- Polypropylene hub
- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Anti-coring stainless steel needle

CODE	DESCRIPTION	BOX/CASE
8881-305117	Filter Needle 18g x 1 1/2"	100/1000
8881-305018	Filter Needle 20g x 1 1/2"	100/1000

MONOJECT™ FILTER ASPIRATOR - Sterile

- Polypropylene hub
- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Specially designed stainless steel blunt cannula

CODE	DESCRIPTION	BOX/CASE
8881-305109	Filter Aspirator 18g x 3"	25/100

ARGYLE™ FILTER CONNECTOR - Sterile

- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Sufficient flow for gravity infusions
- Male-female luer connector
- Latex free

CODE	DESCRIPTION	BOX/CASE
8888-118059	Filter Connector	100/1000

MONOJECT™ MEDICATION TRANSFER NEEDLE - Sterile

- Double-ended stainless steel needle with special anti-coring bevels
- Smooth, rounded aluminum hub
- Short butt end for evacuation
- Latex free

CODE	DESCRIPTION	BOX/CASE
8881-204005	Medication Transfer Needle 20g x 1"	500/1000

**Monoject™ Filter Needle****Argyle™ Filter Connector**

MONOJECT™ SYRINGES - Sterile

- Rigid pack packaging
- Bold graduations
- May be autoclaved or gas sterilised for use in O.R./procedure packs

CODE	DESCRIPTION	BOX/CASE
8881-513918	3mL Monoject Syringe Regular Luer Tip	100/1000
8881-513934	3mL Monoject Syringe Luer Lock Tip	100/1000
8881-516911	6mL Monoject Syringe Regular Luer Tip	50/500
8881-516937	6mL Monoject Syringe Luer Lock Tip	50/500
8881-512860	12mL Monoject Syringe Eccentric Luer Tip	80/480
8881-520665	20mL Monoject Syringe Eccentric Luer Tip	50/300
8881-535788	35mL Monoject Syringe Eccentric Luer Tip	30/180
8881-560125	60mL Monoject Syringe Luer Lock Tip	20/100
8881-560141	60mL Monoject Syringe Catheter Tip	20/100
8881-560182	60mL Monoject Syringe Eccentric Luer Tip	20/100
8881-560224	60mL Monoject Syringe Regular Luer Tip	20/100
8881-501160	1mL Tuberculin Syringe 25g X 5/8"	100/500
8881-501210	1mL Insulin Syringe U100	100/500
8881-501400	1mL Tuberculin Syringe Reg Luer Tip	100/500
8881-682101	Monoject Luer Tip Caps - Flexible/Sterile	25/1000



Monoject™ Syringes

MONOJECT™ RIGID PACK BLUNT CANNULA - Sterile

- Rigid pack packaging
- Aluminum luer lock hub
- May be autoclaved or gas sterilised for use in O.R./procedure packs

CODE	DESCRIPTION	BOX/CASE
8881-202314	Blunt Cannula 15g x 1 1/2"	25/100
8881-202322	Blunt Cannula 16g x 1 1/2"	25/100
8881-202348	Blunt Cannula 18g x 1"	25/100
8881-202355	Blunt Cannula 19g x 1 1/2"	25/100
8881-202363	Blunt Cannula 20g x 1 1/2"	25/100
8881-202371	Blunt Cannula 21g x 1"	25/100
8881-202389	Blunt Cannula 22g x 1 1/2"	25/100
8881-202397	Blunt Cannula 23g x 1"	25/100



Rigid Pack Blunt Cannula

Dental Products



MONOJECT™ PLASTIC HUB DENTAL NEEDLES

- Sterile, disposable
- Rigid autoclavable polypropylene hardpack
- Tamper-evident heat stake
- Ultra-Sharp tri-bevelled anti-coring lancet point
- Anti-coring butt end point design



Plastic Hub Dental Needles

CODE	DESCRIPTION	NEEDLE LENGTH	COLOUR	BOX/CASE
8881-400033	Non Metric (Imperial) Plastic Hub	25 Long (36.5mm)	Red	100/1000
8881-400041	Non Metric (Imperial) Plastic Hub	25 Short (26mm)	Red	100/1000
8881-400058	Non Metric (Imperial) Plastic Hub	27 Long (36.5mm)	Yellow	100/1000
8881-400066	Non Metric (Imperial) Plastic Hub	27 Short (26mm)	Yellow	100/1000
8881-400074	Non Metric (Imperial) Plastic Hub	30 Short (26mm)	Blue	100/1000
8881-400173	Non Metric (Imperial) Plastic Hub	30G Extra Short (11mm)	Blue	100/1000

MONOJECT™ METAL HUB DENTAL NEEDLES

- Ultra-Sharp tri-bevelled lancet point
- Anti-coring heel
- Stainless steel needle
- Rigid polypropylene protective hardpack
- Tamper-evident heat seal
- Sterile, disposable
- Fits all standard thread dental syringes



Metal Hub Dental Needles

CODE	DESCRIPTION	NEEDLE LENGTH	COLOUR	BOX/CASE
8881-401031	Non Metric (Imperial) Metal Hub	25G Long (36.5mm)	Red	100/1000
8881-401049	Non Metric (Imperial) Metal Hub	25G Short (26mm)	Red	100/1000
8881-401056	Non Metric (Imperial) Metal Hub	27G Long (36.5mm)	Yellow	100/1000
8881-401064	Non Metric (Imperial) Metal Hub	27G Short (26mm)	Yellow	100/1000
8881-401072	Non Metric (Imperial) Metal Hub	30G Short (26mm)	Blue	100/1000
8881-401171	Non Metric (Imperial) Metal Hub	30G Extra Short (11mm)	Blue	100/1000

ENDODONTIC IRRIGATION SYRINGE WITH NEEDLE

- Specially designed for safe flushing and irrigation of root canals
- Luer lock design
- Sterile, individually packaged, single patient use product
- Flow control notched tip permits greater irrigation of tooth canal
- Fluid pressure at apex is significantly reduced by notch tip design

ENDODONTIC IRRIGATION NEEDLES

- Rigid polypropylene hardpack, autoclavable
- Specially designed for safe flushing and irrigation of root canals
- Sterile

CODE	DESCRIPTION	NEEDLE GAUGE	COLOUR	BOX/CASE
8881-471232	Endodontic Needle with notched tip	23G (31mm)	Orange	25/100
8881-471273	Endodontic Needle with notched tip	27G (31mm)	Yellow	25/100
8881-513843	3cc Endodontic Syringe with Needle	23G (31mm)	Orange	100/1000
8881-513850	3cc Endodontic Syringe with Needle	27G (31mm)	Yellow	100/1000



Endodontic Irrigation Needles & Syringes

DENTAL INJECTOR

Disposable Injection Syringe with Reusable Plunger

- Rigid polypropylene hardpack protects against accidental product contamination
- Tamper-evident heat seal protects against accidental exposure to contamination
- Ultra-Sharp tri-bevelled, anti-coring lancet point permits greater injection control and patient comfort
- Anti-coring butt end of needle
- Clear view plastic syringe
- Colour coded cap and Needle sheath
- Rigid, shatter-proof construction
- Reduced potential for accidental needle sticks

CODE	DESCRIPTION	NEEDLE LENGTH	COLOUR	BOX/CASE
8881-418027	1.8ml Dental Injector	25G Short (26mm)	Red	100 + 2 plungers / 500 +10 plungers
8881-418019	1.8ml Dental Injector	25G Long (36.5mm)	Red	100 + 2 plungers / 500 +10 plungers
8881-418035	1.8ml Dental Injector	27G Long (36.5mm)	Yellow	100 + 2 plungers / 500 +10 plungers
8881-418043	1.8ml Dental Injector	27G Short (26mm)	Yellow	100 + 2 plungers / 500 +10 plungers
8881-418050	1.8ml Dental Injector	30G Short (26mm)	Blue	100 + 2 plungers / 500 +10 plungers



Dental Injector

ORAL EVACUATOR

- Soft rigid wall design allows oral evacuator to be used straight or shaped to any angle desired
- Rigid wall construction provides extra strength for cheek/tongue retraction
- Smooth rounded edges allow for maximum patient comfort
- Fits all high speed evacuation units
- Single use, disposable
- Non-autoclavable

CODE	DESCRIPTION	QUANTITY
8881-443009	Oral Evacuator	100/500

ASPIRATION TIP

- Pre-bent tapered tip allows access to all areas. Tapered tip permits pinpoint aspiration
- Smooth, rounded comfort tip prevents tissue damage and improves patient comfort
- Rigid wall design provides extra strength for cheek and tongue retraction
- Fits all high speed evacuation units
- Sterile, disposable

CODE	DESCRIPTION	QUANTITY
8881-444015	Aspiration Tip	25/100



Oral Evacuator



Aspiration Tip

SALIVA EJECTOR

- Basket strainer tip permits maximum fluid aspiration without clogging
- Soft contoured design increases patient comfort by eliminating sharp edges
- Permanently attached tip increases safety
- Wire embedded in tubing wall allows tubing to be shaped to satisfy virtually all procedure requirements
- Prevents tubing from kinking and/or losing shape
- Clean, ready to use
- Single use

CODE	DESCRIPTION	QUANTITY
8881-450004	Saliva Ejector	250/1000



Saliva Ejector


CURVED TIP SYRINGE

- Precision curved tapered tip
- Ungraduated
- Ideal for pressure irrigation, perio irrigation or impression materials
- Clean, ready to use

CODE	DESCRIPTION	QUANTITY
8881-412012	Curved Tip Syringe 12cc	50/500



Curved Tip Syringe



Sharps Disposal

PHLEBOTOMY SHARPS CONTAINERS

- A variety of container sizes and openings allow for the disposal of other sharps
- Facilitates the disconnection of multi-sample needles from needle holders
- Injection moulded design provides greater puncture resistance
- Containers lock for final disposal
- Holders and brackets are available

CODE	DESCRIPTION	DIMENSIONS	CASE QTY
8906	Phlebotomy Sharps Container adapts to present blood drawing trays and techniques 1 litre capacity	16 x 11 x 11cm	100
8970Y	Multi-purpose phlebotomy sharps container with rotor opening. 7.6 litre capacity	25 x 18 x 26cm	20
8912SA	Locking Wall Bracket for use with code: 8906	950ml Phlebotomy size bracket	5
8963	Locking Wall Bracket For use with code: 8970Y	7.6L Phlebotomy size bracket	5
8910	Holder for use with code: 8906	—	10



Phlebotomy Sharps Containers

SHARPSTAR IN-ROOM SYSTEM

- Unique counter-balanced lid:
 - Sharps drop into container away from operator
 - Anti-kick back engineering prevents sharps being thrown out of container
 - Shuts automatically when full
- Injection moulded design provides greater puncture resistance
- Nestable (stackable) for easy storage

CODE	DESCRIPTION	DIMENSIONS	CASE QTY
8506Y	Sharps Container with SHARPSTAR* Lid and counter balanced door 4.75 litre capacity	31 x 14 x 27cm	20
8537Y	Sharps Container with SHARPSTAR* Lid and counter balanced door 11.35 litre capacity	42 x 15 x 35cm	10
8518X	Locking Brackets for use with code: 8506Y	4.75L In-Room size bracket	5
8528C	Locking Brackets for use with code: 8537Y	11.30L In-Room size bracket	5
8550LG	Glove Dispenser for SHARPSTAR 4.75L & 11.35L	—	10
8515	Holder for 4.75L SHARPSTAR	—	5
8529	Holder for 11.35L SHARPSTAR	—	5
85161H	Wall enclosure for 4.75L container	—	1
85301H	Wall enclosure for 11.35L container	—	1



Sharpstar In-Room System

MULTI-PURPOSE SHARPS CONTAINERS

- Adjustable rotor opening or hinged lid accommodate a variety of sizes of sharps and provide temporary and permanent closure
- Easy, simple to use assembly
- Containers lock for final disposal
- Nestable (stackable) for easy storage
- Holders and brackets available can secure containers wherever needed
- Plastics and pigments are safe to autoclave or incinerate

CODE	DESCRIPTION	DIMENSIONS	CASE/QTY
8535Y	Multi-purpose container hinge opening lid. with rotor and 7.6 litre capacity	22 x 15 x 31cm	20
8967Y	Multi-purpose container with horizontal drop opening lid - 7.6 litre capacity	32 x 18 x 26cm	20
8990Y	Multi-purpose container with BIOMAX* sealing gasket hinged lid - 7.6 litre capacity	25 x 18 x 26cm	20
8975	Wall/Cart Bracket for use with codes: 8990Y	–	1
8525Y	Multi-purpose container with rotor opening and hinged lid - 11.35 litre capacity	35 x 15 x 35cm	10
8963	Locking Wall Bracket for use with codes: 8535Y, 8967Y, 8990Y	7.6L Multipurpose size bracket	5
8528C	Locking Wall Bracket for use with code: 8525Y	11.35L Multipurpose size bracket	5
31307005	7.6L & 11.35L Wall Enclosure	–	1
31323358	7.6L Sharps Container Yellow	–	12
31317483	11.35L Sharps Container Yellow	–	12
31320008	7.6L & 11.35L Tabletop Holder	–	1
31320016	7.6L & 11.35L Tabletop Holder	–	1
31320032	Gatorguard Glove Dispenser	–	1
31353553	4.75L Wall Enclosure	–	1
31353611	4.75L Sharps Container Yellow	–	14
31369955	4.75L Tabletop Holder	–	1
31369963	4.75L Tabletop Holder	–	1



Multi-Purpose Containers

LARGE VOLUME CONTAINERS WITH SLIDING, SEALING GASKET AND HINGED LIDS

- Adjustable rotor opening and hinged lids accommodate a variety of sizes of sharps and provide temporary closure.
- Easy, simple to use assembly.
- Containers lock for final disposal.
- Nestable containers save valuable storage space.
- Lids can provide temporary closure when container is not in use to reduce exposure to contents.
- Injection moulded design provides greater puncture resistance.

CODE	DESCRIPTION	DIMENSIONS	CASE/QTY
8980Y	Large Volume Sharps Container with hinged lid - 30 litre capacity	45 x 28 x 39cm	10
8935Y	Large Volume Sharps Container with sliding lid - 45 litre capacity	47 x 32 x 46cm	10
8938Y	Large Volume Sharps Container with sliding lid - 68 litre capacity	66 x 32 x 46cm	5
8991Y	Large Volume Sharps Container with hinged lid - 68 litre capacity	66 x 32 x 46cm	5
8938FP	Wire Foot Pedal Sharps Cart	–	1
8940FP	Molded Foot Pedal Sharps Cart	–	1
8981FP	Wire Foot Pedal Sharps Cart	–	1
8984	Locking Wire Bracket for use with 30L	–	1



Large Volume Containers

Needle & Syringe Solutions Index by Product Code

CODE	PRODUCT DESCRIPTION	PAGE	CODE	PRODUCT DESCRIPTION	PAGE
31307005	7.6L & 11.35L Wall Enclosure	26	8881-200011	Monoject Metal Hub Hypodermic Needle 14g x 1.5" IN	14
31323358	7.6L Sharps Container Yellow	26	8881-200029	Monoject Metal Hub Hypodermic Needle 15g x 1.5" B	14
31317483	11.35L Sharps Container Yellow	26	8881-200037	Monoject Metal Hub Hypodermic Needle 16g x 1.5" A	14
31320008	7.6L & 11.35L Tabletop Holder	26	8881-200045	Monoject Metal Hub Hypodermic Needle 16g x 1.5" B	14
31320016	7.6L & 11.35L Tabletop Holder	26	8881-200052	Monoject Metal Hub Hypodermic Needle 16g x 1" A	14
31320032	Gatorguard Glove Dispenser	26	8881-200060	Monoject Metal Hub Hypodermic Needle 16g x 3/4" IN	14
31353553	4.75L Wall Enclosure	26	8881-200078	Monoject Metal Hub Hypodermic Needle 18g x 1.5" A	14
31353611	4.75L Sharps Container Yellow	26	8881-200136	Monoject Metal Hub Hypodermic Needle 19g TW*x 1.5" A	14
31369955	4.75L Tabletop Holder	26	8881-200169	Monoject Metal Hub Hypodermic Needle 19g TW*x 1" A	14
31369963	4.75L Tabletop Holder	26	8881-200185	Monoject Metal Hub Hypodermic Needle 20g x 1.5" A	14
5551-850010	Magellan Safety Hypodermic Needle 20g x 1"	6	8881-200219	Monoject Metal Hub Hypodermic Needle 20g x 1" A	14
5551-850015	Magellan Safety Hypodermic Needle 20g x 1.5"	6	8881-200235	Monoject Metal Hub Hypodermic Needle 20g x 3" A	14
5551-850110	Magellan Safety Hypodermic Needle 21g x 1"	6	8881-200260	Monoject Metal Hub Hypodermic Needle 21g x 1 1/2"	14
5551-850115	Magellan Safety Hypodermic Needle 21g x 1.5"	6	8881-200383	Monoject Metal Hub Hypodermic Needle 23g x 1" A	14
5551-850158	Magellan Safety Hypodermic Needle 21g x 5/8"	6	8881-200441	Monoject Metal Hub Hypodermic Needle 25g x 2" A	14
5551-850210	Magellan Safety Hypodermic Needle 22g x 1"	6	8881-200508	Monoject Metal Hub Hypodermic Needle 27g x 1.25" A	14
5551-850215	Magellan Safety Hypodermic Needle 22g x 1.5"	6	8881-200573	Monoject Metal Hub Hypodermic Needle 14g x 2" IN	14
5551-850310	Magellan Safety Hypodermic Needle 23g x 1"	6	8881-200672	Monoject Metal Hub Hypodermic Needle 18g x 5/8" B	14
5551-850358	Magellan Safety Hypodermic Needle 23g x 5/8"	6	8881-200714	Monoject Metal Hub Hypodermic Needle 18g x 1" A	14
5551-850510	Magellan Safety Hypodermic Needle 25g x 1"	6	8881-200755	Monoject Metal Hub Hypodermic Needle 16g x 5/8" B	14
5551-850558	Magellan Safety Hypodermic Needle 25g x 5/8"	6	8881-200805	Monoject Metal Hub Hypodermic Needle 14g x 1" IN	14
5551-850810	Magellan Safety Hypodermic Needle 18g x 1"	6	8881-202314	Blunt Cannula 15g x 1.5"	16
5551-850815	Magellan Safety Hypodermic Needle 18g x 1.5"	6	8881-202322	Blunt Cannula 16g x 1.5"	16
5551-850910	Magellan Safety Hypodermic Needle 19g x 1"	6	8881-202348	Blunt Cannula 18g x 1"	16
5551-850915	Magellan Safety Hypodermic Needle 19g x 1.5"	6	8881-202355	Blunt Cannula 19g x 1.5"	16
8506Y	Sharpstar Sharps Container - 4.75L	25	8881-202363	Blunt Cannula 20g x 1.5"	16
8515	Holder for 4.75L Sharpstar	25	8881-202371	Blunt Cannula 21g x 1"	16
85161H	Wall Enclosure for 4.75L Container	25	8881-202389	Blunt Cannula 22g x 1.5"	16
8518X	Locking Bracket for 8506Y	25	8881-202397	Blunt Cannula 23g x 1"	16
8525Y	Multi Purpose Container - hinge with rotor 11.35L	26	8881-204005	Monoject Medication Transfer Needle	15
8528C	Locking Bracket for 8537Y	25	8881-225119	Angel Wing Luer Lock Collection Set	11
8528C	11.35L Multi Purpose Bracket	26	8881-225174	Angel Wing Blood Collection Set with Female Luer 19g	10
8529	Holder for 11.35L Sharpstar	25	8881-225182	Angel Wing Blood Collection Set with Female Luer 21g	10
85301H	Wall Enclosure for 11.35L Container	25	8881-225190	Angel Wing Blood Collection Set with Female Luer 23g	10
8535Y	Multi Purpose Container - hinge with rotor 7.6L	26	8881-225208	Angel Wing Blood Collection Set with Female Luer 25g	10
8537Y	Sharpstar Sharps Container - 11.35L	25	8881-225221	Angel Wing Blood Collection Set with Holder 21g	10
8550LG	Glove Dispenser for Sharpstar	14	8881-225232	Angel Wing Tube Holder with Female Adapter	11

CODE	PRODUCT DESCRIPTION	PAGE	CODE	PRODUCT DESCRIPTION	PAGE
8881-225240	Angel Wing Multi-Sample Dome Transfer Set	11	8881-401049	Monoject Metal Hub Dental Needle 25g Short	18
8881-225245	Angel Wing BacT/Alert Safety Collection Device	11	8881-401056	Monoject Metal Hub Dental Needle 27g Long	18
8881-225265	Angel Wing Blood Collection Set with Holder 23g	10	8881-401064	Monoject Metal Hub Dental Needle 27g Short	18
8881-225273	Angel Wing Blood Collection Set with Holder 25g	10	8881-401072	Monoject Metal Hub Dental Needle 30g Short	18
8881-225281	Angel Wing Blood Collection Set w Multi-Sample Adaptor 19g	10	8881-401171	Monoject Metal Hub Dental Needle 30g Extra-Short	18
8881-225299	Angel Wing Blood Collection Set w Multi-Sample Adaptor 21g	10	8881-412012	Curved Tip Syringe 12mL	22
8881-225307	Angel Wing Blood Collection Set w Multi-Sample Adaptor 23g	10	8881-418019	1.8mL Dental Injector 25G Long	20
8881-225315	Angel Wing Blood Collection Set w Multi-Sample Adaptor 25g	10	8881-418027	1.8mL Dental Injector 25G Short	20
8881-225339	Angel Wing Blood Collection Set w Multi-Sample Adaptor 21g	10	8881-418035	1.8mL Dental Injector 27G Long	20
8881-225382	Angel Wing Blood Collection Set w Female Luer 21g	10	8881-418043	1.8mL Dental Injector 27G Short	20
8881-225390	Angel Wing Blood Collection Set w Female Luer 23g	10	8881-418050	1.8mL Dental Injector 30G Short	20
8881-225707	Angel Wing Blood Collection Set w Multi-Sample Adaptor 23g	10	8881-443009	Oral Evacuator	21
8881-227121	Magellan Safety Blood Collector 21 g x 1.25"	7	8881-444015	Aspiration Tip	21
8881-227122	Magellan Safety Blood Collector 22 g x 1.25"	7	8881-450004	Saliva Ejector	22
8881-250016	Monoject Standard Hypodermic Needle 18g x 1.5" A	14	8881-471232	Endodontic Needle with Notched Tip 23g	19
8881-250040	Monoject Standard Hypodermic Needle 18g x 1" A	14	8881-471273	Endodontic Needle with Notched Tip 27g	19
8881-250107	Monoject Standard Hypodermic Needle 20g x 1.5" A	14	8881-501160	1mL Tuberculin Syringe 25g X 5/8"	16
8881-250123	Monoject Standard Hypodermic Needle 20g x 1" A	14	8881-501210	1mL Insulin Syringe U100	16
8881-250206	Monoject Standard Hypodermic Needle 22g x 1.5" A	14	8881-501400	1mL Tuberculin Syringe Reg Luer Tip	16
8881-250222	Monoject Standard Hypodermic Needle 22g x 1" A	14	8881-511110	Monoject Safety Insulin Syringe 29g x 13mm x 1.0mL	8
8881-250230	Monoject Standard Hypodermic Needle 22g x 1" A	14	8881-511136	Monoject Safety Insulin Syringe 29g x 13mm x 0.5mL	8
8881-250248	Monoject Standard Hypodermic Needle 22g x 3/4" A	14	8881-511201	Monoject Safety Tuberculin Syringe 28g x 13mm x 1mL	8
8881-250289	Monoject Standard Hypodermic Needle 23g x 3/4" A	14	8881-511235	Monoject Safety Tuberculin Syringe 25g x 5/8" x 1mL	8
8881-250305	Monoject Standard Hypodermic Needle 25g x 1" A	14	8881-512860	12mL Monoject Syringe Eccentric Luer Tip	16
8881-250313	Monoject Standard Hypodermic Needle 25g x 5/8" A	14	8881-513843	3mL Endodontic Syringe with Needle 23g	19
8881-250321	Monoject Standard Hypodermic Needle 26g x 0.5" IN	14	8881-513850	3mL Endodontic Syringe with Needle 27g	19
8881-250339	Monoject Standard Hypodermic Needle 20g x 0.5" IN	14	8881-513918	3mL Monoject Syringe Regular Luer Tip	16
8881-305018	Monoject Filter Needle 20g x 1.5"	15	8881-513934	3mL Monoject Syringe Luer Lock Tip	16
8881-305109	Monoject Filter Aspirator 18g x 3"	15	8881-516911	6mL Monoject Syringe Regular Luer Tip	16
8881-305117	Monoject Filter Needle 18g x 1.5"	15	8881-516937	6mL Monoject Syringe Luer Lock Tip	16
8881-400033	Monoject Plastic Hub Dental Needle 25g Long	18	8881-520665	20mL Monoject Syringe Eccentric Luer Tip	16
8881-400041	Monoject Plastic Hub Dental Needle 25g Short	18	8881-522000	Monoject Safety Syringe 12mL	9
8881-400058	Monoject Plastic Hub Dental Needle 27g Long	18	8881-535788	35mL Monoject Syringe Eccentric Luer Tip	16
8881-400066	Monoject Plastic Hub Dental Needle 27g Short	18	8881-540101	BlunTip Safety IV Access Cannula	12
8881-400074	Monoject Plastic Hub Dental Needle 30g Short	18	8881-540111	SmarTip Needles Medication Prep Cannula	12
8881-400173	Monoject Plastic Hub Dental Needle 30g Extra-Short	18	8881-540120	BlunTip Safety IV Access Cannula & Vial Access Pin	12
8881-401031	Monoject Metal Hub Dental Needle 25g Long	18	8881-560125	60mL Monoject Syringe Luer Lock Tip	16

Needle & Syringe Solutions **Index by Product Code**

CODE	PRODUCT DESCRIPTION	PAGE
8881-560141	60mL Monoject Syringe Catheter Tip	16
8881-560182	60mL Monoject Syringe Eccentric Luer Tip	16
8881-560224	60mL Monoject Syringe Regular Luer Tip	16
8881-566007	Monoject Safety Syringe 6mL	9
8881-682101	Monoject Luer Tip Caps - Flexible/Sterile	16
8881-882712	Magellan Safety Tuberculin Syringe 27g x 13mm x 1mL	7
8881-882812	Magellan Safety Tuberculin Syringe 28g x 13mm x 1mL	7
8881-892910	Magellan Safety Insulin Syringe 29g x 13mm x 1.0mL	7
8881-892930	Magellan Safety Insulin Syringe 29g x 13mm x 0.3mL	7
8881-892950	Magellan Safety Insulin Syringe 29g x 13mm x 0.5mL	7
8881-893010	Magellan Safety Insulin Syringe 30g x 8mm x 1.0mL	7
8881-893030	Magellan Safety Insulin Syringe 30g x 8mm x 0.3mL	7
8881-893050	Magellan Safety Insulin Syringe 30g x 8mm x 0.5mL	7
8881-901006	Monoject Oral Medication Syringe Amber 1mL	9
8881-901014	Monoject Oral Medication Syringe Clear 1mL	9
8881-903002	Monoject Oral Medication Syringe Clear 3mL	9
8881-903010	Monoject Oral Medication Syringe Amber 3mL	9
8881-906005	Monoject Oral Medication Syringe Amber 6mL	9
8881-906104	Monoject Oral Medication Syringe Clear 6mL	9
8881-907003	Monoject Oral Medication Syringe Amber 10mL	9
8881-907102	Monoject Oral Medication Syringe Clear 10mL	9
8888-118059	Argyle Filter Connector	15
8906	Phlebotomy Sharps Container 1L	24
8910	Holder for use with code 8906	24
8912SA	Locking Bracket for 8906	24
8935Y	Large Volume Sharps Container with Sliding Lid - 40L	27
8938Y	Large Volume Sharps Container with Sliding Lid - 68L	27
8963	7.6L Multi Purpose Bracket	24, 26
8967Y	7.6L Multi Purpose Container	26
8970Y	Multi Purpose Phlebotomy Container 7.6L	24
8975	Wall/Cart Bracket	26
8980Y	Large Volume Sharps Container with Hinged Lid - 30L	27
8990Y	Multi Purpose Container - Sealing Gasket & Hinge 7.6L	26
8991Y	Large Volume Sharps Container with Hinged Lid - 68L	27
8938FP	Wire Foot Pedal Sharps Cart	27
8940FP	Molded Foot Pedal Sharps Cart	27
8981FP	Wire Foot Pedal Sharps Cart	27
8984	Locking Bracket for 30L	27

COVIDIEN, COVIDIEN with logo and TM marked
brands are trademarks of Covidien AG or its affiliate.
© 2011 Covidien AG or its affiliate. All rights reserved.
Needles 047-02-11



COVIDIEN

COVIDIEN PTY LTD
166 EPPING ROAD,
LANE COVE NSW 2066
AUSTRALIA
(T) 1800 252 467

COVIDIEN NEW ZEALAND LTD
GROUND FLOOR, 15B VESTEY DRIVE,
MOUNT WELLINGTON, AUCKLAND
NEW ZEALAND
(T) 0508 489 264

WWW.COVIDIEN.COM