



Disposable Dental Injection Needles

Pre-threaded plastic hub needles. Bevel indicator on the hub for quick orientation of the needle. Color-coded caps indicate gauge size.

FEATURES

Pre-threaded Plastic Hubs

Bevel Indicator

Siliconized Cannula

Color-Coded Caps

BENEFITS

Easy assembly with less risk of damage to syringes

For quick orientation of the needle for different injections

Ensures a smooth insertion

Easy identification

CODE	DESCRIPTION	COLOR	QTY
DP42830	30 Gauge X-Short $\frac{3}{8}$ "	Purple	100/box
DP47830	30 Gauge Short 1"	Blue	100/box
DP47827	27 Gauge Short 1"	Orange	100/box
DP47825	25 Gauge Short 1"	Pink	100/box
DP40827	27 Gauge Long $1\frac{3}{8}$ "	Yellow	100/box
DP40825	25 Gauge Long $1\frac{3}{8}$ "	Red	100/box

Distributed by:



Quala
DENTAL PRODUCTS
www.quala.com

#QSS010R1