

## Redefining Excellence



589 Davies Drive  
York, PA 17402  
toll-free phone 866-854-8300  
toll-free fax 866-854-8400  
phone 717-840-9335  
fax 717-840-9347  
email [dental@miltex.com](mailto:dental@miltex.com)  
[www.miltex.com](http://www.miltex.com)



## INTRODUCTION

The Precision™ by Miltex and Precise™ by Miltex product lines are the same high quality product lines that were once brought to you by the Precise/Precision Dental company. Our product offering includes but is not limited to the following products:

- Cabinet Trays
- Code Rings
- Cotton Roll Holders
- Germicide Soaking Trays
- True-Stop Endo Stop Dispenser
- Mouth Props
- Precision™ Drills
- Precision™ Posts
- Dental Mouth Mirrors
- Gutta Percha
- Paper Points



## TABLE OF CONTENTS

|                                     |       |
|-------------------------------------|-------|
| Amalgam Well .....                  | 2     |
| Cabinet Trays .....                 | 2-3   |
| Instrument Mats .....               | 4     |
| Cotton Roll Holders .....           | 5     |
| Deluxe Cotton Roll Holders .....    | 5     |
| Sponge Dispensers .....             | 6     |
| Germicide Soaking Trays .....       | 7     |
| Code Rings .....                    | 8-9   |
| X-Ray Holders .....                 | 10    |
| Bib Clips .....                     | 10    |
| Mouth Props .....                   | 11    |
| Black Rubber .....                  | 11    |
| Tan Silicone .....                  | 11    |
| Dental Mouth Mirrors .....          | 12    |
| Silicone Endo Stops .....           | 13    |
| True-Stop Endo Stop Dispenser ..... | 13    |
| Precision™ Drills .....             | 14    |
| Precision™ Post Kit .....           | 14    |
| Precision™ Posts .....              | 15    |
| Gutta Percha .....                  | 16-17 |
| Paper Points .....                  | 18    |



### Amalgam Well

- Designed for mixing and dispensing amalgam by guiding carrier to amalgam filling material
- Corrosion resistant stainless steel
- Skid-resistant silicone "o" ring base
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes

| REF             | Description                    |
|-----------------|--------------------------------|
| 700100-0000-001 | Amalgam Well - Stainless Steel |

### Cabinet Trays

- Organize instruments and other dental products in drawers or for procedures
- Heavy gauge plastic prevents warping and resists staining
- Available in 6 different configurations and in multiple colors
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes
- Trays measure 3-3/4" × 7-3/4" (9.5 cm × 19.7 cm)



#### Cabinet Tray - 16

| REF             | Description            |
|-----------------|------------------------|
| 700210-1000-001 | Cabinet Tray 16, Beige |
| 700210-2000-001 | Cabinet Tray 16, Blue  |
| 700210-4500-001 | Cabinet Tray 16, Mauve |
| 700210-8500-001 | Cabinet Tray 16, White |



#### Cabinet Tray - 16A/20A

| REF             | Description                 |
|-----------------|-----------------------------|
| 700215-1000-001 | Cabinet Tray 16A/20A, Beige |
| 700215-2000-001 | Cabinet Tray 16A/20A, Blue  |
| 700215-4500-001 | Cabinet Tray 16A/20A, Mauve |
| 700215-8500-001 | Cabinet Tray 16A/20A, White |

## Cabinet Trays

### Cabinet Tray - 17

| REF             | Description            |
|-----------------|------------------------|
| 700211-1000-001 | Cabinet Tray 17, Beige |
| 700211-2000-001 | Cabinet Tray 17, Blue  |
| 700211-4500-001 | Cabinet Tray 17, Mauve |
| 700211-8500-001 | Cabinet Tray 17, White |



### Cabinet Tray - 18

| REF             | Description            |
|-----------------|------------------------|
| 700212-1000-001 | Cabinet Tray 18, Beige |
| 700212-2000-001 | Cabinet Tray 18, Blue  |
| 700212-4500-001 | Cabinet Tray 18, Mauve |
| 700212-8500-001 | Cabinet Tray 18, White |



### Cabinet Tray - 19

| REF             | Description            |
|-----------------|------------------------|
| 700213-1000-001 | Cabinet Tray 19, Beige |
| 700213-2000-001 | Cabinet Tray 19, Blue  |
| 700213-4500-001 | Cabinet Tray 19, Mauve |
| 700213-8500-001 | Cabinet Tray 19, White |



### Cabinet Tray - 20

| REF             | Description            |
|-----------------|------------------------|
| 700214-1000-001 | Cabinet Tray 20, Beige |
| 700214-2000-001 | Cabinet Tray 20, Blue  |
| 700214-4500-001 | Cabinet Tray 20, Mauve |
| 700214-8500-001 | Cabinet Tray 20, White |





## Instrument Mats

- Hold instruments during chairside use
- Each mat holds up to 12 instruments
- Made from a special polymer compound which resists staining
- Available in a variety of colors, both rigid and flexible styles
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes
- Mats measure 3-5/8" × 7-1/2" (9.2 cm × 19.1 cm)

### Instrument Mat - Flexible

| REF             | Description                    |
|-----------------|--------------------------------|
| 700230-1000-001 | Instrument Mat Flexible, Beige |
| 700230-2000-001 | Instrument Mat Flexible, Blue  |
| 700230-4500-001 | Instrument Mat Flexible, Mauve |
| 700230-8500-001 | Instrument Mat Flexible, White |



### Instrument Mat - Rigid

| REF             | Description                 |
|-----------------|-----------------------------|
| 700235-1000-001 | Instrument Mat Rigid, Beige |
| 700235-2000-001 | Instrument Mat Rigid, Blue  |
| 700235-4500-001 | Instrument Mat Rigid, Mauve |
| 700235-8500-001 | Instrument Mat Rigid, White |

## Cotton Roll Holders

- Hold cotton rolls upright for easy dispensing
- Easy insertion of an entire bundle of cotton rolls
- Hinged lid for easy access to cotton rolls

| REF             | Description |
|-----------------|-------------|
| 700400-1000-001 | Beige       |
| 700400-2000-001 | Blue        |
| 700400-4500-001 | Mauve       |
| 700400-8500-001 | White       |



## Deluxe Cotton Roll Holders

- Store cotton rolls for easy dispensing of 3 cotton rolls at a time
- Easy loading, holds over one pack of rolls
- Lid keeps cotton rolls covered to reduce the risk of airborne contamination
- Weighted, non-skid base for maximum stability

| REF             | Description |
|-----------------|-------------|
| 700410-1000-001 | Beige       |
| 700410-2000-001 | Blue        |
| 700410-3500-001 | Gray        |
| 700410-4500-001 | Mauve       |
| 700410-8500-001 | White       |





## Dispensers

- Hold sponges and gauze for easy dispensing
- Spring loaded from the bottom to push sponges to the top of the dispenser
- Leakproof; allows for disinfecting solutions or alcohol to be added to the gauze or sponges
- Bottom weighted to prevent tipping
- Available in 2 sizes and a variety of colors

### Sponge Dispenser 2" x 2" (5.1 cm x 5.1 cm)

| REF             | Descriptions                       |
|-----------------|------------------------------------|
| 700700-1008-001 | Sponge Dispenser, Beige, Leakproof |
| 700700-2008-001 | Sponge Dispenser, Blue, Leakproof  |
| 700700-4508-001 | Sponge Dispenser, Mauve, Leakproof |
| 700700-8508-001 | Sponge Dispenser, White, Leakproof |



### Sponge Dispenser 4" x 4" (10.2 cm x 10.2 cm)

| REF             | Descriptions                       |
|-----------------|------------------------------------|
| 700701-1010-001 | Sponge Dispenser, Beige, Leakproof |
| 700701-2010-001 | Sponge Dispenser, Blue, Leakproof  |
| 700701-4510-001 | Sponge Dispenser, Mauve, Leakproof |
| 700701-8510-001 | Sponge Dispenser, White, Leakproof |

## Germicide Soaking Trays

- Designed for soaking instruments in a variety of solutions including phenols, glutaraldehydes or other sterilizing and disinfecting solutions
- Perforated inner tray has a heavy duty handle that enables secure handling when transferring instruments
- Inner tray automatically lowers into the solution as the lid is closed and raises to drain when the lid is opened
- Easy pour corners allow for clean disposal of expired solution
- Stainless steel lids and plastic lids are available
- Available in a variety of colors and 2 sizes



### Germicide Tray with Stainless Steel Lid

Tray measures: 5" × 10-3/4" × 3" (12.7 cm × 27.3 cm × 7.6 cm)

| REF             | Description          |
|-----------------|----------------------|
| 700225-1010-001 | Standard Size, Beige |
| 700225-2010-001 | Standard Size, Blue  |
| 700225-3510-001 | Standard Size, Gray  |
| 700225-4510-001 | Standard Size, Mauve |
| 700225-8510-001 | Standard Size, White |

### Germicide Tray with Plastic Lid

Tray measures: 4-1/2" × 10-3/8" × 3" (11.4 cm × 26.4 cm × 7.6 cm)

| REF             | Description          |
|-----------------|----------------------|
| 700220-1010-001 | Standard Size, Beige |
| 700220-2010-001 | Standard Size, Blue  |
| 700220-3510-001 | Standard Size, Gray  |
| 700220-4510-001 | Standard Size, Mauve |
| 700220-8510-001 | Standard Size, White |



### Germicide Tray with Plastic Lid

Tray measures: 3-3/4" × 8-1/8" × 2-3/8" (9.5 cm × 20.6 cm × 6.1 cm)

| REF             | Description       |
|-----------------|-------------------|
| 700226-1008-001 | Small Size, Beige |
| 700226-2008-001 | Small Size, Blue  |
| 700226-3508-001 | Small Size, Gray  |
| 700226-4508-001 | Small Size, Mauve |
| 700226-8508-001 | Small Size, White |





## Code Rings Small

- Used to color code standard hand instrument handles sizes 2-6 (e.g. mirror handles, scalers, curettes, etc.)
- Code rings are round and are made from silicone
- Available in a variety of colors and package configurations
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes

| Color     | Box of 50       | Box of 100      |
|-----------|-----------------|-----------------|
| Baby Blue | 700301-0508-001 | 700303-0508-001 |
| Beige     | 700301-1008-001 | 700303-1008-001 |
| Black     | 700301-1508-001 | 700303-1508-001 |
| Blue      | 700301-2008-001 | 700303-2008-001 |
| Brown     | 700301-2508-001 | 700303-2508-001 |
| Dark Rose | 700301-3008-001 | 700303-3008-001 |
| Gray      | 700301-3508-001 | 700303-3508-001 |
| Green     | 700301-4008-001 | 700303-4008-001 |
| Mauve     | 700301-4508-001 | 700303-4508-001 |
| Orange    | 700301-5008-001 | 700303-5008-001 |
| Pink      | 700301-5508-001 | 700303-5508-001 |
| Red       | 700301-6508-001 | 700303-6508-001 |
| Teal      | 700301-7508-001 | 700303-7508-001 |
| White     | 700301-8508-001 | 700303-8508-001 |
| Yellow    | 700301-9008-001 | 700303-9008-001 |



## Code Rings Small - Assorted

| REF             | Description                             |
|-----------------|-----------------------------------------|
| 700304-9908-001 | Box of 120, (12 Colors x 10 Rings)      |
| 700305-9908-001 | Box of 175, (7 Colors x 25 Code Rings)  |
| 700306-9908-001 | Box of 300, (12 Colors x 25 Code Rings) |
| 700307-9908-001 | Box of 600, (12 Colors x 50 Code Rings) |

## Code Rings Medium

- Used to color code larger hand instruments (e.g. extracting forceps, ortho pliers, etc.)
- Code rings are round and are made from silicone
- Available in a variety of colors and package configurations
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes


| Color     | Box of 50       | Box of 100      |
|-----------|-----------------|-----------------|
| Baby Blue | 700308-0509-001 | 700310-0509-001 |
| Beige     | 700308-1009-001 | 700310-1009-001 |
| Black     | 700308-1509-001 | 700310-1509-001 |
| Blue      | 700308-2009-001 | 700310-2009-001 |
| Brown     | 700308-2509-001 | 700310-2509-001 |
| Dark Rose | 700308-3009-001 | 700310-3009-001 |
| Gray      | 700308-3509-001 | 700310-3509-001 |
| Green     | 700308-4009-001 | 700310-4009-001 |
| Mauve     | 700308-4509-001 | 700310-4509-001 |
| Orange    | 700308-5009-001 | 700310-5009-001 |
| Pink      | 700308-5509-001 | 700310-5509-001 |
| Red       | 700308-6509-001 | 700310-6509-001 |
| Teal      | 700308-7509-001 | 700310-7509-001 |
| White     | 700308-8509-001 | 700310-8509-001 |
| Yellow    | 700308-9009-001 | 700310-9009-001 |



## Code Rings Medium - Assorted

| REF             | Description                             |
|-----------------|-----------------------------------------|
| 700308-9909-001 | Box of 50                               |
| 700311-9909-001 | Box of 120, (12 Colors × 10 Rings)      |
| 700312-9909-001 | Box of 175, (7 Colors × 25 Code Rings)  |
| 700313-9909-001 | Box of 300, (12 Colors × 25 Code Rings) |
| 700314-9909-001 | Box of 600, (12 Colors × 50 Code Rings) |





## X-Ray Film Holders

- Holds x-ray film for anterior and posterior x-rays
- Both ends securely hold film without bending
- Available with an angulated end or with a straight end
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes

## X-Ray Film Holders

| REF             | Description – Pack of 3           |
|-----------------|-----------------------------------|
| 700800-1000-001 | X-Ray Film Holder, Beige, Angled  |
| 700801-2000-001 | X-Ray Film Holder, Blue, Straight |

## Bib Clips

- Helps keep napkin bib in place during dental procedures
- Sold in packs of 3



700902-2000-001

| REF             | Description – Pack of 3 |
|-----------------|-------------------------|
| 700902-1000-001 | Bib Chains - Beige      |
| 700902-2000-001 | Bib Chains - Blue       |
| 700902-3000-001 | Bib Chains - Dusty Rose |
| 700902-6000-001 | Bib Chains - Purple     |
| 700902-4000-001 | Bib Chains - Green      |
| 700902-5000-001 | Bib Chains - White      |
| 700902-7000-001 | Bib Chains - Red        |

## Mouth Props – Black Rubber

- Designed to help position the mouth open during dental procedures
- Made from black synthetic rubber and natural rubber proteins, they are LATEX FREE
- Adult and medium child sizes feature a stainless steel insert for superior protection against deformation and to provide patient comfort
- Available in three different sizes:  
Small Child, Medium Child, Adult (Large)
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes  
(Not recommended to sterilize by dry heat)



| Descriptions                                                   | Single Mouth Prop | Box of 2 Mouth Props | Box of 3 Mouth Props |
|----------------------------------------------------------------|-------------------|----------------------|----------------------|
| Mouth Prop - Rubber, Large                                     | 700600-1510-001   | 700601-1510-001      | 700603-1510-001      |
| Mouth Prop - Rubber, Medium                                    | 700600-1509-001   | 700601-1509-001      | 700603-1509-001      |
| Mouth Prop - Rubber, Small                                     | 700600-1508-001   | 700601-1508-001      | 700603-1508-001      |
| Mouth Prop - Rubber, Assorted,<br>(1 small, 1 Medium, 1 Large) | -                 | -                    | 700603-1589-001      |

## Mouth Props – Tan Silicone

- Designed to help position the mouth open during dental procedures
- Made from silicone, they are LATEX FREE
- Adult and medium child sizes feature a stainless steel insert for superior protection against deformation and to provide patient comfort
- Available in three different sizes:  
Small Child, Medium Child, Adult (Large)
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes



| Descriptions                  | Single Mouth Prop | Box of 2 Mouth Props |
|-------------------------------|-------------------|----------------------|
| Mouth Prop - Silicone, Large  | 700605-7210-001   | 700606-7210-001      |
| Mouth Prop - Silicone, Medium | 700605-7209-001   | 700606-7209-001      |
| Mouth Prop - Silicone, Small  | 700605-7208-001   | 700606-7208-001      |

## Note: mouth props are also available with a chain

| Descriptions                   | REF     |
|--------------------------------|---------|
| Mouth Prop - Large with chain  | PM-0058 |
| Mouth Prop - Medium with chain | PM-0062 |
| Mouth Prop - Small with chain  | PM-0064 |



Rhodium coated Front Surface Mirrors with cone socket





Image Front Surface Mirrors with simple stem



Plane Surface Mirrors - silver backed with simple stem;  
size 6 shown with cone socket

## Dental Mouth Mirrors

- Designed for visualization of hard to see areas of the mouth
- All mirror shells and shanks are made from high quality stainless steel
- The stems are triple welded to the shells which helps ensure a long instrument life
- Available in the following sizes: 3 (21 mm), 4 (23 mm), 5 (25 mm) and 6 (27 mm)
- Available in the following surfaces:
  - Rhodium coated – Front Surface
  - Image – Front Surface
  - Plane (Silver backed) – Rear Surface

### Front Surface Mirror - Rhodium Coating

| Size      | Cone Socket<br>Box of 12 | Simple Stem<br>Box of 12 |
|-----------|--------------------------|--------------------------|
| 3 (21 mm) | 49001021                 | 49001017                 |
| 4 (23 mm) | 49001022                 | 49001018                 |
| 5 (25 mm) | 49001023                 | 49001019                 |
| 6 (27 mm) | 49001024                 | 49001020                 |

### Front Surface Mirror - Image Mirrors

| Size      | Cone Socket<br>Box of 12 | Simple Stem<br>Box of 12 |
|-----------|--------------------------|--------------------------|
| 3 (21 mm) | 49001013                 | 49001009                 |
| 4 (23 mm) | 49001014                 | 49001010                 |
| 5 (25 mm) | 49001015                 | 49001011                 |
| 6 (27 mm) | 49001016                 | 49001012                 |

### Plane Surface Mirrors - Silver Backed

| Size      | Cone Socket<br>Box of 12 | Simple Stem<br>Box of 12 |
|-----------|--------------------------|--------------------------|
| 3 (21 mm) | 49001005                 | 49001001                 |
| 4 (23 mm) | 49001006                 | 49001002                 |
| 5 (25 mm) | 49001007                 | 49001003                 |
| 6 (27 mm) | 49001008                 | 49001004                 |

## Silicone Endo Stops

- Aids in marking measurements on endo files during procedures
- Pre-punched center hole for easy instrument placement
- Available in ISO colors and 2 additional colors

### Endo Stops - Pre-Punched Center Hole

| Color    | Box of 100      | Bag of 100      |
|----------|-----------------|-----------------|
| Assorted | 700502-9900-001 | 700500-9900-001 |
| Black    | 700502-1500-001 | 700500-1500-001 |
| Blue     | 700502-2000-001 | 700500-2000-001 |
| Green    | 700502-4000-001 | 700500-4000-001 |
| Red      | 700502-6500-001 | 700500-6500-001 |
| Rust     | 700502-7000-001 | 700500-7000-001 |
| White    | 700502-8500-001 | 700500-8500-001 |
| Yellow   | 700502-9000-001 | 700500-9000-001 |



## True-Stop Endo Stop Dispenser

- Makes dispensing endo stops easier and provides a convenient storage compartment for endo stops
- Pre-punched center hole for easy instrument placement
- Available in ISO colors and 2 additional colors

| Color  | 1 Dispenser with 250 Silicone Stops |
|--------|-------------------------------------|
| Black  | 700510-1500-001                     |
| Blue   | 700510-2000-001                     |
| Green  | 700510-4000-001                     |
| Red    | 700510-6500-001                     |
| White  | 700510-8500-001                     |
| Yellow | 700510-9000-001                     |
| Rust   | 700510-7000-001                     |





## Precision™ Drills

- Feature a non-aggressive cutting edge to form safe and concentric post space
- Drills are manufactured from stainless steel
- The shank of the drill has a colored silicone ring for easy size identification
- Precision drills are compatible with the ParaPost®\* system
- Sterilizable by standard table top sterilizer:  
Temperature: 250° F (121° C)  
Exposure time: 30 minutes

\*ParaPost® is a registered trademark of Whaledent International



| Size | Pack of 2       |
|------|-----------------|
| #3   | 700005-5001-001 |
| #4   | 700005-5502-001 |
| #4.5 | 700005-3503-001 |
| #5   | 700005-2504-001 |
| #5.5 | 700005-4505-001 |
| #6   | 700005-8506-001 |
| #7   | 700005-2007-001 |

## Precision™ Post Kit

Kit contains

- 7 Drills
- 21 Posts (3 of each size)
- Vial of endo stops



## Precision™ Post Kit

| REF             | Descriptions               |
|-----------------|----------------------------|
| 700001-0000-001 | Post Kit - Stainless Steel |
| 700002-0000-001 | Post Kit - Titanium        |

## Precision™ Posts

- Are designed as a parallel sided/flute style for maximum retention when placed into a root canal
- Vented to reduce pressure and stress during placement
- Available in a medical grade stainless steel or a high strength corrosion resistant titanium
- Precision posts are compatible with the ParaPost®\* system
- Posts are color coded to match the Precision™ drills color coding



\*ParaPost® is a registered trademark of Whaledent International

## Precision™ Stainless Steel Post

| Size                               | Pack of 10      |
|------------------------------------|-----------------|
| #3 (.036 diameter × .9mm long)     | 700003-0001-001 |
| #4 (.040 diameter × 1.0mm long)    | 700003-0002-001 |
| #4.5 (.045 diameter × 1.12mm long) | 700003-0003-001 |
| #5 (.050 diameter × 1.25mm long)   | 700003-0004-001 |
| #5.5 (.055 diameter × 1.37mm long) | 700003-0005-001 |
| #6 (.060 diameter × 1.50mm long)   | 700003-0006-001 |
| #7 (.070 diameter × 1.75mm long)   | 700003-0007-001 |



## Precision™ Titanium Post

| Size                               | Pack of 10      |
|------------------------------------|-----------------|
| #3 (.036 diameter × .9mm long)     | 700001-0001-001 |
| #4 (.040 diameter × 1.0mm long)    | 700001-0002-001 |
| #4.5 (.045 diameter × 1.12mm long) | 700001-0003-001 |
| #5 (.050 diameter × 1.25mm long)   | 700001-0004-001 |
| #5.5 (.055 diameter × 1.37mm long) | 700001-0005-001 |
| #6 (.060 diameter × 1.50mm long)   | 700001-0006-001 |
| #7 (.070 diameter × 1.75mm long)   | 700001-0007-001 |





700-1021

## Gutta Percha

- Used in the obturation stage of a root canal procedure
- ISO tips and tapers available
- Color coding by ISO sizes
- Hand-rolled with a 4 year shelf life from date rolled
- Available in standard and non-standard (accessory) sizes
- Available in a variety of package configurations



700-121



700-321

## Gutta Percha - 0.02 mm Taper

| Size   | Slide Box<br>120/Box<br>Color Coding | Slide Box<br>120/Box<br>Color Coding<br>Length Marked | Slide Box<br>120/Box | Slide Box<br>100/Box<br>Color Coding | Slide Box<br>100/Box | Round Vial<br>120/Box<br>Color Coding | Square Vial<br>120/Box<br>Color Coding | D-4<br>100/Box |
|--------|--------------------------------------|-------------------------------------------------------|----------------------|--------------------------------------|----------------------|---------------------------------------|----------------------------------------|----------------|
| 8      | 700-101                              | 700-301                                               | 700-401              | 700-201                              | -                    | 700-901                               | 700-1001                               | -              |
| 10     | 700-102                              | 700-302                                               | 700-402              | 700-202                              | -                    | 700-902                               | 700-1002                               | -              |
| 15     | 700-103                              | 700-303                                               | 700-403              | 700-203                              | -                    | 700-903                               | 700-1003                               | -              |
| 20     | 700-104                              | 700-304                                               | 700-404              | 700-204                              | -                    | 700-904                               | 700-1004                               | -              |
| 25     | 700-105                              | 700-305                                               | 700-405              | 700-205                              | -                    | 700-905                               | 700-1005                               | -              |
| 30     | 700-106                              | 700-306                                               | 700-406              | 700-206                              | -                    | 700-906                               | 700-1006                               | -              |
| 35     | 700-107                              | 700-307                                               | 700-407              | 700-207                              | -                    | 700-907                               | 700-1007                               | -              |
| 40     | 700-108                              | 700-308                                               | 700-408              | 700-208                              | -                    | 700-908                               | 700-1008                               | -              |
| 45     | 700-109                              | 700-309                                               | 700-409              | 700-209                              | -                    | 700-909                               | 700-1009                               | -              |
| 50     | 700-110                              | 700-310                                               | 700-410              | 700-210                              | -                    | 700-910                               | 700-1010                               | -              |
| 55     | 700-111                              | 700-311                                               | 700-411              | 700-211                              | -                    | 700-911                               | 700-1011                               | -              |
| 60     | 700-112                              | 700-312                                               | 700-412              | 700-212                              | -                    | 700-912                               | 700-1012                               | -              |
| 70     | 700-113                              | 700-313                                               | 700-413              | 700-213                              | -                    | 700-913                               | 700-1013                               | -              |
| 80     | 700-114                              | 700-314                                               | 700-414              | 700-214                              | -                    | 700-914                               | 700-1014                               | -              |
| 90     | 700-115                              | 700-315                                               | 700-415              | 700-215                              | -                    | -                                     | 700-1015                               | -              |
| 100    | 700-116                              | 700-316                                               | 700-416              | 700-216                              | -                    | -                                     | 700-1016                               | -              |
| 110    | 700-117                              | 700-317                                               | 700-417              | 700-217                              | -                    | -                                     | 700-1017                               | -              |
| 120    | 700-118                              | 700-318                                               | 700-418              | 700-218                              | -                    | -                                     | 700-1018                               | -              |
| 130    | 700-119                              | 700-319                                               | 700-419              | 700-219                              | -                    | -                                     | 700-1019                               | -              |
| 140    | 700-120                              | 700-320                                               | 700-420              | 700-220                              | -                    | -                                     | 700-1020                               | -              |
| 15-40  | 700-121                              | 700-321                                               | 700-421              | 700-221                              | -                    | 700-921                               | 700-1021                               | -              |
| 45-80  | 700-122                              | 700-322                                               | 700-422              | 700-222                              | -                    | 700-922                               | 700-1022                               | -              |
| 90-140 | 700-123                              | 700-323                                               | 700-423              | 700-223                              | -                    | -                                     | 700-1023                               | -              |
| XF     | -                                    | -                                                     | 700-424              | -                                    | 700-224              | -                                     | 700-1024                               | 700-1201       |
| F      | -                                    | -                                                     | 700-425              | -                                    | 700-225              | -                                     | 700-1025                               | 700-1202       |
| MF     | -                                    | -                                                     | 700-426              | -                                    | 700-226              | -                                     | 700-1026                               | 700-1203       |
| FF     | -                                    | -                                                     | 700-427              | -                                    | 700-227              | -                                     | 700-1027                               | 700-1204       |
| FM     | -                                    | -                                                     | 700-428              | -                                    | 700-228              | -                                     | 700-1028                               | 700-1205       |
| MF     | -                                    | -                                                     | 700-429              | -                                    | 700-229              | -                                     | 700-1029                               | 700-1206       |
| ML     | -                                    | -                                                     | 700-430              | -                                    | 700-230              | -                                     | 700-1030                               | 700-1207       |
| L      | -                                    | -                                                     | 700-431              | -                                    | 700-231              | -                                     | 700-1031                               | 700-1208       |
| XL     | -                                    | -                                                     | 700-432              | -                                    | 700-232              | -                                     | 700-1032                               | 700-1209       |
| XF-M   | -                                    | -                                                     | 700-433              | -                                    | 700-233              | -                                     | 700-1033                               | -              |

Note: Endoaid also available #15/80 color coded, 240/Box: 700-1101


## Gutta Percha

Color coding by ISO sizes  
for easy reference

Length control marks  
provide quick visual  
depth reference

ISO-sized in three tapers  
(.02/.04-.06)

(Enlarged to show detail)



700-613

## Gutta Percha - 0.04 mm Taper

| Size        | Slide Box<br>60/Box<br>Color Coding | Slide box<br>60/Box | Slide box<br>100/Box<br>Color Coding | Slide box<br>100/Box |
|-------------|-------------------------------------|---------------------|--------------------------------------|----------------------|
| .04, #15    | 700-501                             | 700-801             | 700-601                              | 700-701              |
| .04, #20    | 700-502                             | 700-802             | 700-602                              | 700-702              |
| .04, #25    | 700-503                             | 700-803             | 700-603                              | 700-703              |
| .04, #30    | 700-504                             | 700-804             | 700-604                              | 700-704              |
| .04, #35    | 700-505                             | 700-805             | 700-605                              | 700-705              |
| .04, #40    | 700-506                             | 700-806             | 700-606                              | 700-706              |
| .04, #45    | 700-507                             | 700-807             | 700-607                              | 700-707              |
| .04, #50    | 700-508                             | 700-808             | 700-608                              | 700-708              |
| .04, #55    | 700-509                             | 700-809             | 700-609                              | 700-709              |
| .04, #60    | 700-510                             | 700-810             | 700-610                              | 700-710              |
| .04, #70    | 700-511                             | 700-811             | 700-611                              | 700-711              |
| .04, #80    | 700-512                             | 700-812             | 700-612                              | 700-712              |
| .04, #15-40 | 700-513                             | 700-813             | 700-613                              | 700-713              |
| .04, #45-80 | 700-514                             | 700-814             | 700-614                              | 700-714              |



700-705



700-713

## Gutta Percha - 0.06 mm Taper

| Size        | Slide Box<br>60/Box<br>Color Coding | Slide box<br>60/Box | Slide box<br>100/Box<br>Color Coding | Slide box<br>100/Box |
|-------------|-------------------------------------|---------------------|--------------------------------------|----------------------|
| .06, #15    | 700-515                             | 700-715             | 700-615                              | 700-815              |
| .06, #20    | 700-516                             | 700-716             | 700-616                              | 700-816              |
| .06, #25    | 700-517                             | 700-717             | 700-617                              | 700-817              |
| .06, #30    | 700-518                             | 700-718             | 700-618                              | 700-818              |
| .06, #35    | 700-519                             | 700-719             | 700-619                              | 700-819              |
| .06, #40    | 700-520                             | 700-720             | 700-620                              | 700-820              |
| .06, #45    | 700-521                             | 700-721             | 700-621                              | 700-821              |
| .06, #50    | 700-522                             | 700-722             | 700-622                              | 700-822              |
| .06, #55    | 700-523                             | 700-723             | 700-623                              | 700-823              |
| .06, #60    | 700-524                             | 700-724             | 700-624                              | 700-824              |
| .06, #70    | 700-525                             | 700-725             | 700-625                              | 700-825              |
| .06, #80    | 700-526                             | 700-726             | 700-626                              | 700-826              |
| .06, #15-40 | 700-527                             | 700-727             | 700-627                              | 700-827              |
| .06, #45-80 | 700-528                             | 700-728             | 700-628                              | 700-828              |



700-827



800-129

## Paper Points

- Used for drying a root canal during the procedure
- ISO tips and tapers available
- Color coding by ISO sizes
- Hand-rolled and sold sterile with a 3 year shelf life
- Available in standard and non-standard (accessory) sizes
- Available in a variety of package configurations



800-303



800-508



800-325

| Size    | Cell Pack<br>200/Box<br>Color Coding | Cell Pack,180/Bx<br>180/Box<br>Color Coding | Slide Box<br>200/Box<br>Color Coding | Round Vial<br>200/Box<br>Color Coding | Square Vial<br>200/Box<br>Color Coding | Slide Box<br>200/Box<br>Color Coding<br>Length Marked |
|---------|--------------------------------------|---------------------------------------------|--------------------------------------|---------------------------------------|----------------------------------------|-------------------------------------------------------|
| #08     | 800-101                              | 800-201                                     | 800-301                              | 800-401                               | 800-501                                | -                                                     |
| #10     | 800-102                              | 800-202                                     | 800-302                              | 800-402                               | 800-502                                | -                                                     |
| #15     | 800-103                              | 800-203                                     | 800-303                              | 800-403                               | 800-503                                | 800-601                                               |
| #20     | 800-104                              | 800-204                                     | 800-304                              | 800-404                               | 800-504                                | 800-602                                               |
| #25     | 800-105                              | 800-205                                     | 800-305                              | 800-405                               | 800-505                                | 800-603                                               |
| #30     | 800-106                              | 800-206                                     | 800-306                              | 800-406                               | 800-506                                | 800-604                                               |
| #35     | 800-107                              | 800-207                                     | 800-307                              | 800-407                               | 800-507                                | 800-605                                               |
| #40     | 800-108                              | 800-208                                     | 800-308                              | 800-408                               | 800-508                                | 800-606                                               |
| #45     | 800-109                              | 800-209                                     | 800-309                              | 800-409                               | 800-509                                | 800-607                                               |
| #50     | 800-110                              | 800-210                                     | 800-310                              | 800-410                               | 800-510                                | 800-608                                               |
| #55     | 800-111                              | 800-211                                     | 800-311                              | 800-411                               | 800-511                                | 800-609                                               |
| #60     | 800-112                              | 800-212                                     | 800-312                              | 800-412                               | 800-512                                | 800-610                                               |
| #70     | 800-113                              | 800-213                                     | 800-313                              | 800-413                               | 800-513                                | 800-611                                               |
| #80     | 800-114                              | 800-214                                     | 800-314                              | 800-414                               | 800-514                                | 800-612                                               |
| #90     | 800-115                              | 800-215                                     | 800-315                              | 800-415                               | 800-515                                | -                                                     |
| #100    | 800-116                              | 800-216                                     | 800-316                              | 800-416                               | 800-516                                | -                                                     |
| #110    | 800-117                              | 800-217                                     | 800-317                              | 800-417                               | 800-517                                | -                                                     |
| #120    | 800-118                              | 800-218                                     | 800-318                              | 800-418                               | 800-518                                | -                                                     |
| #130    | 800-119                              | 800-219                                     | 800-319                              | 800-419                               | 800-519                                | -                                                     |
| #140    | 800-120                              | 800-220                                     | 800-320                              | 800-420                               | 800-520                                | -                                                     |
| #15/40  | 800-121                              | 800-221                                     | 800-321                              | 800-421                               | 800-521                                | 800-613                                               |
| #45/80  | 800-122                              | 800-222                                     | 800-322                              | 800-422                               | 800-522                                | 800-614                                               |
| #90/140 | 800-123                              | 800-223                                     | 800-323                              | 800-423                               | 800-523                                | -                                                     |
| #XXF    | 800-124                              | 800-224                                     | 800-324                              | 800-424                               | 800-524                                | -                                                     |
| #XF     | 800-125                              | 800-225                                     | 800-325                              | 800-425                               | 800-525                                | -                                                     |
| #F      | 800-126                              | 800-226                                     | 800-326                              | 800-426                               | 800-526                                | -                                                     |
| #M      | 800-127                              | 800-227                                     | 800-327                              | 800-427                               | 800-527                                | -                                                     |
| #C      | 800-128                              | 800-228                                     | 800-328                              | 800-428                               | 800-528                                | -                                                     |
| #XC     | 800-129                              | 800-229                                     | 800-329                              | 800-429                               | 800-529                                | -                                                     |
| #XXF/XC | 800-130                              | 800-230                                     | 800-330                              | 800-430                               | 800-530                                | -                                                     |
| #FM     | -                                    | -                                           | 800-331                              | -                                     | -                                      | -                                                     |
| #MF     | -                                    | -                                           | 800-332                              | -                                     | -                                      | -                                                     |

Note: Endoaid also available #15/80 without color coding in the following sizes: 800-701 (200/box), 800-702 (360/box), and 800-703 (400/box)

## INDEX

|              |    |              |    |              |    |              |    |
|--------------|----|--------------|----|--------------|----|--------------|----|
| 700-101..... | 16 | 700-232..... | 16 | 700-430..... | 16 | 700-623..... | 17 |
| 700-102..... | 16 | 700-233..... | 16 | 700-431..... | 16 | 700-624..... | 17 |
| 700-103..... | 16 | 700-301..... | 16 | 700-432..... | 16 | 700-625..... | 17 |
| 700-104..... | 16 | 700-302..... | 16 | 700-433..... | 16 | 700-626..... | 17 |
| 700-105..... | 16 | 700-303..... | 16 | 700-501..... | 17 | 700-627..... | 17 |
| 700-106..... | 16 | 700-304..... | 16 | 700-502..... | 17 | 700-628..... | 17 |
| 700-107..... | 16 | 700-305..... | 16 | 700-503..... | 17 | 700-701..... | 17 |
| 700-108..... | 16 | 700-306..... | 16 | 700-504..... | 17 | 700-702..... | 17 |
| 700-109..... | 16 | 700-307..... | 16 | 700-505..... | 17 | 700-703..... | 17 |
| 700-110..... | 16 | 700-308..... | 16 | 700-506..... | 17 | 700-704..... | 17 |
| 700-111..... | 16 | 700-309..... | 16 | 700-507..... | 17 | 700-705..... | 17 |
| 700-112..... | 16 | 700-310..... | 16 | 700-508..... | 17 | 700-706..... | 17 |
| 700-113..... | 16 | 700-311..... | 16 | 700-509..... | 17 | 700-707..... | 17 |
| 700-114..... | 16 | 700-312..... | 16 | 700-510..... | 17 | 700-708..... | 17 |
| 700-115..... | 16 | 700-313..... | 16 | 700-511..... | 17 | 700-709..... | 17 |
| 700-116..... | 16 | 700-314..... | 16 | 700-512..... | 17 | 700-710..... | 17 |
| 700-117..... | 16 | 700-315..... | 16 | 700-513..... | 17 | 700-711..... | 17 |
| 700-118..... | 16 | 700-316..... | 16 | 700-514..... | 17 | 700-712..... | 17 |
| 700-119..... | 16 | 700-317..... | 16 | 700-515..... | 17 | 700-713..... | 17 |
| 700-120..... | 16 | 700-318..... | 16 | 700-516..... | 17 | 700-714..... | 17 |
| 700-121..... | 16 | 700-319..... | 16 | 700-517..... | 17 | 700-715..... | 17 |
| 700-122..... | 16 | 700-320..... | 16 | 700-518..... | 17 | 700-716..... | 17 |
| 700-123..... | 16 | 700-321..... | 16 | 700-519..... | 17 | 700-717..... | 17 |
| 700-201..... | 16 | 700-322..... | 16 | 700-520..... | 17 | 700-718..... | 17 |
| 700-202..... | 16 | 700-323..... | 16 | 700-521..... | 17 | 700-719..... | 17 |
| 700-203..... | 16 | 700-401..... | 16 | 700-522..... | 17 | 700-720..... | 17 |
| 700-204..... | 16 | 700-402..... | 16 | 700-523..... | 17 | 700-721..... | 17 |
| 700-205..... | 16 | 700-403..... | 16 | 700-524..... | 17 | 700-722..... | 17 |
| 700-206..... | 16 | 700-404..... | 16 | 700-525..... | 17 | 700-723..... | 17 |
| 700-207..... | 16 | 700-405..... | 16 | 700-526..... | 17 | 700-724..... | 17 |
| 700-208..... | 16 | 700-406..... | 16 | 700-527..... | 17 | 700-725..... | 17 |
| 700-209..... | 16 | 700-407..... | 16 | 700-528..... | 17 | 700-726..... | 17 |
| 700-210..... | 16 | 700-408..... | 16 | 700-601..... | 17 | 700-727..... | 17 |
| 700-211..... | 16 | 700-409..... | 16 | 700-602..... | 17 | 700-728..... | 17 |
| 700-212..... | 16 | 700-410..... | 16 | 700-603..... | 17 | 700-801..... | 17 |
| 700-213..... | 16 | 700-411..... | 16 | 700-604..... | 17 | 700-802..... | 17 |
| 700-214..... | 16 | 700-412..... | 16 | 700-605..... | 17 | 700-803..... | 17 |
| 700-215..... | 16 | 700-413..... | 16 | 700-606..... | 17 | 700-804..... | 17 |
| 700-216..... | 16 | 700-414..... | 16 | 700-607..... | 17 | 700-805..... | 17 |
| 700-217..... | 16 | 700-415..... | 16 | 700-608..... | 17 | 700-806..... | 17 |
| 700-218..... | 16 | 700-416..... | 16 | 700-609..... | 17 | 700-807..... | 17 |
| 700-219..... | 16 | 700-417..... | 16 | 700-610..... | 17 | 700-808..... | 17 |
| 700-220..... | 16 | 700-418..... | 16 | 700-611..... | 17 | 700-809..... | 17 |
| 700-221..... | 16 | 700-419..... | 16 | 700-612..... | 17 | 700-810..... | 17 |
| 700-222..... | 16 | 700-420..... | 16 | 700-613..... | 17 | 700-811..... | 17 |
| 700-223..... | 16 | 700-421..... | 16 | 700-614..... | 17 | 700-812..... | 17 |
| 700-224..... | 16 | 700-422..... | 16 | 700-615..... | 17 | 700-813..... | 17 |
| 700-225..... | 16 | 700-423..... | 16 | 700-616..... | 17 | 700-814..... | 17 |
| 700-226..... | 16 | 700-424..... | 16 | 700-617..... | 17 | 700-815..... | 17 |
| 700-227..... | 16 | 700-425..... | 16 | 700-618..... | 17 | 700-816..... | 17 |
| 700-228..... | 16 | 700-426..... | 16 | 700-619..... | 17 | 700-817..... | 17 |
| 700-229..... | 16 | 700-427..... | 16 | 700-620..... | 17 | 700-818..... | 17 |
| 700-230..... | 16 | 700-428..... | 16 | 700-621..... | 17 | 700-819..... | 17 |
| 700-231..... | 16 | 700-429..... | 16 | 700-622..... | 17 | 700-820..... | 17 |

|               |    |               |    |              |    |              |    |
|---------------|----|---------------|----|--------------|----|--------------|----|
| 700-821.....  | 17 | 700-1202..... | 16 | 800-222..... | 18 | 800-419..... | 18 |
| 700-822.....  | 17 | 700-1203..... | 16 | 800-223..... | 18 | 800-420..... | 18 |
| 700-823.....  | 17 | 700-1204..... | 16 | 800-224..... | 18 | 800-421..... | 18 |
| 700-824.....  | 17 | 700-1205..... | 16 | 800-225..... | 18 | 800-422..... | 18 |
| 700-825.....  | 17 | 700-1206..... | 16 | 800-226..... | 18 | 800-423..... | 18 |
| 700-826.....  | 17 | 700-1207..... | 16 | 800-227..... | 18 | 800-424..... | 18 |
| 700-827.....  | 17 | 700-1208..... | 16 | 800-228..... | 18 | 800-425..... | 18 |
| 700-828.....  | 17 | 700-1209..... | 16 | 800-229..... | 18 | 800-426..... | 18 |
| 700-901.....  | 16 | 800-101.....  | 18 | 800-230..... | 18 | 800-427..... | 18 |
| 700-902.....  | 16 | 800-102.....  | 18 | 800-301..... | 18 | 800-428..... | 18 |
| 700-903.....  | 16 | 800-103.....  | 18 | 800-302..... | 18 | 800-429..... | 18 |
| 700-904.....  | 16 | 800-104.....  | 18 | 800-303..... | 18 | 800-430..... | 18 |
| 700-905.....  | 16 | 800-105.....  | 18 | 800-304..... | 18 | 800-501..... | 18 |
| 700-906.....  | 16 | 800-106.....  | 18 | 800-305..... | 18 | 800-502..... | 18 |
| 700-907.....  | 16 | 800-107.....  | 18 | 800-306..... | 18 | 800-503..... | 18 |
| 700-908.....  | 16 | 800-108.....  | 18 | 800-307..... | 18 | 800-504..... | 18 |
| 700-909.....  | 16 | 800-109.....  | 18 | 800-308..... | 18 | 800-505..... | 18 |
| 700-910.....  | 16 | 800-110.....  | 18 | 800-309..... | 18 | 800-506..... | 18 |
| 700-911.....  | 16 | 800-111.....  | 18 | 800-310..... | 18 | 800-507..... | 18 |
| 700-912.....  | 16 | 800-112.....  | 18 | 800-311..... | 18 | 800-508..... | 18 |
| 700-913.....  | 16 | 800-113.....  | 18 | 800-312..... | 18 | 800-509..... | 18 |
| 700-914.....  | 16 | 800-114.....  | 18 | 800-313..... | 18 | 800-510..... | 18 |
| 700-921.....  | 16 | 800-115.....  | 18 | 800-314..... | 18 | 800-511..... | 18 |
| 700-922.....  | 16 | 800-116.....  | 18 | 800-315..... | 18 | 800-512..... | 18 |
| 700-1001..... | 16 | 800-117.....  | 18 | 800-316..... | 18 | 800-513..... | 18 |
| 700-1002..... | 16 | 800-118.....  | 18 | 800-317..... | 18 | 800-514..... | 18 |
| 700-1003..... | 16 | 800-119.....  | 18 | 800-318..... | 18 | 800-515..... | 18 |
| 700-1004..... | 16 | 800-120.....  | 18 | 800-319..... | 18 | 800-516..... | 18 |
| 700-1005..... | 16 | 800-121.....  | 18 | 800-320..... | 18 | 800-517..... | 18 |
| 700-1006..... | 16 | 800-122.....  | 18 | 800-321..... | 18 | 800-518..... | 18 |
| 700-1007..... | 16 | 800-123.....  | 18 | 800-322..... | 18 | 800-519..... | 18 |
| 700-1008..... | 16 | 800-124.....  | 18 | 800-323..... | 18 | 800-520..... | 18 |
| 700-1009..... | 16 | 800-125.....  | 18 | 800-324..... | 18 | 800-521..... | 18 |
| 700-1010..... | 16 | 800-126.....  | 18 | 800-325..... | 18 | 800-522..... | 18 |
| 700-1011..... | 16 | 800-127.....  | 18 | 800-326..... | 18 | 800-523..... | 18 |
| 700-1012..... | 16 | 800-128.....  | 18 | 800-327..... | 18 | 800-524..... | 18 |
| 700-1013..... | 16 | 800-129.....  | 18 | 800-328..... | 18 | 800-525..... | 18 |
| 700-1014..... | 16 | 800-130.....  | 18 | 800-329..... | 18 | 800-526..... | 18 |
| 700-1015..... | 16 | 800-201.....  | 18 | 800-330..... | 18 | 800-527..... | 18 |
| 700-1016..... | 16 | 800-202.....  | 18 | 800-331..... | 18 | 800-528..... | 18 |
| 700-1017..... | 16 | 800-203.....  | 18 | 800-332..... | 18 | 800-529..... | 18 |
| 700-1018..... | 16 | 800-204.....  | 18 | 800-401..... | 18 | 800-530..... | 18 |
| 700-1019..... | 16 | 800-205.....  | 18 | 800-402..... | 18 | 800-601..... | 18 |
| 700-1020..... | 16 | 800-206.....  | 18 | 800-403..... | 18 | 800-602..... | 18 |
| 700-1021..... | 16 | 800-207.....  | 18 | 800-404..... | 18 | 800-603..... | 18 |
| 700-1022..... | 16 | 800-208.....  | 18 | 800-405..... | 18 | 800-604..... | 18 |
| 700-1023..... | 16 | 800-209.....  | 18 | 800-406..... | 18 | 800-605..... | 18 |
| 700-1024..... | 16 | 800-210.....  | 18 | 800-407..... | 18 | 800-606..... | 18 |
| 700-1025..... | 16 | 800-211.....  | 18 | 800-408..... | 18 | 800-607..... | 18 |
| 700-1026..... | 16 | 800-212.....  | 18 | 800-409..... | 18 | 800-608..... | 18 |
| 700-1027..... | 16 | 800-213.....  | 18 | 800-410..... | 18 | 800-609..... | 18 |
| 700-1028..... | 16 | 800-214.....  | 18 | 800-411..... | 18 | 800-610..... | 18 |
| 700-1029..... | 16 | 800-215.....  | 18 | 800-412..... | 18 | 800-611..... | 18 |
| 700-1030..... | 16 | 800-216.....  | 18 | 800-413..... | 18 | 800-612..... | 18 |
| 700-1031..... | 16 | 800-217.....  | 18 | 800-414..... | 18 | 800-613..... | 18 |
| 700-1032..... | 16 | 800-218.....  | 18 | 800-415..... | 18 | 800-614..... | 18 |
| 700-1033..... | 16 | 800-219.....  | 18 | 800-416..... | 18 | 800-701..... | 18 |
| 700-1101..... | 16 | 800-220.....  | 18 | 800-417..... | 18 | 800-702..... | 18 |
| 700-1201..... | 16 | 800-221.....  | 18 | 800-418..... | 18 | 800-703..... | 18 |

|                       |    |                       |   |                       |    |                       |    |
|-----------------------|----|-----------------------|---|-----------------------|----|-----------------------|----|
| 700001-0000-001 ..... | 14 | 700226-2008-001 ..... | 7 | 700308-8509-001 ..... | 9  | 700601-1510-001 ..... | 11 |
| 700001-0001-001 ..... | 15 | 700226-3508-001 ..... | 7 | 700308-9009-001 ..... | 9  | 700603-1508-001 ..... | 11 |
| 700001-0002-001 ..... | 15 | 700226-4508-001 ..... | 7 | 700308-9909-001 ..... | 9  | 700603-1509-001 ..... | 11 |
| 700001-0003-001 ..... | 15 | 700226-8508-001 ..... | 7 | 700310-0509-001 ..... | 9  | 700603-1510-001 ..... | 11 |
| 700001-0004-001 ..... | 15 | 700230-1000-001 ..... | 4 | 700310-1009-001 ..... | 9  | 700603-1589-001 ..... | 11 |
| 700001-0005-001 ..... | 15 | 700230-2000-001 ..... | 4 | 700310-1509-001 ..... | 9  | 700605-7208-001 ..... | 11 |
| 700001-0006-001 ..... | 15 | 700230-4500-001 ..... | 4 | 700310-2009-001 ..... | 9  | 700605-7209-001 ..... | 11 |
| 700001-0007-001 ..... | 15 | 700230-8500-001 ..... | 4 | 700310-2509-001 ..... | 9  | 700605-7210-001 ..... | 11 |
| 700002-0000-001 ..... | 14 | 700235-1000-001 ..... | 4 | 700310-3009-001 ..... | 9  | 700606-7208-001 ..... | 11 |
| 700003-0001-001 ..... | 15 | 700235-2000-001 ..... | 4 | 700310-3509-001 ..... | 9  | 700606-7209-001 ..... | 11 |
| 700003-0002-001 ..... | 15 | 700235-4500-001 ..... | 4 | 700310-4009-001 ..... | 9  | 700606-7210-001 ..... | 11 |
| 700003-0003-001 ..... | 15 | 700235-8500-001 ..... | 4 | 700310-4509-001 ..... | 9  | 700700-1008-001 ..... | 6  |
| 700003-0004-001 ..... | 15 | 700301-0508-001 ..... | 8 | 700310-5009-001 ..... | 9  | 700700-2008-001 ..... | 6  |
| 700003-0005-001 ..... | 15 | 700301-1008-001 ..... | 8 | 700310-5509-001 ..... | 9  | 700700-4508-001 ..... | 6  |
| 700003-0006-001 ..... | 15 | 700301-1508-001 ..... | 8 | 700310-6509-001 ..... | 9  | 700700-8508-001 ..... | 6  |
| 700003-0007-001 ..... | 15 | 700301-2008-001 ..... | 8 | 700310-7509-001 ..... | 9  | 700701-1010-001 ..... | 6  |
| 700005-2007-001 ..... | 14 | 700301-2508-001 ..... | 8 | 700310-8509-001 ..... | 9  | 700701-2010-001 ..... | 6  |
| 700005-2504-001 ..... | 14 | 700301-3008-001 ..... | 8 | 700310-9009-001 ..... | 9  | 700701-4510-001 ..... | 6  |
| 700005-3503-001 ..... | 14 | 700301-3508-001 ..... | 8 | 700311-9909-001 ..... | 9  | 700701-8510-001 ..... | 6  |
| 700005-4505-001 ..... | 14 | 700301-4008-001 ..... | 8 | 700312-9909-001 ..... | 9  | 700800-1000-001 ..... | 10 |
| 700005-5001-001 ..... | 14 | 700301-4508-001 ..... | 8 | 700313-9909-001 ..... | 9  | 700801-2000-001 ..... | 10 |
| 700005-5502-001 ..... | 14 | 700301-5008-001 ..... | 8 | 700314-9909-001 ..... | 9  | 700902-2000-001 ..... | 10 |
| 700005-8506-001 ..... | 14 | 700301-5508-001 ..... | 8 | 700400-1000-001 ..... | 5  | 700902-3000-001 ..... | 10 |
| 700100-0000-001 ..... | 2  | 700301-6508-001 ..... | 8 | 700400-2000-001 ..... | 5  | 700902-4000-001 ..... | 10 |
| 700210-1000-001 ..... | 2  | 700301-7508-001 ..... | 8 | 700400-4500-001 ..... | 5  | 700902-5000-001 ..... | 10 |
| 700210-2000-001 ..... | 2  | 700301-8508-001 ..... | 8 | 700400-8500-001 ..... | 5  | 700902-6000-001 ..... | 10 |
| 700210-4500-001 ..... | 2  | 700301-9008-001 ..... | 8 | 700410-1000-001 ..... | 5  | 700902-7000-001 ..... | 10 |
| 700210-8500-001 ..... | 2  | 700303-0508-001 ..... | 8 | 700410-2000-001 ..... | 5  | 49001001 .....        | 12 |
| 700211-1000-001 ..... | 3  | 700303-1008-001 ..... | 8 | 700410-3500-001 ..... | 5  | 49001002 .....        | 12 |
| 700211-2000-001 ..... | 3  | 700303-1508-001 ..... | 8 | 700410-4500-001 ..... | 5  | 49001003 .....        | 12 |
| 700211-4500-001 ..... | 3  | 700303-2008-001 ..... | 8 | 700410-8500-001 ..... | 5  | 49001004 .....        | 12 |
| 700211-8500-001 ..... | 3  | 700303-2508-001 ..... | 8 | 700500-1500-001 ..... | 13 | 49001005 .....        | 12 |
| 700212-1000-001 ..... | 3  | 700303-3008-001 ..... | 8 | 700500-2000-001 ..... | 13 | 49001006 .....        | 12 |
| 700212-2000-001 ..... | 3  | 700303-3508-001 ..... | 8 | 700500-4000-001 ..... | 13 | 49001007 .....        | 12 |
| 700212-4500-001 ..... | 3  | 700303-4008-001 ..... | 8 | 700500-6500-001 ..... | 13 | 49001008 .....        | 12 |
| 700212-8500-001 ..... | 3  | 700303-4508-001 ..... | 8 | 700500-7000-001 ..... | 13 | 49001009 .....        | 12 |
| 700213-1000-001 ..... | 3  | 700303-5008-001 ..... | 8 | 700500-8500-001 ..... | 13 | 49001010 .....        | 12 |
| 700213-2000-001 ..... | 3  | 700303-5508-001 ..... | 8 | 700500-9000-001 ..... | 13 | 49001011 .....        | 12 |
| 700213-4500-001 ..... | 3  | 700303-6508-001 ..... | 8 | 700500-9900-001 ..... | 13 | 49001012 .....        | 12 |
| 700213-8500-001 ..... | 3  | 700303-7508-001 ..... | 8 | 700502-1500-001 ..... | 13 | 49001013 .....        | 12 |
| 700214-1000-001 ..... | 3  | 700303-8508-001 ..... | 8 | 700502-2000-001 ..... | 13 | 49001014 .....        | 12 |
| 700214-2000-001 ..... | 3  | 700303-9008-001 ..... | 8 | 700502-4000-001 ..... | 13 | 49001015 .....        | 12 |
| 700214-4500-001 ..... | 3  | 700304-9908-001 ..... | 8 | 700502-6500-001 ..... | 13 | 49001016 .....        | 12 |
| 700214-8500-001 ..... | 3  | 700305-9908-001 ..... | 8 | 700502-7000-001 ..... | 13 | 49001017 .....        | 12 |
| 700215-1000-001 ..... | 2  | 700306-9908-001 ..... | 8 | 700502-8500-001 ..... | 13 | 49001018 .....        | 12 |
| 700215-2000-001 ..... | 2  | 700307-9908-001 ..... | 8 | 700502-9000-001 ..... | 13 | 49001019 .....        | 12 |
| 700215-4500-001 ..... | 2  | 700308-0509-001 ..... | 9 | 700502-9900-001 ..... | 13 | 49001020 .....        | 12 |
| 700215-8500-001 ..... | 2  | 700308-1009-001 ..... | 9 | 700510-1500-001 ..... | 13 | 49001021 .....        | 12 |
| 700220-1010-001 ..... | 7  | 700308-1509-001 ..... | 9 | 700510-2000-001 ..... | 13 | 49001022 .....        | 12 |
| 700220-2010-001 ..... | 7  | 700308-2009-001 ..... | 9 | 700510-4000-001 ..... | 13 | 49001023 .....        | 12 |
| 700220-3510-001 ..... | 7  | 700308-2509-001 ..... | 9 | 700510-6500-001 ..... | 13 | 49001024 .....        | 12 |
| 700220-4510-001 ..... | 7  | 700308-3009-001 ..... | 9 | 700510-7000-001 ..... | 13 |                       |    |
| 700220-8510-001 ..... | 7  | 700308-3509-001 ..... | 9 | 700510-8500-001 ..... | 13 |                       |    |
| 700225-1010-001 ..... | 7  | 700308-4009-001 ..... | 9 | 700510-9000-001 ..... | 13 |                       |    |
| 700225-2010-001 ..... | 7  | 700308-4509-001 ..... | 9 | 700600-1508-001 ..... | 11 |                       |    |
| 700225-3510-001 ..... | 7  | 700308-5009-001 ..... | 9 | 700600-1509-001 ..... | 11 |                       |    |
| 700225-4510-001 ..... | 7  | 700308-5509-001 ..... | 9 | 700600-1510-001 ..... | 11 |                       |    |
| 700225-8510-001 ..... | 7  | 700308-6509-001 ..... | 9 | 700601-1508-001 ..... | 11 |                       |    |
| 700226-1008-001 ..... | 7  | 700308-7509-001 ..... | 9 | 700601-1509-001 ..... | 11 |                       |    |





## MILTEX® WARRANTY

### INSTRUMENTS AND PRODUCTS

Miltex hereby warrants all Miltex® products purchased directly from Miltex or from an Authorized Dealer in accordance with the following terms and conditions:

All Miltex® instruments and products are guaranteed to be free from defects in material and workmanship when used for their intended surgical or dental purpose. Any Miltex® product that proves defective in workmanship or material will either be repaired or replaced without charge, at Miltex's discretion. Products not purchased from an authorized source may not be authentic or may not meet Miltex's quality standards and may have been tampered with.

#### Other Instrument and Product Warranties:

- . Union Broach® Chiropractic Mallet – 90 days
- . Union Broach® Electrical Devices – 1 year except for Sprint™ TCM – 2 years
- . Cassettes - 4-year warranty on all moving parts and lifetime guarantee against corrosion

### REPAIRS

Instruments and products repaired by the Miltex® Instrument Repair Service are guaranteed to be free from defects in material and workmanship for 90 days when used for their intended surgical or dental purpose. Any Miltex repair that proves defective in workmanship or material within this 90-day period will either be repaired again or replaced without charge, at Miltex's discretion. Instruments must be cleaned and sterilized prior to returning to Miltex for repair.

### WARRANTY EXCEPTIONS

This warranty is void for Miltex® instruments and products serviced by any person or by any facility other than the Miltex® Instrument Repair Services. Any custom made, custom etched or neutralized products are not returnable. Warranties are not valid for Miltex® instruments and products that prove defective as a result of improper care and cleaning or misuse. Instruments or products that are damaged in fire, flood or other acts or disasters are also excluded from warranty. Except as set forth in this warranty, Miltex makes no other representations and warranties regarding its instruments or products, either expressed or implied. In addition, in no event will Miltex be liable for any incidental, special or exemplary damages or loss of profits in connection with the use of Miltex® instruments and products.



[www.miltex.com](http://www.miltex.com)  
toll-free phone 866-854-8300  
toll-free fax 866-854-8400  
phone 717-840-9335  
fax 717-840-9347

**Please call Customer  
Service today for additional  
information, an in-office  
demonstration, to place an  
order, or to receive other  
Miltex catalogs.**

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation.  
Precise and Precision are trademarks of Miltex, Inc.  
The Integra wave logo is a trademark of Integra LifeSciences Corporation.  
Copyright© 2009 Integra LifeSciences Corporation

