

There Are Two Ways To Reduce Your Exposure...

One Actually Works.



Ultra Safety Plus XL Safety Syringe

Ultra Safety Plus XL safety syringes feature a patented system with a sliding sheath that locks into place for disposing of contaminated needles. This system allows for use of multiple cartridges and has been proven in dental clinics around the world.

Hope will not protect you and your staff, but Ultra Safety Plus XL will reduce your risk. Your team deserves the best needle safety system available today, exclusively from Septodont.



Now Available!
Pre-sterilized Single Use Plastic Handles

*Some lead...
others follow.*



Features and Benefits

Bevel Indicator - Orient your needle quickly with our new bevel indicator.

Transparent Barrel - Provides easy inspection of solution during aspiration.

Unique Sliding Protective Sheath - No longer necessary to unthread contaminated needle, just slide the protective sheath to cover the needle.

Reloadable - Accepts multiple cartridges for single patient use.

Re-engineered locking mechanism - Improves stability and safety.

Autoclavable Handle - Reusable handle - no additional waste.

Pre-sterilized Single Use Handle - individually packed - protection against cross-contamination

Indications

Ultra Safety Plus XL is used in the routine administration of dental local anesthetics.

Order Information

Ultra Safety Plus XL	Item#
25 gauge long (red), 100 injection units & 1 syringe handle	01-N4252
27 gauge short (orange), 100 injection units & 1 syringe handle	01-N4271
27 gauge long (yellow), 100 injection units & 1 syringe handle	01-N4272
30 gauge extra short (purple), 100 injection units & 1 syringe handle	01-N4300
30 gauge short (blue), 100 injection units & 1 syringe handle	01-N4301
Safety Plus XL Syringe Handle Plastic	01-N3000
Safety Plus XL Syringe Handle Pre-sterilized single use, 50 handles, individually packed	01-N3400

To order Septodont Ultra Safety Plus XL, contact your dental dealer or call Septodont at 800-872-8305.

