

Easy to handle pen
style design

Only .06ml of local
anesthetic with
each "click"

Uses standard disposable
dental needles

Compact. Sleek. Discreet.



Paroject® Intraligamental Syringe

Paroject is a intraligamental syringe used for "PDL" injections. This pen style syringe is compact, durable and most of all easy to use. Its unique design makes it versatile and less intimidating to the patient.

Paroject:

- Permits slow and controlled injection
- Releases calibrated doses - injects only .06 ml of local anesthetic with each "click" of the activator level
- Design reduces risk of cartridge breakage

*Some lead...
others follow.*



From the manufacturers of Septocaine®

800-872-8305 • www.septodontusa.com

Features and Benefits

“Pen” style activator

Injects calibrated .06ml dose of local anesthetic with each “click” of the activator level

Permanent hub design

Indications

Paroject is used to perform “PDL” injections.

Order Information

Description	Item#
Paroject Intraligamental Syringe	01N2040

To order Septodont Paroject syringe, contact your dental dealer or call Septodont at 800-872-8305.

