



BLACK LINE ENDO COLLECTION

DESIGNED FOR PREMIUM PERFORMANCE

How the best perform



HU-FRIEDY'S ENDODONTIC INSTRUMENTATION

Hu-Friedy's comprehensive line of endodontic instruments and accessories give you versatile options to meet all of your clinical needs. If you are looking for exceptional instrument design, quality and reliability, and unmatched customer service, choose Hu-Friedy's Endodontic Instruments.

The Challenge

The sensitive nature of endodontic procedures necessitates instruments with exceptional response, integrity and durability to accommodate a wide variety of endodontic procedures – from routine root canals to complex microsurgery.

The Solution

Hu-Friedy instruments are individually handcrafted to precise specifications from quality materials to ensure superior performance, longevity and reliability. Instrument design combined with proprietary materials and incomparable workmanship result in instruments with unique balance, working leverage and control. They function as a seamless extension of the practitioner's skill.

The Result

High-quality Hu-Friedy instrumentation provides unsurpassed clinical efficiency, so practitioners can deliver superior outcomes and focus on the patient. The result is outstanding patient care.



BLACK LINE FOR ENDO

Featuring a performance engineered coating to enhance lubricity, reduce light reflection under harsh lighting, and optimize edge retention, the Black Line Endo instruments are designed to provide exceptional clinical outcomes and improved efficiency during endodontic procedures.

AT INSTRUMENT TIP

- Aluminum Titanium Nitride (AlTiN) coating offers an exceptional non-stick surface when performing vertical or lateral condensation
- Reduced light reflection and enhanced contrast afforded by a matte black finish

AT HANDLE

- The unique smooth, large, lightweight handle reduces fatigue and offers an increased level of control



EXPLORERS

Hu-Friedy's endodontic explorers are made from a special blend of alloy which provides the flexibility, strength, tactile sensitivity and resilience clinicians need when exploring and assessing the treatment site. They are available in several double-ended designs and handles.




DG16
| EXDG16X

DG16/17
| EXDG16/17X

DG16/23
| EXDG16/23X

EXCAVATORS

Hu-Friedy's endodontic excavators feature extended-length shanks and are meticulously handcrafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin. The result is a superior clinical outcome and improved efficiency.



31L
| EXC31LX


31W
Oval Spoon
| EXC31WX

32L
| EXC32LX

33L
| EXC33LX



SPREADERS

Hu-Friedy manufactures a wide variety of spreaders, all handcrafted from the finest stainless steel to precise specifications. Our spreaders help the clinician gain access into the most difficult canals for optimal compaction of filling material during lateral condensation.



PLUGGERS

Hu-Friedy's pluggers aid clinicians in performing vertical condensation as they are handcrafted to exact specifications using the finest materials.



BLACK LINE RUBBER DAM CLAMPS

Hu-Friedy rubber dam clamps are made in the USA from proprietary Satin Steel for corrosion resistance, superior strength, and a fine balance between flexibility and memory. The black matte finish on all 48 styles absorbs and dissipates light instead of reflecting it back, so you can experience a comfortable field of vision that is easier on the eyes – to help you perform at your best.



POINTS OF PERFORMANCE

- Reduced light reflection and enhanced contrast
- Improved visual acuity
- Corrosion resistant

PART CODE	DESCRIPTION
RDCM12AX	Rubber Dam Clamp #12A, Black Line
RDCM13AX	Rubber Dam Clamp #13A, Black Line
RDCM14AX	Rubber Dam Clamp #14A, Black Line
RDCM2AX	Rubber Dam Clamp #2A, Black Line
RDCM14X	Rubber Dam Clamp #14, Black Line
RDCM9X	Rubber Dam Clamp #9, Black Line
RDCM2X	Rubber Dam Clamp #2, Black Line
RDCM00X	Rubber Dam Clamp #00, Black Line
RDCM7X	Rubber Dam Clamp #7, Black Line
RDCM212X	Rubber Dam Clamp #212, Black Line
RDCM8AX	Rubber Dam Clamp #8A, Black Line
RDCMW8AX	Rubber Dam Clamp #W8A, Black Line
RDCM1X	Rubber Dam Clamp #1, Black Line
RDCM0X	Rubber Dam Clamp #0, Black Line
RDCM8X	Rubber Dam Clamp #8, Black Line
RDCM27NX	Rubber Dam Clamp #27N, Black Line
RDCMW3X	Rubber Dam Clamp #W3, Black Line
RDCM3X	Rubber Dam Clamp #3, Black Line
RDCM27X	Rubber Dam Clamp #27, Black Line
RDCMW14AX	Rubber Dam Clamp #W14A, Black Line
RDCM26NX	Rubber Dam Clamp #26N, Black Line
RDCM4X	Rubber Dam Clamp #4, Black Line
RDCM6X	Rubber Dam Clamp #6, Black Line
RDCMW2X	Rubber Dam Clamp #W2, Black Line

PART CODE	DESCRIPTION
RDCM2ASX	Rubber Dam Clamp #2AS, Black Line
RDCM56X	Rubber Dam Clamp #56, Black Line
RDCM201X	Rubber Dam Clamp #201, Black Line
RDCM18X	Rubber Dam Clamp #18, Black Line
RDCM7AX	Rubber Dam Clamp #7A, Black Line
RDCM56SX	Rubber Dam Clamp #56S, Black Line
RDCM205X	Rubber Dam Clamp #205, Black Line
RDCM5X	Rubber Dam Clamp #5, Black Line
RDCM8ADX	Rubber Dam Clamp #8AD, Black Line
RDCM9SX	Rubber Dam Clamp #9S, Black Line
RDCM212SAX	Rubber Dam Clamp #212SA, Black Line
RDCM1AX	Rubber Dam Clamp #1A, Black Line
RDCMW7X	Rubber Dam Clamp #W7, Black Line
RDCM209X	Rubber Dam Clamp #209, Black Line
RDCMW4X	Rubber Dam Clamp #W4, Black Line
RDCM26X	Rubber Dam Clamp #26, Black Line
RDCM24X	Rubber Dam Clamp #24, Black Line
RDCM25X	Rubber Dam Clamp #25, Black Line
RDCM29X	Rubber Dam Clamp #29, Black Line
RDCM10X	Rubber Dam Clamp #10, Black Line
RDCM28X	Rubber Dam Clamp #28, Black Line
RDCM11X	Rubber Dam Clamp #11, Black Line
RDCM31X	Rubber Dam Clamp #31, Black Line
RDCM30X	Rubber Dam Clamp #30, Black Line

PERIODONTAL



RESTORATIVE



ORTHODONTIC



INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



ULTRASONIC
SCALING



For more information visit
Hu-Friedy.com/EndoEssentials

Hu-Friedy Mfg. Co., LLC: 3232 N. Rockwell St. | Chicago, IL 60618 | USA
Telephone: 1-800-Hu-Friedy (1-773-975-6100) | Website: Hu-Friedy.com

©2017 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-656/0317

Hu-Friedy