

Medicut™ Scalpels and Blades



WE CARE LIKE FAMILY



Medicut™ Make The Cut

Medicut™ Scalpels (Disposable & Sterile)

- Disposable scalpels...for maximum safety, efficiency, and economy
- Sterile...individual peel open pouch
- Sheffield stainless steel blade
- Length of the handle is 4.9" (125 mm)

Medicut™ Scalpels

Reorder No.	Size	Qty
4110	#10	10/bx
4111	#11	10/bx
4112	#12	10/bx
4115	#15	10/bx
4120	#20	10/bx
4121	#21	10/bx
4122	#22	10/bx
4123	#23	10/bx

Medicut™ Blades (Surgical Disposable & Sterile)

- For guaranteed sharpness and exacting performance
- Sterile...individually foil packaged
- Sheffield stainless steel blade
- Finely honed cutting edge
- Gamma Sterilized

Medicut™ Blades

Reorder No.	Size	Qty
4130	#10	100/bx
4131	#11	100/bx
4132	#12	100/bx
4135	#15	100/bx
4140	#20	100/bx
4141	#21	100/bx
4142	#22	100/bx
4143	#23	100/bx
4144	#24	100/bx

Medicut™ Safety Scalpels (Sterile)

- Satisfies regulatory guidelines for safe handling of sharp edges
- Blade retracts with one simple motion and can be permanently locked inside the handle for safe disposal
- Sheffield stainless steel blade for exacting performance
- Packaged in individual pouches
- Gamma Sterilized
- Length of the handle is 4.9" (125 mm)

Medicut™ Safety Scalpels

Reorder No.	Size	Qty
4160	#10	10/bx
4161	#11	10/bx
4165	#15	10/bx



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com