

3M Glass Ionomer Solutions

Which is best for your practice?

The unique properties of 3M glass ionomers can fit a wide range of patient situations where simplicity, speed, fluoride release and radiopacity are especially important.

Ketac™ Universal Aplicap™ Glass Ionomer Restorative, our latest addition to the 3M family of glass ionomer restoratives, is fast, strong and simple.

To try it for yourself ... ask your 3M representative for a sample and demonstration.

Ketac™ Universal Aplicap™ Glass Ionomer Restorative

- No need for conditioner or coating
- Bulk fill—no need for layering
- Enables fast, one-step, stress-bearing restorative solutions
- Extended indications for use compared to other glass ionomer materials
- Low stickiness allows an easy placement
- Self adhesive
- Available in 6 shades

RotoMix™ Capsule Mixing Unit,
Aplicap™ Activator, Aplicap™ Applier



Ketac™ Nano Light-Curing Glass Ionomer Restorative

- The esthetic 3M glass ionomer restorative
- Nano technology inspired by Filtek™ Supreme Ultra Universal Restorative
- Low wear and high initial polish in comparison with other glass ionomer restoratives
- Significantly higher gloss reflectance after polishing
- Faster, easier dispensing with the innovative Quick Mix Capsule or Clicker™ Dispenser
- Available in 6 shades

Elipar™ DeepCure-S LED Curing Light,
Aplicap™ Applier



Ordering Information

Ketac™ Universal Aplicap™ Glass Ionomer Restorative

Item #	Shade/Product Description	Item #	Shade/Product Description
	Refills Includes: 50 Aplicap™ capsules		
61085	White, A2, A3	61097	Giant Pack A3 Includes: 300 Aplicap™ capsules
61086	White	73040	Aplicap™ Activator
61087	A1	73050	Aplicap™ Applier
61088	A2	76310	RotoMix™ Capsule Mixing Unit
61089	A3		
61095	A3.5		
61096	A4		

Ketac™ Nano Light-Curing Glass Ionomer Restorative

Item #	Shade/Product Description	Item #	Shade/Product Description
	Refills Includes: 20 Quick Mix capsules		
3305A1	A1	3305TK	Trial Kit—Quick Mix Capsule Kit Includes:
3305A2	A2		10–A2 Quick Mix capsules;
3305A3	A3		10–A3 Quick Mix capsules;
3305A3.5	A3.5		Primer
3305A4	A4		
3305B2	B2		

Saving time and money is important for every practice.

Which 3M glass ionomer restorative is more economical, covering more indications?

	Indications										
											
Product	Class I	Class II	Class III	Class V	Core Build-Up	Base/Liner	Fissure Sealant	Delivery Options	Conditioner/ Primer Required	Glaze Required	Curing
Conventional											
Ketac™ Universal Aplicap™ Glass Ionomer Restorative											Self Cure
Ketac™ Molar Quick Aplicap™ Glass Ionomer Restorative										Optional	Self Cure
Ketac™ Fil Plus Aplicap™ Glass Ionomer Restorative											Self Cure
Metal Reinforced											
Ketac™ Silver Glass Ionomer Restorative								 			Self Cure
Resin Modified											
Ketac™ Nano Light-Curing Glass Ionomer Restorative								 			Light Cure
Photac™ Fil Quick Aplicap™ Light-Curing Glass Ionomer Restorative											Light Cure
Vitremer™ Core Buildup/Restorative											Tri Cure
 Indicated  Indicated with restriction  Quick Mix Capsule  Capsule  Clicker™ Dispenser  Handmix  C P Conditioner/Primer  Glaze											

● Indicated

● Indicated with restriction

 Quick Mix Capsule

 Capsule

 Clicker™ Dispenser

 Handmix

  Conditioner/Primer

 Glaze

Refer to Instructions for Use for complete product information.

For more information on the product you are most interested in, visit

www.3M.com/KetacUniversal



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, Aplicap, Clicker, Ketac, Photac, Vitremer and RotoMix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are not trademarks of 3M. Please recycle. Printed in USA. © 3M 2016. All rights reserved. 70-2013-0772-8