

# Versa-Temp<sup>2</sup>

TEMPORARY CROWN AND BRIDGE RESIN

Up to 70% stronger  
than leading competitors

Up to 42%  
faster set time  
than leading competitors

Easy to adjust  
for beautiful results

Versa-Temp<sup>2</sup>

TEMPORARY CROWN AND BRIDGE RESIN

 **Sultan**Healthcare  
sultanhealthcare.com

The anatomy of a great temporary.


# Versa-Temp<sup>2</sup>

## TEMPORARY CROWN AND BRIDGE RESIN

The temporary is an essential part of the restorative process... important to your patient, and a reflection of the work you and your practice produce. **Versa-Temp<sup>2</sup>** combines the superb performance and natural-looking aesthetics you want—for temporaries that set quickly, handle easily, look great and stay strong until the final crown or bridge arrives.

UP TO  
**70%**  
STRONGER\*  
THAN LEADING COMPETITORS

At 450 MPa compressive strength, **Versa-Temp<sup>2</sup>** is formulated to last—to help ensure your patients don't have to come back before their crown or bridge is ready.



UP TO  
**42%**  
FASTER\*\* SET TIME  
THAN LEADING COMPETITORS

That means less chair time, less expense, and happier patients.

4 minute  
total set time!



**EASY** TO  
**SHAPE**  
WITH BEAUTIFUL RESULTS

**Versa-Temp<sup>2</sup>** makes the temporary process smooth and easy, for aesthetically pleasing temporaries in minutes!

- ✓ Simple to adjust and trim
- ✓ Natural fluorescence
- ✓ Easy to polish with an alcohol wipe (e.g. Isopropyl)



## VERSA-TEMP<sup>2</sup> CASE STUDY

A 45-year-old female patient reported sensitivity on #14 and dental hygiene difficulties on #12 and #13. A temporary matrix was made using **Genie<sup>®</sup> Putty** and **Genie<sup>®</sup> VPS**. A failed restoration was removed and #14 was endodontically treated. A temporary was created using **Versa-Temp<sup>2</sup>** Temporary Crown and Bridge Resin and **SensiTemp<sup>®</sup>** Resin Non-Eugenol Temporary Cement.



Failed restoration on #12 and #13 with endodontic repair needed on #14.

Three-unit temporary made using **Versa-Temp<sup>2</sup>**, Shade A2.

### A great temporary begins with a great matrix

Sultan Healthcare offers a variety of materials that can be used for a temporary matrix, all of which work hand-in-hand with **Versa-Temp<sup>2</sup>**. A great impression means less time adjusting and polishing... and a high-quality temporary.



**SILGIMIX<sup>®</sup>**



**GENIE<sup>®</sup> PUTTY**



**GENIE<sup>®</sup> VPS**

TEMPORARY MATRIX MATERIAL	Alginate Replacement	Putty+VPS	VPS
DESCRIPTION	An alginate replacement combines the low cost of an alginate with the detail and versatility of a VPS.	Using a putty/VPS matrix eliminates the need for an impression tray, saving on cost while still benefiting from high quality details of a VPS.	A dual-phase technique leverages the high quality of a final impression, leading to unmatched detail.
MATERIALS USED	<i>Silgimix<sup>®</sup>; Genie<sup>®</sup> VPS Tray Adhesive; DenTray<sup>®</sup> II</i>	<i>Genie<sup>®</sup> Putty; Genie<sup>®</sup> VPS Light Body</i>	<i>Genie<sup>®</sup> VPS Heavy Body; Genie<sup>®</sup> VPS Light Body; Genie<sup>®</sup> VPS Tray Adhesive; DenTray<sup>®</sup> II</i>

To learn more about **Silgimix<sup>®</sup>** and **Genie<sup>®</sup>**, contact your Sultan Healthcare sales representative or visit [sultanhc.com](http://sultanhc.com)

### Versa-Temp<sup>2</sup> Technical Specifications

3-Point Flexural Strength	141 MPa
3-Point Flexural Strength (10 min)	100 MPa
Compressive Strength	450 MPa
Flexural Modulus	2500 MPa
Edge Stability	1000 MPa
Setting Temperature	Max 3.6 F above mouth temperature

### Versa-Temp<sup>2</sup> Temporary Crown & Bridge Resin

REF #	50ml Cartridge, 15 Mixing Tips	REF #	3 x 50ml Cartridge
<b>72041</b>	Shade BL		
<b>72051</b>	Shade A1	<b>72066</b>	Shade A1
<b>72056</b>	Shade A2	<b>72071</b>	Shade A2
<b>72061</b>	Shade A3	<b>72076</b>	Shade A3



**Sultan** Healthcare

\*Based on compressive strength of leading competitors shown. Data on file.

\*\*Based on manufacturers' directions for use.

Protent<sup>TM</sup> and Luxatemp<sup>TM</sup> are not registered trademarks of Sultan Healthcare.