


# Material Safety Data Sheet

Dr.Jean Bausch KG  
Oskar-Schindler-Str. 4  
50769 Cologne - Germany  
Phone: ++49-221-70936-0  
Fax: ++49-221-70936-66  
Internet:www.bauschdental.de  
eMail: info@bauschdental.de

**US-Sales Office:**  
**Bausch Articulating Papers Inc.**  
**12 Murphy Drive,Unit 4**  
**Nashua, NH 03062 - U.S.A**  
**Phone & Fax: (603) 883-2155**  
**Internet: www.bauschdental.com**  
**eMail: info@bauschdental.com**



## Identity (As Used on Label and List)

**Trade Name:** Bausch Articulating Paper 100  $\mu$  blue with progressive Color transfer

**Code:** BK 51 (precut strips), BK 53 horseshoeshape, BK 57 (precut strips)

**Chemical** natural vegetable waxes, paraffin oil (Pharmaceutical Standard) DAB 6, cellulose (Paper oxygen

**Description:** bleached), C.I.Pigment blue 27

**Product Use:** Dental Material, for Marking of Occlusal Surface

## Section I

|                                                     |                                  |
|-----------------------------------------------------|----------------------------------|
| Manufacturer's Name                                 | Emergency Telephone Number       |
| Dr.Jean Bausch KG                                   | ++49-221-70936-0, (603)883-2155  |
| Address (Number, Street, City, State, and Zip Code) | Telephone Number for Information |
| Oskar-Schindler-Str. 4                              | ++49-221-70936-0, (603)883-2155  |
| 50769 Cologne-Germany                               | Date Prepared                    |
| Bausch Articulating Papers Inc.                     | 11.02.1999                       |
| 12 Murphy Drive, Unit 4                             | Signature of Preparer (optional) |
| Nashua NH, 03062 -U.S.A                             | Peter Bausch                     |

## Section II - Hazardous Ingredients / Identity Information

| Hazardous Components (Specific Chemical Identity; Common Name(s)) | OSHA PEL | ACGIH TLV | Other Limits Recommended | % (optional) |
|-------------------------------------------------------------------|----------|-----------|--------------------------|--------------|
| none                                                              |          |           |                          |              |

## Section III - Physical / Chemical Characteristics

|                         |                                         |                                      |
|-------------------------|-----------------------------------------|--------------------------------------|
| Boiling Point           | Specific Gravity (H <sub>2</sub> O = 1) | Melting Point                        |
| n.a.                    | 0.91 (D)                                |                                      |
| Vapor Pressure (mm Hg.) | Vapor Density (AIR=1)                   | Evaporation Rate (Butyl Acetate = 1) |
| n.a.                    | n.a.                                    | n                                    |
| Solubility in Water     | Appearance and Odor                     |                                      |
| none                    | blue strips of paper                    |                                      |

## Section IV - Fire and Explosion Hazard Data

|                              |                  |      |     |
|------------------------------|------------------|------|-----|
| Flash Point (Method Used)    | Flammable Limits | LEL  | UEL |
| > 300°C (572°F) (Literature) |                  | n.a. | n.a |

Extinguishing Media

Use dry chemicals or carbon dioxide; water may be ineffective, but should be used to keep fire exposed containers cool.

Special Fire Fighting Procedures

not required

Unusual Fire and Explosion Hazards

Firefighters should wear self-contained breathing apparatus in the positive pressure mode with a full facepiece when there is a possibility of exposure to smoke, fumes, or hazardous decomposition products.