

elementsTMfree

OBTURATION SYSTEM


n., Powered obturation with
no strings attached.
(In other words...it's cord-free!)



Dentistry *unleashed.*

Visit elementsfreeendo.com or call 800.KERR.123 for more information.

n. , Powered obturation with
no strings attached.



Downpack



Ready to go when you're ready to go...


- Downpack unit heats to 200°C in less than 0.5 seconds and with an adjustable heating range of 100°-400°C
- 360° activation ring to improve your experience and clinical outcome
- Digital temperature control for procedural accuracy
- Perfect for Continuous Wave™ and single motion downpack obturation

Backfill



Powered obturation, no strings attached...

- Swiss motor and gear drive for exceptional performance
- Motorized extrusion for less hand fatigue
- Single-use cartridges with built-in finger grips for ease of use
- Aerogel insulation surrounds the heating element for your comfort and your patient's safety



Protect your Patients...

Thermal images of three backfill devices, after being heated to 200°C for 5 minutes.

Data on file at SybronEndo

B&L®, Calamus®, and Obtura® are trademarks or registered trademarks of their respective owners.

Ready to give elementsfree a try? Go to...

elementsfreeendo.com

elementsTMfree

n., Powered obturation with
no strings attached



Dual Charging Base



Single Charging Base

PART #	UNITS
973-0500-TypeB	Elements Free System Type B with Dual Charging Base
973-0502-TypeB	Elements Free Downpack Type B with Single Charging Base
973-0504-TypeB	Elements Free Backfill Type B with Single Charging Base

PART #	REPLACEMENT PARTS
973-0516-TypeB	Elements Free Cord Type B
973-0520	Elements Free, Battery
973-0515	Elements Free Transformer
973-0337	Elements Cartridge Bender Tool

PART #	ACCESSORIES
973-0510	Elements Free, Dual Charging Base
973-0512	Elements Free, Single Charging Base
972-2500	Elements Gutta Percha Cartridge 23GA
972-2501	Elements Gutta Percha Cartridge Light 25GA - Pack 10
972-2502	Elements Gutta Percha Cartridge Heavy 23GA - Pack 10
952-0003	Buchanan Heat Pluggers (XF, F, FM, M, ML)
952-0004	Buchanan Heat Plugger Fine .06 Taper
952-0005	Buchanan Heat Plugger Fine Medium .08 Taper
952-0006	Buchanan Heat Plugger Medium .10 Taper
952-0007	Buchanan Heat Plugger Medium-Large .12 Taper
952-0031	Buchanan Heat Plugger Extra Fine .04 Taper
952-0017	System B Tip Stand
974-0058	Buchanan Hand Plugger 1
974-0059	Buchanan Hand Plugger 2
974-0060	Buchanan Hand Plugger 0

Dentistry *unleashed.*

Visit elementsfreeendo.com or call 800.KERR.123 for more information.