



Rotary when you **want it**.
Reciprocation when you **need it**.

Introducing the *NEW* Elements Motor.

Adaptive Motion -

The new Elements Motor adapts the motion based on the amount of pressure on the file. It may or may not reciprocate and does not pull into the canal. Reciprocating angles vary: No load 600° forward and 0° backwards, loaded 370° forward and up to 50° backwards, based on file load.

6 presets -

This motor provides recommended settings for TF Adaptive, K3/K3XF, TF, Lightspeed, M4 and custom settings.

Intuitive User Interface -

With an easy to read LCD screen, the Elements Motor is easy to learn how to use.

8:1 Contra Angle -

Smooth operation. Handpieces are available in latch-pin and push-button.

Fully Autoclavable motor/cable -

Easily detaches from the console.



PART #	DESCRIPTION
815-1500	Elements Motor - includes motor console, motor handpiece, foot pedal, 8:1 push button contra-angle, power cords for N.America & Aus/NZ, IFU.
815-1510	Elements Motor Handpiece
815-1512	Elements Motor Foot Pedal
815-1515	Elements Motor Power Cord

Visit KerrDental.com or call 800.KERR.123 for more information.

ELEMENTS MOTOR



Rotary when you **want it**.
Reciprocation when you **need it**.

Introducing the *NEW* Elements Motor.

Adaptive Motion -

The new Elements Motor adapts the motion based on the amount of pressure on the file. It may or may not reciprocate and does not pull into the canal. Reciprocating angles vary: No load 600° forward and 0° backwards, loaded 370° forward and up to 50° backwards, based on file load.

6 presets -

This motor provides recommended settings for TF Adaptive, K3/K3XF, TF, Lightspeed, M4 and custom settings.

Intuitive User Interface -

With an easy to read LCD screen, the Elements Motor is easy to learn how to use.

8:1 Contra Angle -

Smooth operation. Handpieces are available in latch-pin and push-button.

Fully Autoclavable motor/cable -

Easily detaches from the console.



PART #	DESCRIPTION
815-1500	Elements Motor - includes motor console, motor handpiece, foot pedal, 8:1 push button contra-angle, power cords for N.America & Aus/NZ, IFU.
815-1510	Elements Motor Handpiece
815-1512	Elements Motor Foot Pedal
815-1515	Elements Motor Power Cord

Visit KerrDental.com or call 800.KERR.123 for more information.