

A combo of flexibility and durability resulting in efficiency

TF Adaptive gives you rotary when you want it. Reciprocation when you need it.
And peace of mind throughout.




FLEXIBILITY

Designed with R-Phase technology, TF Adaptive files are up to 70% more flexible than other files on the market. This allows you to comfortably navigate curved canals and decreases the risk of canal transportation when compared to other file systems.

LESS force to bend.

TF Adaptive Files are more flexible than comparable files.
Over 80% less force required to bend the TFA file vs. ProTaper®




Force Required to Bend 45 degrees (g*cm)

*Testing data available upon request

MORE resistant to fracture.

TF Adaptive Files took over 2x more the number of cycles to fracture vs. WaveOne Gold.



DURABILITY

TF Adaptive files are the only endodontic files that are twisted, rather than ground, into their shape. Grinding causes microfractures which increases the likelihood of file separation. Because of their unique design, TF Adaptive files are up to 58% more resistant to fracture than other leading file systems.

TF ADAPTIVE STARTER KIT #815-1535

1



**Elements Motor
with 8:1 Contra Angle
Handpiece**

2



**6 Packs of Assorted
23mm NiTi files:
3-Small &
3 Medium/Large**

3



**4 Packs of Assorted
27mm NiTi Files:
2-Small &
2-Medium/Large**

4



**2 Packs Each of
#8, #10 & #15 K-File
Stainless Steel
Hand Files**

5



LA Axxess™ 2.0 Bur Kit

6



EndoRing

PART # DESCRIPTION

815-1535	TF Adaptive Starter Kit
815-1500	Elements Motor
815-1510	Elements Motor Handpiece
815-1512	Elements Motor Foot Pedal
815-1515	Elements Motor Power Cord
817-3003	TF Adaptive Small Procedure Pack 23mm
817-3007	TF Adaptive Small Procedure Pack 27mm
817-4003	TF Adaptive Med/Lrg Procedure Pack 23mm
817-4007	TF Adaptive Med/Lrg Procedure Pack 27mm
817-4203	TF Adaptive File SM1 23mm
817-4207	TF Adaptive File SM1 27mm
817-6253	TF Adaptive File SM2 23mm
817-6257	TF Adaptive File SM2 27mm
817-4353	TF Adaptive File SM3 23mm
817-4357	TF Adaptive File SM3 27mm
817-8253	TF Adaptive File ML1 23mm
817-8257	TF Adaptive File ML1 27mm
817-6353	TF Adaptive File ML2 23mm
817-6357	TF Adaptive File ML2 27mm
817-4503	TF Adaptive File ML3 23mm
817-4507	TF Adaptive File ML3 27mm
815-1540	TF Adaptive Gutta Percha Small Yellow Pack - SM2
815-1541	TF Adaptive Gutta Percha Small Red Pack - SM3
815-1542	TF Adaptive Gutta Percha Small Asst Pack - SM2, SM3
815-1543	TF Adaptive Gutta Percha Medium/Large Green Pack - ML1
815-1544	TF Adaptive Gutta Percha Medium/Large Yellow Pack - ML2
815-1545	TF Adaptive Gutta Percha Medium/Large Red Pack - ML3
815-1546	TF Adaptive Gutta Percha Medium/Large Asst Pack - ML1, ML2, ML3
815-1550	TF Adaptive Obturator Small Yellow Pack - SM2
815-1551	TF Adaptive Obturator Small Red Pack - SM3
815-1554	TF Adaptive Obturator Medium/Large Green Pack - ML1
815-1555	TF Adaptive Obturator Medium/Large Yellow Pack - ML2
815-1556	TF Adaptive Obturator Medium/Large Red Pack - ML3
815-1560	TF Adaptive Bulk Obturator Small Yellow Pack - SM2
815-1561	TF Adaptive Bulk Obturator Small Red Pack - SM3
815-1562	TF Adaptive Bulk Obturator Medium/Large Green Pack - ML1
815-1563	TF Adaptive Bulk Obturator Medium/Large Yellow Pack - ML2
815-1564	TF Adaptive Bulk Obturator Medium/Large Red Pack - ML3
815-1552	TF Adaptive Verifier Small Yellow Pack - SM2
815-1553	TF Adaptive Verifier Small Red Pack - SM3
815-1557	TF Adaptive Verifier Medium/Large Green Pack - ML1
815-1558	TF Adaptive Verifier Medium/Large Yellow Pack - ML2
815-1559	TF Adaptive Verifier Medium/Large Red Pack - ML3
815-1570	TF Adaptive Paper Points Small Yellow Pack - SM2
815-1571	TF Adaptive Paper Points Small Red Pack - SM3
815-1572	TF Adaptive Paper Points Medium/Large Green Pack - ML1
815-1573	TF Adaptive Paper Points Medium/Large Yellow Pack - ML2
815-1574	TF Adaptive Paper Points Medium/Large Red Pack - ML3
815-1575	TF Adaptive Paper Points Small Asst Pack - SM2, SM3
815-1576	TF Adaptive Paper Points Medium/Large Asst Pack - ML1, ML2, ML3
077-3142	TF Adaptive Technique Card

For more information visit TFAdaptive.com or call us at 800-KERR-123