

TF™ Adaptive

A complete system from Shape to Fill. TF Adaptive offers integrated simplicity and peace of mind.

TF ADAPTIVE



SHAPE

Smart, intuitive instruments—designed with you in mind.

You're in control—TF Adaptive offers an intuitive color-coded system to help you through your procedures with ease. With the safety features of Elements™ Adaptive Motion Technology, you now have both safety and efficiency in your hands.



TF ADAPTIVE



FILL

TF Adaptive Obturators—a perfect match.

Utilizing the same color-coded system as the files, obturating with TF Adaptive couldn't be simpler. The obturators are sized to precisely match the file that you finish with, giving you predictable results every time.



Kerr™
ENDODONTICS

TF ADAPTIVE OBTURATOR KIT

PART # 815-1590



1

Sealapex Xpress

Non-eugenol, calcium hydroxide root canal sealer in automix syringes—no manual mixing.

2

Soft-Core Heater

Specially designed to eliminate the risk of overheating. Now your obturators will be ready when you are.

3

5 Packs of Verifiers

Each pack contains six size verifiers to help ensure a proper fit every time.

4

5 Packs of Obturators

Each pack contains six obturators that provide an exact match to the file system.

PART #	DESCRIPTION	QUANTITY
972-3002	Soft-Core Oven	1
815-1550	TF Adaptive Obturator SM2 (6pk)	1
815-1551	TF Adaptive Obturator SM3 (6pk)	1
815-1554	TF Adaptive Obturator ML1 (6pk)	1
815-1555	TF Adaptive Obturator ML2 (6pk)	1
815-1556	TF Adaptive Obturator ML3 (6pk)	1

PART #	DESCRIPTION	QUANTITY
815-1552	TF Adaptive Verifier SM2 (6pk)	1
815-1553	TF Adaptive Verifier SM3 (6pk)	1
815-1557	TF Adaptive Verifier ML1 (6pk)	1
815-1558	TF Adaptive Verifier ML2 (6pk)	1
815-1559	TF Adaptive Verifier ML3 (6pk)	1
33639	Sealapex Xpress (2 10.5 syringes)	1

For more information visit TFAdaptive.com or call us at 800-Kerr-123