Safety Data Sheet

According to Hazard Communication Standard (29 CFR 1910.1200)

Arti-Spot Frühkontaktindikator BK 87

Issue date: 02/06/2015 Version 1.0 Revision date: 02/06/2015

1. Identification

Product name Arti-Spot Frühkontaktindikator BK 87

Synonyms -

CAS# See section 3

Product code -

Product use Indicator solution. Paint.

Manufacturer/Supplier

Supplier(Manufacturer): Dr. Jean Bausch GmbH & Co. KG

Address: Oskar-Schindler-Str. 4, D-50769 Köln

 Contact person(E-mail):
 info@BauschDental.de

 Telephone:
 +49 (0)221-70936-0

 Fax:
 +49 (0)221-70936-66

Emergency telephone Number: +49 30 19240 (D-13437 Berlin, 24 hour)

2. Hazard(s) identification

GHS classification

Physical hazardsFlammable liquidsCategory 1Health hazardsSerious eye damage/eye irritationCategory 2A

Specific target organ toxicity after single

exposure

Category 3

Environmental hazards Not classified

GHS label elements

Hazard Pictograms



Signal word Danger

Hazard statement Extremely flammable liquid and vapor.

Causes serious eye irritation.

May cause drowsiness or dizziness.

Precautionary statement

Prevention Keep away from heat/sparks/open flames/hot surfaces. — No smoking.

Keep container tightly closed.

Ground/bond container and receiving equipment.

Use explosion-proof electrical/ventilating/lighting/equipment.

Use only non-sparking tools.

Take precautionary measures against static discharge.

Avoid breathing dust/fume/gas/mist/vapors/spray.

Wash thoroughly after handling.

Use only outdoors or in a well-ventilated area.

Wear protective gloves/protective clothing/eye protection/face protection.

Response If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with

water/shower.

If inhaled: Remove person to fresh air and keep comfortable for breathing.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing.

Call a poison center /doctor if you feel unwell.

If eye irritation persists: Get medical advice/attention.

In case of fire: Use CO2, extinction powder, water jet spray, alcohol resistant foam for

extinction.

Storage Store in a well-ventilated place. Keep container tightly closed. Keep cool.

Disposal Dispose of contents/container in corroding with local regulation.

3. Composition / information on ingredients

Components	CAS#	Percent
ethanol	64-17-5	30-40%
ethyl acetate	141-78-6	20-30%
diethyl ether	60-29-7	10-<25%
butanone	78-93-3	0-1%

4. First-aid Measures

First aid procedures

Eye contact Remove contact lenses. Wash thoroughly for several minutes using copious water.

Seek medical help if necessary.

Skin contact Remove polluted, soaked clothing immediately, wash thoroughly with plenty of water

and soap, in case of irritation of the skin (flare), consult a doctor. Unsuitable cleaning

product: Solvent, thinners.

Inhalation Remove person from danger area. Supply person with fresh air and consult doctor

according to symptoms. If the person is unconscious, place in a stable side position

and consult a doctor.

Ingestion Rinse the mouth thoroughly with water. Do not induce vomiting - give copious water to

drink. Consult doctor immediately. Danger of aspiration.

Notes to physician Treat symptoms.

5. Fire-fighting measures

Flammable properties Extremely flammable.

Extinguishing media

Suitable extinguishing media Carbon dioxide, extinction powder, water jet spray, alcohol resistant foam.

Unsuitable extinguishing media High volume water jet.

Firefighting equipment/instructions In case of fire and/or explosion do not breathe fumes. Protective respirator with

independent air supply. According to size of fire Full protection, if necessary. Dispose

of contaminated extinction water according to official regulations.

Hazardous combustion products

Oxides of carbon, oxides of nitrogen, toxic pyrolysis products, explosive vapour/air

Material name: Arti-Spot Frühkontaktindikator BK 87

SDS US

mixture, dangerous vapours heavier than air.

6. Accidental release measures

Personal precautions Remove possible causes of ignition - do not smoke. Ensure sufficient supply of air.

Avoid inhalation, and contact with eyes or skin. If applicable, caution - risk of slipping.

For personal protection see section 8.

Environmental precautions If leakage occurs, dam up. Resolve leaks if this possible without risk. Prevent surface

and ground-water infiltration, as well as ground penetration. Prevent from entering drainage system. If accidental entry into drainage system occurs, inform responsible

authorities.

Methods for cleaning up

Soak up with absorbent material (e.g. universal binding agent, sand, diatomaceous

earth) and dispose of according to Section 13. Use no flammable substances. Fill the absorbed material into lockable containers. Keep moist. Do not let the solution dry up.

Flush residue using copious water.

7. Handling and storage

Handling Avoid inhalation of the vapours. Ensure good ventilation. Keep away from sources of

ignition - Do not smoke. Take measures against electrostatic charging, if appropriate.

Avoid contact with eyes or skin. Eating, drinking, smoking, as well as food-storage, is prohibited in work-room. Observe directions on label and instructions for use. Use

working methods according to operating instructions. For precautions see section 2.2.

Keep out of access to unauthorised individuals. Not to be stored in gangways or stair

wells. Store product closed and only in original packing. Do not store with flammable or

self-igniting materials. Store in a well-ventilated place. Protect from direct sunlight and

warming. Store at room temperature. Do not store over 30°C.

8. Exposure controls / personal protection

Control parameters:

Storage

OCCUPATIONAL EXPOSURE LIMITS (OEL)

INGREDIENT DATA:

Source	Ingredient	TWA	STEL	Peak
US OSHA Permissible	ethanol	1900 mg/m3 /1000 ppm	Not Available	Not Available
Exposure Levels (PELs)				
- Table Z1				
US ACGIH Threshold	ethanol	Not Available	1000 ppm	Not Available
Limit Values (TLV)				
US NIOSH	ethanol	1900 mg/m3 /1000 ppm	Not Available	Not Available
Recommended				
Exposure Limits (RELs)				
US OSHA Permissible	ethyl acetate	1400 mg/m3 /400 ppm	Not Available	Not Available
Exposure Levels (PELs)				
- Table Z1				
US ACGIH Threshold	ethyl acetate	400 ppm	Not Available	Not Available
Limit Values (TLV)				

US NIOSH	ethyl acetate	1400 mg/m3 /400 ppm	Not Available	Not Available
Recommended				
Exposure Limits (RELs)				
US OSHA Permissible	diethyl ether	1200 mg/m3 /400 ppm	Not Available	Not Available
Exposure Levels (PELs)				
- Table Z1				
US ACGIH Threshold	diethyl ether	400 ppm	500 ppm	Not Available
Limit Values (TLV)				
US NIOSH	diethyl ether	Not Available	Not Available	Not Available
Recommended				
Exposure Limits (RELs)				
US OSHA Permissible	butanone	590 mg/m3 /200 ppm	Not Available	Not Available
Exposure Levels (PELs)				
- Table Z1				
US ACGIH Threshold	butanone	200 ppm	300 ppm	Not Available
Limit Values (TLV)				
US NIOSH	butanone	590 mg/m3 /200 ppm	885 mg/m3 /300 ppm	Not Available
Recommended				
Exposure Limits (RELs)				

EMERGENCY LIMITS:

Ingredient	TEEL-1	TEEL-2	TEEL-3
ethanol	Not Available Not Available		Not Available
ethyl acetate	400 ppm	400 ppm	10000 ppm
diethyl ether	500 ppm	ppm 500 ppm 190	
butanone	Not Available	Not Available	Not Available

Ingredient	Original IDLH	Revised IDLH	
ethanol	15,000 ppm	3,300 [LEL] ppm	
ethyl acetate	10,000 ppm	2,000 [LEL] ppm	
diethyl ether	19,000 [LEL] ppm	1,900 [LEL] ppm	
butanone	3,000 ppm	3,000 [Unch] ppm	

Exposure controls:

Appropriate engineering controls:

Ensure good ventilation. This can be achieved by local suction or general air extraction. If this is insufficient to maintain the concentration under the WEL or AGW values, suitable breathing protection should be worn. Applies only if maximum permissible exposure values are listed here.

Individual protection measures, such as personal protective equipment:

Eye / face protection Tight fitting protective goggles (EN 166).

Skin protection Solvent resistant protective gloves (EN 374). If applicable, safety gloves made of butyl (EN

374), protective Neoprene® / polychloroprene gloves (EN 374), protective nitrile gloves (EN 374), protective hand cream recommended. Protective working garments (e.g. safety shoes

EN ISO 20345, long-sleeved protective working garments).

Respiratory protection Normally not necessary. If OES or MEL is exceeded. Gas mask filter A (EN 14387), code colour

brown. Observe wearing time limitations for respiratory protection equipment.

General hygiene Wash hands, forearms and face thoroughly after handling chemical products, before eating, considerations smoking and using the lavatory and at the end of the working period. Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing.

9. Physical and chemical properties

Appearance

Physical stateLiquidFormLiquidColorBlue

Odor Alcoholic, Characteristic

Odor threshold

PH

Not available

Vapor pressure

Not available

Melting point/Freezing point

Not available

initial boiling point and boiling range

Not available

Flash point -28,5 °C (ISO 1523 (Rapid Equilibrium, closed cup, RECC))

Evaporation rate Not available

Flammability (solid, gas) 200 °C (Ignition temperature)

Explosion limits Lower: 2,1 Vol-% upper: 13,5 Vol-%

Vapor density Not available Relative density Not available Solubility (water) Insoluble Not available Partition coefficient Not available Auto-ignition temperature **Decomposition temperature** Not available Not available Specific gravity Density 0,902 g/ml Flammability limits in air, upper, %by volume Not available Flammability limits in air, lower, % by volume Not available VOC Not available Percent volatile Not available

Other data

Viscosity Not available

10. Stability and reactivity

Reactivity Can form explosive peroxides.

Chemical stability Explosive when dry.

Conditions to avoid Incompatible materials. Heating, open flame, ignition sources, electrostatic charge.

Protect from direct sunlight. Product is light sensitive.

Incompatible materials Avoid contact with strong oxidizing agents, strong acids, alkali metals.

Hazardous decomposition products

Oxides of carbon, oxides of nitrogen, toxic pyrolysis products, explosive vapour/air

mixture, dangerous vapours heavier than air.

Possibility of hazardous reactions

Possible build up of explosive/highly flammable vapour/air mixture.

11. Toxicological information

Toxicokinetics, metabolism and distribution:

Non-human toxicological data: Not available

Information on toxicological effects:

Acute toxicity:

ethanol (CAS#: 64-17-5)

 LD50(Oral, Rat):
 10470 mg/kg

 LD50(Dermal, Rabbit):
 >2000 mg/kg

 LC50(Inhalation, Rat):
 117-125 mg/l/4h

Acute toxicity:

ethyl acetate (CAS#:141-78-6)

 LD50(Oral, Rat):
 5620 mg/kg

 LD50(Dermal, Rabbit):
 >18000 mg/kg

 LC50(Inhalation, Rat):
 >28,6 mg/l/4h

Acute toxicity:

diethyl ether (CAS#: 60-29-7)

LD50(Oral, Rat):1215 mg/kgLD50(Dermal, Rabbit):>20000 mg/kgLC50(Inhalation, Rat):>20 mg/l/4hSkin corrosion/Irritation:Not classified.

Serious eye damage/irritation: Causes serious eye irritation.

Respiratory or skin sensitization:Not classifiedGerm cell mutagenicity:Not classifiedCarcinogenicity:Not classifiedReproductive toxicity:Not classified

STOT- single exposure:May cause drowsiness or dizziness.

STOT-repeated exposure: Not classified
Aspiration hazard: Not classified

12. Ecological information

Toxicity:

ethanol (CAS#: 64-17-5)

Acute to	xicity	Time	Species	Method	Evaluation	Remarks
LC50	13000 mg/l	96h	Fish	OECD 203	N/A	N/A
EC50	12340 mg/l	48h	Daphnia	OECD 202	N/A	N/A
EC50	275 mg/l	72h	Algae	OECD 201	N/A	N/A

Toxicity:

ethyl acetate (CAS#:141-78-6)

Acute to	xicity	Time	Species	Method	Evaluation	Remarks
LC50	230 mg/m3	96h	Fish	OECD 203	N/A	N/A
EC50	610 mg/l	48h	Daphnia	OECD 202	N/A	N/A

EC50	N/A	72h	Algae	OECD 201	N/A	N/A	
------	-----	-----	-------	----------	-----	-----	--

Toxicity:

diethyl ether (CAS#: 60-29-7)

Acute to	xicity	Time	Species	Method	Evaluation	Remarks
LC50	2600 mg/l	96h	Fish	OECD 203	N/A	N/A
EC50	165	24h	Daphnia	OECD 202	N/A	N/A
EC50	>100 mg/l	72h	Algae	OECD 201	N/A	N/A

Persistence and degradability: Not available.

Bioaccumulative potential: Not available.

Mobility in soil: The product is Insoluble in water.

Results of PBT&vPvB assessment: Not available.

Other adverse effects: No known significant effects or critical hazards.

13. Disposal considerations

Disposal instructions

Dispose of contents/container in accordance with local/regional/national/international

regulations.

Contaminated packaging Since emptied containers may retain product residue, follow label warnings even after

container is emptied.

14. Transport information

DOT

Basic shipping requirements:

UN number UN1263

Proper shipping name Paint (including paint, lacquer, enamel, stain, shellac, varnish, polish, liquid filler and

liquid lacquer base)

Hazard class 3
Packing group I
Environmental hazards No

IATA

UN number UN1263

UN proper shipping name Paint (including paint, lacquer, enamel, stain, shellac, varnish, polish, liquid filler and

liquid lacquer base)

Transport hazard class(es) 3
Packing group I
Environmental hazards No

IMDG

UN number UN1263

UN proper shipping name Paint (including paint, lacquer, enamel, stain, shellac, varnish, polish, liquid filler and

liquid lacquer base)

Transport hazard class(es) 3
Packing group I
Environmental hazards No

15. Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture:

ethanol (64-17-5) is found on the following	"US - Hawaii Air Contaminant Limits" List.
regulatory lists	"US -Idaho - Limits for Air Contaminants" List.
	"US - Alaska Limits for Air Contaminants" List.
	"US Spacecraft Maximum Allowable Concentrations (SMACs) for Airborne
	Contaminants" List.
	"US Toxic Substances Control Act (TSCA) - Chemical Substance Inventory"
	List.
	"US – Wyoming Toxic and Hazardous Substances Table Z1 Limits for Air
	Contaminants" List.
ethyl acetate (141-78-6) is found on the following	"US - Hawaii Air Contaminant Limits" List.
regulatory lists	"US -Idaho - Limits for Air Contaminants" List.
	"US - Alaska Limits for Air Contaminants" List.
	"US Toxic Substances Control Act (TSCA) - Chemical Substance Inventory"
	List.
	"US – Wyoming Toxic and Hazardous Substances Table Z1 Limits for Air
	Contaminants" List.
diethyl ether(60-29-7) is found on the following	"US - Hawaii Air Contaminant Limits" List.
regulatory lists	"US -Idaho - Limits for Air Contaminants" List.
	"US - Alaska Limits for Air Contaminants" List.
	"US Toxic Substances Control Act (TSCA) - Chemical Substance Inventory"
	List.
	"US – Wyoming Toxic and Hazardous Substances Table Z1 Limits for Air
	Contaminants" List.
butanone (78-93-3) is found on the following	"US - Hawaii Air Contaminant Limits" List.
regulatory lists	"US -Idaho - Limits for Air Contaminants" List.
	"US - Alaska Limits for Air Contaminants" List.
	"US - Washington Toxic air pollutants and their ASIL, SQER and de minimis
	emission values" List.
	"US Spacecraft Maximum Allowable Concentrations (SMACs) for Airborne
	Contaminants" List.
	"US Toxic Substances Control Act (TSCA) - Chemical Substance Inventory"
	List.
	"US – Wyoming Toxic and Hazardous Substances Table Z1 Limits for Air
	Contaminants" List.

16. Other information, including date of preparation or last revision

HMIS®ratings Health: 2

Flammability: 3 Physical hazard: 0

NFPA ratings Health: 2

Flammability: 3

Material name: Arti-Spot Frühkontaktindikator BK 87

Instability: 0

Disclaimer The information in the sheet was written based on the best knowledge and experience

currently available.

Issue date 02-06-2015

Material name: Arti-Spot Frühkontaktindikator BK 87 Version #:1.0 Revision date: 02-06-2015. Issue date: 02-06-2015.