

# **ETHICON PRODUCT CATALOG**

## **Sutures/Adhesives/Drains/Hernia Repair**



# Table of Contents

|   |     |
|---|-----|
| <b>General information</b> .....        | 1   |
| <b>Needle Index</b> .....               | 11  |
| <b>True to Scale Needle Chart</b> ..... | 46  |
| <b>Sutures</b> .....                    | 48  |
| <b>Adhesives</b> .....                  | 235 |
| <b>Adhesives</b> .....                  | 236 |
| <b>Adhesives</b> .....                  | 237 |
| <b>Drains</b> .....                     | 238 |
| <b>Hernia Repair</b> .....              | 241 |
| <b>Miscellaneous</b> .....              | 257 |

## **Explanation of a Typical Selection And Page**

### **ATRALOC Needle**

The needle type and size/shape are shown only for needle/suture combinations. A repeated needle type (eg RB-1, RB-1) indicates a "double-armed" suture, that is, a suture with a needle attached to each end.

### **Needle Type/Needle Shape Symbol**

Needle type and needle shape are shown in this subhead.

### **Material**

Suture materials are designated by color and construction characteristics (eg Undyed Braided, Violet Braided, etc.)

### **Length**

Suture length is shown in both inches and the approximate metric equivalent. Multi-strand means 3 or 8 strands per packet. These products are packaged 12, 24 or 36 packets per box.

The CR mark refers to control release needle sutures.

### **Product Name**

Sutures and other ETHICON products are designated by trademarked and generic names.

### **Size**

### **About ATRALOC Needles**

ATRALOC Needles are supplied in a variety of sizes, shapes and chord lengths. Surgeons should

Suture size, denoting the diameter of the material, is in accordance with the United States Pharmacopeia (USP).

### **Needle Silhouette**

Needle Silhouettes are shown in ACTUAL SIZE. The designated number above the needle, within the Needle Section, is the actual length measured in millimeters.

### **Code**

Product codes ending in suffix G are packaged dozen packets to a box. Product codes ending in suffix H are packaged 3 dozen to a box. Product codes ending in suffix T are packaged 2 dozen packets to a box. In the few instances where no alphabetic suffix is in the product code, the number of packets per box is identified by a footnote. The number of strands per Multi-Strand packet appears on the left of suture length.

### **Needle Length**

Needle length is the actual length measured in millimeters. Needle lengths can be found within the Needle Section and throughout the catalog to the right of the needle code. As a user you will find that many of these needle lengths have been rounded off. Recognize that the needles are not changing, but for consistency ETHICON has applied a rounding convention to all needle lengths. Due to timing differences, the box and label dimensions may not match the catalog exactly, but will be updated accordingly in the near future.

select the type of needle that, in their experience, is appropriate for the specific procedure and

tissue.

Needle shapes are generally classified according to the degree of curvature of the body — 5/8, 1/2, 3/8, or 1/4

circle and straight — with taper, cutting, blunt

points or TAPERCUT\* Surgical Needle tips.

In general, the same size of needle may be made from finer gauge wire for use in soft or delicate tissues and from heavier gauge wire for use in tough or fibrosed tissues (the surgeon's choice).

### **The Major Characteristics of ATRALOC Needles are:**

- They must be made from high quality stainless steel.
- They resist bending but are processed so that they will tend to bend before breaking.
- Taper points must be sharp and contoured for easy passage into tissues.
- Cutting points or edges must be sharp and free of burrs.
- On most ATRALOC Needles, a Super-Smooth finish is provided that permits the needle to penetrate and pass through with minimal resistance or drag.
- Ribbed needles — longitudinal ribs are provided on many ATRALOC Needles to increase the stability of the needle in the needleholder.
- The attachment of the needle to the suture material must be secure so that the needle will not separate from the suture material under normal use. An exception is the CONTROL RELEASE Needle suture principle, which facilitates fast separation of needle from suture when desirable.

## **MULTIPASS\* Needles**

### Elements of Needle Design


Needle design involves analyzing a surgical procedure and the density of the tissue involved in great detail. ETHICON engineers work continuously to improve upon their needle line.

The anatomy of the ideal surgical needle has three key factors that make up the ideal needle. - Alloy - Geometry - tip and body - Coating


The combination of these attributes is ETHICON MULTIPASS Needle Technology.

MULTIPASS Needles are the highest performing needles that ETHICON Products offers; a combination [or a balance] of three proprietary technologies ensure excellent strength, tissue penetration and control; pass after pass.

## **CONTROL RELEASE Needle Sutures**



- The needle is held securely in holder. Suture is grasped securely just below needle, pulling strand taut.



- The needle is released with a straight tug of the needleholder.

## Selection and Use of Needleholders

- The needleholder must be an appropriate size instrument for the size needle selected.
- It must be made from good quality steel with a secure jaw design.



- Needles should be grasped in an area about 1/3 to 1/2 of the distance from the swaged area to the point. Avoid placement on or near the swaged area.




- The needle should be placed securely in the tip of the needleholder jaws.



- Do not damage taper points or cutting edges when using the needleholder to pull the needle out through the tissue. Grasp as far back as possible.
- When placing the needle in tissue, any force applied should be in the direction following the curve of the needle.
- Do not take excessively large bites of tissue with small needles.



- Do not force or twist the needle in an effort to bring the point out through the tissue – withdraw the needle and replace tissue.
- Do not force a dull needle through tissue – obtain a new one.
- Avoid using the needle to bridge or approximate tissues for suturing.
- If the needle is held too tightly in a sharp or hard jawed or defective needleholder, the needle may be damaged or notched in such a manner that it will have more of a tendency to bend or break on successive passes through tissue.



- When the needle holder with needle is passed to the surgeon, the needle should be pointing in the right direction in which it will start to be used without need for readjustments.
- In some patients, the tissues may be tougher or fibrosed more than normal and require the use of a heavier gauge needle.
- In deep confined area, ideal positioning of the needle may not be possible. The surgeon should proceed with caution or use a heavier gauge needle

An ATRALOC Needle is an instrument that complements the surgeon's skills. The variety of ETHICON needles offers the surgeon a choice of needles which are dependable for consistent quality and for practically all surgical applications.

## Steps to Create Performance

### Controlled Manufacturing

ETHICON makes its own needle/sutures combinations. Each step, from raw material to finished needle, is carefully controlled for uniformly superior performance.



Conventional needle



Super smooth

### Strong, Sharp Micrograin Stainless Steel

ETHICON makes its needles from premium grade

stainless steel formulated specifically for surgical use.

The ETHICON process controls the orientation of the micrograins so that ductility and strength are closely controlled – for strong needles that start sharp, and stay sharp, thus assuring easy passage through tissue.

### Super-Smooth

Virtually all ATRALOC Needles have a micro-thin Super-Smooth finish that permits the needle to penetrate and pass through the toughest tissue with an almost fluid-like flow, with minimal resistance and virtually no drag, giving the surgeon a feeling of control.


### Needle Swaging

Needles are permanently attached to the suture materials, eliminating the need for threading and providing a new, sharp needle each time. ATRALOC Needle swaging gives the surgeon a needle with an unbroken surface from end to end; without double strands, it draws through tissue easier and with less tissue trauma.

### Ribbed for Security and Control

Longitudinal ribbing on the inside and outside curvatures of most needles has a cross-locking action in the needleholder, minimizing turning, rocking or twisting and thus providing unsurpassed needle control.


## **Shapes**



## Types


### Taper Point Needle

For soft, easily penetrated tissues.




### ETHIGUARD® Blunt Point Needle

Taper body. For blunt dissection and suturing friable tissue.




### Conventional Cutting Needle

Two opposing cutting edges, with a third on inside curve. Change in cross-section from a triangle cutting tip to a flattened body.




### Reverse Cutting Needle

Cutting edge on outer curve. For tough, difficult-to-penetrate tissues.




### TAPERCUT® Needle

Cutting tip, taper body. For tough tissue, like two needles in one.



### Precision Point Needle


For delicate plastic or cosmetic surgery. Cutting tip electropolished for added sharpness.



### Precision Cosmetic-Conventional Cutting

#### PRIME Needle

For delicate plastic or cosmetic surgery. Conventional cutting tip and PRIME geometry for increased sharpness.



### MICRO-POINT® Surgical Needle (Spatula)

Thin, flat profile. Specially designed for ophthalmic anterior segment surgery.



### Conventional Spatula Needle

Visibility of the point at bottom provides control of depth penetration.

**SABRELOC\* Needle (Spatula)**


Side-cutting spatula shaped edges. For layers of scleral or corneal tissue. Spatulated needle point is centered for maximum needle stability in thin sclera. Four equidistant and properly defined edges provide greater control.

**ULTIMA\* Needle (Spatula)**

Reduced edge-angles provide better penetration.  
Readily facilitates knot rotation in ophthalmic surgery.

**VISI-BLACK\* Surgical Needle**

Slim Taper point needles with a black finish for improved visibility and penetration.

**MICRO-POINT\* Reverse Cutting Needle**

Cutting edge on outer curve. Extremely smooth. Extremely sharp for ophthalmic surgery.



## Trademarks

|   |  |
|---|--|
| ATRALOC Surgical Needle   | PERMA-HAND Silk Suture   |
| CONTROL RELEASE Needle or Needle Suture  | PROCEED Surgical Mesh  |
| CS ULTIMA Ophthalmic Needles  | PROLENE Polypropylene Hernia System  |
| DERMABOND ProPen Topical Skin Adhesive  | PROLENE Polypropylene Mesh   |
| DERMABOND ProPen XL   | PROLENE Polypropylene Suture   |
| DERMABOND Topical Skin Adhesive   | PROLENE 3D Patch Polypropylene Mesh  |
| ETHALLOY Needle Alloy   | PROLENE Soft Polypropylene Mesh  |
| ETHIBOND Extra Polyester Suture   | PRONOVA Poly (Hexafluoropropylene – VDF)   |
| ETHIBOND EXCEL Polyester Suture   | Suture   |
| ETHIGUARD Blunt Point Needle  | PROXI-STRIP Skin Closure   |
| ETHILON Nylon Suture  | PS PRIME Needle  |
| ETHI-PACK Pre-cut Sutures   | PVP™ Device  |
| HEMO-SEAL Needle Suture                  | RELAY Suture Delivery System   |
| J-VAC Closed Wound Drainage System  | RELAY Suture Package   |
| J-VAC Bulb Suction Reservoir  | SABRELOC Spatula Needle  |
| J-VAC Reservoir   | SUTUPAK Pre-cut Sterile Sutures  |
| J-VAC Drainage Adapter  | TAPERCUT Surgical Needle   |
| J-VAC Silicone Drains   | ULTRAPRO Partially Absorbable Lightweight Mesh   |
| LAPRA-TY Suture Clip  | ULTRAPRO Hernia System   |
| LIGAPAK Dispensing Reel   | VICRYL (polyglactin 910) Knitted Mesh  |
| LIGAPAK Ligature  | VICRYL (polyglactin 910) Woven Mesh  |
| MERSILENE Polyester Fiber Suture  | VICRYL (polyglactin 910) Mesh Bags   |
| MERSILENE Polyester Fiber Mesh  | VICRYL (polyglactin 910) Perirenal Mesh Bag  |
| MERSILENE Polyester Fiber Strip   | VICRYL (polyglactin 910) Perisplenic Mesh Bag  |
| MICRO-POINT Surgical Needle   | VICRYL (polyglactin 910) Perihepatic Mesh Bags   |
| MONOCRYL (poliglecaprone 25) Suture   | VICRYL (polyglactin 910) Perihepatic Right Lobe  |
| MONOCRYL Plus Antibacterial (poliglecaprone 25) Suture  | Mesh Bag   |
| MULTIPLASS Needle   | VICRYL (polyglactin 910) Perihepatic Left Lobe   |
| NUROLON Nylon Suture  | Mesh Bag   |
| PC PRIME Needle   | Coated VICRYL (polyglactin 910) Suture   |
| PDS (polydioxanone) Suture  | Coated VICRYL Plus Antibacterial (polyglactin 910) Suture  |
| PDS II (polydioxanone) Suture   | VICRYL RAPIDE (polyglactin 910) Suture   |
| PDS Plus Antibacterial (polydioxanone) Suture   | VISI-BLACK Surgical Needle  |

# Blunt Point

## 1/2 Circle

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>BP</b> 80mm<br><br>Surgical Gut Suture - Chromic                      |    |
| <b>BP</b> 80mm<br><b>BP</b><br><br>MERSILENE Suture                      |    |
| <b>BP-1</b> 65mm<br><br>PDS Plus Suture<br>Surgical Gut Suture - Chromic |   |
| <b>BP-1</b> 65mm<br><b>BP-1</b><br><br>MERSILENE Suture                  |  |
| <b>TPB-1</b> 65mm<br><br>Coated VICRYL Plus Suture<br>PDS Plus Suture    |  |

# Conventional Cutting - Sternum

## 1/2 Circle

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>CCS</b> 48mm<br><br>ETHIBOND EXCEL Suture<br>Surgical Stainless Steel Suture |  |

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>CCS-1</b> 60mm<br><br>Coated VICRYL Suture<br>PDS II Suture<br>Surgical Stainless Steel Suture |  |

## Endoscopic Needles

### Straight Taper Point

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>ST-4</b> 19mm<br><b>ST-4</b><br><br>PDS II Suture<br>PDS Plus Suture |  |
| <b>ST100-3</b><br><br>ETHILON Suture                                    |  |

## ETHIGUARD - Safety needle

### 1/2 Circle

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>CTB</b> 40mm<br><br>Coated VICRYL Plus Suture<br>PDS Plus Suture   |  |
| <b>CTB-1</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture |  |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>CTB-2</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PDS Plus Suture   |   |
| <b>CTX-B</b> 48mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture |   |
| <b>SHB</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS Plus Suture              |  |

## EVERPOINT Cardiovascular Needles

### 1/2 Circle

| Needles/<br>Avialable Products                        | Needle image  |
|---|---|
| <b>RB-2</b> 13mm<br><b>RB-2</b><br><br>PROLENE Suture |  |

### 3/8 Circle

| Needles/<br>Avialable Products                    | Needle image  |
|---|---|
| <b>BV</b> 11mm<br><b>BV</b><br><br>PROLENE Suture |  |

| Needles/<br>Avialable Products                               | Needle image  |
|--|---|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br><br>PROLENE Suture       | <br><br>                                       |
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b><br><br>PROLENE Suture | <br><br>EVERPOINT™<br>Cardiovascular Needles   |
| <b>BV175-6</b> 8.0mm<br><b>BV175-6</b><br><br>PROLENE Suture | <br><br>EVERPOINT™<br>Cardiovascular Needles   |
| <b>BV175-7</b> 8.0mm<br><b>BV175-7</b><br><br>PROLENE Suture | <br><br>  |
| <b>BV175-8</b> 9.3mm<br><b>BV175-8</b><br><br>PROLENE Suture | <br><br>EVERPOINT™<br>Cardiovascular Needles   |
| <b>C-1</b> 13mm<br><b>C-1</b><br><br>PROLENE Suture          | <br><br>EVERPOINT™<br>Cardiovascular Needles |
| <b>C1</b> 13mm<br><br>PROLENE Suture                         | <br><br>EVERPOINT™<br>Cardiovascular Needles |
| <b>CC</b> 9.3mm<br><b>CC</b><br><br>PROLENE Suture           | <br><br>                                     |
| <b>CC-1</b> 13mm<br><b>CC-1</b><br><br>PROLENE Suture        | <br><br>                                     |
| <b>CC175-8</b> 9.3mm<br><b>CC175-8</b><br><br>PROLENE Suture | <br><br>EVERPOINT™<br>Cardiovascular Needles |

## MICROPOINT - Spatula

## 1/2 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>TG160-4</b> 5.5mm<br><b>TG160-4</b><br><br>ETHILON Suture<br>MERSILENE Suture       |    |
| <b>TG160-4-3M</b> 4.3mm<br><b>TG160-4-3M</b><br><br>ETHILON Suture<br>MERSILENE Suture |    |
| <b>TG160-6</b> 5.5mm<br><b>TG160-6</b><br><br>ETHILON Suture                           |    |
| <b>TG160-8</b> 5.5mm<br><b>TG160-8</b><br><br>Coated VICRYL Suture                     |    |
| <b>TG175-8</b> 7.0mm<br><b>TG175-8</b><br><br>ETHILON Suture                           |  |

## 1/4 Circle

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>TG100-8</b> 6.5mm<br><b>TG100-8</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>Surgical Gut Suture - Chromic |  |

## 3/8 Circle

| Needles/<br>Avialable Products | Needle image |
|--------------------------------|--------------|
|                                |              |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>TG140-6</b> 6.5mm<br><b>TG140-6</b><br><br>Coated VICRYL Suture<br>ETHILON Suture  |  |
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |  |
| <b>VAS100-4</b> 5.1mm<br><br>ETHILON Suture   |  |

## MICROPOINT-Reverse Cutting

### 1/2 Circle

| Needles/<br>Avialable Products                                      | Needle image  |
|---|---|
| <b>G-2</b> 8.0mm<br><b>G-2</b><br><br>Surgical Gut Suture - Chromic |  |
| <b>G-7</b> 8.0mm<br><b>G-7</b><br><br>PERMA-HAND Suture             |  |

### 3/8 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>G-1</b> 11mm<br><br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |  |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>G-1</b> 11mm<br><b>G-1</b><br><br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain                          |  |
| <b>G-3</b> 13mm<br><br>PERMA-HAND Suture  |  |
| <b>G-3</b> 13mm<br><b>G-3</b><br><br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain                          |  |
| <b>G-6</b> 8.0mm<br><b>G-6</b><br><br>Coated VICRYL Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |  |

## Precision Cosmetic - Conventional Cutting PRIME

### 1/2 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>PC-12</b> 16mm<br><br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Suture |  |

### 3/8 Circle

| Needles/<br>Avialable Products | Needle image |
|--------------------------------|--------------|
|                                |              |

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>PC-1</b> 13mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |    |
| <b>PC-3</b> 16mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>VICRYL RAPIDE Suture                             |    |
| <b>PC-5</b> 19mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture   |  |

## Precision Point - Reverse Cutting

### 1/2 Circle

| Needles/<br>Avialable Products | Needle image |
|--------------------------------|--------------|
|                                |              |

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>P-2</b> 8.0mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain   |   |
| <b>PS-4</b> 16mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain | <br><b>MultIPass<sup>®</sup></b><br><small>NEEDLES</small> |
| <b>PS-5</b> 13mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic  |    |
| <b>PS-6</b> 11mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic   |   |

## 1/2 Circle Compound Curve

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>PS-4C</b> 16mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>Surgical Gut Suture - Chromic |  |

## 3/8 Circle

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>P-1</b> 11mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture                                  |    |
| <b>P-3</b> 13mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |     |
| <b>P-6</b> 8.0mm<br><br>ETHILON Suture<br>PROLENE Suture  |  |
| <b>PS</b> 26mm<br><br>ETHILON Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic  |  |

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>PS-1</b> 24mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture   |    |
| <b>PS-2</b> 19mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |    |
| <b>PS-3</b> 16mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain  |  |
| <b>PSL</b> 30mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic   |  |
| <b>PSLX</b> 40mm<br><br>ETHILON Suture<br>PERMA-HAND Suture  |  |

# Reverse Cutting

## 1/2 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>CP</b> 40mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>Surgical Gut Suture - Chromic  |     |
| <b>CP-1</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture |     |
| <b>CP-2</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>Surgical Gut Suture - Chromic   |  |
| <b>CPX</b> 48mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PDS II Suture<br>PDS Plus Suture   |   |
| <b>J-1</b> 17mm<br><br>Coated VICRYL Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic  |  |
| <b>M-1</b> 13mm<br><br>PERMA-HAND Suture<br>Surgical Gut Suture - Plain  |  |

## Needle Index

| <b>Needles/<br/>Avialable Products</b>   | <b>Needle image</b>   |
|--|---|
| <b>M-1</b> 13mm<br><b>M-1</b><br><br>Surgical Gut Suture - Chromic   |    |
| <b>OS-2</b> 17mm<br><br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture  |    |
| <b>OS-4</b> 22mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture                           |    |
| <b>OS-6</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture                           |     |
| <b>OS-8</b> 40mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>NUROLON Suture   |   |
| <b>X-1</b> 22mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>NUROLON Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Plain |  |

## 3/8 Circle

| <b>Needles/<br/>Avialable Products</b> | <b>Needle image</b>   |
|--|---|
| <b>C-2</b> 13mm<br><br>ETHILON Suture  |  |

## Needle Index

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>C-2</b> 13mm<br><b>C-2</b><br><br>PERMA-HAND Suture   |    |
| <b>C-3</b> 13mm<br><br>Coated VICRYL Suture<br>ETHILON Suture<br>PERMA-HAND Suture   |    |
| <b>C-3</b> 13mm<br><b>C-3</b><br><br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic  |    |
| <b>FS</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>NUROLON Suture<br>PERMA-HAND Suture<br>PROLENE Suture  |    |
| <b>FS-1</b> 24mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture |  |

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>FS-2</b> 19mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture |    |
| <b>FSL</b> 30mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>NUROLON Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Stainless Steel Suture  |     |
| <b>FSLX</b> 40mm<br><br>ETHILON Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Stainless Steel Suture   |   |
| <b>LR</b> 75mm<br><br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>PROLENE Suture  |   |
| <b>LR</b> 75mm<br><b>LR</b><br><br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>PERMA-HAND Suture<br>Surgical Stainless Steel Suture   |   |
| <b>LS-1</b> 45mm<br><br>Surgical Stainless Steel Suture  |   |
| <b>M-2</b> 13mm<br><br>Surgical Gut Suture - Chromic   |  |

| Needles/<br>Avialable Products                                     | Needle image  |
|--|---|
| <b>M-2</b> 13mm<br><b>M-2</b><br><br>Surgical Gut Suture - Chromic |  |

## SABRELOC - Center Point Spatula

### 1/4 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>RD-1</b> 8.0mm<br><b>RD-1</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>MERSILENE Suture |  |

## SABRELOC - Conventional Spatula

### 1/2 Circle

| Needles/<br>Avialable Products                               | Needle image  |
|--|---|
| <b>S-28</b> 8.0mm<br><b>S-28</b><br><br>Coated VICRYL Suture |  |

### 1/4 Circle

| Needles/<br>Avialable Products                               | Needle image  |
|--|---|
| <b>S-29</b> 8.0mm<br><b>S-29</b><br><br>Coated VICRYL Suture |  |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>S-4</b> 8.0mm<br><b>S-4</b><br><br>MERSILENE Suture<br>Surgical Gut Suture - Chromic |  |

## SABRELOC - Spatula

### 1/2 Circle

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>S-2</b> 8.0mm<br><b>S-2</b><br><br>MERSILENE Suture<br>Surgical Gut Suture - Chromic |    |
| <b>S-22</b> 8.0mm<br><b>S-22</b><br><br>ETHILON Suture                                  |  |

### 1/4 Circle

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br><br>Coated VICRYL Suture<br>MERSILENE Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |  |
| <b>S-24</b> 8.0mm<br><b>S-24</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>MERSILENE Suture  |  |

## 1/8 Circle

| Needles/<br>Avialable Products                     | Needle image  |
|--|---|
| <b>SM-1</b> 5.5mm<br><b>SM-1</b><br>PROLENE Suture |  |

## SABRELOC - Straight Transchamber

### Straight Straight

| Needles/<br>Avialable Products                           | Needle image   |
|--|--|
| <b>STC-6</b> 16mm<br><b>STC-6</b> 16mm<br>PROLENE Suture |  |

## Straight Cutting Needles

### Straight

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>KS</b> 60mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture |  |

| Needles/<br>Avialable Products                                     | Needle image  |
|--|---|
| <b>SC-1</b> 13mm<br><br>Surgical Gut Suture - Plain                |  |
| <b>SC-1</b> 13mm<br><b>SC-1</b><br><br>Surgical Gut Suture - Plain |  |
| <b>TS</b> 48mm<br><br>PERMA-HAND Suture                            |   |
| <b>TS</b> 48mm<br><b>TS</b><br><br>PROLENE Suture                  |   |

## Straight Taper Point

### Straight

| Needles/<br>Avialable Products            | Needle image   |
|---|--|
| <b>ST-1</b> 45mm<br><br>PERMA-HAND Suture |  |

### Straight Straight

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>ST-1</b> 45mm<br><b>SH</b> 26mm<br><br>Surgical Gut Suture - Chromic |  |

## Straight TAPERCUT

### Straight

| Needles/<br>Avialable Products                    | Needle image   |
|---|--|
| <b>V-47</b> 40mm<br><b>V-47</b><br>PROLENE Suture |  |

## Taper Point

### 1/2 Circle

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>BVH100-3</b> 4.8mm<br><br>ETHILON Suture  |   |
| <b>CT</b> 40mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |    |
| <b>CT</b> 40mm<br><b>CT</b><br><br>Surgical Gut Suture - Chromic   |  |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>CT-1</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MERSILENE Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture<br>VICRYL RAPIDE Suture |     |
| <b>CT-2</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic   |    |
| <b>CT-3</b> 22mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain  |  |

| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>CTX</b> 48mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>ETHILON Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>Surgical Stainless Steel Suture<br>VICRYL RAPIDE Suture |    |
| <b>CTX</b> 48mm<br><b>CTX</b><br><br>MERSILENE Suture   |    |
| <b>MH</b> 36mm<br><br>PERMA-HAND Suture   |   |
| <b>MH</b> 36mm<br><b>MH</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture<br>PRONOVA Suture  |  |
| <b>MO-2</b> 40mm<br><br>Surgical Gut Suture - Chromic   |  |
| <b>MO-4</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MERSILENE Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic   |  |
| <b>MO-5</b> 31mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic   |  |

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>MO-6</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>NUROLON Suture<br>PERMA-HAND Suture<br>PROLENE Suture  |    |
| <b>MO-7</b> 22mm<br><br>ETHIBOND EXCEL Suture<br>NUROLON Suture<br>PERMA-HAND Suture<br>PROLENE Suture   |    |
| <b>RB-1</b> 17mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |    |
| <b>RB-1</b> 17mm<br><b>RB-1</b><br><br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>PRONOVA Suture  |  |
| <b>RB-2</b> 13mm<br><b>RB-2</b><br><br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>PRONOVA Suture  |  |
| <b>RB-3</b> 11mm<br><b>RB-3</b><br><br>PROLENE Suture  |  |

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>SH</b> 26mm<br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MERSILENE Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |    |
| <b>SH</b> 26mm<br><b>SH</b><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MONOCRYL Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>PRONOVA Suture<br>Surgical Gut Suture - Chromic   |    |
| <b>SH-1</b> 22mm<br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic   |  |
| <b>SH-1</b> 22mm<br><b>SH-1</b><br>ETHIBOND EXCEL Suture<br>PROLENE Suture<br>PRONOVA Suture   |  |
| <b>SH-2</b> 20mm<br><b>SH-2</b><br>ETHIBOND EXCEL Suture   |  |

## Needle Index

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>TF</b> 13mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>NUROLON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture                          |    |
| <b>TF</b> 13mm<br><b>TF</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture<br>PRONOVA Suture   |    |
| <b>TF-4</b> 8.0mm<br><br>PERMA-HAND Suture   |    |
| <b>TP-1</b> 65mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHILON Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic             |    |
| <b>UCL</b> 14mm<br><br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture   |  |
| <b>XLH</b> 70mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |   |
| <b>XLH</b> 70mm<br><b>XLH</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture   |   |

## 3/8 Circle

| Needles/<br>Available Products   | Needle image  |
|--|---|
| <b>BB</b> 17mm<br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic  |    |
| <b>BB</b> 17mm<br><b>BB</b><br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>PRONOVA Suture  |    |
| <b>BV</b> 11mm<br><b>BV</b><br><br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture<br>PRONOVA Suture  |    |
| <b>BV-1</b> 9.3mm<br><br>PERMA-HAND Suture   |  |
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br><br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>PRONOVA Suture |  |
| <b>BV100-3</b> 3.9mm<br><br>Coated VICRYL Suture<br>ETHILON Suture<br>PROLENE Suture   |  |
| <b>BV100-4</b> 5.1mm<br><br>ETHILON Suture<br>PROLENE Suture   |  |
| <b>BV130-3</b> 3.8mm<br><br>ETHILON Suture   |  |

## Needle Index

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>BV130-4</b> 5.0mm<br><br>Coated VICRYL Suture<br>ETHILON Suture   |                                 |
| <b>BV130-5</b> 6.5mm<br><br>Coated VICRYL Suture<br>ETHILON Suture<br>PROLENE Suture                                 | <br><b>Multipass® NEEDLES</b>    |
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b><br><br>ETHILON Suture<br>PROLENE Suture<br>PRONOVA Suture                     | <br><b>Multipass® NEEDLES</b>   |
| <b>BV175-6</b> 8.0mm<br><b>BV175-6</b><br><br>PROLENE Suture<br>PRONOVA Suture                                       | <br><b>Multipass® NEEDLES</b>   |
| <b>BV175-7</b> 8.0mm<br><b>BV175-7</b><br><br>PROLENE Suture<br>PRONOVA Suture                                       | <br><b>Multipass® NEEDLES</b>  |
| <b>BV175-8</b> 9.3mm<br><b>BV175-8</b><br><br>PROLENE Suture<br>PRONOVA Suture                                       | <br><b>Multipass® NEEDLES</b> |
| <b>BV50-3</b> 4.1mm<br><br>ETHILON Suture  |                               |
| <b>BV75-3</b> 3.8mm<br><br>ETHILON Suture<br>PROLENE Suture  |                               |
| <b>BV75-4</b> 5.1mm<br><br>ETHILON Suture  |                               |
| <b>C-1</b> 13mm<br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic |                               |

| Needles/<br>Avialable Products  | Needle image  |
|---|---|
| <b>C-1</b> 13mm<br><b>C-1</b><br><br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PERMA-HAND Suture<br>PROLENE Suture<br>PRONOVA Suture |  |

## 5/8 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>UR-4</b> 40mm<br><br>Coated VICRYL Suture<br>Surgical Gut Suture - Chromic  |    |
| <b>UR-5</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic                                |   |
| <b>UR-6</b> 26mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Plus Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain |  |

## Taper Point - Reverse Cutting


### 1/2 Circle 1/2 Circle

| Needles/<br>Avialable Products | Needle image |
|--------------------------------|--------------|
|--------------------------------|--------------|


| Needles/<br>Avialable Products  | Needle image   |
|---|--|
| <b>CT</b> 40mm<br><b>FSL</b> 30mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>Surgical Gut Suture - Chromic   |  |
| <b>CT-1</b> 36mm<br><b>CP-1</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>MONOCRYL Suture<br>Surgical Gut Suture - Chromic<br>VICRYL RAPIDE Suture |  |

## TAPERCUT - Straight Cutting

### 3/8 Circle

| Needles/<br>Avialable Products                          | Needle image  |
|---|---|
| <b>V-26</b> 26mm<br><br>Surgical Stainless Steel Suture | <br>Polypropylene |

### TS Straight

| Needles/<br>Avialable Products                           | Needle image   |
|--|--|
| <b>TS_TS</b> 48mm<br><br>Surgical Stainless Steel Suture |  |

## TAPERCUT

**1/2 Circle**

| Needles/<br>Available Products  | Needle image  |
|---|---|
| <b>V-34</b> 36mm<br><br>Coated VICRYL Plus Suture<br>Coated VICRYL Suture<br>ETHIBOND EXCEL Suture<br>MONOCRYL Suture<br>PROLENE Suture<br>Surgical Gut Suture - Chromic<br>Surgical Gut Suture - Plain<br>VICRYL RAPIDE Suture |     |
| <b>V-37</b> 40mm<br><br>ETHIBOND EXCEL Suture<br>Surgical Stainless Steel Suture  |     |
| <b>V-40</b> 48mm<br><br>ETHIBOND EXCEL Suture<br>Surgical Stainless Steel Suture  |     |
| <b>V-5</b> 17mm<br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Chromic  |  |
| <b>V-5</b> 17mm<br><b>V-5</b><br><br>ETHIBOND EXCEL Suture<br>PDS II Suture<br>PDS Plus Suture<br>PROLENE Suture  |  |
| <b>V-56</b> 60mm<br><br>Surgical Stainless Steel Suture   |   |
| <b>V-7</b> 26mm<br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture<br>Surgical Gut Suture - Plain  |  |
| <b>V-7</b> 26mm<br><b>V-7</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture<br>PRONOVA Suture  |  |

## 1/4 Circle

| Needles/<br>Avialable Products                          | Needle image  |
|---|---|
| <b>CIF-4</b> 13mm<br><b>CIF-4</b><br><br>PROLENE Suture |    |
| <b>V100-3</b> 3.5mm<br><br>ETHILON Suture               |    |
| <b>V130-3</b> 3.8mm<br><br>ETHILON Suture               |    |
| <b>V130-4</b> 4.2mm<br><br>Coated VICRYL Suture         |   |
| <b>V130-5</b> 6.5mm<br><br>Coated VICRYL Suture         |   |
| <b>V75-3</b> 3.5mm<br><br>ETHILON Suture                |  |

## 3/8 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>CC</b> 9.3mm<br><b>CC</b><br><br>PROLENE Suture<br>PRONOVA Suture                             |  |
| <b>CC-1</b> 13mm<br><br>PERMA-HAND Suture  |  |
| <b>CC-1</b> 13mm<br><b>CC-1</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture<br>PRONOVA Suture |  |

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>V-4</b> 17mm<br><br>ETHIBOND EXCEL Suture<br>PERMA-HAND Suture            |  |
| <b>V-4</b> 17mm<br><b>V-4</b><br><br>ETHIBOND EXCEL Suture<br>PROLENE Suture |  |

## Temporary Pacing Wire

### 1/2 Circle Straight Cutting

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>RB-1</b> 17mm<br><b>SKS</b> 60mm<br><br>Surgical Stainless Steel Suture |   |
| <b>SH</b> 26mm<br><b>SKS-3</b> 90mm<br><br>Surgical Stainless Steel Suture |  |

### 3/8 Circle 3/8 Circle

| Needles/<br>Avialable Products   | Needle image   |
|--|--|
| <b>BB</b> 17mm<br><b>SKS</b> 60mm<br><br>Surgical Stainless Steel Suture     |  |
| <b>BB-1</b> 22mm<br><b>SKS-3</b> 90mm<br><br>Surgical Stainless Steel Suture |  |
| <b>TE</b> 31mm<br><b>SKS-3</b> 90mm<br><br>Surgical Stainless Steel Suture   |  |

# Trocar Point

## Straight

| Needles/<br>Avialable Products   | Needle image |
|--|--------------|
| <b>STP-10</b> 254mm<br><b>STP-10</b><br><br>ETHIBOND EXCEL Suture<br>PDS II Suture |              |

# ULTIMA - Spatula

## 1/2 Circle

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>CS160-6</b> 5.5mm<br><br>ETHILON Suture   |  |
| <b>CS160-6</b> 5.5mm<br><b>CS160-6</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>MERSILENE Suture<br>PROLENE Suture |  |
| <b>CS160-8</b> 5.5mm<br><b>CS160-8</b><br><br>Coated VICRYL Suture<br>ETHILON Suture<br>PROLENE Suture                     |  |
| <b>CS175-6</b> 7.0mm<br><b>CS175-6</b><br><br>ETHILON Suture   |  |
| <b>CS175-8</b> 7.0mm<br><b>CS175-8</b><br><br>ETHILON Suture<br>Surgical Gut Suture - Chromic                              |  |

## 1/2 Circle 1/2 Circle

| Needles/<br>Avialable Products                                  | Needle image  |
|---|---|
| <b>CS160-6</b> 5.5mm<br><b>CTC-6</b> 13mm<br><br>PROLENE Suture |  |

## 1/4 Circle

| Needles/<br>Avialable Products                            | Needle image  |
|---|---|
| <b>CS90-6</b> 7.0mm<br><br>Coated VICRYL Suture           |  |
| <b>CTC-6L</b> 13mm<br><b>CTC-6L</b><br><br>PROLENE Suture |  |

## 1/8 Circle

| Needles/<br>Avialable Products            | Needle image  |
|---|---|
| <b>CS44-6</b> 5.5mm<br><br>ETHILON Suture |  |

## 3/8 Circle

| Needles/<br>Avialable Products                                     | Needle image  |
|--|---|
| <b>CS140-6</b> 6.5mm<br><br>Coated VICRYL Suture<br>ETHILON Suture |  |

| Needles/<br>Avialable Products   | Needle image  |
|--|---|
| <b>CS140-6</b> 6.5mm<br><b>CS140-6</b><br><br>ETHILON Suture<br>MERSILENE Suture<br>PROLENE Suture |  |
| <b>CS140-8</b> 6.5mm<br><b>CS140-8</b><br><br>ETHILON Suture                                       |  |


## Bi-Curve

| Needles/<br>Avialable Products                             | Needle image  |
|--|---|
| <b>CSB-6</b> 5.0mm<br><br>Coated VICRYL Suture             |    |
| <b>CSB-6</b> 5.0mm<br><b>CSB-6</b><br><br>MERSILENE Suture |  |

 **MultiPass® NEEDLES** MULTIPASS Needles

# ETHICON NEEDLE SALES TYPES

## CUTTING EDGE



\*TRADEMARK

# **ETHICON NEEDLE SALES TYPES**


## TAPER POINT

|            |            |            |
|------------|------------|------------|
| 1/4 CIRCLE | 5/8 CIRCLE | 1/2 CIRCLE |
| BV30-2     | TF-5       | BVH100-3   |
| 3/8 CIRCLE | MF-1       | TF-6       |
| BV50-3     | MF         | TF-4       |
| BV75-3     | UR-6       | TF-1       |
| BV75-4     | UR-5       | RB-4       |
| BV100-3    | UR-4       | RB-3       |
| BV100-4    |            | TF         |
| BV130-3    |            | UCL        |
| BV130-4    |            | RB-2       |
| BV130-5    |            | UCLX       |
| BV175-5    |            | RB-1       |
| BV175-6    |            | SH-2       |
| BV175-7    |            | RS         |
| BV175-8    |            | SH-1       |
| BV-1       |            | SH         |
| BV         |            | MH-1       |
| C-1        |            | MH         |
| BB         |            | LH         |
| BB-1       |            | XLH        |
| TE-2       |            | XLH-1      |
| TE         |            | XXLH       |
| TE-1       |            |            |

## TAPERCUT\*


| 3/8 CIRCLE | 1/4 CIRCLE |
|------------|------------|
| V130-3     | V75-3      |
| V130-4     | V100-3     |
| V130-5     | RV-2       |
| CC         | CIF-4      |
| CC-1       |            |
| V-18       |            |
| V-4        |            |
| V-3        |            |
| V-26       |            |

## TROCAR




STRAIGHT  
V-47


## BLUNT POINT




ETHIGUARD®



SKI




## STRAIGHT




ENDOSCOPY

TAPERPOINT



TAPERCUT\*



\*TRADEMARK

MW v.4

# Coated VICRYL Plus Suture

## Blunt Point

| Needle   | Material       | Length       | 0 | 1        |
|--|----------------|--------------|---|----------|
| <b>TPB-1</b> 65mm<br>1/2 Circle<br> | Violet Braided | 4-27" ~ 70cm |   | VCPB650G |

## ETHIGUARD - Safety needle


| Needle   | Material       | Length       | 0  | 1  | 2-0  | 3-0                  |
|--|----------------|--------------|--|--|--|----------------------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Undyed Braided | 8-18" ~ 45cm | VCPB840D<br> VCPB840Z | VCPB841D<br> VCPB841Z | VCPB839D<br> VCPB839Z |                      |
|  | Undyed Braided | 27" ~ 70cm   | VCPB260H<br>VCPB260W   |  | VCPB259H<br>VCPB259W   |                      |
|  | Undyed Braided | 8-27" ~ 70cm | VCPB41D<br>         |  |  |                      |
|  | Undyed Braided | 36" ~ 90cm   | VCPB946H<br>VCPB946W   | VCPB947H<br>VCPB947W   | VCPB945H   | VCPB944H<br>VCPB944W |
|  | Violet Braided | 8-18" ~ 45cm | VCPB740D<br>          | VCPB741D<br>          | VCPB739D<br>          |                      |
|  | Violet Braided | 27" ~ 70cm   | VCPB340H<br>VCPB340W   | VCPB341H<br>VCPB341W   |  |                      |
|  | Violet Braided | 8-27" ~ 70cm | VCPB341D<br>          |  |  |                      |
|  | Violet Braided | 36" ~ 90cm   | VCPB346H   |  |  |                      |


| Needle | Material | Length | 0 | 1 |
|--------|----------|--------|---|---|
|        |          |        |   |   |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>  |
|--|-----------------|---------------|---|---|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  | VCPB724D<br><br>VCPB724Z | VCPB725D<br><br>VCPB725Z |
|  | Undyed Braided  | 36" ~ 90cm    | VCPB978H<br>VCPB978W  | VCPB977H <sup>3</sup><br>VCPB977W   |
|  | Violet Braided  | 8-18" ~ 45cm  |   | VCPB765D<br>             |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>  |
|--|-----------------|---------------|------------|---|
| <b>SHB</b> 26mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  |            | VCPB864D<br><br>VCPB864Z |
|  | Undyed Braided  | 27" ~ 70cm    | VCPB417H   | VCPB416H  |
|  | Violet Braided  | 8-18" ~ 45cm  |            | VCPB774D<br>             |
|  | Violet Braided  | 27" ~ 70cm    |            | VCPB316H  |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>              | <b>2-0</b> |
|--|-----------------|---------------|-----------------------|------------|
| <b>CTB</b> 40mm<br>1/2 Circle<br> | Violet Braided  | 36" ~ 90cm    | VCPB330H <sup>3</sup> | VCPB357H   |
|  | Undyed Braided  | 36" ~ 90cm    | VCPB958H              |            |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>2-0</b>  |
|--|-----------------|---------------|---|---|
| <b>CTB-2</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCPB727D<br> | VCPB726D<br> |
|  | Undyed Braided  | 27" ~ 70cm    |   | VCPB269H  |

## **LIGAPAK Ligating Reel**

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b>           | <b>4-0</b>           |
|--|-----------------|---------------|----------|----------|------------|----------------------|----------------------|
| <b>LIGAPAK Ligating Reel</b><br>N/A<br> | Undyed Braided  | 144" ~ 366cm  |          |          |            | VCP885G <sup>1</sup> |                      |
|  | Undyed Braided  | 54" ~ 135cm   | VCP287G  |          | VCP286G    | VCP285G              | VCP284G <sup>1</sup> |
|  | Violet Braided  | 54" ~ 135cm   | VCP207G  | VCP208G  |            | VCP205G              |                      |

## **Precision Cosmetic - Conventional Cutting PRIME**

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP835G    | VCP834G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>           | <b>5-0</b> |
|---|-----------------|---------------|----------------------|------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP845G<br>VCP845G15 | VCP844G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b>         | <b>5-0</b> |
|---|-----------------|---------------|------------|--------------------|------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP824G    | VCP823G<br>VCP823H | VCP822G    |

## Precision Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP504G    | VCP503G    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br><br> | Undyed Braided  | 18" ~ 45cm    | VCP507G    | VCP506G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>           | <b>5-0</b> |
|---|-----------------|---------------|----------------------|------------|
| <b>PS-5</b> 13mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP594G <sup>1</sup> | VCP593G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-6</b> 11mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP511G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-4C</b> 16mm<br>1/2 Circle Compound Curve<br> | Undyed Braided  | 18" ~ 45cm    | VCP656G    |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP490G    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b>           | <b>5-0</b>                      |
|--|-----------------|---------------|----------------------|---------------------------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | VCP464G <sup>1</sup> | VCP463G                         |
|  | Undyed Braided  | 18" ~ 45cm    | VCP494G<br>VCP494H   | VCP493G<br>VCP493H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>         | <b>4-0</b>         |
|---|-----------------|---------------|--------------------|--------------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP683G<br>VCP683H | VCP682G<br>VCP682H |
|   | Undyed Braided  | 27" ~ 70cm    | VCP936H            | VCP935H            |


| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>                      | <b>4-0</b>                                   | <b>5-0</b>                                   |
|---|-----------------|---------------|------------|---------------------------------|--|--|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    |            | VCP497G<br>VCP497H <sup>3</sup> | VCP496G <sup>1</sup><br>VCP496H <sup>3</sup> | VCP495G <sup>1</sup><br>VCP495H <sup>3</sup> |
|   | Undyed Braided  | 27" ~ 70cm    | VCP428H    | VCP427H <sup>3</sup>            | VCP426H <sup>3</sup>                         |  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b>         |
|---|-----------------|---------------|--------------------|
| <b>PS-3</b> 16mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | VCP500G<br>VCP500H |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>PSL</b> 30mm<br>3/8 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP596H    | VCP502H    |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b>             | <b>2-0</b> | <b>2</b> |
|---|-----------------|---------------|----------|----------------------|------------|----------|
| <b>CP</b> 40mm<br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP197H  | VCP196H <sup>1</sup> |            | VCP195H  |
|   | Violet Braided  | 27" ~ 70cm    |          | VCP480H              | VCP478H    |          |
|   | Violet Braided  | 36" ~ 90cm    |          | VCP486H              |            |          |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|----------|----------|------------|------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP267H  | VCP268H  | VCP266H    | VCP265H    |
|   | Violet Braided  | 27" ~ 70cm    | VCP467H  | VCP468H  | VCP466H    |            |
|   | Violet Braided  | 36" ~ 90cm    |          | VCP474H  | VCP472H    |            |
|   | Violet Braided  | 54" ~ 135cm   |          |          |            | VCP531H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b>   | <b>3-0</b>   |
|---|-----------------|---------------|----------|----------|--|--|
| <b>CP-2</b> 26mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  |          |          | VCP762D<br> | VCP761D<br> |
|   | Undyed Braided  | 27" ~ 70cm    | VCP870H  | VCP871H  | VCP869H <sup>3</sup>   | VCP868H  |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> |
|--|-----------------|---------------|----------|----------|
| <b>CPX 48mm</b><br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP568H  | VCP569H  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   |
|---|-----------------|---------------|---|--|--|
| <b>OS-4 22mm</b><br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | VCP707T<br> | VCP708T<br> | VCP706T<br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | VCP749T<br> |  | VCP748T<br> |
|   | Undyed Braided  | 27" ~ 70cm    | VCP694H   | VCP695H  |  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>             | <b>2-0</b>   |
|---|-----------------|---------------|---|----------------------|--|
| <b>OS-6 36mm</b><br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | VCP711T<br> |                      | VCP710T<br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | VCP754T<br> |                      | VCP750T<br> |
|   | Undyed Braided  | 27" ~ 70cm    | VCP534H   | VCP535H<br>VCP535H05 | VCP533H <sup>3</sup>   |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   | <b>2</b>   |
|---|-----------------|---------------|---|--|--|--|
| <b>OS-8 40mm</b><br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | VCP717T<br> | VCP718T<br> | VCP716T<br> | VCP719T<br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | VCP756T<br> | VCP757T<br> |  |  |
|   | Undyed Braided  | 27" ~ 70cm    | VCP698H   | VCP699H  |  |  |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>  |
|--|-----------------|---------------|------------|---|
| <b>X-1</b> 22mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |            | VCP790D<br><br>VCP790Z |
|  | Violet Braided  | 27" ~ 70cm    |            | VCP460H   |
|  | Undyed Braided  | 27" ~ 70cm    | VCP459H    | VCP458H   |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>   |
|---|-----------------|---------------|--|
| <b>FS</b> 26mm<br>3/8 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP787D<br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>           | <b>4-0</b> |
|---|-----------------|---------------|------------|----------------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP443H    | VCP442H <sup>3</sup> | VCP441H    |
|   | Violet Braided  | 27" ~ 70cm    |            | VCP452H              | VCP451H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b>           | <b>5-0</b> |
|---|-----------------|---------------|------------|----------------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | VCP393H    | VCP392H <sup>3</sup> | VCP391H    |
|   | Undyed Braided  | 27" ~ 70cm    | VCP423H    | VCP422H <sup>3</sup> | VCP421H    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>FSL</b> 30mm<br>3/8 Circle<br> | Undyed Braided  | 27" ~ 70cm    | VCP589H    | VCP588H    |

## Standard and Short Length Sutures

| Needle   | Material       | Length        | 0       | 1                    | 2-0     | 2                    | 3-0     |
|--|----------------|---------------|---------|----------------------|---------|----------------------|---------|
| <b>Standard and Short Length Sutures</b><br>N/A<br> | Violet Braided | 3-18" ~ 45cm  | VCP636T |                      | VCP635T |                      | VCP634T |
|  | Violet Braided | 54" ~ 135cm   | VCP616H | VCP617H <sup>3</sup> | VCP615H | VCP618H <sup>3</sup> |         |
|  | Undyed Braided | 3-18" ~ 45cm  | VCP646T |                      | VCP645T |                      | VCP644H |
|  | Undyed Braided | 12-18" ~ 45cm |         |                      |         |                      | VCP110G |
|  | Undyed Braided | 54" ~ 135cm   | VCP608H |                      | VCP607H |                      |         |

## Straight Cutting Needles

| Needle  | Material       | Length     | 3-0     | 4-0     |
|---|----------------|------------|---------|---------|
| <b>KS 60mm</b><br>Straight<br> | Undyed Braided | 27" ~ 70cm | VCP663H | VCP662H |

## SUTUPAK Pre-Cut Sutures

| Needle   | Material       | Length        | 0       | 1       | 2-0     | 3-0     | 4-0     |
|--|----------------|---------------|---------|---------|---------|---------|---------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Violet Braided | 6-18" ~ 45cm  | VCP106G | VCP107G |         |         |         |
|  | Violet Braided | 12-18" ~ 45cm |         |         | VCP105G | VCP104G | VCP103G |
|  | Undyed Braided | 6-18" ~ 45cm  | VCP112G |         |         |         |         |
|  | Undyed Braided | 12-18" ~ 45cm |         |         | VCP111G |         | VCP109G |

## Taper Point

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   | <b>2-0</b>   | <b>3-0</b> |
|---|-----------------|---------------|--|--|--|------------|
| <b>CT 40mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP752D<br> | VCP753D<br> | VCP751D<br> |            |
|   | Violet Braided  | 27" ~ 70cm    | VCP352H  | VCP353H  | VCP351H  |            |
|   | Violet Braided  | 36" ~ 90cm    | VCP358H  | VCP359H  | VCP357H  | VCP356H    |
|   | Undyed Braided  | 27" ~ 70cm    | VCP280H  | VCP281H  | VCP275H  |            |
|   | Undyed Braided  | 36" ~ 90cm    | VCP958H  | VCP959H  | VCP957H  | VCP956H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   | <b>2-0</b>   | <b>3-0</b>   |
|---|-----------------|---------------|--|--|--|--|
| <b>CT-1 36mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP740D<br> VCP740Z | VCP741D <sup>13</sup><br> VCP741Z | VCP739D<br> | VCP738D<br> |
|   | Violet Braided  | 27" ~ 70cm    | VCP340H<br>VCP340W   | VCP341H<br>VCP341W   | VCP339H<br>VCP339W   | VCP338H<br>VCP338W   |
|   | Violet Braided  | 8-27" ~ 70cm  | VCPP31D<br>         |  |  |  |
|   | Violet Braided  | 36" ~ 90cm    | VCP346H  | VCP347H <sup>3</sup>   | VCP345H  | VCP344H  |
|   | Undyed Braided  | 8-18" ~ 45cm  | VCP840D<br>         | VCP841D<br>                       | VCP839D<br> | VCP838D<br> |
|   | Undyed Braided  | 27" ~ 70cm    | VCP260H<br>VCP946H   | VCP261H  | VCP259H  | VCP258H  |
|   | Undyed Braided  | 8-27" ~ 70cm  | VCPP41D<br>         | VCPP40D<br>                       | VCPP42D<br> |  |
|   | Undyed Braided  | 36" ~ 90cm    |  | VCP947H  | VCP945H  | VCP944H  |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|----------|----------|------------|------------|
|               |                 |               |          |          |            |            |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>             | <b>2-0</b>  | <b>3-0</b> |
|---|-----------------|---------------|---|----------------------|---|------------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP727D<br> |                      | VCP726D<br><br>VCP726Z |            |
|   | Violet Braided  | 27" ~ 70cm    | VCP334H   | VCP335H <sup>3</sup> | VCP333H<br>VCP333W  | VCP332H    |
|   | Undyed Braided  | 27" ~ 70cm    | VCP270H   |                      | VCP269H   | VCP232H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|----------|------------|------------|
| <b>CT-3</b> 22mm<br>1/2 Circle<br> | Violet Braided  | 27" ~ 70cm    | VCP329H  | VCP328H    | VCP327H    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   | <b>3-0</b>           |
|--|-----------------|---------------|---|--|--|----------------------|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  | VCP724D<br> |  | VCP723D<br> |                      |
|  | Undyed Braided  | 36" ~ 90cm    | VCP978H   | VCP977H  | VCP979H  | VCP980H <sup>3</sup> |
|  | Violet Braided  | 8-18" ~ 45cm  | VCP764D<br> | VCP765D<br> |  |                      |
|  | Violet Braided  | 27" ~ 70cm    | VCP364H   | VCP365H  | VCP363H  |                      |
|  | Violet Braided  | 36" ~ 90cm    | VCP370H<br>VCP370H05  | VCP371H<br>VCP371H05   | VCP369H  |                      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   |
|---|-----------------|---------------|---|--|--|
| <b>MO-4</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP701D<br> | VCP702D<br> | VCP700D<br> |
|   | Violet Braided  | 27" ~ 70cm    | VCP436H   | VCP437H  |  |
|   | Violet Braided  | 8-27" ~ 70cm  | VCPP71D<br> |  |  |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-5</b> 31mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP801D<br> |  |
|   | Violet Braided  | 8-27" ~ 70cm  | VCPP81D<br> | VCPP80D<br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |  |  |
|---|-----------------|---------------|------------|--|--|
| <b>MO-6</b> 26mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |            | VCP789D<br> |  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>   | <b>4-0</b>   | <b>5-0</b> |
|---|-----------------|---------------|------------|--|--|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |            | VCP713D<br> VCP713Z |  |            |
|   | Violet Braided  | 27" ~ 70cm    | VCP306H    | VCP305H<br>VCP305W   | VCP304H<br>VCP304W   | VCP303H    |
|   | Undyed Braided  | 8-18" ~ 45cm  |            |  | VCP714D<br> VCP714Z |            |
|   | Undyed Braided  | 27" ~ 70cm    |            | VCP215H  | VCP214H  | VCP213H    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|------------|
|               |                 |               |          |            |            |            |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b>  | <b>3-0</b>  | <b>4-0</b>  |
|---|-----------------|---------------|----------|---|---|---|
| <b>SH 26mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |          | VCP775D<br><br>VCP775Z | VCP774D<br><br>VCP774Z | VCP773D<br><br>VCP773Z |
|   | Violet Braided  | 27" ~ 70cm    |          | VCP317H<br>VCP317W  | VCP316H <sup>3</sup><br>VCP316W   | VCP315H   |
|   | Violet Braided  | 8-27" ~ 70cm  |          | VCP785D<br>            | VCP784D<br>            |   |
|   | Undyed Braided  | 8-18" ~ 45cm  |          |   | VCP864D<br>            |   |
|   | Undyed Braided  | 27" ~ 70cm    | VCP418H  | VCP417H   | VCP416H   | VCP415H   |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|
| <b>SH 26mm</b><br><b>SH</b><br>1/2 Circle<br> | Violet Braided  | 36" ~ 90cm    | VCP527H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>  | <b>4-0</b>  |
|---|-----------------|---------------|---|---|
| <b>SH-1 22mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP772D<br><br>VCP772Z | VCP771D<br><br>VCP771Z |
|   | Violet Braided  | 27" ~ 70cm    | VCP311H<br>VCP311W  | VCP310H<br>VCP310W  |
|   | Undyed Braided  | 27" ~ 70cm    | VCP219H   | VCP218H   |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>  | <b>5-0</b> |
|---|-----------------|---------------|---|------------|
| <b>TF 13mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | VCP743D<br><br>VCP743Z |            |
|   | Undyed Braided  | 27" ~ 70cm    | VCP434H   | VCP433H    |

## Sutures : Absorbable / Coated VICRYL Plus Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>1</b> | <b>2</b> |
|---|-----------------|---------------|----------|----------|
| <b>TP-1</b> 65mm<br>1/2 Circle<br> | Violet Braided  | 4-27" ~ 70cm  | VCP650G  | VCP649G  |
|   | Undyed Braided  | 2-27" ~ 70cm  |          | VCP849G  |
|   | Undyed Braided  | 54" ~ 135cm   |          | VCP880T  |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>1</b> | <b>2-0</b> |
|--|-----------------|---------------|----------|------------|
| <b>XLH</b> 70mm<br>1/2 Circle<br> | Violet Braided  | 27" ~ 70cm    | VCP583G  | VCP581G    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|-----------------|---------------|----------|------------|
| <b>UR-5</b> 36mm<br>5/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | VCP376H  | VCP375H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>             | <b>2-0</b>                      |
|---|-----------------|---------------|----------------------|---------------------------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | VCP603H <sup>3</sup> | VCP602H <sup>3</sup><br>VCP602W |

## Taper Point - Reverse Cutting

| Needle   | Material       | Length      | 3-0                  |
|--|----------------|-------------|----------------------|
| <b>CT</b> 40mm<br><b>FSL</b> 30mm<br>1/2 Circle<br> | Undyed Braided | 54" ~ 135cm | VCP940H <sup>3</sup> |

| Needle  | Material       | Length      | 3-0                  |
|---|----------------|-------------|----------------------|
| <b>CT-1</b> 36mm<br><b>CP-1</b> 36mm<br>1/2 Circle<br> | Undyed Braided | 54" ~ 135cm | VCP931H <sup>3</sup> |

## TAPERCUT

| Needle  | Material       | Length     | 0       | 1                    | 2-0     | 3-0     |
|---|----------------|------------|---------|----------------------|---------|---------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Violet Braided | 36" ~ 90cm | VCP518H | VCP519H <sup>3</sup> | VCP517H | VCP516H |
|   | Undyed Braided | 36" ~ 90cm |         |                      | VCP917H | VCP916H |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>13</sup> Safety organizer tray

 **MULTIPASS** NEEDLES

 CONTROL RELEASE

# MONOCRYL Plus Suture

## ETHIGUARD - Safety needle

| Needle   | Material            | Length     | 0                     | 2-0      | 3-0                   |
|--|---------------------|------------|-----------------------|----------|-----------------------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm | MCPB346H <sup>3</sup> | MCPB345H |                       |
|  | Undyed Monofilament | 36" ~ 90cm | MCPB946H              | MCPB945H | MCPB944H <sup>3</sup> |

| Needle  | Material            | Length     | 3-0                   |
|---|---------------------|------------|-----------------------|
| <b>SHB</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | MCPB316H <sup>3</sup> |
|   | Undyed Monofilament | 27" ~ 70cm | MCPB416H              |

## Precision Cosmetic - Conventional Cutting PRIME

| Needle   | Material            | Length     | 4-0     | 5-0     |
|--|---------------------|------------|---------|---------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br><br> | Undyed Monofilament | 18" ~ 45cm | MCP835G | MCP834G |

| Needle   | Material            | Length     | 4-0     | 5-0     |
|--|---------------------|------------|---------|---------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br><br> | Undyed Monofilament | 18" ~ 45cm | MCP845G | MCP844G |

## Sutures : Absorbable / MONOCRYL Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | MCP823G    | MCP822G    |

## Precision Point - Reverse Cutting

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|--|---------------------|---------------|------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | MCP490G    | MCP489G    |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|---------------------|---------------|------------|------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | MCP494G    | MCP493G    | MCP492G    |
|  | Violet Monofilament | 18" ~ 45cm    |            | MCP463G    |            |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    |            | MCP682H    |
|   | Undyed Monofilament | 27" ~ 70cm    | MCP936H    | MCP935H    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|
|               |                 |               |            |            |            |

## Sutures : Absorbable / MONOCRYL Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | MCP497G    | MCP496G    | MCP495G    |
|   | Undyed Monofilament | 27" ~ 70cm    | MCP427H    | MCP426H    |            |
|   | Violet Monofilament | 18" ~ 45cm    |            | MCP513G    |            |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> |
|---|---------------------|---------------|------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | MCP266H    |

## Straight Cutting Needles

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> |
|---|---------------------|---------------|------------|
| <b>KS</b> 60mm<br>Straight<br> | Undyed Monofilament | 27" ~ 70cm    | MCP523H    |

## Taper Point

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> |
|---|---------------------|---------------|----------|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | MCP358H  |

## Sutures : Absorbable / MONOCRYL Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---|---------------------|---------------|----------|----------|------------|------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | MCP340H  |          | MCP339H    |            |
|   | Violet Monofilament | 36" ~ 90cm    | MCP346H  | MCP347H  | MCP345H    | MCP344H    |
|   | Undyed Monofilament | 36" ~ 90cm    | MCP946H  | MCP947H  | MCP945H    | MCP944H    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> |
|---|---------------------|---------------|------------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | MCP333H    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | MCP214H    | MCP213H    |
|   | Violet Monofilament | 27" ~ 70cm    | MCP304H    |            |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | MCP317H    | MCP316H    | MCP315H    |
|   | Undyed Monofilament | 27" ~ 70cm    | MCP417H    | MCP416H    | MCP415H    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Absorbable / MONOCRYL Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> |
|---|---------------------|---------------|------------|
| <b>UR-5</b> 36mm<br>5/8 Circle<br> | Violet Monofilament | 27" ~ 70cm    | MCP377H    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Violet Monofilament | 27" ~ 70cm    | MCP605H    | MCP604H    |

<sup>3</sup> Three dozen packets per box

 **Multipass<sup>®</sup>** MULTIPASS Needles


# PDS Plus Suture

## Blunt Point

| Needle  | Material                 | Length      | 0        | 1                        |
|---|--------------------------|-------------|----------|--------------------------|
| <b>BP-1</b> 65mm<br>1/2 Circle<br> | Violet Monofilament Loop | 48" ~ 122cm |          | PDPB880G <sup>1,15</sup> |
|   | Violet Monofilament Loop | 60" ~ 150cm | PDPB991G |                          |

| Needle   | Material            | Length      | 1                     |
|--|---------------------|-------------|-----------------------|
| <b>TPB-1</b> 65mm<br>1/2 Circle<br> | Violet Monofilament | 54" ~ 135cm | PDPB879G <sup>1</sup> |

## Endoscopic Needles

| Needle   | Material            | Length     | 4-0                  |
|--|---------------------|------------|----------------------|
| <b>ST-4</b> 19mm<br><b>ST-4</b><br>Straight Taper Point<br> | Violet Monofilament | 14" ~ 35cm | PDP420G <sup>1</sup> |

## ETHIGUARD - Safety needle

| Needle | Material | Length | 0 | 1 | 2-0 |
|--------|----------|--------|---|---|-----|
|        |          |        |   |   |     |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>              | <b>2-0</b>            |
|--|---------------------|---------------|----------------------|-----------------------|-----------------------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDPB340 <sup>3</sup> | PDPB341H <sup>3</sup> | PDPB339H <sup>3</sup> |
|  | Violet Monofilament | 36" ~ 90cm    | PDPB346              | PDPB347 <sup>3</sup>  | PDPB345 <sup>3</sup>  |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             |
|--|---------------------|---------------|----------------------|----------------------|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | PDPB370 <sup>6</sup> | PDPB371 <sup>6</sup> |
|  | Violet Monofilament | 60" ~ 150cm   | PDPB990G             |                      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>2-0</b>            | <b>3-0</b>            |
|--|---------------------|---------------|-----------------------|-----------------------|
| <b>SHB</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDPB317H <sup>3</sup> | PDPB316H <sup>3</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b>              | <b>1</b>              |
|--|---------------------|---------------|-----------------------|-----------------------|
| <b>CTB</b> 40mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | PDPB358T <sup>6</sup> | PDPB359T <sup>6</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b>              | <b>2-0</b>            | <b>3-0</b>            |
|--|---------------------|---------------|-----------------------|-----------------------|-----------------------|
| <b>CTB-2</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDPB334H <sup>3</sup> | PDPB333H <sup>3</sup> | PDPB332H <sup>3</sup> |


## Precision Cosmetic - Conventional Cutting PRIME

| Needle  | Material            | Length     | 6-0                  |
|---|---------------------|------------|----------------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm | PDP833G <sup>1</sup> |

| Needle  | Material            | Length     | 4-0                  | 5-0                  |
|---|---------------------|------------|----------------------|----------------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm | PDP845G <sup>1</sup> | PDP844G <sup>1</sup> |

| Needle  | Material            | Length     | 3-0                  | 4-0                  |
|---|---------------------|------------|----------------------|----------------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm | PDP824G <sup>1</sup> | PDP823G <sup>1</sup> |


## Precision Point - Reverse Cutting

| Needle  | Material            | Length     | 4-0                  | 5-0                  |
|---|---------------------|------------|----------------------|----------------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm | PDP504G <sup>1</sup> | PDP503G <sup>1</sup> |

| Needle | Material | Length | 4-0 |
|--------|----------|--------|-----|
|        |          |        |     |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b>           |
|---|---------------------|---------------|----------------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | PDP507G <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>6-0</b>           |
|--|---------------------|---------------|----------------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Violet Monofilament | 18" ~ 45cm    | PDP487G <sup>1</sup> |
|  | Undyed Monofilament | 18" ~ 45cm    | PDP489G              |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>4-0</b>           | <b>5-0</b> |
|--|---------------------|---------------|----------------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | PDP494G <sup>1</sup> | PDP493G    |
|  | Violet Monofilament | 18" ~ 45cm    |                      | PDP463G    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | PDP683G <sup>1</sup> | PDP682G <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           | <b>5-0</b> |
|---|---------------------|---------------|----------------------|----------------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | PDP497G <sup>1</sup> | PDP496G <sup>1</sup> | PDP495G    |
|   | Violet Monofilament | 18" ~ 45cm    |                      | PDP513G              |            |

## Reverse Cutting

| Needle  | Material            | Length     | 0                    | 1                    | 2-0                  |
|---|---------------------|------------|----------------------|----------------------|----------------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | PDP467H <sup>3</sup> | PDP468H <sup>3</sup> | PDP466H <sup>3</sup> |

| Needle  | Material            | Length     | 0                    | 2-0                  | 3-0                  |
|---|---------------------|------------|----------------------|----------------------|----------------------|
| <b>CP-2</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | PDP970H <sup>3</sup> | PDP969H <sup>3</sup> | PDP968H <sup>3</sup> |


| Needle   | Material            | Length     | 1                    |
|--|---------------------|------------|----------------------|
| <b>CPX</b> 48mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | PDP569T <sup>6</sup> |

| Needle  | Material            | Length     | 1                    |
|---|---------------------|------------|----------------------|
| <b>OS-4</b> 22mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | PDP695H <sup>3</sup> |

| Needle | Material | Length | 0 | 1 |
|--------|----------|--------|---|---|
|        |          |        |   |   |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             |
|---|---------------------|---------------|----------------------|----------------------|
| <b>OS-6</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP534T <sup>6</sup> | PDP535H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b>           | <b>3-0</b>           | <b>4-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|----------------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | PDP443H <sup>3</sup> | PDP442H <sup>3</sup> | PDP441H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           | <b>5-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|----------------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | PDP423H <sup>3</sup> | PDP422H <sup>3</sup> | PDP421H <sup>3</sup> |

## Straight Cutting Needles

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|
| <b>KS</b> 60mm<br>Straight<br> | Undyed Monofilament | 27" ~ 70cm    | PDP663H <sup>3</sup> | PDP662H <sup>3</sup> |

## Taper Point

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|---------------|-----------------|---------------|----------|----------|------------|
|               |                 |               |          |          |            |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             | <b>2-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|----------------------|
| <b>CT 40mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP352H <sup>3</sup> | PDP353H <sup>3</sup> |                      |
|   | Violet Monofilament | 36" ~ 90cm    | PDP358T              | PDP359T <sup>6</sup> | PDP357H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             | <b>2-0</b> | <b>3-0</b> |
|---|---------------------|---------------|----------------------|----------------------|------------|------------|
| <b>CT-1 36mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP340H <sup>3</sup> | PDP341H <sup>3</sup> | PDP339H    | PDP338H    |
|   | Violet Monofilament | 36" ~ 90cm    | PDP346H              | PDP347H <sup>3</sup> |            |            |
|   | Undyed Monofilament | 27" ~ 70cm    |                      |                      | PDP259H    | PDP258H    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>2-0</b>           | <b>3-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|----------------------|
| <b>CT-2 26mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP334H <sup>3</sup> | PDP333H <sup>3</sup> | PDP332H <sup>3</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             |
|--|---------------------|---------------|----------------------|----------------------|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | PDP370T <sup>6</sup> | PDP371T <sup>6</sup> |
|  | Violet Monofilament | 60" ~ 150cm   | PDP990G              |                      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           | <b>5-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|----------------------|
| <b>RB-1 17mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP305H <sup>3</sup> | PDP304H <sup>3</sup> | PDP303H <sup>3</sup> |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b>           |
|--|---------------------|---------------|----------------------|
| <b>RB-1</b> 17mm<br><b>RB-1</b><br>1/2 Circle<br> | Violet Monofilament | 30" ~ 75cm    | PDP320H <sup>3</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b>           | <b>6-0</b>           |
|--|---------------------|---------------|----------------------|----------------------|
| <b>RB-2</b> 13mm<br><b>RB-2</b><br>1/2 Circle<br> | Violet Monofilament | 30" ~ 75cm    | PDP148H <sup>3</sup> | PDP149H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b>           | <b>3-0</b>           | <b>4-0</b> |
|---|---------------------|---------------|----------------------|----------------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP317H <sup>3</sup> | PDP316H <sup>3</sup> | PDP315H    |
|   | Undyed Monofilament | 27" ~ 70cm    |                      | PDP416H              |            |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>           | <b>4-0</b>           |
|---|---------------------|---------------|----------------------|----------------------|
| <b>SH-1</b> 22mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP311H <sup>3</sup> | PDP310H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>6-0</b>           |
|---|---------------------|---------------|----------------------|
| <b>TF</b> 13mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | PDP432H <sup>3</sup> |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>                  | <b>Material</b>          | <b>Length</b> | <b>0</b> | <b>1</b> |
|--------------------------------|--------------------------|---------------|----------|----------|
| <b>TP-1</b> 65mm<br>1/2 Circle | Violet Monofilament      | 2-27" ~ 70cm  |          | PDP690G  |
|                                | Violet Monofilament      | 54" ~ 135cm   |          | PDP879G  |
|                                | Violet Monofilament Loop | 60" ~ 150cm   | PDP991G  |          |
|                                | Violet Monofilament Loop | 96" ~ 244cm   |          | PDP880G  |

| <b>Needle</b>                 | <b>Material</b>     | <b>Length</b> | <b>2-0</b>           |
|-------------------------------|---------------------|---------------|----------------------|
| <b>UCL</b> 14mm<br>1/2 Circle | Violet Monofilament | 27" ~ 70cm    | PDP114H <sup>3</sup> |

| <b>Needle</b>                 | <b>Material</b>     | <b>Length</b> | <b>0</b>             | <b>1</b>             |
|-------------------------------|---------------------|---------------|----------------------|----------------------|
| <b>XLH</b> 70mm<br>1/2 Circle | Violet Monofilament | 27" ~ 70cm    | PDP582G <sup>1</sup> | PDP583G <sup>1</sup> |
|                               | Violet Monofilament | 54" ~ 135cm   |                      | PDP877G              |
|                               | Violet Monofilament | 96" ~ 244cm   |                      | PDP881G              |

| <b>Needle</b>                             | <b>Material</b>     | <b>Length</b> | <b>6-0</b>           |
|---|---------------------|---------------|----------------------|
| <b>BV</b> 11mm<br><b>BV</b><br>3/8 Circle | Violet Monofilament | 30" ~ 75cm    | PDP123H <sup>3</sup> |

## Sutures : Absorbable / PDS Plus Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>6-0</b>           |
|---|---------------------|---------------|----------------------|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Violet Monofilament | 30" ~ 75cm    | PDP117H <sup>3</sup> |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b>           | <b>6-0</b>           |
|--|---------------------|---------------|----------------------|----------------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Violet Monofilament | 30" ~ 75cm    | PDP126H <sup>3</sup> | PDP127H <sup>3</sup> |

## TAPERCUT

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>4-0</b>           |
|--|---------------------|---------------|----------------------|
| <b>V-5</b> 17mm<br><b>V-5</b><br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | PDP137H <sup>3</sup> |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>6</sup> Two dozen packets per box

<sup>15</sup> This product has a loop length that is half of the length specified here.

 **Multipass® NEEDLES**


## VICRYL RAPIDE Suture

**Precision Cosmetic - Conventional Cutting PRIME**

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>         | <b>5-0</b>          |
|---|----------------------|---------------|--------------------|---------------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 18" ~ 45cm    | VR835 <sup>1</sup> | VR834 <sup>01</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>         | <b>5-0</b>         |
|---|----------------------|---------------|--------------------|--------------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 18" ~ 45cm    | VR845 <sup>1</sup> | VR844 <sup>1</sup> |

**Precision Point - Reverse Cutting**

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>5-0</b>         |
|--|----------------------|---------------|--------------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 18" ~ 45cm    | VR490 <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b>         | <b>5-0</b>         |
|--|----------------------|---------------|--------------------|--------------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 18" ~ 45cm    | VR494 <sup>1</sup> | VR493 <sup>1</sup> |

## Sutures : Absorbable / VICRYL RAPIDE Suture

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b>         |
|---|----------------------|---------------|--------------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 27" ~ 70cm    | VR935 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b>         | <b>4-0</b>         |
|---|----------------------|---------------|--------------------|--------------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed<br>Braid/coat | 18" ~ 45cm    | VR497 <sup>1</sup> | VR496 <sup>1</sup> |
|   | Undyed<br>Braid/coat | 27" ~ 70cm    |                    | VR426              |

## Taper Point

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>            |
|---|----------------------|---------------|---------------------|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Undyed<br>Braid/coat | 36" ~ 90cm    | VR9580 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b>         | <b>3-0</b>         |
|---|----------------------|---------------|--------------------|--------------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Undyed<br>Braid/coat | 36" ~ 90cm    | VR945 <sup>1</sup> | VR944 <sup>1</sup> |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Absorbable / VICRYL RAPIDE Suture

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>0</b>            |
|--|----------------------|---------------|---------------------|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Undyed<br>Braid/coat | 36" ~ 90cm    | VR1695 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>         |
|---|----------------------|---------------|--------------------|
| <b>RB-1 17mm</b><br>1/2 Circle<br> | Undyed<br>Braid/coat | 27" ~ 70cm    | VR214 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b>         | <b>3-0</b>         |
|---|----------------------|---------------|--------------------|--------------------|
| <b>SH 26mm</b><br>1/2 Circle<br> | Undyed<br>Braid/coat | 27" ~ 70cm    | VR417 <sup>1</sup> | VR416 <sup>1</sup> |

## Taper Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b>         | <b>3-0</b>         |
|---|----------------------|---------------|--------------------|--------------------|
| <b>CT-1 36mm</b><br><b>CP-1 36mm</b><br>1/2 Circle<br> | Undyed<br>Braid/coat | 54" ~ 135cm   | VR932 <sup>1</sup> | VR931 <sup>1</sup> |

## TAPERCUT

## Sutures : Absorbable / VICRYL RAPIDE Suture

| Needle  | Material             | Length     | 0                   | 2-0                | 3-0                |
|---|----------------------|------------|---------------------|--------------------|--------------------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Undyed<br>Braid/coat | 36" ~ 90cm | VR2218 <sup>1</sup> | VR917 <sup>1</sup> | VR916 <sup>1</sup> |

Nonsterile

<sup>1</sup> One dozen packets per box

 **Multipass**<sup>®</sup> MULTIPASS Needles

# Coated VICRYL Suture

## Conventional Cutting - Sternum

| Needle   | Material       | Length     | 3     |
|--|----------------|------------|-------|
| <b>CCS-1</b> 60mm<br>1/2 Circle<br> | Undyed Braided | 27" ~ 70cm | J808T |

## ETHIGUARD - Safety needle

| Needle   | Material       | Length       | 0   | 1  | 2-0  | 3-0                |
|--|----------------|--------------|---|--|--|--------------------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Undyed Braided | 8-18" ~ 45cm | JB840<br> | JB841 <sup>1-13</sup><br> | JB839<br> |                    |
|  | Undyed Braided | 27" ~ 70cm   | JB260   |  | JB259  |                    |
|  | Undyed Braided | 36" ~ 90cm   | JB946   | JB947 <sup>3</sup>   | JB945 <sup>3</sup>   | JB944 <sup>3</sup> |
|  | Violet Braided | 27" ~ 70cm   | JB340   | JB341 <sup>3</sup>   |  |                    |

| Needle   | Material       | Length       | 0   | 1   |
|--|----------------|--------------|---|---|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Undyed Braided | 8-18" ~ 45cm | JB724 <sup>1-13</sup><br><br>JB72433 | JB725 <sup>1-13</sup><br><br>JB72533 |
|  | Undyed Braided | 36" ~ 90cm   | JB978<br>JB97833  | JB977 <sup>3</sup><br>JB97730   |

| Needle | Material | Length | 2-0 | 3-0 |
|--------|----------|--------|-----|-----|
|        |          |        |     |     |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>   |
|--|-----------------|---------------|------------|--|
| <b>SHB</b> 26mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  |            | JB864<br> |
|  | Undyed Braided  | 27" ~ 70cm    | JB417      | JB416 <sup>3</sup>   |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b>         | <b>3-0</b>         |
|--|-----------------|---------------|--------------------|--------------------|
| <b>CTB-2</b> 36mm<br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | JB269 <sup>3</sup> | JB232 <sup>3</sup> |

## LIGAPAK Ligating Reel

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|----------|----------|------------|------------|------------|------------|
| <b>LIGAPAK Ligating Reel</b><br>N/A<br> | Undyed Braided  | 144" ~ 366cm  |          |          |            | J885G      |            |            |
|  | Undyed Braided  | 54" ~ 135cm   | J287G    |          | J286G      | J285G      | J284G      | J283G      |
|  | Violet Braided  | 54" ~ 135cm   | J207G    | J208G    | J206G      | J205G      | J204G      |            |

## MICROPOINT - Spatula

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|---|-----------------|---------------|------------|
| <b>TG160-6</b> 5.5mm<br><b>TG160-6</b><br>1/2 Circle<br> | Violet Braided  | 12" ~ 30cm    | J975G33    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>TG160-8</b> 5.5mm<br><b>TG160-8</b><br>1/2 Circle<br> | Violet Braided  | 12" ~ 30cm    | J575G      |            |
|   | Violet Braided  | 18" ~ 45cm    | J576G      |            |
|   | Violet Braided  | 8" ~ 20cm     |            | J574G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>TG100-8</b> 6.5mm<br><b>TG100-8</b><br>1/4 Circle<br> | Violet Braided  | 18" ~ 45cm    | J544G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|---|-----------------|---------------|------------|
| <b>TG140-6</b> 6.5mm<br><b>TG140-6</b><br>3/8 Circle<br> | Violet Braided  | 12" ~ 30cm    | J974G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b>                  | <b>9-0</b> |
|---|-----------------|---------------|------------|-----------------------------|------------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Violet Braided  | 12" ~ 30cm    | J566G      | J548G<br>J599G <sup>7</sup> | V549G      |
|   | Violet Braided  | 18" ~ 45cm    | J546G      |                             |            |
|   | Violet Braided  | 8" ~ 20cm     |            | J547G                       |            |

## **MICROPOINT-Reverse Cutting**

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>G-6</b> 8.0mm<br><b>G-6</b><br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | J551G      |

## Precision Cosmetic - Conventional Cutting PRIME

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b>          | <b>5-0</b>          |
|--|-----------------|---------------|---------------------|---------------------|
| <b>PC-12</b> 16mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J814G <sup>65</sup> | J813G <sup>65</sup> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b>          | <b>5-0</b>          | <b>6-0</b>          |
|--|-----------------|---------------|---------------------|---------------------|---------------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br><br> | Undyed Braided  | 18" ~ 45cm    | J835G <sup>65</sup> | J834G <sup>65</sup> | J833G <sup>65</sup> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br><br> | Undyed Braided  | 18" ~ 45cm    | J845G      | J844G      | J843G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b>     | <b>5-0</b> |
|--|-----------------|---------------|------------|----------------|------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br><br> | Undyed Braided  | 18" ~ 45cm    | J824G      | J823G<br>J823H | J822G      |

## Precision Point - Reverse Cutting

| <b>Needle</b>                      | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|------------------------------------|-----------------|---------------|------------|------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J504G      | J503G      |

| <b>Needle</b>                          | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br><br> | Undyed Braided  | 18" ~ 45cm    | J507G      | J506G      |

| <b>Needle</b>                      | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|------------------------------------|-----------------|---------------|------------|------------|
| <b>PS-5</b> 13mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J594G      | J593G      |

| <b>Needle</b>                      | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|------------------------------------|-----------------|---------------|------------|------------|
| <b>PS-6</b> 11mm<br>1/2 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J511G      | J510G      |

| <b>Needle</b>                                      | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>PS-4C</b> 16mm<br>1/2 Circle Compound Curve<br> | Undyed Braided  | 18" ~ 45cm    | J657G      | J656G      |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> | <b>7-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J490G      | J489G      | J488G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b>     | <b>5-0</b>     | <b>6-0</b>     |
|--|-----------------|---------------|----------------|----------------|----------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | J464G          | J463G          |                |
|  | Undyed Braided  | 18" ~ 45cm    | J494G<br>J494H | J493G<br>J493H | J492G<br>J492H |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>       | <b>4-0</b>     |
|---|-----------------|---------------|------------------|----------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J683G<br>J683H   | J682G<br>J682H |
|   | Undyed Braided  | 27" ~ 70cm    | J936H<br>J936H05 | J935H          |


| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>     | <b>4-0</b>     | <b>5-0</b>     |
|---|-----------------|---------------|------------|----------------|----------------|----------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    |            | J497G<br>J497H | J496G<br>J496H | J495G<br>J495H |
|   | Undyed Braided  | 27" ~ 70cm    | J428H      | J427H          | J426H          |                |


| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b>     | <b>6-0</b> |
|---|-----------------|---------------|----------------|------------|
| <b>PS-3</b> 16mm<br>3/8 Circle<br> | Undyed Braided  | 18" ~ 45cm    | J500G<br>J500H | J499G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|----------|------------|------------|
| <b>PSL 30mm</b><br>3/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | J610H    |            |            |
|  | Undyed Braided  | 27" ~ 70cm    |          | J596H      | J502H      |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>           | <b>1</b> | <b>2-0</b> | <b>2</b> |
|---|-----------------|---------------|--------------------|----------|------------|----------|
| <b>CP 40mm</b><br>1/2 Circle<br> | Violet Braided  | 27" ~ 70cm    | J479H <sup>3</sup> | J480H    | J478H      |          |
|   | Violet Braided  | 36" ~ 90cm    | J485H              | J486H    |            |          |
|   | Undyed Braided  | 27" ~ 70cm    |                    | J196H    |            | J195H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b>           | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|----------|--------------------|------------|------------|
| <b>CP-1 36mm</b><br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | J267H    | J268H              | J266H      |            |
|   | Violet Braided  | 27" ~ 70cm    | J467H    | J468H              | J466H      |            |
|   | Violet Braided  | 36" ~ 90cm    | J473H    | J474H <sup>3</sup> | J472H      | J471H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b>   | <b>3-0</b>   |
|---|-----------------|---------------|----------|----------|--|--|
| <b>CP-2 26mm</b><br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  |          |          | J762D <sup>13</sup><br> | J761D <sup>13</sup><br> |
|   | Undyed Braided  | 27" ~ 70cm    | J870H    | J871H    | J869H  | J868H  |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> |
|--|-----------------|---------------|----------|----------|
| <b>CPX 48mm</b><br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | J568H    | J569H    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>J-1 17mm</b><br>1/2 Circle<br> | Undyed Braided  | 27" ~ 70cm    | J455H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  |
|---|-----------------|---------------|---|
| <b>OS-2 17mm</b><br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | J704T <sup>4</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>  |
|---|-----------------|---------------|---|--|---|
| <b>OS-4 22mm</b><br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | J707T <sup>4</sup><br> | J708T<br> | J706T <sup>4</sup><br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | J749T<br>              |  | J748T<br>              |
|   | Undyed Braided  | 27" ~ 70cm    | J694H   | J695H  | J693H <sup>4</sup>  |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|---------------|-----------------|---------------|----------|----------|------------|
|---------------|-----------------|---------------|----------|----------|------------|

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b> | <b>2-0</b>   |
|---|-----------------|---------------|--|----------|--|
| <b>OS-6</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | J711T <sup>4</sup><br> |          | J710T<br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | J754T<br>              |          | J750T<br> |
|   | Undyed Braided  | 27" ~ 70cm    | J534H  | J535H    | J533H  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>  | <b>2-0</b>  | <b>2</b>   |
|---|-----------------|---------------|--|---|---|--|
| <b>OS-8</b> 40mm<br>1/2 Circle<br> | Violet Braided  | 3-18" ~ 45cm  | J717T <sup>4</sup><br> | J718T<br> | J716T <sup>4</sup><br> | J719T<br> |
|   | Undyed Braided  | 3-18" ~ 45cm  | J756T<br>              | J757T<br> |   |  |
|   | Undyed Braided  | 27" ~ 70cm    | J698H  | J699H   |   |  |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>   |
|--|-----------------|---------------|------------|--|
| <b>X-1</b> 22mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |            | J790D <sup>13</sup><br> |
|  | Violet Braided  | 27" ~ 70cm    | J461H      | J460H  |
|  | Undyed Braided  | 27" ~ 70cm    | J459H      | J458H  |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>C-3</b> 13mm<br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | J386H      | J385H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>   |
|---|-----------------|---------------|--|
| <b>FS 26mm</b><br>3/8 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J787D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|------------|------------|
| <b>FS-1 24mm</b><br>3/8 Circle<br> | Undyed Braided  | 27" ~ 70cm    | J443H      | J442H      | J441H      |
|   | Violet Braided  | 27" ~ 70cm    |            |            | J451H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|------------|
| <b>FS-2 19mm</b><br>3/8 Circle<br> | Violet Braided  | 18" ~ 45cm    | J393H      | J392H      | J391H      |
|   | Undyed Braided  | 27" ~ 70cm    | J423H      | J422H      | J421H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>FSL 30mm</b><br>3/8 Circle<br> | Undyed Braided  | 27" ~ 70cm    | J589H      | J588H      |

## SABRELOC - Center Point Spatula

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>RD-1 8.0mm</b><br><b>RD-1</b><br>1/4 Circle<br> | Violet Braided  | 12" ~ 30cm    | J446G      |

## **SABRELOC - Conventional Spatula**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>S-28</b> 8.0mm<br><b>S-28</b><br>1/2 Circle<br> | Violet Braided  | 18" ~ 45cm    | J562G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>S-29</b> 8.0mm<br><b>S-29</b><br>1/4 Circle<br> | Violet Braided  | 12" ~ 30cm    | J556G      |
|   | Violet Braided  | 18" ~ 45cm    | J555G      |

## **SABRELOC - Spatula**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b>       | <b>6-0</b> |
|---|-----------------|---------------|------------------|------------|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br>1/4 Circle<br> | Violet Braided  | 18" ~ 45cm    | J571G            | J570G      |
|   | Violet Braided  | 8" ~ 20cm     | J591G<br>J591G05 | J590G      |
|   | Undyed Braided  | 18" ~ 45cm    | J671G            | J670G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>S-24</b> 8.0mm<br><b>S-24</b><br>1/4 Circle<br> | Violet Braided  | 12" ~ 30cm    | J553G      | J552G      |
|   | Violet Braided  | 8" ~ 20cm     | J579G      | J578G      |

## **Standard and Short Length Sutures**

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b>           | <b>2</b>           | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|--|-----------------|---------------|----------|--------------------|--------------------|------------|------------|------------|
| <b>Standard and Short Length Sutures</b><br>N/A<br> | Violet Braided  | 3-18" ~ 45cm  | J636H    |                    |                    | J635H      | J634H      | J633H      |
|  | Violet Braided  | 54" ~ 135cm   | J616H    | J617H <sup>3</sup> | J618H <sup>3</sup> | J615H      | J614H      |            |
|  | Undyed Braided  | 3-18" ~ 45cm  | J646H    |                    |                    | J645H      | J644H      | J643H      |
|  | Undyed Braided  | 54" ~ 135cm   | J608H    |                    |                    | J607H      | J606H      |            |

## **Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>KS 60mm</b><br>Straight<br> | Violet Braided  | 27" ~ 70cm    | J523H      | J522H      |
|   | Undyed Braided  | 27" ~ 70cm    | J663H      | J662H      |

## **SUTUPAK Pre-Cut Sutures**

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>                    | <b>1</b>                    | <b>2-0</b>                           | <b>3-0</b>                           | <b>4-0</b>                  |
|--|-----------------|---------------|-----------------------------|-----------------------------|--------------------------------------|--------------------------------------|-----------------------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Violet Braided  | 6-18" ~ 45cm  | J106G<br>J106T <sup>5</sup> |                             |                                      |                                      |                             |
|  | Violet Braided  | 12-18" ~ 45cm | J906G                       | J907G                       | J105G<br>J105T <sup>3</sup><br>J905T | J104G<br>J104T <sup>3</sup><br>J904T | J103G<br>J103T <sup>5</sup> |
|  | Undyed Braided  | 6-18" ~ 45cm  | J112G<br>J112T <sup>5</sup> | J113G<br>J113T <sup>5</sup> |                                      |                                      |                             |
|  | Undyed Braided  | 12-18" ~ 45cm | J912G                       | J913G                       | J111G<br>J111T <sup>3</sup><br>J911T | J110G<br>J110T <sup>3</sup><br>J910T | J109G<br>J109T <sup>5</sup> |

## Taper Point

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>  | <b>2-0</b>  | <b>3-0</b> |
|---|-----------------|---------------|---|---|---|------------|
| <b>CT 40mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J752D <sup>13</sup><br> | J753D <sup>13</sup><br> | J751D <sup>13</sup><br> |            |
|   | Violet Braided  | 27" ~ 70cm    | J352H   | J353H   | J351H   |            |
|   | Violet Braided  | 36" ~ 90cm    | J358H   | J359H   | J357H   | J356H      |
|   | Undyed Braided  | 27" ~ 70cm    | J280H   | J281H   | J275H<br>J275H33  |            |
|   | Undyed Braided  | 36" ~ 90cm    | J958H<br>J958H33  | J959H<br>J959H33  | J957H   | J956H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>  | <b>2-0</b>  | <b>3-0</b>  | <b>4-0</b> |
|---|-----------------|---------------|---|---|---|---|------------|
| <b>CT-1 36mm</b><br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J740D <sup>13</sup><br> | J741D<br> | J739D<br> | J738D<br> |            |
|   | Violet Braided  | 27" ~ 70cm    | J340H   | J341H   | J339H   | J338H   | J337H      |
|   | Violet Braided  | 3-27" ~ 70cm  | J430T   |   |   |   |            |
|   | Violet Braided  | 8-27" ~ 70cm  | JJ31G <sup>4</sup><br>  |   |   |   |            |
|   | Violet Braided  | 36" ~ 90cm    | J346H   | J347H   | J345H   | J344H   |            |
|   | Undyed Braided  | 8-18" ~ 45cm  | J840D<br>               | J841D<br> | J839D<br> | J838D<br> |            |
|   | Undyed Braided  | 27" ~ 70cm    | J260H   | J261H   | J259H   | J258H   |            |
|   | Undyed Braided  | 8-27" ~ 70cm  | JJ41G<br>               | JJ40G<br> | JJ42G<br> |   |            |
|   | Undyed Braided  | 36" ~ 90cm    | J946H<br>J946H33  | J947H   | J945H   | J944H   |            |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b> | <b>2-0</b>   | <b>3-0</b>       |
|---|-----------------|---------------|---|----------|--|------------------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J727D <sup>13</sup><br> |          | J726D <sup>13</sup><br> |                  |
|   | Violet Braided  | 27" ~ 70cm    | J334H   | J335H    | J333H  | J332H            |
|   | Undyed Braided  | 27" ~ 70cm    | J270H   |          | J269H  | J232H<br>J232H33 |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|----------|------------|------------|
| <b>CT-3</b> 22mm<br>1/2 Circle<br> | Violet Braided  | 27" ~ 70cm    | J329H    | J328H      | J327H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   | <b>3-0</b> |
|--|-----------------|---------------|---|--|--|------------|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Undyed Braided  | 8-18" ~ 45cm  | J724D <sup>13</sup><br> |  | J723D<br> |            |
|  | Undyed Braided  | 36" ~ 90cm    | J978H   | J977H  | J979H  | J980H      |
|  | Violet Braided  | 8-18" ~ 45cm  | J764D<br>               | J765D <sup>13</sup><br> |  |            |
|  | Violet Braided  | 27" ~ 70cm    | J364H   | J365H  | J363H  |            |
|  | Violet Braided  | 36" ~ 90cm    | J370H   | J371H  | J369H  |            |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b>   |
|---|-----------------|---------------|---|--|--|
| <b>MO-4</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J701D <sup>13</sup><br> | J702D <sup>13</sup><br> | J700D <sup>13</sup><br> |
|   | Violet Braided  | 27" ~ 70cm    | J436H   | J437H  | J435H <sup>4</sup>   |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>  | <b>2-0</b>  |
|---|-----------------|---------------|--|---|---|
| <b>MO-5</b> 31mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J801D <sup>13</sup><br> |   |   |
|   | Violet Braided  | 8-27" ~ 70cm  | JJ81G<br>               | JJ80G <sup>4</sup><br> | JJ82G <sup>4</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b>   |  |  |  |
|---|-----------------|---------------|--|--|--|--|
| <b>MO-6</b> 26mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J789D <sup>13</sup><br> |  |  |  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>   | <b>4-0</b>   | <b>5-0</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|--|--|------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |            | J713D <sup>13</sup><br> |  |            |            |
|   | Violet Braided  | 27" ~ 70cm    | J306H      | J305H  | J304H  | J303H      |            |
|   | Undyed Braided  | 8-18" ~ 45cm  |            |  | J714D<br> |            |            |
|   | Undyed Braided  | 27" ~ 70cm    |            | J215H  | J214H  | J213H      | J212H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b>   | <b>3-0</b>   | <b>4-0</b>   |
|---|-----------------|---------------|----------|--|--|--|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  |          | J775D <sup>13</sup><br> | J774D <sup>13</sup><br> | J773D<br> |
|   | Violet Braided  | 27" ~ 70cm    |          | J317H  | J316H  | J315H  |
|   | Violet Braided  | 8-27" ~ 70cm  |          | J785G<br>               | J784G<br>               | J783G<br> |
|   | Undyed Braided  | 8-18" ~ 45cm  |          |  | J864D<br>               | J864D33  |
|   | Undyed Braided  | 27" ~ 70cm    | J418H    | J417H  | J416H  | J415H  |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|
| <b>SH</b> 26mm<br><b>SH</b><br>1/2 Circle<br> | Violet Braided  | 36" ~ 90cm    | J527H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>   | <b>4-0</b>   |
|---|-----------------|---------------|--|--|
| <b>SH-1</b> 22mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J772D <sup>13</sup><br> | J771D <sup>13</sup><br> |
|   | Violet Braided  | 27" ~ 70cm    | J311H  | J310H  |
|   | Undyed Braided  | 27" ~ 70cm    | J219H  | J218H  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>   | <b>5-0</b> |
|---|-----------------|---------------|--|------------|
| <b>TF</b> 13mm<br>1/2 Circle<br> | Violet Braided  | 8-18" ~ 45cm  | J743D <sup>13</sup><br> |            |
|   | Undyed Braided  | 27" ~ 70cm    | J434H  | J433H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>1</b> | <b>2</b>         |
|---|-----------------|---------------|----------|------------------|
| <b>TP-1</b> 65mm<br>1/2 Circle<br> | Violet Braided  | 4-27" ~ 70cm  | J650G    | J649G            |
|   | Undyed Braided  | 2-27" ~ 70cm  |          | J849G<br>J849G05 |
|   | Undyed Braided  | 54" ~ 135cm   |          | J880T            |

## Sutures : Absorbable / Coated VICRYL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|--|-----------------|---------------|----------|----------|------------|
| <b>XLH 70mm</b><br>1/2 Circle<br> | Violet Braided  | 27" ~ 70cm    | J582G    | J583G    | J581G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|---|-----------------|---------------|------------|
| <b>BV130-4</b> 5.0mm<br>3/8 Circle<br> | Violet Braided  | 5" ~ 13cm     | J405G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|--|-----------------|---------------|------------|
| <b>BV130-5</b> 6.5mm<br>3/8 Circle<br><br> | Violet Braided  | 5" ~ 13cm     | J401G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---|-----------------|---------------|----------|
| <b>UR-4</b> 40mm<br>5/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | J381H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|-----------------|---------------|----------|------------|
| <b>UR-5</b> 36mm<br>5/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | J376H    | J375H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|-----------------|---------------|----------|------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Violet Braided  | 27" ~ 70cm    | J603H    | J602H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>9-0</b> |
|---|---------------------|---------------|------------|
| <b>BV100-3</b> 3.9mm<br>3/8 Circle<br> | Violet Monofilament | 5" ~ 13cm     | V402G      |

## Taper Point - Reverse Cutting

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|
| <b>CT</b> 40mm<br><b>FSL</b> 30mm<br>1/2 Circle<br> | Undyed Braided  | 54" ~ 135cm   | J940H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>CT-1</b> 36mm<br><b>CP-1</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 54" ~ 135cm   | J532H      | J531H      |
|   | Undyed Braided  | 54" ~ 135cm   |            | J931H      |

## **TAPERCUT**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|----------|----------|------------|------------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Violet Braided  | 36" ~ 90cm    | J518H    | J519H    | J517H      | J516H      |
|   | Undyed Braided  | 36" ~ 90cm    |          |          | J917H      | J916H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|--|-----------------|---------------|------------|
| <b>V130-4</b> 4.2mm<br>1/4 Circle<br> | Violet Braided  | 9" ~ 23cm     | J408G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|--|-----------------|---------------|------------|
| <b>V130-5</b> 6.5mm<br>1/4 Circle<br> | Violet Braided  | 9" ~ 23cm     | J409G      |

## **ULTIMA - Spatula**

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>10-0</b>      |
|---|---------------------|---------------|------------------|
| <b>CS160-6</b> 5.5mm<br><b>CS160-6</b><br>1/2 Circle<br> | Violet Monofilament | 12" ~ 30cm    | V448G<br>V448G33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---------------|-----------------|---------------|-------------|------------|
|               |                 |               |             |            |

## Sutures : Absorbable / Coated VICRYL Suture


| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---|---------------------|---------------|-------------|------------|
| <b>CS160-8</b> 5.5mm<br><b>CS160-8</b><br>1/2 Circle<br> | Violet Monofilament | 12" ~ 30cm    | V450G       | V449G      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>10-0</b>      |
|--|---------------------|---------------|------------------|
| <b>CS90-6</b> 7.0mm<br>1/4 Circle<br> | Violet Monofilament | 4" ~ 10cm     | V980G<br>V980G33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>10-0</b>      | <b>9-0</b>       |
|---|---------------------|---------------|------------------|------------------|
| <b>CS140-6</b> 6.5mm<br>3/8 Circle<br> | Violet Monofilament | 4" ~ 10cm     | V960G<br>V960G05 |                  |
|   | Violet Braided      | 4" ~ 10cm     |                  | V950G<br>V950G33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>10-0</b> |
|---|---------------------|---------------|-------------|
| <b>CSB-6</b> 5.0mm<br>Bi-Curve<br> | Violet Monofilament | 4" ~ 10cm     | V966G       |

- <sup>1</sup> One dozen packets per box
- <sup>3</sup> Three dozen packets per box
- <sup>4</sup> Single strand delivery organizer
- <sup>5</sup> Single strand delivery
- <sup>7</sup> Beaded retraction sutures
- <sup>13</sup> Safety organizer tray
- <sup>65</sup> Pliabilized

 **Multipass<sup>®</sup>** NEEDLES MULTIPASS Needles

 **CR** CONTROL RELEASE

# MONOCRYL Suture

## ETHIGUARD - Safety needle

| Needle   | Material            | Length     | 0     | 1     | 2-0   | 3-0   |
|--|---------------------|------------|-------|-------|-------|-------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | YB258 |       | YB257 |       |
|  | Undyed Monofilament | 36" ~ 90cm | YB946 | YB947 |       | YB944 |

| Needle  | Material            | Length     | 0                |
|---|---------------------|------------|------------------|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Undyed Monofilament | 36" ~ 90cm | YB978<br>YB97833 |

| Needle   | Material            | Length     | 3-0                |
|--|---------------------|------------|--------------------|
| <b>SHB</b> 26mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm | YB416 <sup>3</sup> |

## Precision Cosmetic - Conventional Cutting PRIME

| Needle   | Material            | Length     | 4-0   |
|--|---------------------|------------|-------|
| <b>PC-12</b> 16mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm | Y814G |


## Sutures : Absorbable / MONOCRYL Suture


| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b>       |
|---|---------------------|---------------|------------|------------|------------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y835G      | Y834G      | Y833G<br>Y833G33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b>       | <b>5-0</b> |
|---|---------------------|---------------|------------------|------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y845G<br>Y845G05 | Y844G      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y823G      | Y822G      |

## Precision Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>6-0</b>       |
|---|---------------------|---------------|------------------|
| <b>PS-6</b> 11mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y510G<br>Y510G33 |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b>       | <b>6-0</b>       |
|--|---------------------|---------------|------------------|------------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y490G<br>Y490G05 | Y489G<br>Y489G33 |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b>       | <b>6-0</b>       |
|--|---------------------|---------------|------------|------------------|------------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Violet Monofilament | 18" ~ 45cm    | Y464G      | Y463G            |                  |
|  | Undyed Monofilament | 18" ~ 45cm    | Y494G      | Y493G<br>Y493G33 | Y492G<br>Y492G33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    |            | Y682H      |
|   | Undyed Monofilament | 27" ~ 70cm    | Y936H      | Y935H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b>       | <b>5-0</b> |
|---|---------------------|---------------|------------|------------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Y497G      | Y496G            | Y495G      |
|   | Undyed Monofilament | 27" ~ 70cm    | Y427H      | Y426H<br>Y426H33 |            |
|   | Violet Monofilament | 18" ~ 45cm    |            | Y513G            |            |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|---------------------|---------------|----------|------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Y267H    | Y266H      |

## Straight Cutting Needles

| Needle  | Material            | Length     | 3-0   |
|---|---------------------|------------|-------|
| <b>KS</b> 60mm<br>Straight<br> | Undyed Monofilament | 27" ~ 70cm | Y523H |

## SUTUPAK Pre-Cut Sutures

| Needle   | Material            | Length       | 0                  | 1                  | 2-0                | 3-0                | 4-0                |
|--|---------------------|--------------|--------------------|--------------------|--------------------|--------------------|--------------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Violet Monofilament | 6-18" ~ 45cm | Y910G <sup>5</sup> | Y911G <sup>5</sup> | Y909G <sup>5</sup> | Y908G <sup>5</sup> | Y907G <sup>5</sup> |

## Taper Point

| Needle  | Material            | Length       | 0   | 1     | 2-0   | 3-0   |
|---|---------------------|--------------|---|-------|-------|-------|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm | Y752D<br> |       |       |       |
|   | Violet Monofilament | 27" ~ 70cm   | Y352H   | Y353H | Y351H | Y350H |
|   | Violet Monofilament | 36" ~ 90cm   | Y358H   | Y359H | Y357H | Y356H |

| Needle | Material | Length | 0 | 1 | 2-0 | 3-0 |
|--------|----------|--------|---|---|-----|-----|
|        |          |        |   |   |     |     |

## Sutures : Absorbable / MONOCRYL Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>   | <b>1</b>   | <b>2-0</b>   | <b>3-0</b>   |
|---|---------------------|---------------|--|--|--|--|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm  | Y740D<br> |  | Y739D<br> | Y738D<br> |
|   | Violet Monofilament | 27" ~ 70cm    | Y340H  | Y341H  | Y339H  | Y338H  |
|   | Violet Monofilament | 4-27" ~ 70cm  | YY31G<br> | YY30G<br> | YY32G<br> |  |
|   | Violet Monofilament | 36" ~ 90cm    | Y346H  | Y347H  | Y345H  | Y344H  |
|   | Undyed Monofilament | 36" ~ 90cm    | Y946H  | Y947H  | Y945H  | Y944H  |


| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>1</b>         | <b>2-0</b> |
|--|---------------------|---------------|----------|------------------|------------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Y334H    | Y335H<br>Y335H33 | Y333H      |
|  | Violet Monofilament | 36" ~ 90cm    |          |                  | Y762H      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>1</b>         | <b>2-0</b>       |
|--|---------------------|---------------|----------|------------------|------------------|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | Y398H    | Y399H<br>Y399H33 | Y397H<br>Y397H33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|---|---------------------|---------------|--|--|
| <b>MO-4</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm  | Y701D<br> |  |
|   | Violet Monofilament | 27" ~ 70cm    | Y436H  | Y437H  |
|   | Violet Monofilament | 4-27" ~ 70cm  | YY71G<br> | YY70G<br> |

## Sutures : Absorbable / MONOCRYL Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>   |
|---|---------------------|---------------|--|
| <b>MO-5</b> 31mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm  | Y801D<br> |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Y215H      | Y214H      | Y213H      |
|   | Violet Monofilament | 27" ~ 70cm    | Y305H      | Y304H      | Y303H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|----------|------------|------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Y318H    | Y317H      | Y316H      | Y315H      |
|   | Violet Monofilament | 36" ~ 90cm    |          | Y732H      |            |            |
|   | Undyed Monofilament | 27" ~ 70cm    |          | Y417H      | Y416H      | Y415H      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>3-0</b> |
|--|---------------------|---------------|------------|
| <b>SH</b> 26mm<br><b>SH</b><br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm    | Y527H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>   | <b>4-0</b> |
|---|---------------------|---------------|--|------------|
| <b>SH-1</b> 22mm<br>1/2 Circle<br> | Undyed Monofilament | 8-18" ~ 45cm  | Y872D<br> |            |
|   | Undyed Monofilament | 27" ~ 70cm    | Y219H  | Y218H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>TF</b> 13mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Y434H      | Y433H      | Y432H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|---------------------|---------------|----------|------------|
| <b>UR-5</b> 36mm<br>5/8 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Y378H    | Y377H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---|---------------------|---------------|----------|------------|------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Y606H    | Y605H      | Y604H      |

## Taper Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>CT-1</b> 36mm<br><b>CP-1</b> 36mm<br>1/2 Circle<br> | Undyed Monofilament | 54" ~ 135cm   | Y932H      | Y931H      |

**TAPERCUT**

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> |
|---|---------------------|---------------|------------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Y817H      |
|   | Undyed Monofilament | 36" ~ 90cm    | Y917H      |

<sup>3</sup> Three dozen packets per box<sup>5</sup> Single strand delivery
 **Multipass<sup>®</sup>** NEEDLES

 **CONTROL RELEASE**

# PDS II Suture

## Conventional Cutting - Sternum

| Needle   | Material            | Length     | 2     |
|--|---------------------|------------|-------|
| <b>CCS-1</b> 60mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | Z807T |

## Endoscopic Needles

| Needle   | Material            | Length     | 4-0              |
|--|---------------------|------------|------------------|
| <b>ST-4</b> 19mm<br><b>ST-4</b><br>Straight Taper Point<br> | Violet Monofilament | 14" ~ 35cm | Z420G<br>Z420G33 |

## ETHIGUARD - Safety needle

| Needle   | Material            | Length       | 0   | 1     | 2-0                |
|--|---------------------|--------------|---|-------|--------------------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm | ZB740 <sup>113</sup><br> |       |                    |
|  | Violet Monofilament | 27" ~ 70cm   | ZB340   |       |                    |
|  | Violet Monofilament | 36" ~ 90cm   | ZB346   | ZB347 | ZB345 <sup>3</sup> |

| Needle | Material | Length | 0 | 1 |
|--------|----------|--------|---|---|
|        |          |        |   |   |

## Sutures : Absorbable / PDS II Suture

| Needle   | Material            | Length     | 0                | 1                |
|--|---------------------|------------|------------------|------------------|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Violet Monofilament | 36" ~ 90cm | ZB370<br>ZB37033 | ZB371<br>ZB37133 |

## Precision Cosmetic - Conventional Cutting PRIME

| Needle   | Material            | Length     | 6-0   |
|--|---------------------|------------|-------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br><br> | Undyed Monofilament | 18" ~ 45cm | Z833G |


| Needle   | Material            | Length     | 4-0              | 5-0              |
|--|---------------------|------------|------------------|------------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br><br> | Undyed Monofilament | 18" ~ 45cm | Z845G<br>Z845G31 | Z844G<br>Z844G31 |

| Needle   | Material            | Length     | 3-0   | 4-0   |
|--|---------------------|------------|-------|-------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br><br> | Undyed Monofilament | 18" ~ 45cm | Z824G | Z823G |

## Precision Point - Reverse Cutting

| Needle | Material | Length | 4-0 | 5-0 |
|--------|----------|--------|-----|-----|
|        |          |        |     |     |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b>       | <b>5-0</b>       |
|---|---------------------|---------------|------------------|------------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Z504G<br>Z504G33 | Z503G<br>Z503G33 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Z507G      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>6-0</b>       |
|--|---------------------|---------------|------------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Violet Monofilament | 18" ~ 45cm    | Z487G            |
|  | Undyed Monofilament | 18" ~ 45cm    | Z489G<br>Z489G31 |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>4-0</b> | <b>5-0</b>       |
|--|---------------------|---------------|------------|------------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Z494G      | Z493G<br>Z493G31 |
|  | Violet Monofilament | 18" ~ 45cm    |            | Z463G            |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>       | <b>4-0</b> |
|---|---------------------|---------------|------------------|------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Z683G<br>Z683G31 | Z682G      |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b>       | <b>4-0</b>       | <b>5-0</b> |
|---|---------------------|---------------|------------------|------------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | Z497G<br>Z497G31 | Z496G<br>Z496G31 | Z495G      |
|   | Violet Monofilament | 18" ~ 45cm    |                  | Z513G            |            |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2</b>         |
|---|---------------------|---------------|------------------|
| <b>CP</b> 40mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z195T<br>Z195T31 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b>       |
|---|---------------------|---------------|----------|----------|------------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z467H    | Z468H    | Z466H<br>Z466H31 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b> | <b>2-0</b>   | <b>3-0</b>   |
|---|---------------------|---------------|----------|--|--|
| <b>CP-2</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm  |          | Z762D <sup>13</sup><br> | Z761D <sup>13</sup><br> |
|   | Violet Monofilament | 27" ~ 70cm    | Z970H    | Z969H  | Z968H<br>Z968H33   |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>1</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>1</b>         |
|--|---------------------|---------------|------------------|
| <b>CPX</b> 48mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z569T<br>Z569T31 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>1</b>         |
|---|---------------------|---------------|------------------|
| <b>OS-4</b> 22mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z695H<br>Z695H31 |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>0</b>         | <b>1</b> |
|---|---------------------|---------------|------------------|----------|
| <b>OS-6</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z534T<br>Z534T31 | Z535H    |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Z443H      | Z442H      | Z441H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Undyed Monofilament | 27" ~ 70cm    | Z423H      | Z422H      | Z421H      |

## Straight Cutting Needles

| Needle  | Material            | Length     | 3-0   | 4-0   |
|---|---------------------|------------|-------|-------|
| <b>KS</b> 60mm<br>Straight<br> | Undyed Monofilament | 27" ~ 70cm | Z663H | Z662H |

## Taper Point

| Needle  | Material            | Length     | 0                | none    | 1                | 2-0              |
|---|---------------------|------------|------------------|---------|------------------|------------------|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | Z352H<br>Z352H31 |         | Z353H<br>Z353H31 |                  |
|   | Violet Monofilament | 36" ~ 90cm | Z358T            |         | Z359T<br>Z359T31 | Z357H<br>Z357H33 |
|   |                     | 36" ~ 90cm |                  | Z358T05 |                  |                  |

| Needle  | Material            | Length       | 0  | 1  | 2-0  | 3-0              |
|---|---------------------|--------------|--|--|--|------------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm | Z740D <sup>13</sup><br> | Z741D <sup>13</sup><br> | Z739D <sup>13</sup><br> |                  |
|   | Violet Monofilament | 27" ~ 70cm   | Z340H<br>Z340H31   | Z341H<br>Z341H31   | Z339H<br>Z339H31   | Z338H<br>Z338H31 |
|   | Violet Monofilament | 36" ~ 90cm   | Z346H<br>Z346H31   | Z347H<br>Z347H31   |  |                  |
|   | Undyed Monofilament | 27" ~ 70cm   |  |  | Z259H  | Z258H            |

| Needle  | Material            | Length     | 0                | 2-0              | 3-0   |
|---|---------------------|------------|------------------|------------------|-------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm | Z334H<br>Z334H31 | Z333H<br>Z333H31 | Z332H |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|--|--------------------------|---------------|--|--|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Violet Monofilament      | 8-18" ~ 45cm  | Z776D <sup>13</sup><br> | Z765D <sup>13</sup><br> |
|  | Violet Monofilament      | 36" ~ 90cm    | Z370T  | Z371T<br>Z371T31   |
|  | Violet Monofilament      | 60" ~ 150cm   | Z990G05<br>Z990G31   |  |
|  | Violet Monofilament Loop | 60" ~ 150cm   | Z990G  |  |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|---------------------|---------------|------------|------------|------------|
| <b>RB-1 17mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z305H      | Z304H      | Z303H      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b> |
|--|---------------------|---------------|------------|
| <b>RB-1 17mm</b><br><b>RB-1</b><br>1/2 Circle<br> | Violet Monofilament | 30" ~ 75cm    | Z320H      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>5-0</b> | <b>6-0</b>       |
|--|---------------------|---------------|------------|------------------|
| <b>RB-2 13mm</b><br><b>RB-2</b><br>1/2 Circle<br> | Violet Monofilament | 30" ~ 75cm    | Z148H      | Z149H<br>Z149H33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|
|               |                 |               |            |            |            |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>2-0</b>  | <b>3-0</b>   | <b>4-0</b>   |
|---|---------------------|---------------|---|--|--|
| <b>SH 26mm</b><br>1/2 Circle<br> | Violet Monofilament | 8-18" ~ 45cm  | Z775D <sup>13</sup><br> | Z774D <sup>13</sup><br> | Z773D <sup>13</sup><br> |
|   | Violet Monofilament | 27" ~ 70cm    | Z317H   | Z316H  | Z315H  |
|   | Undyed Monofilament | 27" ~ 70cm    |   | Z416H  |  |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>SH-1 22mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z311H      | Z310H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>6-0</b> |  |  |
|---|---------------------|---------------|------------|--|--|
| <b>TF 13mm</b><br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z432H      |  |  |

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>none</b> | <b>0</b> | <b>1</b> |
|---|--------------------------|---------------|-------------|----------|----------|
| <b>TP-1 65mm</b><br>1/2 Circle<br> |                          | 2-27" ~ 70cm  | Z690G31     |          |          |
|   | Violet Monofilament      | 2-27" ~ 70cm  |             |          | Z690G    |
|   | Violet Monofilament      | 48" ~ 122cm   |             |          | Z880G31  |
|   | Violet Monofilament      | 54" ~ 135cm   |             |          | Z879G    |
|   | Violet Monofilament      | 60" ~ 150cm   |             | Z991G33  |          |
|   | Violet Monofilament Loop | 60" ~ 150cm   |             | Z991G    |          |
|   | Violet Monofilament Loop | 96" ~ 244cm   |             |          | Z880G    |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>2-0</b>       |
|--|---------------------|---------------|------------------|
| <b>UCL</b> 14mm<br>1/2 Circle<br> | Violet Monofilament | 27" ~ 70cm    | Z114H<br>Z114H33 |

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>0</b> | <b>1</b> |
|--|--------------------------|---------------|----------|----------|
| <b>XLH</b> 70mm<br>1/2 Circle<br> | Violet Monofilament      | 27" ~ 70cm    | Z582G    | Z583G    |
|  | Violet Monofilament      | 54" ~ 135cm   |          | Z877G    |
|  | Violet Monofilament Loop | 96" ~ 244cm   |          | Z881G    |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>6-0</b> |
|--|---------------------|---------------|------------|
| <b>BV</b> 11mm<br><b>BV</b><br>3/8 Circle<br> | Violet Monofilament | 30" ~ 75cm    | Z123H      |

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|---|---------------------|---------------|------------|------------|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Violet Monofilament | 24" ~ 60cm    |            | Z155H      |
|   | Violet Monofilament | 30" ~ 75cm    | Z117H      | Z135H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Absorbable / PDS II Suture

| <b>Needle</b>  | <b>Material</b>        | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|--|------------------------|---------------|------------|------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Violet<br>Monofilament | 30" ~ 75cm    | Z126H      | Z127H      |

## TAPERCUT

| <b>Needle</b>  | <b>Material</b>        | <b>Length</b> | <b>4-0</b> |
|--|------------------------|---------------|------------|
| <b>V-5</b> 17mm<br><b>V-5</b><br>1/2 Circle<br> | Violet<br>Monofilament | 36" ~ 90cm    | Z137H      |

## Trocar Point

| <b>Needle</b>   | <b>Material</b>        | <b>Length</b> | <b>2-0</b>         |
|---|------------------------|---------------|--------------------|
| <b>STP-10</b> 254mm<br><b>STP-10</b><br>Straight<br> | Violet<br>Monofilament | 27" ~ 70cm    | Z997G <sup>1</sup> |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>13</sup> Safety organizer tray

 **MULTIPASS** Needles

 CONTROL RELEASE

# Surgical Gut Suture - Chromic

## Blunt Point

| Needle  | Material             | Length     | 1                |
|---|----------------------|------------|------------------|
| <b>BP</b> 80mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm | 48G <sup>1</sup> |

| Needle  | Material             | Length     | 0   |
|---|----------------------|------------|-----|
| <b>BP-1</b> 65mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm | 47T |

## LIGAPAK Dispensing Reel

| Needle   | Material             | Length      | 0     | 1     | 2-0   | 3-0   | 4-0   |
|--|----------------------|-------------|-------|-------|-------|-------|-------|
| <b>LIGAPAK Dispensing Reel</b><br>N/A<br> | Chromic Monofilament | 54" ~ 135cm | L114G | L115G | L113G | L112G | L111G |

## MICROPOINT - Spatula

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>7-0</b> |
|---|----------------------|---------------|------------|
| <b>TG100-8</b> 6.5mm<br><b>TG100-8</b><br>1/4 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1744G      |

| <b>Needle</b>   | <b>Material</b>           | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|---|---------------------------|---------------|------------|------------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Blue Chromic Monofilament | 18" ~ 45cm    | 1731G      | 1797G      |

## MICROPOINT-Reverse Cutting

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>none</b> | <b>4-0</b> |
|---|----------------------|---------------|-------------|------------|
| <b>G-2</b> 8.0mm<br><b>G-2</b><br>1/2 Circle<br> |                      | 18" ~ 45cm    | 798G34      |            |
|   | Chromic Monofilament | 18" ~ 45cm    |             | 798G       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>6-0</b>     |
|--|----------------------|---------------|----------------|
| <b>G-1</b> 11mm<br>3/8 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 796G<br>796G33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>6-0</b>     |
|--|----------------------|---------------|----------------|
| <b>G-1</b> 11mm<br><b>G-1</b><br>3/8 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 790G<br>790G33 |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|----------------------|---------------|------------|------------|
| <b>G-3</b> 13mm<br><b>G-3</b><br>3/8 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 793G       | 792G       |

| <b>Needle</b>   | <b>Material</b>           | <b>Length</b> | <b>6-0</b> |
|---|---------------------------|---------------|------------|
| <b>G-6</b> 8.0mm<br><b>G-6</b><br>3/8 Circle<br> | Chromic Monofilament      | 18" ~ 45cm    | 794G       |
|   | Blue Chromic Monofilament | 18" ~ 45cm    | 795G       |

## Precision Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|----------------------|---------------|------------|------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1659G      | 1658G      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|
|               |                 |               |            |            |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|----------------------|---------------|------------|------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br><br> | Chromic Monofilament | 18" ~ 45cm    | 1643G      | 1642G      | 1641G      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>       | <b>5-0</b> |
|---|----------------------|---------------|------------------|------------|
| <b>PS-5</b> 13mm<br>1/2 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1653G<br>1653G33 | 1652G      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>6-0</b> |
|---|----------------------|---------------|------------|
| <b>PS-6</b> 11mm<br>1/2 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1816G      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|----------------------|---------------|------------|------------|
| <b>PS-4C</b> 16mm<br>1/2 Circle Compound Curve<br> | Chromic Monofilament | 18" ~ 45cm    | 2643G      | 2644G      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b> | <b>5-0</b>     |
|---|----------------------|---------------|------------|----------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br><br> | Chromic Monofilament | 18" ~ 45cm    | 1654G      | 687G<br>687G33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> |
|---|----------------------|---------------|----------|
| <b>PS 26mm</b><br>3/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 598H     |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|----------------------|---------------|------------|------------|------------|
| <b>PS-2 19mm</b><br>3/8 Circle<br> | Chromic Monofilament | 18" ~ 45cm    |            |            | 1637G      |
|   | Chromic Monofilament | 27" ~ 70cm    | 1621H      | 1638H      |            |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---|----------------------|---------------|------------|------------|
| <b>PS-3 16mm</b><br>3/8 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1636G      | 1635G      |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>2-0</b> |
|--|----------------------|---------------|------------|
| <b>PSL 30mm</b><br>3/8 Circle<br> | Chromic Monofilament | 36" ~ 90cm    | 929H       |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>1</b>                    | <b>2-0</b>        | <b>3-0</b> |
|---|----------------------|---------------|-----------------------------|-------------------|------------|
| <b>CP 40mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 805H <sup>3</sup><br>805H33 | 833H <sup>3</sup> | 832H       |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

---

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|---|----------------------|---------------|----------|----------|------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 817H     |          | 816H       |
|   | Chromic Monofilament | 36" ~ 90cm    |          | 975H     |            |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>2-0</b>     | <b>3-0</b> |
|---|----------------------|---------------|----------|----------------|------------|
| <b>CP-2</b> 26mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 888H     | 887H<br>887H33 | 886H       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b>       |
|--|----------------------|---------------|------------------|
| <b>J-1</b> 17mm<br>1/2 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 724G             |
|  | Chromic Monofilament | 27" ~ 70cm    | U213H<br>U213H33 |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b> |
|--|----------------------|---------------|------------|
| <b>M-1</b> 13mm<br><b>M-1</b><br>1/2 Circle<br> | Chromic Monofilament | 12" ~ 30cm    | 744G       |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>5-0</b> |
|---|----------------------|---------------|------------|
| <b>C-3</b> 13mm<br>C-3<br><br>3/8 Circle<br> | Chromic Monofilament | 12" ~ 30cm    | 723G       |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b>     |
|---|----------------------|---------------|----------------|
| <b>FS-1</b> 24mm<br><br>3/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 637H<br>637H33 |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b>     | <b>4-0</b>     | <b>5-0</b> |
|---|----------------------|---------------|----------------|----------------|------------|
| <b>FS-2</b> 19mm<br><br>3/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 636H<br>636H33 | 635H<br>635H33 | 634G       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>3-0</b> |
|--|----------------------|---------------|------------|
| <b>FSL</b> 30mm<br><br>3/8 Circle<br> | Chromic Monofilament | 36" ~ 90cm    | 932H       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b>     |
|--|----------------------|---------------|----------------|
| <b>M-2</b> 13mm<br><br>3/8 Circle<br> | Chromic Monofilament | 12" ~ 30cm    | 752G<br>752G33 |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b> |
|--|----------------------|---------------|------------|
| <b>M-2</b> 13mm<br><b>M-2</b><br>3/8 Circle<br> | Chromic Monofilament | 12" ~ 30cm    | 751G       |

## SABRELOC - Conventional Spatula

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>       |
|---|----------------------|---------------|------------------|
| <b>S-4</b> 8.0mm<br><b>S-4</b><br>1/4 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1793G<br>1793G33 |

## SABRELOC - Spatula

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b> |
|---|----------------------|---------------|------------|
| <b>S-2</b> 8.0mm<br><b>S-2</b><br>1/2 Circle<br> | Chromic Monofilament | 18" ~ 45cm    | 1798G      |

| <b>Needle</b>   | <b>Material</b>           | <b>Length</b> | <b>5-0</b>       | <b>6-0</b> |
|---|---------------------------|---------------|------------------|------------|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br>1/4 Circle<br> | Chromic Monofilament      | 12" ~ 30cm    | 1766G            |            |
|   | Chromic Monofilament      | 18" ~ 45cm    | 1792G<br>1792G33 | 1791G33    |
|   | Blue Chromic Monofilament | 18" ~ 45cm    |                  | 1791G      |

## **Standard and Short Length Sutures**

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>2</b> | <b>3-0</b> | <b>3</b> | <b>4-0</b> |
|--|----------------------|---------------|----------|----------|------------|----------|------------|----------|------------|
| <b>Standard and Short Length Sutures</b><br>N/A<br> | Chromic Monofilament | 54" ~ 135cm   | S114H    | S115H    | S113H      | S116H    | S112H      | S117H    | S111H      |

## **Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b>     |
|---|----------------------|---------------|----------------|
| <b>KS</b> 60mm<br>Straight<br> | Chromic Monofilament | 27" ~ 70cm    | 654H<br>654H33 |

## **Straight Taper Point**

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b> |
|---|----------------------|---------------|------------|
| <b>ST-1</b> 45mm<br><b>SH</b> 26mm<br>Straight<br> | Chromic Monofilament | 36" ~ 90cm    | G152H      |

## **SUTUPAK Pre-Cut Sutures**

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|--|----------------------|---------------|----------|----------|------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Chromic Monofilament | 12-18" ~ 45cm | SG14T    | SG15T    | SG13T      |

## Taper Point

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>       | <b>1</b>       | <b>2-0</b>     | <b>3-0</b> |
|---|----------------------|---------------|----------------|----------------|----------------|------------|
| <b>CT 40mm</b><br><br>1/2 Circle | Chromic Monofilament | 3-27" ~ 70cm  | 3802T          | 3803T          |                |            |
|   | Chromic Monofilament | 27" ~ 70cm    | 802H<br>802H33 | 803H<br>803H33 | 801H<br>801H33 | 800H       |
|   | Chromic Monofilament | 36" ~ 90cm    | 914H<br>914H33 | 915H<br>915H33 | 913H<br>913H33 | 912H       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|----------------------|---------------|------------|------------|
| <b>CT 40mm</b><br><b>CT</b><br><br>1/2 Circle | Chromic Monofilament | 54" ~ 135cm   | 991H       | 990H       |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>   | <b>1</b>   | <b>2-0</b>     | <b>2</b>       | <b>3-0</b>     | <b>4-0</b> |
|---|----------------------|---------------|--|--|----------------|----------------|----------------|------------|
| <b>CT-1 36mm</b><br><br>1/2 Circle | Chromic Monofilament | 8-18" ~ 45cm  | CC41G <sup>13</sup><br> | CC40G <sup>13</sup><br> |                |                |                |            |
|   | Chromic Monofilament | 3-27" ~ 70cm  | 3812T  | 3813G  | 3811T          |                |                |            |
|   | Chromic Monofilament | 27" ~ 70cm    | 812H<br>812H33   | 813H<br>813H33   | 811H<br>811H33 | 814H<br>814H33 | 810H<br>810H33 |            |
|   | Chromic Monofilament | 36" ~ 90cm    | 924H   | 925H<br>925H33   | 923H           |                | 922H           | 921H       |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|----------|----------|------------|------------|
|               |                 |               |          |          |            |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>       | <b>1</b>       | <b>2-0</b>     | <b>3-0</b> |
|---|----------------------|---------------|----------------|----------------|----------------|------------|
| <b>CT-2 26mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 884H<br>884H33 | 885H<br>885H33 | 883H<br>883H33 | 882H       |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>2-0</b>     | <b>3-0</b>     |
|---|----------------------|---------------|----------|----------------|----------------|
| <b>CT-3 22mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 894H     | 893H<br>893H33 | 892H<br>892H33 |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>0</b>       | <b>1</b>       | <b>2-0</b>     | <b>3-0</b> |
|--|----------------------|---------------|----------------|----------------|----------------|------------|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 864H           | 865H<br>865H33 | 863H<br>863H33 | 862H       |
|  | Chromic Monofilament | 36" ~ 90cm    | 904H<br>904H33 | 905H<br>905H33 | 903H           |            |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>1</b>         |
|---|----------------------|---------------|------------------|
| <b>MO-2 40mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 4915H<br>4915H33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|---------------|-----------------|---------------|----------|----------|------------|
|               |                 |               |          |          |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>  | <b>1</b>   | <b>2-0</b> |
|---|----------------------|---------------|---|--|------------|
| <b>MO-4</b> 36mm<br>1/2 Circle<br> | Chromic Monofilament | 8-18" ~ 45cm  | CC02G <sup>13</sup><br> | CC01G <sup>13</sup><br> |            |
|   | Chromic Monofilament | 27" ~ 70cm    | 4924H<br>4924H33  | 4925H<br>4925H33   | 4923H      |
|   | Chromic Monofilament | 3-27" ~ 70cm  |   | 5589T  |            |
|   | Chromic Monofilament | 8-27" ~ 70cm  |   | CC38G <sup>13</sup><br> |            |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|---|----------------------|---------------|--|--|
| <b>MO-5</b> 31mm<br>1/2 Circle<br> | Chromic Monofilament | 8-18" ~ 45cm  | CC04G <sup>13</sup><br> | CC03G <sup>13</sup><br> |
|   | Chromic Monofilament | 27" ~ 70cm    | 4994H  | 4995H<br>4995H33   |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|----------------------|---------------|------------|------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | U204H      | U203H      | U202H      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>2-0</b>         | <b>3-0</b> | <b>4-0</b> |
|---|----------------------|---------------|----------|--------------------|------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | G124H    | G123H <sup>3</sup> | G122H      | G121H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>2-0</b>         | <b>3-0</b> |
|--|----------------------|---------------|--------------------|------------|
| <b>SH 26mm</b><br><b>SH</b><br>1/2 Circle<br> | Chromic Monofilament | 36" ~ 90cm    | G173H <sup>3</sup> | G172H      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b>         |
|---|----------------------|---------------|------------|------------|--------------------|
| <b>SH-1 22mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | G182H      | G181H      | G180H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>1</b> | <b>2</b> |
|---|----------------------|---------------|----------|----------|
| <b>TP-1 65mm</b><br>1/2 Circle<br> | Chromic Monofilament | 4-27" ~ 70cm  | M760G    |          |
|   | Chromic Monofilament | 60" ~ 150cm   | GL30G    | GL31G    |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>0</b> |
|--|----------------------|---------------|----------|
| <b>XLH 70mm</b><br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 44T      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Absorbable / Surgical Gut Suture - Chromic

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>4-0</b>         |
|---|----------------------|---------------|--------------------|
| <b>BB</b> 17mm<br>3/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 4887G <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>5-0</b> |
|--|----------------------|---------------|------------|
| <b>C-1</b> 13mm<br>3/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | K895H      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---|----------------------|---------------|----------|------------|
| <b>UR-4</b> 40mm<br>5/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | U256H    | U255H      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|---|----------------------|---------------|----------|----------|------------|
| <b>UR-5</b> 36mm<br>5/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | U246H    | U247H    | U245H      |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---|----------------------|---------------|----------|------------|------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | N879H    | N878H      | N877H      |

## **Taper Point - Reverse Cutting**

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>2-0</b>        |
|--|----------------------|---------------|-------------------|
| <b>CT</b> 40mm<br><b>FSL</b> 30mm<br>1/2 Circle<br> | Chromic Monofilament | 54" ~ 135cm   | 977H <sup>3</sup> |

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|----------------------|---------------|------------|------------|
| <b>CT-1</b> 36mm<br><b>CP-1</b> 36mm<br>1/2 Circle<br> | Chromic Monofilament | 54" ~ 135cm   | 983H       | 982H       |

## **TAPERCUT**

| <b>Needle</b>   | <b>Material</b>      | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|---|----------------------|---------------|----------|----------|------------|------------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | 838H     | 839H     | 837H       |            |
|   | Chromic Monofilament | 36" ~ 90cm    | 944H     | 945H     | 943H       | 942H       |

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>4-0</b> |
|--|----------------------|---------------|------------|
| <b>V-5</b> 17mm<br>1/2 Circle<br> | Chromic Monofilament | 27" ~ 70cm    | U233H      |

## **ULTIMA - Spatula**

| <b>Needle</b>   | <b>Material</b>         | <b>Length</b> | <b>7-0</b>       |
|---|-------------------------|---------------|------------------|
| <b>CS175-8</b> 7.0mm<br><b>CS175-8</b><br>1/2 Circle<br> | Chromic<br>Monofilament | 18" ~ 45cm    | 1745G<br>1745G34 |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>13</sup> Safety organizer tray

 **Multipass<sup>®</sup>** NEEDLES MULTIPASS Needles

 **CR/** CONTROL RELEASE

# Surgical Gut Suture - Plain

## LIGAPAK Dispensing Reel

| Needle   | Material           | Length      | 0     | 2-0   | 3-0   | 4-0   |
|--|--------------------|-------------|-------|-------|-------|-------|
| <b>LIGAPAK Dispensing Reel</b><br>N/A<br> | Plain Monofilament | 54" ~ 135cm | L104G | L103G | L102G | L101G |

## MICROPOINT - Spatula

| Needle  | Material           | Length     | 6-0   |
|---|--------------------|------------|-------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm | 1735G |

## MICROPOINT-Reverse Cutting

| Needle   | Material           | Length     | 6-0            |
|--|--------------------|------------|----------------|
| <b>G-1</b> 11mm<br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm | 774G<br>774G33 |

| Needle   | Material           | Length     | 6-0            |
|--|--------------------|------------|----------------|
| <b>G-1</b> 11mm<br><b>G-1</b><br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm | 770G<br>770G33 |

## Sutures : Absorbable / Surgical Gut Suture - Plain

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|--------------------|---------------|------------|------------|
| <b>G-3</b> 13mm<br><b>G-3</b><br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 773G       | 772G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>6-0</b> |
|---|--------------------|---------------|------------|
| <b>G-6</b> 8.0mm<br><b>G-6</b><br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 775G       |

## Precision Cosmetic - Conventional Cutting PRIME

| <b>Needle</b>  | <b>Material</b>      | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|--|----------------------|---------------|------------|------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br><br> | Plain Fast Absorbing | 18" ~ 45cm    | 1915G      | 1916G      |
|  | Plain Monofilament   | 18" ~ 45cm    |            | 1916G33    |

## Precision Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>5-0</b> |
|---|--------------------|---------------|------------|
| <b>P-2</b> 8.0mm<br>1/2 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 658G       |

## Sutures : Absorbable / Surgical Gut Suture - Plain

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|--------------------|---------------|------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br><br><b>Multipass NEEDLES</b> | Plain Monofilament | 18" ~ 45cm    | 1633G      | 1632G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>6-0</b> |
|--|--------------------|---------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br><br><b>Multipass NEEDLES</b> | Plain Monofilament | 18" ~ 45cm    | 1657G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>4-0</b>       | <b>5-0</b>     |
|--|--------------------|---------------|------------------|----------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br><br><b>Multipass NEEDLES</b> | Plain Monofilament | 18" ~ 45cm    | 1644G<br>1644G33 | 686G<br>686G33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|---|--------------------|---------------|------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br><br><b>Multipass NEEDLES</b> | Plain Monofilament | 27" ~ 70cm    | 1640H      |


| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b> | <b>4-0</b>       |
|---|--------------------|---------------|------------|------------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br><br><b>Multipass NEEDLES</b> | Plain Monofilament | 18" ~ 45cm    |            | 1627H<br>1627H33 |
|   | Plain Monofilament | 27" ~ 70cm    | 1630H      |                  |

## Sutures : Absorbable / Surgical Gut Suture - Plain


| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>5-0</b>       |
|---|--------------------|---------------|------------------|
| <b>PS-3</b> 16mm<br>3/8 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 1626G<br>1626G33 |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b>     |
|---|--------------------|---------------|----------------|
| <b>CP-1</b> 36mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | 898H<br>898H33 |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>4-0</b>     |
|--|--------------------|---------------|----------------|
| <b>M-1</b> 13mm<br>1/2 Circle<br> | Plain Monofilament | 12" ~ 30cm    | 742G<br>742G33 |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|--|--------------------|---------------|------------|
| <b>X-1</b> 22mm<br>1/2 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 612G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b>       | <b>3-0</b>       |
|---|--------------------|---------------|------------------|------------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Plain Monofilament | 27" ~ 70cm    | H823H<br>H823H33 | H810H<br>H810H33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b>       | <b>4-0</b>       | <b>5-0</b>       |
|---|--------------------|---------------|------------------|------------------|------------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Plain Monofilament | 27" ~ 70cm    | H822H<br>H822H33 | H821H<br>H821H33 | H820G<br>H820G33 |

## **SABRELOC - Spatula**

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>5-0</b>       |
|---|--------------------|---------------|------------------|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br>1/4 Circle<br> | Plain Monofilament | 12" ~ 30cm    | 1775G            |
|   | Plain Monofilament | 18" ~ 45cm    | 1772G<br>1772G33 |

## **Standard and Short Length Sutures**

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> |
|--|--------------------|---------------|----------|----------|------------|------------|
| <b>Standard and Short Length Sutures</b><br>N/A<br> | Plain Monofilament | 54" ~ 135cm   | S104H    | S105H    | S103H      | S102H      |

## **Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|---|--------------------|---------------|------------|
| <b>KS</b> 60mm<br>Straight<br> | Plain Monofilament | 27" ~ 70cm    | 665H       |

## Sutures : Absorbable / Surgical Gut Suture - Plain

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>4-0</b>       |
|---|--------------------|---------------|------------------|
| <b>SC-1</b> 13mm<br>Straight<br> | Plain Monofilament | 18" ~ 45cm    | 1824H<br>1824H33 |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>4-0</b> |
|--|--------------------|---------------|------------|
| <b>SC-1</b> 13mm<br><b>SC-1</b><br>Straight<br> | Plain Monofilament | 18" ~ 45cm    | 1828H      |

## SUTUPAK Pre-Cut Sutures

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|--|--------------------|---------------|------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Plain Monofilament | 12-18" ~ 45cm | SG2T       |

## Taper Point

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b>     | <b>3-0</b>   |
|---|--------------------|---------------|----------------|--|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Plain Monofilament | 8-18" ~ 45cm  |                | CP53G <sup>13</sup><br> |
|   | Plain Monofilament | 27" ~ 70cm    | 853H<br>853H33 | 852H   |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|
|               |                 |               |          |            |            |

## Sutures : Absorbable / Surgical Gut Suture - Plain

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>0</b>       | <b>2-0</b>     | <b>3-0</b>     |
|---|--------------------|---------------|----------------|----------------|----------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Plain Monofilament | 3-18" ~ 45cm  |                |                | 5503T          |
|   | Plain Monofilament | 27" ~ 70cm    | 844H<br>844H33 | 843H<br>843H33 | 842H<br>842H33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>0</b> | <b>2-0</b>       | <b>3-0</b>       |
|---|--------------------|---------------|----------|------------------|------------------|
| <b>CT-3</b> 22mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | N864H    | N863H<br>N863H33 | N862H<br>N862H33 |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>2-0</b>     | <b>3-0</b> |
|--|--------------------|---------------|----------------|------------|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | 872H<br>872H33 | 873H       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>4-0</b> |
|---|--------------------|---------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | U207H      |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---|--------------------|---------------|------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | G322H      | G321H      |

## Sutures : Absorbable / Surgical Gut Suture - Plain

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|--------------------|---------------|------------|------------|
| <b>XLH</b> 70mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | 53T        | 52T        |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b> |
|---|--------------------|---------------|------------|
| <b>UR-6</b> 26mm<br>5/8 Circle<br> | Plain Monofilament | 27" ~ 70cm    | N873H      |

## TAPERCUT

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|---|--------------------|---------------|------------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Plain Monofilament | 27" ~ 70cm    | 846H       |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|--|--------------------|---------------|------------|
| <b>V-7</b> 26mm<br>1/2 Circle<br> | Plain Monofilament | 18" ~ 45cm    | 849G       |

<sup>13</sup> Safety organizer tray



# PRONOVA Suture

## Taper Point

| Needle   | Material          | Length     | 3-0   |
|--|-------------------|------------|-------|
| <b>MH</b> 36mm<br><b>MH</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm | 3842H |

| Needle   | Material          | Length     | 3-0   | 4-0   | 5-0   |
|--|-------------------|------------|-------|-------|-------|
| <b>RB-1</b> 17mm<br><b>RB-1</b><br>1/2 Circle<br> | Blue Monofilament | 24" ~ 60cm |       |       | 3555H |
|  | Blue Monofilament | 36" ~ 90cm | 3558H | 3557H | 3556H |

| Needle   | Material          | Length     | 5-0   | 6-0   |
|--|-------------------|------------|-------|-------|
| <b>RB-2</b> 13mm<br><b>RB-2</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm | 3710H | 3711H |

| Needle   | Material          | Length     | 2-0   | 3-0   | 4-0   |
|--|-------------------|------------|-------|-------|-------|
| <b>SH</b> 26mm<br><b>SH</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm | 3523H | 3522H | 3521H |

## Sutures : Non Absorbable / PRONOVA Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>SH-1</b> 22mm<br><b>SH-1</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 3762H      | 3526H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>5-0</b>         |
|--|-------------------|---------------|--------------------|
| <b>TF</b> 13mm<br><b>TF</b><br>1/2 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 3205H<br><b>BV</b> |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-------------------|---------------|------------|------------|------------|
| <b>BB</b> 17mm<br><b>BB</b><br>3/8 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 3862H      | 3861H      |            |
|  | Blue Monofilament | 36" ~ 90cm    |            | 3581H      | 3580H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>BV</b> 11mm<br><b>BV</b><br>3/8 Circle<br> | Blue Monofilament | 4-24" ~ 60cm  |            | M3600      |
|  | Blue Monofilament | 30" ~ 75cm    | 3776H      |            |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Non Absorbable / PRONOVA Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b>   | <b>7-0</b>   |
|---|-------------------|---------------|--|--|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 3806H  | 3701H  |
|   | Blue Monofilament | 4-18" ~ 45cm  |  | M3701  |
|   | Blue Monofilament | 24" ~ 60cm    | 3305H<br> 3805H | 3304H<br> 3702H |
|   | Blue Monofilament | 4-24" ~ 60cm  | M3805  | M3304<br> M3702 |
|   | Blue Monofilament | 30" ~ 75cm    | 3709H  | 3703H  |
|   | Blue Monofilament | 4-30" ~ 75cm  | M3709  | M3703  |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>8-0</b>       |
|--|-------------------|---------------|------------------|
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b><br>3/8 Circle<br><br> | Blue Monofilament | 18" ~ 45cm    | 3730H<br>3730H33 |
|  | Blue Monofilament | 24" ~ 60cm    | 3732H<br>3732H33 |
|  | Blue Monofilament | 2-24" ~ 60cm  | M3739            |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>BV175-6</b> 8.0mm<br><b>BV175-6</b><br>3/8 Circle<br><br> | Blue Monofilament | 24" ~ 60cm    | 3735H      | 3741H      |
|  | Blue Monofilament | 4-24" ~ 60cm  | M3735      |            |
|  | Blue Monofilament | 2-24" ~ 60cm  | M3737      | M3743      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Non Absorbable / PRONOVA Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---|-------------------|---------------|------------|------------|
| <b>BV175-8</b> 9.3mm<br><b>BV175-8</b><br>3/8 Circle<br> | Blue Monofilament | 2-24" ~ 60cm  | M3747      | M3753      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>5-0</b>  | <b>6-0</b>  |
|--|-------------------|---------------|---|---|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 3717H   |   |
|  | Blue Monofilament | 24" ~ 60cm    | 3725H   | 3726H   |
|  | Blue Monofilament | 4-24" ~ 60cm  | M3725   | M3726   |
|  | Blue Monofilament | 30" ~ 75cm    |   | 3307H<br><br>3706H   |
|  | Blue Monofilament | 4-30" ~ 75cm  |   | M3307<br><br>M3706 |
|  | Blue Monofilament | 36" ~ 90cm    | 3321H<br><br>3720H |   |
|  | Blue Monofilament | 4-36" ~ 90cm  | M3720   |   |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---|-------------------|---------------|------------|------------|
| <b>BV175-7</b> 8.0mm<br><b>BV175-7</b><br>3/8 Circle<br> | Blue Monofilament | 2-18" ~ 45cm  | M3768      | M3769      |
|   | Blue Monofilament | 2-24" ~ 60cm  | M3766      | M3767      |
|   | Blue Monofilament | 4-30" ~ 75cm  | M3770      |            |

## TAPERCUT

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>V-7</b> 26mm<br><b>V-7</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 3976H      | 3975H      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|---|-------------------|---------------|------------|------------|
| <b>CC</b> 9.3mm<br><b>CC</b><br>3/8 Circle<br> | Blue Monofilament | 4-24" ~ 60cm  | M3807      |            |
|   | Blue Monofilament | 24" ~ 60cm    |            | 3704H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b> |
|--|-------------------|---------------|------------|
| <b>CC-1</b> 13mm<br><b>CC-1</b><br>3/8 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 3707H      |

 **Multipass<sup>®</sup>** NEEDLES MULTIPASS Needles  
**VB** VISI-BLACK

# ETHIBOND EXCEL Suture

## Conventional Cutting - Sternum

| Needle  | Material      | Length       | 5     |
|---|---------------|--------------|-------|
| CCS 48mm<br>1/2 Circle<br> | Green Braided | 4-30" ~ 75cm | MB47G |

## Precision Point - Reverse Cutting

| Needle  | Material      | Length     | 4-0   |
|---|---------------|------------|-------|
| PS-4 16mm<br>1/2 Circle<br><br> | White Braided | 18" ~ 45cm | X695G |

| Needle  | Material      | Length     | 6-0   |
|---|---------------|------------|-------|
| P-1 11mm<br>3/8 Circle<br> | Green Braided | 18" ~ 45cm | X697G |

| Needle   | Material      | Length     | 5-0   |
|--|---------------|------------|-------|
| P-3 13mm<br>3/8 Circle<br><br> | Green Braided | 18" ~ 45cm | X698G |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> | White Braided   | 18" ~ 45cm    | X692G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-3</b> 16mm<br>3/8 Circle<br> | Green Braided   | 18" ~ 45cm    | X681G      |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>1</b>  |
|---|-----------------|---------------|---|---|
| <b>OS-2</b> 17mm<br>1/2 Circle<br> | Green Braided   | 3-18" ~ 45cm  | CX74T <sup>5</sup><br> | CX75T <sup>5</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2</b> |
|---|-----------------|---------------|----------|----------|----------|
| <b>OS-4</b> 22mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X517H    | X518H    | X519H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>1</b> |
|---|-----------------|---------------|----------|
| <b>OS-6</b> 36mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X538H    |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

---

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>1</b> |
|---|-----------------|---------------|----------|
| <b>OS-8</b> 40mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X548H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>FS</b> 26mm<br>3/8 Circle<br> | Green Braided   | 18" ~ 45cm    | X664H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Green Braided   | 18" ~ 45cm    | X663H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Green Braided   | 18" ~ 45cm    | X683G      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2</b> | <b>5</b> |
|---------------|-----------------|---------------|----------|----------|
|               |                 |               |          |          |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2</b> | <b>5</b> |
|---|-----------------|---------------|----------|----------|
| <b>LR 75mm</b><br>3/8 Circle<br> | Green Braided   | 20" ~ 50cm    | X406T    | B409T    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2</b> | <b>5</b> |
|--|-----------------|---------------|----------|----------|
| <b>LR 75mm</b><br><b>LR</b><br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | X496T    | B499T    |

## Straight Cutting Needles

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|
| <b>KS 60mm</b><br>Straight<br> | Green Braided   | 30" ~ 75cm    | X622H      |

## SUTUPAK Pre-Cut Sutures

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|------------|------------|
|---------------|-----------------|---------------|----------|------------|------------|------------|------------|

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|----------|------------|------------|------------|------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Green Braided   | 6-18" ~ 45cm  | X186H    |            |            |            |            |
|  | Green Braided   | 12-18" ~ 45cm |          | X185H      | X184H      | X183H      | X182H      |
|  | Green Braided   | 6-30" ~ 75cm  | X306H    |            |            |            |            |
|  | Green Braided   | 12-30" ~ 75cm |          | X305H      | X304H      | X303H      |            |

## Taper Point

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b> | <b>2-0</b>   |
|--|-----------------|---------------|--|----------|--|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | CX21D <sup>13</sup><br> |          | CX22D <sup>13</sup><br> |
|  | Green Braided   | 30" ~ 75cm    | X424H  | X425H    | X423H  |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | CX27D <sup>13</sup><br> | CX26D <sup>13</sup><br> |
|   | Green Braided   | 30" ~ 75cm    | X412H  | X411H  |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|--|-----------------|---------------|--|--|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | CX31D <sup>13</sup><br> | CX30D <sup>13</sup><br> |
|  | Green Braided   | 30" ~ 75cm    | X864H  | X865H  |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b>     |
|--|-----------------|---------------|----------|----------------|
| <b>MH</b> 36mm<br><b>MH</b><br>1/2 Circle<br> | Green Braided   | 36" ~ 90cm    | X844H    | X843H<br>X843H |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-6</b> 26mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | CX45D <sup>13</sup><br> | CX46D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-7</b> 22mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | CX41D <sup>13</sup><br> | CX42D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|------------|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X873H      | X872H      | X871H      | X870H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|------------|
|               |                 |               |            |            |            |            |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>2-0</b>                 | <b>3-0</b> | <b>4-0</b>         | <b>5-0</b> |
|--|-----------------------|---------------|----------------------------|------------|--------------------|------------|
| <b>RB-1</b> 17mm<br><b>RB-1</b><br>1/2 Circle<br> | Green & White Braided | 10-30" ~ 75cm | 10X42 SX43 <sup>1,22</sup> |            |                    |            |
|  | Green & White Braided | 8-30" ~ 75cm  | PXX41 PXX43 <sup>20</sup>  |            |                    |            |
|  | White Braided         | 4-30" ~ 75cm  | MX523                      |            |                    |            |
|  | Green Braided         | 4-30" ~ 75cm  | MX553                      | MX552      | MX551 <sup>1</sup> |            |
|  | Green Braided         | 30" ~ 75cm    |                            |            | X551H              | X550H      |
|  | Green Braided         | 36" ~ 90cm    |                            | X558H      | X557H              | X556H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b>              | <b>3-0</b>              | <b>4-0</b>              |
|---|-----------------|---------------|----------|-------------------------|-------------------------|-------------------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  |          | CX12D <sup>13</sup> CR/ | CX13D <sup>13</sup> CR/ | CX14D <sup>13</sup> CR/ |
|   | Green Braided   | 30" ~ 75cm    | X834H    | X833H                   | X832H                   | X831H                   |
|   | Green Braided   | 4-30" ~ 75cm  |          | MX833                   |                         |                         |
|   | White Braided   | 4-30" ~ 75cm  |          | MX823                   |                         |                         |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|------------|
|               |                 |               |          |            |            |            |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>0</b>              | <b>2-0</b>   | <b>3-0</b> | <b>4-0</b> |
|--|-----------------------|---------------|-----------------------|--|------------|------------|
| <b>SH 26mm SH</b><br>1/2 Circle<br> | Green & White Braided | 10-30" ~ 75cm | PXX87 <sup>1 21</sup> | PXX76 <sup>1 23</sup><br>SXP76 <sup>1 21</sup><br>SXX86 <sup>1 20</sup><br>10X82 |            |            |
|  | Green & White Braided | 10-36" ~ 90cm |                       | 10X86<br>PXX86 <sup>1 21</sup><br>SXP86  |            |            |
|  | Green & White Braided | 6-36" ~ 90cm  |                       | PXX85  |            |            |
|  | Green Braided         | 4-30" ~ 75cm  |                       | MX563  |            |            |
|  | Green Braided         | 30" ~ 75cm    |                       | X563H  | X562H      |            |
|  | Green Braided         | 36" ~ 90cm    | X524H                 | SX14H <sup>20</sup><br>PX514<br>SXP02<br>X523H                                   | X522H      |            |
|  | White Braided         | 36" ~ 90cm    |                       | X513H  |            |            |
|  | Green                 | 36" ~ 90cm    |                       |  |            | X521H      |

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>2-0</b>         | <b>3-0</b> | <b>4-0</b> |
|--|-----------------------|---------------|--------------------|------------|------------|
| <b>SH-1 22mm SH-1</b><br>1/2 Circle<br> | Green Braided         | 4-30" ~ 75cm  | MX763 <sup>1</sup> |            |            |
|  | Green Braided         | 36" ~ 90cm    | X763H              | X762H      |            |
|  | Green & White Braided | 10-30" ~ 75cm | PXX92              |            |            |
|  | Green                 | 36" ~ 90cm    |                    |            | X761H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>2-0</b>   |
|--|-----------------------|---------------|--|
| <b>SH-2 20mm<br/>SH-2</b><br><br>1/2 Circle<br> | Green Braided         | 30" ~ 75cm    | PX82H <sup>21</sup><br>X582H<br>PX82H <sup>21</sup><br>X582H |
|  | Green Braided         | 36" ~ 90cm    | X559H<br>X559H   |
|  | White Braided         | 30" ~ 75cm    | PX83H<br>PX83H   |
|  | Green & White Braided | 10-30" ~ 75cm | PXX80  |
|  | White & Green Braided | 10-30" ~ 75cm | PXX82<br>SXP82 <sup>123</sup>                                |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>TF 13mm<br/>TF</b><br><br>1/2 Circle<br> | Green Braided   | 24" ~ 60cm    | X204H      | X205H      |


| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|--|-----------------|---------------|------------|
| <b>UCL 14mm</b><br><br>1/2 Circle<br> | Green Braided   | 27" ~ 70cm    | X114H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>1</b> |
|---------------|-----------------|---------------|----------|
|---------------|-----------------|---------------|----------|

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>1</b> |
|--|-----------------|---------------|----------|
| <b>XLH</b> 70mm<br><b>XLH</b><br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X190G    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>BB</b> 17mm<br><b>BB</b><br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | X881H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>BB</b> 17mm<br><b>BB</b><br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    |            | X581H      | X580H      |
|  | Green Braided   | 36" ~ 90cm    | X588H      |            |            |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> | <b>7-0</b>         |
|---|-----------------|---------------|------------|--------------------|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Green Braided   | 10-18" ~ 45cm |            | 10X07 <sup>1</sup> |
|   | Green Braided   | 18" ~ 45cm    |            | X701H              |
|   | Green Braided   | 24" ~ 60cm    | X805H      |                    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>C-1</b> 13mm<br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | X890H      | X889H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Green Braided   | 18" ~ 45cm    |            | X717H      | X706H      |
|  | Green Braided   | 30" ~ 75cm    | X726H      | X727H      |            |
|  | Green Braided   | 36" ~ 90cm    |            |            | X719H      |

## TAPERCUT

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---|-----------------|---------------|----------|
| <b>V-34</b> 36mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X444H    |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2</b> | <b>5</b> |
|---|-----------------|---------------|----------|----------|
| <b>V-37</b> 40mm<br>1/2 Circle<br> | Green Braided   | 4-30" ~ 75cm  | MX69G    | MB66G    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5</b>           |
|---|-----------------|---------------|--------------------|
| <b>V-40</b> 48mm<br>1/2 Circle<br> | Green Braided   | 4-30" ~ 75cm  | MB46G <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|
| <b>V-5</b> 17mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X951H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>V-7</b> 26mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | X964H      | X963H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|------------|
|               |                 |               |            |            |            |            |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>2-0</b>   | <b>3-0</b>          | <b>4-0</b> | <b>5-0</b> |
|--|-----------------------|---------------|--|---------------------|------------|------------|
| <b>V-5</b> 17mm<br><b>V-5</b><br>1/2 Circle<br> | White & Green Braided | 10-30" ~ 75cm | 10X52 <sup>1</sup><br>PXX52 <sup>1 21</sup><br>SXP52 <sup>1 23</sup><br>SXP54 <sup>1 22</sup><br>SXX54 <sup>1 20</sup> | 10X53               |            |            |
|  | White & Green Braided | 8-30" ~ 75cm  | PXX50<br>SXX50 <sup>20</sup>   | SXX53 <sup>20</sup> |            |            |
|  | White Braided         | 4-30" ~ 75cm  | MPX17  |                     |            |            |
|  | White Braided         | 30" ~ 75cm    | X917H  | X916H               |            |            |
|  | White Braided         | 36" ~ 90cm    | PX17H<br>SX17H   |                     |            |            |
|  | Green Braided         | 4-30" ~ 75cm  | MPX37  |                     |            |            |
|  | Green Braided         | 30" ~ 75cm    | PX937 <sup>21 3</sup><br>X937H   | X936H               | X935H      | X934H      |
|  | Green Braided         | 36" ~ 90cm    | PX54H<br>SX54H<br>X807H  | X806H               |            |            |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|------------|
|---------------|-----------------|---------------|----------|------------|------------|------------|

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>0</b> | <b>2-0</b>   | <b>3-0</b> | <b>4-0</b> |
|--|-----------------------|---------------|----------|--|------------|------------|
| <b>V-7</b> 26mm<br><b>V-7</b><br>1/2 Circle<br> | White & Green Braided | 10-30" ~ 75cm |          | 10X72 <sup>1</sup><br>PXX77 <sup>1 21</sup><br>SXP77 |            |            |
|  | White & Green Braided | 8-30" ~ 75cm  |          | PXX72  |            |            |
|  | White & Green Braided | 10-36" ~ 90cm |          | 10X66  |            |            |
|  | White & Green Braided | 8-36" ~ 90cm  |          | PXX79  |            |            |
|  | Green Braided         | 4-30" ~ 75cm  |          | MPX77<br>MX977 <sup>1</sup>                          |            |            |
|  | Green Braided         | 30" ~ 75cm    |          | PX977 <sup>21 3</sup><br>SXP07<br>X977H              | X976H      | X975H      |
|  | Green Braided         | 36" ~ 90cm    | X905H    | PX62H<br>X966H                                       |            |            |
|  | White Braided         | 30" ~ 75cm    |          | X987H  |            |            |
|  | White Braided         | 36" ~ 90cm    | PX65H    | PX64H  |            |            |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>CC-1</b> 13mm<br><b>CC-1</b><br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | X728H      |            |
|  | Green Braided   | 36" ~ 90cm    |            | X720H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>V-4</b> 17mm<br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | X942H      |

## Sutures : Non Absorbable / ETHIBOND EXCEL Suture

| <b>Needle</b>  | <b>Material</b>       | <b>Length</b> | <b>2-0</b>                                  | <b>3-0</b> | <b>4-0</b> |
|--|-----------------------|---------------|---|------------|------------|
| <b>V-4</b> 17mm<br><b>V-4</b><br>3/8 Circle<br> | Green Braided         | 4-30" ~ 75cm  | MPX27 <sup>1 21</sup><br>MX927 <sup>1</sup> |            |            |
|  | Green Braided         | 30" ~ 75cm    | X927H                                       | X926H      | X925H      |
|  | White Braided         | 4-30" ~ 75cm  | MX944                                       |            |            |
|  | White & Green Braided | 10-30" ~ 75cm | PXX22<br>SXP22 <sup>1 23</sup>              |            |            |

## Trocar Point

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>STP-10</b> 254mm<br><b>STP-10</b><br>Straight<br> | Green Braided   | 27" ~ 70cm    | X997G      |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>5</sup> Single strand delivery

<sup>13</sup> Safety organizer tray

<sup>20</sup> 1/8 x 1/8 x 1/16 (3.0mm x 3.0mm x 1.5 mm) Firm TFE Polymer Pledgets

<sup>21</sup> 1/4 x 1/8 x 1/16 (7.0mm x 3.0mm x 1.5 mm) Firm TFE Polymer Pledgets

<sup>22</sup> 1/8 x 1/8 x 1/16 (3.0mm x 3.0mm x 1.5 mm) Soft TFE Polymer Pledgets

<sup>23</sup> 1/4 x 1/8 x 1/16 (7.0mm x 3.0mm x 1.5 mm) Soft TFE Polymer Pledgets

 **Multipass<sup>®</sup>** NEEDLES

 CONTROL RELEASE

# ETHILON Suture

## Endoscopic Needles

| Needle  | Material           | Length    | 10-0  |
|---|--------------------|-----------|-------|
| <b>ST100-3</b><br>Straight Taper Point<br> | Black Monofilament | 5" ~ 13cm | 2849G |

## MICROPOINT - Spatula

| Needle  | Material           | Length     | 10-0  |
|---|--------------------|------------|-------|
| <b>TG160-4</b> 5.5mm<br><b>TG160-4</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm | 7707G |

| Needle  | Material           | Length    | 10-0  |
|---|--------------------|-----------|-------|
| <b>TG160-4-3M</b> 4.3mm<br><b>TG160-4-3M</b><br>1/2 Circle<br> | Black Monofilament | 8" ~ 20cm | 7711G |

| Needle  | Material           | Length     | 10-0             | 9-0   |
|---|--------------------|------------|------------------|-------|
| <b>TG160-6</b> 5.5mm<br><b>TG160-6</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm | 7756G            | 7760G |
|   | Black Monofilament | 6" ~ 15cm  | 7757G<br>7757G33 |       |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>8-0</b>                  |
|---|--------------------|---------------|-----------------------------|
| <b>TG175-8</b> 7.0mm<br><b>TG175-8</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm    | 1716G<br>1716G05<br>1716G33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>8-0</b>       |
|---|--------------------|---------------|------------------|
| <b>TG100-8</b> 6.5mm<br><b>TG100-8</b><br>1/4 Circle<br> | Black Monofilament | 12" ~ 30cm    | 1714G<br>1714G33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>TG140-6</b> 6.5mm<br><b>TG140-6</b><br>3/8 Circle<br> | Black Monofilament | 12" ~ 30cm    | 7770G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---|--------------------|---------------|-------------|------------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Black Monofilament | 12" ~ 30cm    | 7718G       | 7717G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>9-0</b>       |
|--|--------------------|---------------|------------------|
| <b>VAS100-4</b> 5.1mm<br>3/8 Circle<br> | Black Monofilament | 6" ~ 15cm     | 2890G<br>2890G33 |

## Precision Cosmetic - Conventional Cutting PRIME

| Needle   | Material           | Length     | 5-0                                       |
|--|--------------------|------------|---|
| <b>PC-12</b> 16mm<br>1/2 Circle<br> | Black Monofilament | 18" ~ 45cm | 1845G <sup>1</sup><br>1945G <sup>65</sup> |

| Needle  | Material           | Length     | 4-0                          | 5-0                          | 6-0                          |
|---|--------------------|------------|------------------------------|------------------------------|------------------------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm | 1854G<br>1954G <sup>65</sup> | 1855G<br>1955G <sup>65</sup> | 1856G<br>1956G <sup>65</sup> |

| Needle  | Material           | Length     | 4-0            | 5-0                          | 6-0                       |
|---|--------------------|------------|----------------|------------------------------|---------------------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm | 1864G<br>1964G | 1865G<br>1965G <sup>65</sup> | 1866G<br>1966G<br>1966G33 |


| Needle  | Material           | Length     | 3-0  | 4-0   | 5-0                                       |
|---|--------------------|------------|--|---|---|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm | 1893G <sup>1</sup><br>1993G <sup>1</sup> <sub>65</sub><br>1893G33<br>1993G33 | 1894G <sup>1</sup><br>1994G <sup>1</sup> <sub>65</sub><br>1994G33 | 1895G <sup>1</sup><br>1995G <sup>65</sup> |

## Precision Point - Reverse Cutting

| Needle | Material | Length | 4-0 | 5-0 | 6-0 |
|--------|----------|--------|-----|-----|-----|
|        |          |        |     |     |     |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b>       | <b>5-0</b> | <b>6-0</b> |
|---|---------------------|---------------|------------------|------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | 1603G<br>1603G33 |            |            |
|   | Black Monofilament  | 18" ~ 45cm    | 1662G            | 1661G      | 1660G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>none</b> | <b>5-0</b> | <b>6-0</b>   | <b>7-0</b> |
|--|--------------------|---------------|-------------|------------|--------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> |                    | 18" ~ 45cm    | 1696G34     |            |              |            |
|  | Green Monofilament | 18" ~ 45cm    |             | G695G      |              |            |
|  | Undyed             | 18" ~ 45cm    |             |            | 689G         |            |
|  | Black Monofilament | 18" ~ 45cm    |             |            | 697G<br>697H | 1696G      |
|  | Green              | 18" ~ 45cm    |             |            | G697G        |            |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>none</b> | <b>4-0</b>   | <b>5-0</b>   | <b>6-0</b> |
|--|---------------------|---------------|-------------|--------------|--------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> |                     | 18" ~ 45cm    | 698G34      |              |              |            |
|  | Undyed Monofilament | 18" ~ 45cm    |             | 691G         | 690G         |            |
|  | Black Monofilament  | 18" ~ 45cm    |             | 699G<br>699H | 698G<br>698H | 1698G      |
|  | Green               | 18" ~ 45cm    |             |              | G698G        | G696G      |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>7-0</b> |
|---|--------------------|---------------|------------|
| <b>P-6</b> 8.0mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 1647G      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b> |
|---|--------------------|---------------|------------|
| <b>PS 26mm</b><br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 585H       |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b>                |
|--|--------------------|---------------|---------------------------|
| <b>PS-1 24mm</b><br>3/8 Circle<br><br> | Black Monofilament | 18" ~ 45cm    | 1663G<br>1663H<br>1663H33 |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b>       | <b>5-0</b>     |
|--|---------------------|---------------|------------|------------|------------------|----------------|
| <b>PS-2 19mm</b><br>3/8 Circle<br><br> | Black Monofilament  | 18" ~ 45cm    | 593H       | 1669H      | 1667G<br>1667H   | 1666G<br>1666H |
|  | Undyed Monofilament | 18" ~ 45cm    |            |            | 1611G<br>1611G33 |                |
|  | Green               | 18" ~ 45cm    |            |            | G667G            | G666G          |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>5-0</b>     | <b>6-0</b> |
|--|--------------------|---------------|----------------|------------|
| <b>PS-3 16mm</b><br>3/8 Circle<br><br> | Black Monofilament | 18" ~ 45cm    | 1668G<br>1668H | 1665G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|--------------------|---------------|------------|------------|
| <b>PSL 30mm</b><br>3/8 Circle<br> | Black Monofilament | 30" ~ 75cm    | 1691H      | 562H       |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|--------------------|---------------|------------|------------|
| <b>PSLX 40mm</b><br>3/8 Circle<br> | Black Monofilament | 30" ~ 75cm    | 1697H      | 1683H      |

## Reverse Cutting

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> |
|--|--------------------|---------------|------------|
| <b>X-1 22mm</b><br>1/2 Circle<br> | Black Monofilament | 18" ~ 45cm    | 642G       |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>6-0</b> |
|--|--------------------|---------------|------------|
| <b>C-2 13mm</b><br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 667G       |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>5-0</b> |
|--|--------------------|---------------|------------|
| <b>C-3 13mm</b><br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 668G       |


| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b>   |
|---|--------------------|---------------|--------------|
| <b>FS 26mm</b><br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 664G<br>664H |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>3-0</b>   | <b>4-0</b> |
|---|--------------------|---------------|--------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 663G<br>663H | 1629H      |
|   | Black Monofilament | 30" ~ 75cm    | 669H         |            |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>4-0</b>   | <b>5-0</b>   |
|---|--------------------|---------------|--------------|--------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Black Monofilament | 18" ~ 45cm    | 662G<br>662H | 661G<br>661H |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>3-0</b> | <b>4-0</b>       |
|--|--------------------|---------------|------------|------------------|
| <b>FSL</b> 30mm<br>3/8 Circle<br> | Black Monofilament | 30" ~ 75cm    | 1671H      | 1670H<br>1670H33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b> | <b>3-0</b>       |
|---|--------------------|---------------|------------|------------------|
| <b>FSLX</b> 40mm<br>3/8 Circle<br> | Black Monofilament | 30" ~ 75cm    | 1674H      | 1673H<br>1673H33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2</b>                             |
|---|--------------------|---------------|--------------------------------------|
| <b>LR</b> 75mm<br>3/8 Circle<br> | Black Monofilament | 20" ~ 50cm    | 460T<br>470G <sup>40</sup><br>460T33 |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>1</b>       | <b>2</b>       |
|--|--------------------|---------------|----------------|----------------|
| <b>LR</b> 75mm<br><b>LR</b><br>3/8 Circle<br> | Black Monofilament | 30" ~ 75cm    | 489T<br>489T05 | 490T<br>490T05 |

## SABRELOC - Center Point Spatula

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>none</b> | <b>5-0</b>     |
|---|--------------------|---------------|-------------|----------------|
| <b>RD-1</b> 8.0mm<br><b>RD-1</b><br>1/4 Circle<br> |                    | 18" ~ 45cm    | 749G34      |                |
|   | Black Monofilament | 18" ~ 45cm    |             | 749G<br>749G33 |

## SABRELOC - Spatula

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>5-0</b> |
|---|--------------------|---------------|------------|
| <b>S-22</b> 8.0mm<br><b>S-22</b><br>1/2 Circle<br> | Black Monofilament | 18" ~ 45cm    | 1740G      |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>5-0</b> |
|---|--------------------|---------------|------------|
| <b>S-24</b> 8.0mm<br><b>S-24</b><br>1/4 Circle<br> | Black Monofilament | 18" ~ 45cm    | 7731G      |

## **Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|--------------------|---------------|------------|------------|------------|
| <b>KS 60mm</b><br>Straight<br> | Black Monofilament | 30" ~ 75cm    | 628H       | 627H       | 626H       |

## **Taper Point**

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>0</b>            |
|---|--------------------|---------------|---------------------|
| <b>CT 40mm</b><br>1/2 Circle<br> | Black Monofilament | 96" ~ 244cm   | L886T <sup>15</sup> |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>1</b> |
|--|--------------------|---------------|----------|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Black Monofilament | 30" ~ 75cm    | 830H     |

| <b>Needle</b>   | <b>Material</b>         | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2</b>       |
|---|-------------------------|---------------|----------|----------|----------------|
| <b>TP-1 65mm</b><br>1/2 Circle<br> | Black Monofilament      | 48" ~ 122cm   | L880G33  |          |                |
|   | Black Monofilament      | 60" ~ 150cm   |          | 824G     | 825G<br>825G33 |
|   | Black Monofilament Loop | 96" ~ 244cm   | L880G    |          |                |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>                      | <b>Material</b>    | <b>Length</b> | <b>10-0</b> | <b>8-0</b> | <b>9-0</b> |
|------------------------------------|--------------------|---------------|-------------|------------|------------|
| <b>BV130-3</b> 3.8mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2820G       | 2822G      | 2819G      |

| <b>Needle</b>                      | <b>Material</b>    | <b>Length</b> | <b>8-0</b> | <b>9-0</b> |
|------------------------------------|--------------------|---------------|------------|------------|
| <b>BV130-4</b> 5.0mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2815G      | 2813G      |

| <b>Needle</b>                      | <b>Material</b>    | <b>Length</b> | <b>10-0</b> | <b>8-0</b> | <b>9-0</b> |
|------------------------------------|--------------------|---------------|-------------|------------|------------|
| <b>BV130-5</b> 6.5mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2810G       | 2808G      | 2809G      |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>9-0</b> |
|--|--------------------|---------------|------------|
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b><br>3/8 Circle | Black Monofilament | 10" ~ 25cm    | 2800G      |

| <b>Needle</b>                     | <b>Material</b>    | <b>Length</b> | <b>11-0</b> |
|-----------------------------------|--------------------|---------------|-------------|
| <b>BV50-3</b> 4.1mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2881G       |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>                       | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|-------------------------------------|--------------------|---------------|-------------|
| <b>BVH100-3</b> 4.8mm<br>1/2 Circle | Black Monofilament | 5" ~ 13cm     | 2860G       |

| <b>Needle</b>                      | <b>Material</b>    | <b>Length</b> | <b>10-0</b>          |
|------------------------------------|--------------------|---------------|----------------------|
| <b>BV100-3</b> 3.9mm<br>3/8 Circle | Black Monofilament | 9-5" ~ 13cm   | 2889 <sup>0.14</sup> |

| <b>Needle</b>                      | <b>Material</b>    | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|------------------------------------|--------------------|---------------|-------------|------------|
| <b>BV100-4</b> 5.1mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2830G       | 2829G      |

| <b>Needle</b>                     | <b>Material</b>    | <b>Length</b> | <b>10-0</b>          |
|-----------------------------------|--------------------|---------------|----------------------|
| <b>BV75-3</b> 3.8mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2870G                |
|                                   | Black Monofilament | 9-5" ~ 13cm   | 2888 <sup>0.14</sup> |

| <b>Needle</b>                     | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|-----------------------------------|--------------------|---------------|-------------|
| <b>BV75-4</b> 5.1mm<br>3/8 Circle | Black Monofilament | 5" ~ 13cm     | 2850G       |

**TAPERCUT**

| <b>Needle</b>                     | <b>Material</b>       | <b>Length</b> | <b>9-0</b> |
|-----------------------------------|-----------------------|---------------|------------|
| <b>V100-3</b> 3.5mm<br>1/4 Circle | Black<br>Monofilament | 5" ~ 13cm     | 2894G      |

| <b>Needle</b>                     | <b>Material</b>       | <b>Length</b> | <b>8-0</b> | <b>9-0</b> |
|-----------------------------------|-----------------------|---------------|------------|------------|
| <b>V130-3</b> 3.8mm<br>1/4 Circle | Black<br>Monofilament | 5" ~ 13cm     | 2897G      | 2898G      |

| <b>Needle</b>                    | <b>Material</b>       | <b>Length</b> | <b>10-0</b> |
|----------------------------------|-----------------------|---------------|-------------|
| <b>V75-3</b> 3.5mm<br>1/4 Circle | Black<br>Monofilament | 5" ~ 13cm     | 2892G       |

**ULTIMA - Spatula**

| <b>Needle</b>                      | <b>Material</b>       | <b>Length</b> | <b>10-0</b>   |
|------------------------------------|-----------------------|---------------|---|
| <b>CS160-6</b> 5.5mm<br>1/2 Circle | Black<br>Monofilament | 4" ~ 10cm     | 9061G<br><br>9061G33 |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>10-0</b> | <b>11-0</b> | <b>9-0</b> |
|---------------|-----------------|---------------|-------------|-------------|------------|
|               |                 |               |             |             |            |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b>      | <b>11-0</b>                 | <b>9-0</b> |
|---|--------------------|---------------|------------------|-----------------------------|------------|
| <b>CS160-6</b> 5.5mm<br><b>CS160-6</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9000G            |                             | 9011G      |
|   | Black Monofilament | 6" ~ 15cm     | 9001G<br>9001G05 |                             |            |
|   | Black Monofilament | 8" ~ 20cm     |                  | 9016G<br>9016G05<br>9016G33 |            |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>CS160-8</b> 5.5mm<br><b>CS160-8</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9030G       |
|   | Black Monofilament | 6" ~ 15cm     | 9031G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b>      |
|---|--------------------|---------------|------------------|
| <b>CS175-6</b> 7.0mm<br><b>CS175-6</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9002G<br>9002G33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> | <b>8-0</b> | <b>9-0</b>       |
|---|--------------------|---------------|-------------|------------|------------------|
| <b>CS175-8</b> 7.0mm<br><b>CS175-8</b><br>1/2 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9032G       | 9038G      | 9036G<br>9036G33 |

| <b>Needle</b>  | <b>Material</b>    | <b>Length</b> | <b>10-0</b>  |
|--|--------------------|---------------|--|
| <b>CS44-6</b> 5.5mm<br>1/8 Circle<br> | Black Monofilament | 4" ~ 10cm     | 9065G<br> |

## Sutures : Non Absorbable / ETHILON Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b>  |
|---|--------------------|---------------|--|
| <b>CS140-6</b> 6.5mm<br>3/8 Circle<br> | Black Monofilament | 4" ~ 10cm     | 9060G<br><br>9060G05<br>9060G33 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b>      | <b>9-0</b>       |
|---|--------------------|---------------|------------------|------------------|
| <b>CS140-6</b> 6.5mm<br><b>CS140-6</b><br>3/8 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9003G<br>9003G05 | 9013G<br>9013G05 |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>CS140-8</b> 6.5mm<br><b>CS140-8</b><br>3/8 Circle<br> | Black Monofilament | 12" ~ 30cm    | 9033G       |

Nonsterile

<sup>1</sup> One dozen packets per box

<sup>14</sup> 8 packets (72 strands) per box

<sup>15</sup> This product has a loop length that is half of the length specified here.

<sup>40</sup> Lead shot and bolster

<sup>65</sup> Pliabilized

 **Multipass<sup>®</sup>** MULTIPASS Needles

 **HEMO-SEAL**

# MERSILENE Suture

## Blunt Point

| Needle   | Material    | Length     | 5MM                |
|--|-------------|------------|--------------------|
| <b>BP</b> 80mm<br><b>BP</b><br>1/2 Circle<br> | White Woven | 12" ~ 30cm | RS20 <sup>10</sup> |

| Needle   | Material    | Length     | none   | 5MM  |
|--|-------------|------------|--------|------|
| <b>BP-1</b> 65mm<br><b>BP-1</b><br>1/2 Circle<br> | White Woven | 12" ~ 30cm | RS2134 |      |
|  |             | 12" ~ 30cm |        | RS21 |

## MICROPOINT - Spatula

| Needle  | Material           | Length     | 11-0             |
|---|--------------------|------------|------------------|
| <b>TG160-4</b> 5.5mm<br><b>TG160-4</b><br>1/2 Circle<br> | Green Monofilament | 12" ~ 30cm | R750G<br>R750G33 |

## Sutures : Non Absorbable / MERSILENE Suture

| <b>Needle</b>   | <b>Material</b>       | <b>Length</b> | <b>11-0</b> |
|---|-----------------------|---------------|-------------|
| <b>TG160-4-3M</b> 4.3mm<br><b>TG160-4-3M</b><br>1/2 Circle<br> | Green<br>Monofilament | 8" ~ 20cm     | R749G       |

## Precision Point - Reverse Cutting

| <b>Needle</b>  | <b>Material</b>  | <b>Length</b> | <b>6-0</b> |
|--|------------------|---------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | White<br>Braided | 18" ~ 45cm    | R689G      |

| <b>Needle</b>   | <b>Material</b>  | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|------------------|---------------|------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br><br> | White<br>Braided | 18" ~ 45cm    | R691G      | R690G      |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>  | <b>Length</b> | <b>none</b> | <b>2-0</b> |
|---|------------------|---------------|-------------|------------|
| <b>FS</b> 26mm<br>3/8 Circle<br> |                  | 18" ~ 45cm    | R665H25     |            |
|   | White<br>Braided | 18" ~ 45cm    |             | R665H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Non Absorbable / MERSILENE Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | White Braided   | 18" ~ 45cm    | R647H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | White Braided   | 18" ~ 45cm    | R633H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>1</b> |
|--|-----------------|---------------|----------|
| <b>LR</b> 75mm<br><b>LR</b><br>3/8 Circle<br> | Green Braided   | 30" ~ 75cm    | R495T    |

## SABRELOC - Center Point Spatula

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b>     |
|---|-----------------|---------------|----------------|
| <b>RD-1</b> 8.0mm<br><b>RD-1</b><br>1/4 Circle<br> | Green Braided   | 18" ~ 45cm    | 740G           |
|   | White Braided   | 18" ~ 45cm    | 746G<br>746G33 |

## SABRELOC - Conventional Spatula

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Non Absorbable / MERSILENE Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>       |
|---|-----------------|---------------|------------------|
| <b>S-4</b> 8.0mm<br><b>S-4</b><br>1/4 Circle<br> | White Braided   | 18" ~ 45cm    | 1778G<br>1778G05 |

## **SABRELOC - Spatula**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>                  |
|---|-----------------|---------------|-----------------------------|
| <b>S-2</b> 8.0mm<br><b>S-2</b><br>1/2 Circle<br> | White Braided   | 18" ~ 45cm    | 1779G<br>1779G05<br>1779G33 |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br>1/4 Circle<br> | White Braided   | 18" ~ 45cm    | 1760G      | 1722G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b>       |
|---|-----------------|---------------|------------|------------------|
| <b>S-24</b> 8.0mm<br><b>S-24</b><br>1/4 Circle<br> | Green Braided   | 18" ~ 45cm    | 1762G      | 1761G<br>1761G05 |
|   | White Braided   | 18" ~ 45cm    |            | 1764G            |

## **Taper Point**

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Non Absorbable / MERSILENE Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> |
|---|-----------------|---------------|----------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Green Braided   | 8-18" ~ 45cm  | MR21T    |
|   | Green Braided   | 30" ~ 75cm    | R424H    |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>5MM</b>                   |
|--|-----------------|---------------|------------------------------|
| <b>CTX</b> 48mm<br><b>CTX</b><br>1/2 Circle<br> | White Woven     | 12" ~ 30cm    | RS22 <sup>10</sup><br>RS2233 |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5MM</b>         |
|---|-----------------|---------------|--------------------|
| <b>MO-4</b> 36mm<br>1/2 Circle<br> | White Woven     | 12" ~ 30cm    | RS23 <sup>10</sup> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Green Braided   | 30" ~ 75cm    | R833H      | R832H      |

## ULTIMA - Spatula

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>10-0</b> |
|---------------|-----------------|---------------|-------------|
|               |                 |               |             |

## Sutures : Non Absorbable / MERSILENE Suture

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>CS160-6</b> 5.5mm<br><b>CS160-6</b><br>1/2 Circle<br> | Green Monofilament | 12" ~ 30cm    | R756G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>CS140-6</b> 6.5mm<br><b>CS140-6</b><br>3/8 Circle<br> | Green Monofilament | 12" ~ 30cm    | R770G       |

| <b>Needle</b>   | <b>Material</b>    | <b>Length</b> | <b>10-0</b> |
|---|--------------------|---------------|-------------|
| <b>CSB-6</b> 5.0mm<br><b>CSB-6</b><br>Bi-Curve<br> | Green Monofilament | 8" ~ 20cm     | R721G       |

<sup>10</sup> 6 strips per box

 **Multipass<sup>®</sup>** MULTIPASS Needles

# NUROLON Suture

## Precision Cosmetic - Conventional Cutting PRIME

| Needle  | Material      | Length     | 5-0   | 6-0   |
|---|---------------|------------|-------|-------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br> | Black Braided | 18" ~ 45cm | 5665G | 5666G |

## Reverse Cutting

| Needle  | Material      | Length       | 2-0  |
|---|---------------|--------------|--|
| <b>CP</b> 40mm<br>1/2 Circle<br> | Black Braided | 8-18" ~ 45cm | C572D <sup>13</sup><br> |

| Needle  | Material      | Length       | 2-0  | 3-0  |
|---|---------------|--------------|--|--|
| <b>CP-2</b> 26mm<br>1/2 Circle<br> | Black Braided | 8-18" ~ 45cm | C575D <sup>13</sup><br> | C576D <sup>13</sup><br> |

| Needle  | Material      | Length       | 1   |
|---|---------------|--------------|---|
| <b>OS-8</b> 40mm<br>1/2 Circle<br> | Black Braided | 3-18" ~ 45cm | C595T <sup>5</sup><br> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b>   |
|--|-----------------|---------------|--|
| <b>X-1</b> 22mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C563D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>   |
|---|-----------------|---------------|--|
| <b>FS</b> 26mm<br>3/8 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C570D <sup>13</sup><br> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b>   |
|--|-----------------|---------------|--|
| <b>FSL</b> 30mm<br>3/8 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C566D <sup>13</sup><br> |

## Taper Point

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|--|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C521D <sup>13</sup><br> | C520D <sup>13</sup><br> | C522D <sup>13</sup><br> |
|   | Black Braided   | 30" ~ 75cm    | 5424H  | 5425H  |  |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|---------------|-----------------|---------------|----------|------------|
|               |                 |               |          |            |

## Sutures : Non Absorbable / NUROLON Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C527D <sup>13</sup><br> | C526D <sup>13</sup><br> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>1</b>   |
|--|-----------------|---------------|--|--|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C551D <sup>13</sup><br> | C550D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-6</b> 26mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C545D <sup>13</sup><br> | C546D <sup>13</sup><br> |


| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-7</b> 22mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C541D <sup>13</sup><br> | C542D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>   | <b>4-0</b>  |
|---|-----------------|---------------|--|---|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C553D <sup>13</sup><br> | C554D <sup>13</sup><br><br>N124T |

## Sutures : Non Absorbable / NUROLON Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b>   | <b>3-0</b>   | <b>4-0</b>   |
|---|-----------------|---------------|--|--|--|
| <b>SH 26mm</b><br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C512D <sup>13</sup><br> | C513D <sup>13</sup><br> | C514D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b>   | <b>3-0</b>   |
|---|-----------------|---------------|--|--|
| <b>SH-1 22mm</b><br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C502D <sup>13</sup><br> | C503D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>   |
|---|-----------------|---------------|--|
| <b>TF 13mm</b><br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C584D <sup>13</sup><br><br>N104T |

<sup>5</sup> Single strand delivery

<sup>13</sup> Safety organizer tray

 **MULTIPASS** NEEDLES

 CONTROL RELEASE

# PERMA-HAND Suture

## LIGAPAK Dispensing Reel

| Needle   | Material      | Length       | 2-0   | 3-0   | 4-0   |
|--|---------------|--------------|-------|-------|-------|
| <b>LIGAPAK Dispensing Reel</b><br>N/A<br> | Black Braided | 144" ~ 366cm | LA55G | LA54G | LA53G |

## MICROPOINT - Spatula

| Needle  | Material      | Length     | 6-0   | 7-0   |
|---|---------------|------------|-------|-------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Black Braided | 18" ~ 45cm | 1732G | 7733G |

## MICROPOINT-Reverse Cutting

| Needle  | Material      | Length     | 6-0  | 7-0  |
|---|---------------|------------|------|------|
| <b>G-7</b> 8.0mm<br><b>G-7</b><br>1/2 Circle<br> | Black Braided | 18" ~ 45cm | 765G | 766G |

| Needle   | Material      | Length     | 6-0            |
|--|---------------|------------|----------------|
| <b>G-1</b> 11mm<br>3/8 Circle<br> | Black Braided | 18" ~ 45cm | 786G<br>786G05 |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|
| <b>G-1</b> 11mm<br><b>G-1</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 780G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>G-3</b> 13mm<br><b>G-3</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 789G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>G-3</b> 13mm<br><b>G-3</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 783G       |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b>     | <b>7-0</b>     |
|---|-----------------|---------------|----------------|----------------|
| <b>G-6</b> 8.0mm<br><b>G-6</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 769G<br>769G33 | 768G<br>768G05 |

## Precision Cosmetic - Conventional Cutting PRIME

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |


## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 1925G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 1984G      |

## Precision Point - Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br> | Black Braided   | 18" ~ 45cm    | 1687G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 639G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 641G       | 640G       | 1639G      |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>PS 26mm</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 1588H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b> |
|---|-----------------|---------------|------------|
| <b>PS-1 24mm</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 1684G      |


| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|------------|------------|
| <b>PS-2 19mm</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 583H       | 1679H      | 1677G      | 1676G      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b>       |
|---|-----------------|---------------|------------|------------------|
| <b>PS-3 16mm</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 1678H      | 1675H<br>1675H33 |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b>     | <b>3-0</b> |
|--|-----------------|---------------|----------|----------------|------------|
| <b>PSL 30mm</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 580H     | 673H           |            |
|  | Black Braided   | 30" ~ 75cm    | 590H     | 486H<br>486H33 | 1690H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>PSLX 40mm</b><br>3/8 Circle<br> | Black Braided   | 30" ~ 75cm    | 591H       |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>CP-1 36mm</b><br>1/2 Circle<br> | Black Braided   | 30" ~ 75cm    | 443H       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>J-1 17mm</b><br>1/2 Circle<br> | Black Braided   | 18" ~ 45cm    | 734G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>M-1 13mm</b><br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | M154T      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>  |
|--|-----------------|---------------|------------|---|
| <b>X-1</b> 22mm<br>1/2 Circle<br> | Black Braided   | 18" ~ 45cm    | 737G       | 632G  |
|  | Black Braided   | 8-18" ~ 45cm  |            | C063D <sup>13</sup><br><br>C063D33 |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>6-0</b>     |
|--|-----------------|---------------|----------------|
| <b>C-2</b> 13mm<br><b>C-2</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 714G<br>714G33 |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>C-3</b> 13mm<br><b>C-3</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 736G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>C-3</b> 13mm<br><b>C-3</b><br>3/8 Circle<br> | Black Braided   | 12" ~ 30cm    | 735G       |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b>   |
|---|-----------------|---------------|--------------|
| <b>FS</b> 26mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 685G<br>685H |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>   | <b>4-0</b> |
|---|-----------------|---------------|--------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 684G<br>684H | 629G       |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 683G       | 682G       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b>     | <b>2-0</b> | <b>3-0</b> |
|--|-----------------|---------------|--------------|------------|------------|
| <b>FSL</b> 30mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 678G<br>678H | 677G       |            |
|  | Black Braided   | 30" ~ 75cm    | 680H         | 679H       | 675H       |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---|-----------------|---------------|------------|
| <b>FSLX</b> 40mm<br>3/8 Circle<br> | Black Braided   | 30" ~ 75cm    | 1689H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2</b> | <b>5</b> |
|--|-----------------|---------------|----------|----------|
| <b>LR</b> 75mm<br><b>LR</b><br>3/8 Circle<br> | Black Braided   | 30" ~ 75cm    | 482T     | 485T     |

## **SABRELOC - Spatula**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b>       |
|---|-----------------|---------------|------------------|
| <b>S-14</b> 8.0mm<br><b>S-14</b><br>1/4 Circle<br> | Black Braided   | 18" ~ 45cm    | 1780G<br>1780G05 |

## **Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b>        | <b>4-0</b> |
|---|-----------------|---------------|----------|------------|-------------------|------------|
| <b>KS</b> 60mm<br>Straight<br> | Black Braided   | 30" ~ 75cm    | 624H     | 623H       | 622H              | 621H       |
|   | Black Braided   | 3-30" ~ 75cm  |          |            | 3622 <sup>3</sup> |            |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---|-----------------|---------------|------------|
| <b>TS</b> 48mm<br>Straight<br> | Black Braided   | 30" ~ 75cm    | 619H       |

## **Straight Taper Point**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>none</b> | <b>3-0</b> |
|---|-----------------|---------------|-------------|------------|
| <b>ST-1</b> 45mm<br>Straight<br> |                 | 30" ~ 75cm    | K852H34     |            |
|   | Black Braided   | 30" ~ 75cm    |             | K852H      |

## **SUTUPAK Pre-Cut Sutures**

| <b>Needle</b>                         | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>2</b> | <b>3-0</b> | <b>3</b> | <b>4-0</b> | <b>4</b> | <b>5</b> |
|---------------------------------------|-----------------|---------------|----------|----------|------------|----------|------------|----------|------------|----------|----------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A | Black Braided   | 17-18" ~ 45cm | SA66G    |          | SA65H      |          | SA64H      |          | SA63H      |          |          |
|                                       | Black Braided   | 13-24" ~ 60cm | SA76G    | SA77G    | SA75H      |          | SA74H      |          | SA73H      |          |          |
|                                       | Black Braided   | 10-30" ~ 75cm | SA86G    | SA87G    | SA85H      |          | SA84H      |          | SA83H      |          |          |
|                                       | Black Braided   | 2-60" ~ 150cm | SA6H     | SA7H     | SA5H       | SA8H     |            | SA9G     |            | SA10G    | SA11G    |

## **SUTUPAK Pre-Cut Sutures in Labyrinth Package**

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|----------|----------|------------|------------|------------|------------|
| <b>SUTUPAK Pre-Cut Sutures in Labyrinth Package</b><br>N/A | Black Braided   | 6-18" ~ 45cm  | A186H    | A187H    |            |            |            |            |
|  | Black Braided   | 12-18" ~ 45cm |          |          | A185H      | A184H      | A183H      | A182H      |
|  | Black Braided   | 6-30" ~ 75cm  | A306H    | A307H    |            |            |            |            |
|  | Black Braided   | 12-30" ~ 75cm |          |          | A305H      | A304H      | A303H      | A302H      |

## **Taper Point**

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>  | <b>2-0</b>   |
|---|-----------------|---------------|---|--|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C021D <sup>13</sup><br><br>C021D33 | C022D <sup>13</sup><br> |
|   | Black Braided   | 30" ~ 75cm    | 424H  | 423H   |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C027D <sup>13</sup><br> | C026D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b>     |
|---|-----------------|---------------|----------------|
| <b>CT-3</b> 22mm<br>1/2 Circle<br> | Black Braided   | 20" ~ 50cm    | 413H<br>413H33 |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2</b> | <b>2-0</b> |
|---|-----------------|---------------|----------|----------|----------|------------|
| <b>MH</b> 36mm<br>1/2 Circle<br> | Black Braided   | 8-24" ~ 60cm  |          |          |          | M72T       |
|   | Black Braided   | 30" ~ 75cm    | K844H    | K845H    | K846H    | K843H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-6</b> 26mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C045D <sup>13</sup><br> | C046D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b>   | <b>2-0</b>   |
|---|-----------------|---------------|--|--|
| <b>MO-7</b> 22mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C041D <sup>13</sup><br> | C042D <sup>13</sup><br> |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b>   | <b>4-0</b>   | <b>5-0</b> |
|---|-----------------|---------------|------------|--|--|------------|
| <b>RB-1</b> 17mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  |            | C053D <sup>13</sup><br> | C054D <sup>13</sup><br> |            |
|   | Black Braided   | 30" ~ 75cm    | K873H      | K872H  | K871H  | K870H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>RB-1</b> 17mm<br><b>RB-1</b><br>1/2 Circle<br> | Black Braided   | 24" ~ 60cm    | K572H      | K571H      | K570H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b>   | <b>3-0</b>   | <b>4-0</b>   |
|---|-----------------|---------------|----------|----------|--|--|--|
| <b>SH</b> 26mm<br>1/2 Circle<br> | Black Braided   | 5-18" ~ 45cm  |          |          | C0125 <sup>5,6</sup><br>         | C0135 <sup>5,6</sup><br>         | C0145 <sup>5,6</sup><br>         |
|   | Black Braided   | 8-18" ~ 45cm  |          |          | C012D <sup>13</sup><br><br>M12T | C013D <sup>13</sup><br><br>M13T | C014D <sup>13</sup><br><br>M14T |
|   | Black Braided   | 30" ~ 75cm    | K834H    | K835H    | K833H  | K832H  | K831H  |
|   | Black Braided   | 4-30" ~ 75cm  | M834G    |          |  |  |  |
|   | Black Braided   | 8-30" ~ 75cm  |          |          | C016D<br>                       | C017D<br>                       | C018D<br>                       |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> |
|--|-----------------|---------------|----------|------------|
| <b>SH</b> 26mm<br><b>SH</b><br>1/2 Circle<br> | Black Braided   | 24" ~ 60cm    | K534H    | K533H      |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>3-0</b>  | <b>4-0</b>   |
|---|-----------------|---------------|---|--|
| <b>SH-1</b> 22mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C003D <sup>13</sup><br><br>M3T | C004D <sup>13</sup><br> |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>4-0</b>  | <b>5-0</b> |
|---|-----------------|---------------|---|------------|
| <b>TF</b> 13mm<br>1/2 Circle<br> | Black Braided   | 8-18" ~ 45cm  | C084D <sup>13</sup><br><br>M104T |            |
|   | Black Braided   | 18" ~ 45cm    | N267H   | N266H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>TF-4</b> 8.0mm<br>1/2 Circle<br> | Black Braided   | 18" ~ 45cm    | N272H      | N271H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-----------------|---------------|------------|------------|------------|------------|
| <b>BB</b> 17mm<br>3/8 Circle<br> | Black Braided   | 30" ~ 75cm    | K883H      | K882H      | K881H      | K880H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-----------------|---------------|------------|------------|
| <b>BB</b> 17mm<br><b>BB</b><br>3/8 Circle<br> | Black Braided   | 24" ~ 60cm    | K581H      | K580H      |

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|
| <b>BV-1</b> 9.3mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | K802H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> | <b>7-0</b>                                 |
|---|-----------------|---------------|------------|--|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | K801H      | K809H<br>Q809H <sup>b</sup>                |
|   | Black Braided   | 10-18" ~ 45cm |            | 10Q07 <sup>18</sup><br>10S07 <sup>18</sup> |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|-----------------|---------------|------------|------------|------------|
| <b>C-1</b> 13mm<br>3/8 Circle<br> | Black Braided   | 12" ~ 30cm    | 7734G      |            |            |
|  | Black Braided   | 18" ~ 45cm    |            |            | 707G       |
|  | Black Braided   | 30" ~ 75cm    | K891H      | K890H      | K889H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>6-0</b>                 |
|--|-----------------|---------------|------------|----------------------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | 709H       | 706G<br>Q706H <sup>b</sup> |
|  | Black Braided   | 10-18" ~ 45cm |            | 10Q06 <sup>18</sup>        |

## TAPERCUT

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / PERMA-HAND Suture

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>V-5</b> 17mm<br>1/2 Circle<br> | Black Braided   | 30" ~ 75cm    | K952H      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|--|-----------------|---------------|------------|
| <b>V-7</b> 26mm<br>1/2 Circle<br> | Black Braided   | 30" ~ 75cm    | K964H      |

| <b>Needle</b>   | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---|-----------------|---------------|------------|
| <b>CC-1</b> 13mm<br>3/8 Circle<br> | Black Braided   | 12" ~ 30cm    | 7729G      |

| <b>Needle</b>  | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|--|-----------------|---------------|------------|
| <b>V-4</b> 17mm<br>3/8 Circle<br> | Black Braided   | 18" ~ 45cm    | K945G      |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>5</sup> Single strand delivery

<sup>6</sup> Two dozen packets per box

<sup>8</sup> Siliconized silk

<sup>13</sup> Safety organizer tray

 **MULTIPASS** Needles

 CONTROL RELEASE

# PROLENE Suture

## ETHIGUARD - Safety needle

| Needle   | Material          | Length     | none    | 0     | 1     |
|--|-------------------|------------|---------|-------|-------|
| <b>CTB-1</b> 36mm<br>1/2 Circle<br> |                   | 30" ~ 75cm | PB20H34 |       |       |
|  | Blue Monofilament | 30" ~ 75cm |         | PB21H | PB20H |

| Needle  | Material          | Length     | 0                | 1                |
|---|-------------------|------------|------------------|------------------|
| <b>CTX-B</b> 48mm<br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm | PB31H<br>PB31H33 | PB30H<br>PB30H30 |

## EVERPOINT Cardiovascular Needles

| Needle | Material | Length | 7-0 | 8-0 |
|--------|----------|--------|-----|-----|
|        |          |        |     |     |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b>           | <b>8-0</b>           |
|--|-------------------|---------------|----------------------|----------------------|
| <b>BV175-6</b> 8.0mm<br><b>BV175-6</b>   | Blue Monofilament | 4-18" ~ 45cm  | EPM8734 <sup>1</sup> |                      |
| 3/8 Circle<br><br>EVERPOINT™ Cardiovascular Needles | Blue Monofilament | 2-18" ~ 45cm  | EPM8736 <sup>1</sup> | EPM8742 <sup>1</sup> |
|  | Blue Monofilament | 18" ~ 45cm    | EP8734H <sup>3</sup> | EP8740H <sup>3</sup> |
|  | Blue Monofilament | 4-24" ~ 60cm  | EPM8735              | EPM8741 <sup>1</sup> |
|  | Blue Monofilament | 2-24" ~ 60cm  | EPM8737 <sup>1</sup> | EPM8743 <sup>1</sup> |
|  | Blue Monofilament | 24" ~ 60cm    | EP8735H <sup>3</sup> | EP8741H <sup>3</sup> |
|  | Blue Monofilament | 2-30" ~ 75cm  | EPM8733              |                      |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8731 <sup>1</sup> |                      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b>             |
|--|-------------------|---------------|------------------------|
| <b>CC175-8</b> 9.3mm<br><b>CC175-8</b>   | Blue Monofilament | 4-24" ~ 60cm  | EPM8704SL <sup>1</sup> |
| 3/8 Circle<br><br>EVERPOINT™ Cardiovascular Needles | Blue Monofilament | 4-30" ~ 75cm  | EPM8812SL              |
|  | Blue Monofilament | 4-24" ~ 60cm  | EPM8704SL              |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8812SL              |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>7-0</b>                                   | <b>8-0</b>           |
|---|-------------------|---------------|--|----------------------|
| <b>BV175-8</b> 9.3mm<br><b>BV175-8</b>  | Blue Monofilament | 2-18" ~ 45cm  | EPM8746 <sup>1</sup><br>EPM8746 <sup>1</sup> | EPM8752 <sup>1</sup> |
| 3/8 Circle<br><br><br>EVERPOINT™<br>Cardiovascular Needles | Blue Monofilament | 18" ~ 45cm    | EP8746H <sup>3</sup>                         | EP8752H <sup>3</sup> |
|   | Blue Monofilament | 4-24" ~ 60cm  | EPM8745                                      |                      |
|   | Blue Monofilament | 2-24" ~ 60cm  | EPM8747 <sup>1</sup>                         | EPM8753 <sup>1</sup> |
|   | Blue Monofilament | 24" ~ 60cm    | EP8747H <sup>3</sup>                         | EP8753H <sup>3</sup> |
|   | Blue Monofilament | 2-30" ~ 75cm  | EPM8755                                      | EPM8754 <sup>1</sup> |
|   | Blue Monofilament | 30" ~ 75cm    | EPM8751<br>EP8755H <sup>3</sup>              | EP8754H <sup>3</sup> |
|   | Blue Monofilament | 4-30" ~ 75cm  | EPM8751 <sup>1</sup>                         |                      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>8-0</b>           |
|---|-------------------|---------------|----------------------|
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b>  | Blue Monofilament | 4-18" ~ 45cm  | EPM8730 <sup>1</sup> |
| 3/8 Circle<br><br><br>EVERPOINT™<br>Cardiovascular Needles | Blue Monofilament | 2-18" ~ 45cm  | EP8738 <sup>1</sup>  |
|   | Blue Monofilament | 18" ~ 45cm    | EP8730H <sup>1</sup> |
|   | Blue Monofilament | 4-24" ~ 60cm  | EP8732               |
|   | Blue Monofilament | 2-24" ~ 60cm  | EP8739 <sup>1</sup>  |
|   | Blue Monofilament | 24" ~ 60cm    | EP8732H <sup>1</sup> |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>6-0</b> | <b>7-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           | <b>7-0</b>           |
|--|-------------------|---------------|----------------------|----------------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br><br>EVERPOINT™<br>Cardiovascular Needles | Blue Monofilament | 18" ~ 45cm    | EP8718H <sup>1</sup> | EP8800H <sup>1</sup> |
|  | Blue Monofilament | 4-18" ~ 45cm  | EPM8718 <sup>1</sup> | EPM8800 <sup>1</sup> |
|  | Blue Monofilament | 4-24" ~ 60cm  | EPM8726              | EPM8727 <sup>1</sup> |
|  | Blue Monofilament | 24" ~ 60cm    | EP8726H <sup>1</sup> |                      |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8706              |                      |
|  | Blue Monofilament | 30" ~ 75cm    | EP8706H <sup>1</sup> | EP8708H <sup>1</sup> |
|  | Blue Monofilament | 2-30" ~ 75cm  | EPM8206 <sup>1</sup> |                      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           |
|---|-------------------|---------------|----------------------|
| <b>C1</b> 13mm<br>3/8 Circle<br><br>EVERPOINT™<br>Cardiovascular Needles | Blue Monofilament | 30" ~ 75cm    | EP8889H <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           |
|--|-------------------|---------------|----------------------|
| <b>CC-1</b> 13mm<br><b>CC-1</b><br>3/8 Circle<br><br>EVERPOINT™<br>Cardiovascular Needles | Blue Monofilament | 4-18" ~ 45cm  | EPM8719 <sup>1</sup> |
|  | Blue Monofilament | 18" ~ 45cm    | EP8719H              |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8707              |
|  | Blue Monofilament | 30" ~ 75cm    | EP8707H              |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b> |
|--|-------------------|---------------|------------|
| <b>RB-2</b> 13mm<br><b>RB-2</b><br>1/2 Circle<br> | Blue Monofilament | 18" ~ 45cm    | EP8714H    |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8711    |
|  | Blue Monofilament | 30" ~ 75cm    | EP8711H    |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           |
|--|-------------------|---------------|----------------------|
| <b>BV</b> 11mm<br><b>BV</b><br>3/8 Circle<br> | Blue Monofilament | 4-24" ~ 60cm  | EPM8610 <sup>1</sup> |
|  | Blue Monofilament | 24" ~ 60cm    | EP8610H              |
|  | Blue Monofilament | 30" ~ 75cm    | EP8776H              |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8776 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           |
|---|-------------------|---------------|----------------------|
| <b>CC</b> 9.3mm<br><b>CC</b><br>3/8 Circle<br> | Blue Monofilament | 24" ~ 60cm    | EP8807H              |
|   | Blue Monofilament | 4-24" ~ 60cm  | EPM8807 <sup>1</sup> |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>BV175-7</b> 8.0mm<br><b>BV175-7</b>   | Blue Monofilament | 2-18" ~ 45cm  | EPM8765    |            |
| 3/8 Circle<br><br><b>Multipass® NEEDLES</b> | Blue Monofilament | 4-24" ~ 60cm  | EPM8775    |            |
|  | Blue Monofilament | 2-24" ~ 60cm  | EPM8766    | EPM8767    |
|  | Blue Monofilament | 24" ~ 60cm    | EP8766H    |            |
|  | Blue Monofilament | 4-30" ~ 75cm  | EPM8785    |            |
|  | Blue Monofilament | 2-30" ~ 75cm  | EPM8768    |            |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b>           |
|---|-------------------|---------------|----------------------|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b>  | Blue Monofilament | 4-18" ~ 45cm  | EPM8806 <sup>1</sup> |
| 3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | EP8806H              |
|   | Blue Monofilament | 4-24" ~ 60cm  | EPM8805              |
|   | Blue Monofilament | 2-24" ~ 60cm  | EPM8205              |
|   | Blue Monofilament | 24" ~ 60cm    | EP8805H              |
|   | Blue Monofilament | 4-30" ~ 75cm  | EPM8709              |
|   | Blue Monofilament | 30" ~ 75cm    | EP8709H              |

## **MICROPOINT - Spatula**

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---------------|-----------------|---------------|-------------|------------|
|               |                 |               |             |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---|-------------------|---------------|-------------|------------|
| <b>TG140-8</b> 6.5mm<br><b>TG140-8</b><br>3/8 Circle<br> | Blue Monofilament | 12" ~ 30cm    | 1771G       |            |
|   | Blue Monofilament | 6" ~ 15cm     | 1757G       | 1754G      |

## Precision Cosmetic - Conventional Cutting PRIME

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|-------------------|---------------|------------|------------|------------|
| <b>PC-1</b> 13mm<br>3/8 Circle<br><br> | Blue Monofilament | 18" ~ 45cm    | 8619G      | 8618G      | 8617G      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b>       |
|--|-------------------|---------------|-------------|------------|------------|------------------|
| <b>PC-3</b> 16mm<br>3/8 Circle<br><br> |                   | 18" ~ 45cm    | 8634G34     |            |            |                  |
|  | Blue Monofilament | 18" ~ 45cm    |             | 8634G      | 8635G      | 8636G<br>8636G33 |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b>       | <b>5-0</b>       |
|--|-------------------|---------------|------------|------------------|------------------|
| <b>PC-5</b> 19mm<br>3/8 Circle<br><br> | Blue Monofilament | 18" ~ 45cm    | 8632G      | 8631G<br>8631G33 | 8630G<br>8630G33 |


## Precision Point - Reverse Cutting

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>     | <b>Length</b> | <b>4-0</b> |
|---|---------------------|---------------|------------|
| <b>PS-4</b> 16mm<br>1/2 Circle<br> | Undyed Monofilament | 18" ~ 45cm    | 8603G      |

| <b>Needle</b>                  | <b>Material</b>   | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--------------------------------|-------------------|---------------|------------|------------|
| <b>PS-5</b> 13mm<br>1/2 Circle | Blue Monofilament | 18" ~ 45cm    | 8656G      | 8655G      |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>none</b> | <b>6-0</b> | <b>7-0</b> |
|--|---------------------|---------------|-------------|------------|------------|
| <b>P-1</b> 11mm<br>3/8 Circle<br> | Blue Monofilament   | 18" ~ 45cm    | 8697G33     | 8697G      | 8696G      |
|  | Clear Monofilament  | 18" ~ 45cm    |             | 8606G      |            |
|  | Undyed Monofilament | 18" ~ 45cm    |             | 8606G33    |            |

| <b>Needle</b>  | <b>Material</b>     | <b>Length</b> | <b>none</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> |
|--|---------------------|---------------|-------------|------------|------------|------------|
| <b>P-3</b> 13mm<br>3/8 Circle<br> | Blue                | 18" ~ 45cm    | 8698G33     |            |            |            |
|  | Undyed Monofilament | 18" ~ 45cm    |             | 8604G      | 8605G      |            |
|  | Blue Monofilament   | 18" ~ 45cm    |             | 8699G      | 8698G      | 8695G      |

| <b>Needle</b>                  | <b>Material</b>   | <b>Length</b> | <b>7-0</b> |
|--------------------------------|-------------------|---------------|------------|
| <b>P-6</b> 8.0mm<br>3/8 Circle | Blue Monofilament | 18" ~ 45cm    | 8648G      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>3-0</b> |
|---|-------------------|---------------|-------------|------------|
| <b>PS-1</b> 24mm<br>3/8 Circle<br> |                   | 18" ~ 45cm    | 8663G33     |            |
|   | Blue Monofilament | 18" ~ 45cm    |             | 8663G      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>3-0</b> | <b>4-0</b>     | <b>5-0</b> |
|---|-------------------|---------------|-------------|------------|----------------|------------|
| <b>PS-2</b> 19mm<br>3/8 Circle<br> |                   | 18" ~ 45cm    | 8682H33     |            |                |            |
|   | Blue Monofilament | 18" ~ 45cm    |             | 8687H      | 8682G<br>8682H | 8686G      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---|-------------------|---------------|------------|------------|
| <b>PS-3</b> 16mm<br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 8681G      | 8680G      |

## Reverse Cutting

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>2-0</b> |
|---|-------------------|---------------|------------|
| <b>FS</b> 26mm<br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 8685H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|------------|------------|
|               |                 |               |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>3-0</b>       | <b>4-0</b> |
|---|-------------------|---------------|------------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 8684G<br>8684G33 | 8629H      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-------------------|---------------|-------------|------------|------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> |                   | 18" ~ 45cm    | 8683G34     |            |            |            |
|   | Blue Monofilament | 18" ~ 45cm    |             | 8665G      | 8683G      | 8661G      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> |  |  |
|--|-------------------|---------------|------------|--|--|
| <b>FSL</b> 30mm<br>3/8 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8675H      |  |  |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|---|-------------------|---------------|----------|------------|------------|
| <b>FSLX</b> 40mm<br>3/8 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8690H    | 8689H      | 8649H      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>2</b> |
|---|-------------------|---------------|----------|
| <b>LR</b> 75mm<br>3/8 Circle<br> | Blue Monofilament | 3-20" ~ 50cm  | 3846T    |

**SABRELOC - Spatula**

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>5-0</b> |
|---|-------------------|---------------|-------------|------------|
| <b>SM-1</b> 5.5mm<br><b>SM-1</b><br>1/8 Circle<br> |                   | 12" ~ 30cm    | 7740G34     |            |
|   | Blue Monofilament | 12" ~ 30cm    |             | 7740G      |

**SABRELOC - Straight Transchamber**

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>10-0</b> |
|---|-------------------|---------------|-------------|-------------|
| <b>STC-6</b> 16mm<br><b>STC-6</b> 16mm<br>Straight<br> |                   | 8" ~ 20cm     | 1713G34     |             |
|   | Blue Monofilament | 8" ~ 20cm     |             | 1713G       |

**Straight Cutting Needles**

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|-------------------|---------------|------------|------------|------------|
| <b>KS</b> 60mm<br>Straight<br> | Blue Monofilament | 30" ~ 75cm    | 8623H      | 8622H      | 8621H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>2-0</b> |
|--|-------------------|---------------|------------|
| <b>TS</b> 48mm<br><b>TS</b><br>Straight<br> | Blue Monofilament | 14" ~ 35cm    | 8571G      |

## **Straight TAPERCUT**

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>none</b> | <b>4-0</b> |
|--|-------------------|---------------|-------------|------------|
| <b>V-47</b> 40mm<br><b>V-47</b><br>Straight<br> |                   | 14" ~ 35cm    | 8592G34     |            |
|  | Blue Monofilament | 14" ~ 35cm    |             | 8592G      |

## **Taper Point**

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>0</b> | <b>1</b> | <b>2-0</b> |
|--|-------------------|---------------|----------|----------|------------|
| <b>CT</b> 40mm<br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8434H    | 8435H    | 8433H      |


| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b>   | <b>1</b> | <b>2-0</b>   | <b>3-0</b>       |
|---|-------------------|---------------|--|----------|--|------------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Blue Monofilament | 8-18" ~ 45cm  | C821G<br> |          | C822G<br> |                  |
|   | Blue Monofilament | 30" ~ 75cm    | 8424H  | 8425H    | 8423H  | 8422H<br>8422H33 |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b>   | <b>2-0</b>   | <b>3-0</b> |
|---|-------------------|---------------|--|--|------------|
| <b>CT-2</b> 26mm<br>1/2 Circle<br> | Blue Monofilament | 8-18" ~ 45cm  | C827G<br> | C826G<br> |            |
|   | Blue Monofilament | 30" ~ 75cm    | 8412H  | 8411H  | 8410H      |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>0</b>   | <b>1</b> |
|--|-------------------|---------------|--|----------|
| <b>CTX 48mm</b><br>1/2 Circle<br> | Blue Monofilament | 8-18" ~ 45cm  | C811G<br> |          |
|  | Blue Monofilament | 30" ~ 75cm    | 8454H  | 8455H    |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>2-0</b> | <b>3-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>MH 36mm</b><br><b>MH</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 8843H      | 8842H      |
|  | Blue Monofilament | 48" ~ 122cm   | 8853H      |            |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b>   | <b>2-0</b> |
|---|-------------------|---------------|--|------------|
| <b>MO-6 26mm</b><br>1/2 Circle<br> | Blue Monofilament | 8-18" ~ 45cm  | C845G<br> |            |
|   | Blue Monofilament | 30" ~ 75cm    | 8418H  | 8417H      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b>   |
|---|-------------------|---------------|--|
| <b>MO-7 22mm</b><br>1/2 Circle<br> | Blue Monofilament | 8-18" ~ 45cm  | C841G<br> |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---|-------------------|---------------|------------|------------|------------|
| <b>RB-1 17mm</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8872H      | 8871H      | 8870H      |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>2-0</b> | <b>3-0</b>  | <b>4-0</b>  | <b>5-0</b>  | <b>none</b> |
|--|-------------------|---------------|------------|---|---|---|-------------|
| <b>RB-1</b> 17mm<br><b>RB-1</b><br>1/2 Circle<br> | Blue Monofilament | 18" ~ 45cm    |            |   | 8757H   | 8756H   |             |
|  | Blue Monofilament | 24" ~ 60cm    |            |   |   | 8555H   |             |
|  | Blue Monofilament | 36" ~ 90cm    | 8559H      | 8558H<br>9556H<br> | 8357H<br><br>8557H | 8356H<br><br>8556H                           |             |
|  | Blue Monofilament | 4-36" ~ 90cm  |            |   | M8557 <sup>1</sup>  | M8356 <sup>1</sup><br><br>M8556 <sup>1</sup> |             |
|  | Blue              | 24" ~ 60cm    |            |   |   |   | 8555H33     |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>5-0</b>                    | <b>6-0</b>         |
|--|-------------------|---------------|-------------------------------|--------------------|
| <b>RB-2</b> 13mm<br><b>RB-2</b><br>1/2 Circle<br> | Blue Monofilament | 18" ~ 45cm    | 8713H                         | 8714H              |
|  | Blue Monofilament | 24" ~ 60cm    | 8554H                         |                    |
|  | Blue Monofilament | 30" ~ 75cm    | 8710H                         | 8711H              |
|  | Blue Monofilament | 4-30" ~ 75cm  | M8710 <sup>1</sup><br>M871033 | M8711 <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b> |  |
|--|-------------------|---------------|------------|--|
| <b>RB-3</b> 11mm<br><b>RB-3</b><br>1/2 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 8206H      |  |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---------------|-----------------|---------------|----------|------------|------------|------------|
|               |                 |               |          |            |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|---|-------------------|---------------|----------|------------|------------|------------|
| <b>SH 26mm</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8834H    | 8833H      | 8832H      | 8831H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>2-0</b> | <b>3-0</b>                     | <b>4-0</b> |
|--|-------------------|---------------|------------|--------------------------------|------------|
| <b>SH 26mm</b><br><b>SH</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 8523H      | 8522H<br>P8522 <sup>21,3</sup> | 8521H      |
|  | Blue Monofilament | 48" ~ 122cm   | 8533H      | 8534H                          |            |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>SH-1 22mm</b><br><b>SH-1</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8762H      | 8526H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>4-0</b> | <b>5-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>TF 13mm</b><br><b>TF</b><br>1/2 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 8204H      | 8205H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>1</b> | <b>2</b> |
|---------------|-----------------|---------------|----------|----------|
|               |                 |               |          |          |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>1</b> | <b>2</b> |
|---|-------------------|---------------|----------|----------|
| <b>TP-1</b> 65mm<br>1/2 Circle<br> | Blue Monofilament | 60" ~ 150cm   | 8824G    | 8825G    |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>0</b> | <b>1</b> |
|--|-------------------|---------------|----------|----------|
| <b>XLH</b> 70mm<br>1/2 Circle<br> | Blue Monofilament | 3-24" ~ 60cm  | 3844G    |          |
|  | Blue Monofilament | 60" ~ 150cm   |          | 8845G    |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>1</b> |
|--|-------------------|---------------|----------|
| <b>XLH</b> 70mm<br><b>XLH</b><br>1/2 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8190G    |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>4-0</b> |
|---|-------------------|---------------|------------|
| <b>BB</b> 17mm<br>3/8 Circle<br> | Blue Monofilament | 30" ~ 75cm    | 8881H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>3-0</b> | <b>4-0</b> | <b>5-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|
|               |                 |               |            |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b>                  | <b>5-0</b> |
|--|-------------------|---------------|------------|-----------------------------|------------|
| <b>BB</b> 17mm<br><b>BB</b><br>3/8 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 8862H      | 8861H                       | 8860H      |
|  | Blue Monofilament | 36" ~ 90cm    |            | 8381H<br><b>VB</b><br>8581H | 8580H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>6-0</b>         | <b>7-0</b>         |
|--|-------------------|---------------|--------------------|--------------------|
| <b>BV</b> 11mm<br><b>BV</b><br>3/8 Circle<br> | Blue Monofilament | 24" ~ 60cm    | 8610H              |                    |
|  | Blue Monofilament | 4-24" ~ 60cm  | M8610 <sup>1</sup> | M8600 <sup>1</sup> |
|  | Blue Monofilament | 30" ~ 75cm    | 8776H              |                    |
|  | Blue Monofilament | 4-30" ~ 75cm  | M8776 <sup>1</sup> |                    |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> | <b>7-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|
|               |                 |               |            |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>5-0</b>   | <b>6-0</b>  | <b>7-0</b>   |
|---|-------------------|---------------|--|---|--|
| <b>BV-1</b> 9.3mm<br><b>BV-1</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    |  | 8806H   | 8301H<br><br>8701H                                    |
|   | Blue Monofilament | 4-18" ~ 45cm  |  | M8806 <sup>1</sup>  | M8301 <sup>1</sup><br><br>M8701 <sup>1</sup>          |
|   | Blue Monofilament | 24" ~ 60cm    | 9702H<br> | 8305H <sup>5,9</sup><br><br>8805H                                  | 8304H<br><br>8702H                                    |
|   | Blue Monofilament | 2-24" ~ 60cm  |  | M8205   |  |
|   | Blue Monofilament | 4-24" ~ 60cm  |  | M8305 <sup>1</sup><br><br>M8805 <sup>1</sup><br>M8805 <sup>1</sup> | M8304 <sup>1</sup><br><br>M8702 <sup>1</sup>          |
|   | Blue Monofilament | 30" ~ 75cm    |  | 8709H   | 8703H  |
|   | Blue Monofilament | 4-30" ~ 75cm  |  | M8709 <sup>1</sup>  | M8303 <sup>1</sup><br><br>M8703 <sup>1</sup><br>M8703 |
|   | Blue Monofilament | 2-30" ~ 75cm  |  |   | M8203  |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>8-0</b> |
|--|-------------------|---------------|------------|
| <b>BV130-5</b> 6.5mm<br>3/8 Circle<br><br> | Blue Monofilament | 5" ~ 13cm     | 2775G      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|
|               |                 |               |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>8-0</b>                               |
|---|-------------------|---------------|--|
| <b>BV130-5</b> 6.5mm<br><b>BV130-5</b>  | Blue Monofilament | 18" ~ 45cm    | 8730H<br>8730H                           |
| 3/8 Circle<br> | Blue Monofilament | 4-18" ~ 45cm  | M8730 <sup>1</sup><br>M8730 <sup>1</sup> |
|   | Blue Monofilament | 2-18" ~ 45cm  | M8738 <sup>1</sup><br>M8738 <sup>1</sup> |
|   | Blue Monofilament | 24" ~ 60cm    | 8732H<br>8732H                           |
|   | Blue Monofilament | 4-24" ~ 60cm  | M8732 <sup>1</sup><br>M8732 <sup>1</sup> |
|   | Blue Monofilament | 2-24" ~ 60cm  | M8739 <sup>1</sup><br>M8739 <sup>1</sup> |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>7-0</b>                    | <b>8-0</b>                    |
|---|-------------------|---------------|-------------------------------|-------------------------------|
| <b>BV175-6</b> 8.0mm<br><b>BV175-6</b>  | Blue Monofilament | 18" ~ 45cm    | 8734H<br>8734H01              | 8740H                         |
| 3/8 Circle<br> | Blue Monofilament | 4-18" ~ 45cm  | M8734 <sup>1</sup><br>M873401 |                               |
|   | Blue Monofilament | 2-18" ~ 45cm  | M8736 <sup>1</sup><br>M873601 | M8742 <sup>1</sup><br>M874201 |
|   | Blue Monofilament | 24" ~ 60cm    | 8735H                         | 8741H <sup>1</sup>            |
|   | Blue Monofilament | 4-24" ~ 60cm  | M8735 <sup>1</sup><br>M873533 | M8741 <sup>1</sup>            |
|   | Blue Monofilament | 2-24" ~ 60cm  | M8737 <sup>1</sup><br>M873701 | M8743<br>M874301              |
|   | Blue Monofilament | 4-30" ~ 75cm  | M8731<br>M873101              |                               |
|   | Blue Monofilament | 2-30" ~ 75cm  | M8733 <sup>1</sup><br>M873301 |                               |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>7-0</b> | <b>8-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>7-0</b>                  | <b>8-0</b>                  |
|---|-------------------|---------------|-----------------------------|-----------------------------|
| <b>BV175-8</b> 9.3mm<br><b>BV175-8</b>  | Blue Monofilament | 18" ~ 45cm    | 8746H<br>8746H              | 8752H<br>8752H              |
| 3/8 Circle<br> | Blue Monofilament | 2-18" ~ 45cm  | M8746 <sup>1</sup><br>M8746 | M8752 <sup>1</sup><br>M8752 |
|   | Blue Monofilament | 24" ~ 60cm    | 8747H<br>8747H              | 8753H<br>8753H              |
|   | Blue Monofilament | 4-24" ~ 60cm  | M8745 <sup>1</sup><br>M8745 |                             |
|   | Blue Monofilament | 2-24" ~ 60cm  | M8747 <sup>1</sup><br>M8747 | M8753 <sup>1</sup><br>M8753 |
|   | Blue Monofilament | 30" ~ 75cm    | 8755H<br>8755H              | 8754H<br>8754H              |
|   | Blue Monofilament | 2-30" ~ 75cm  | M8755 <sup>1</sup><br>M8755 | M8754 <sup>1</sup><br>M8754 |

| <b>Needle</b>                 | <b>Material</b>   | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|-------------------------------|-------------------|---------------|------------|------------|
| <b>C-1</b> 13mm<br>3/8 Circle | Blue Monofilament | 30" ~ 75cm    | 8890H      | 8889H      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>4-0</b> | <b>5-0</b> | <b>6-0</b> | <b>7-0</b> |
|---------------|-----------------|---------------|------------|------------|------------|------------|
|               |                 |               |            |            |            |            |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>4-0</b>         | <b>5-0</b>                  | <b>6-0</b>  | <b>7-0</b> |
|--|-------------------|---------------|--------------------|-----------------------------|---|------------|
| <b>C-1</b> 13mm<br><b>C-1</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    |                    | 8717H                       | 8318H<br><b>VB</b><br>8718H                       | 8800H      |
|  | Blue Monofilament | 4-18" ~ 45cm  |                    | M8717                       | M8318<br><b>VB</b><br>M8718                       | M8800      |
|  | Blue Monofilament | 24" ~ 60cm    |                    | 8325H<br><b>VB</b><br>8725H | 8726H   |            |
|  | Blue Monofilament | 4-24" ~ 60cm  |                    | M8325<br><b>VB</b><br>M8725 | M8726   | M8727      |
|  | Blue Monofilament | 30" ~ 75cm    | 9706H<br><b>HS</b> |                             | 8307H<br><b>VB</b><br>8706H                       | 8708H      |
|  | Blue Monofilament | 2-30" ~ 75cm  |                    |                             | M8206   |            |
|  | Blue Monofilament | 4-30" ~ 75cm  |                    |                             | M8307<br><b>VB</b><br>M8706<br>M8706 <sup>1</sup> |            |
|  | Blue Monofilament | 36" ~ 90cm    | 8519H              | 8321H<br><b>VB</b><br>8720H |   |            |
|  | Blue Monofilament | 2-36" ~ 90cm  |                    | M8220                       |   |            |
|  | Blue Monofilament | 4-36" ~ 90cm  |                    | M8321<br><b>VB</b><br>M8720 |   |            |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>10-0</b> | <b>9-0</b> |
|---|-------------------|---------------|-------------|------------|
| <b>BV100-3</b> 3.9mm<br>3/8 Circle<br> | Blue Monofilament | 5" ~ 13cm     | 2794G       | 2793G      |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>10-0</b> |
|---------------|-----------------|---------------|-------------|
|---------------|-----------------|---------------|-------------|

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>                      | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|------------------------------------|-------------------|---------------|-------------|
| <b>BV100-4</b> 5.1mm<br>3/8 Circle | Blue Monofilament | 2" ~ 5cm      | 2795G       |

| <b>Needle</b>                     | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|-----------------------------------|-------------------|---------------|-------------|
| <b>BV75-3</b> 3.8mm<br>3/8 Circle | Blue Monofilament | 5" ~ 13cm     | 2790G       |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>7-0</b>       | <b>8-0</b>        |
|--|-------------------|---------------|------------------|-------------------|
| <b>BV175-7</b> 8.0mm<br><b>BV175-7</b><br>3/8 Circle | Blue Monofilament | 4-18" ~ 45cm  | M8764<br>M876433 |                   |
|  | Blue Monofilament | 2-18" ~ 45cm  | M8765<br>M876533 |                   |
|  | Blue Monofilament | 24" ~ 60cm    | 8766H<br>8766H33 | 8767 <sup>1</sup> |
|  | Blue Monofilament | 2-24" ~ 60cm  | M8766<br>M876633 | M8767<br>M876733  |
|  | Blue Monofilament | 4-24" ~ 60cm  | M8775            |                   |
|  | Blue Monofilament | 2-30" ~ 75cm  | M8768            |                   |
|  | Blue Monofilament | 4-30" ~ 75cm  | M8785            |                   |

## TAPERCUT

| <b>Needle</b>                  | <b>Material</b>   | <b>Length</b> | <b>0</b> |
|--------------------------------|-------------------|---------------|----------|
| <b>V-34</b> 36mm<br>1/2 Circle | Blue Monofilament | 30" ~ 75cm    | 8444H    |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b>                    | <b>5-0</b> |
|--|-------------------|---------------|------------|------------|-------------------------------|------------|
| <b>V-5</b> 17mm<br><b>V-5</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 8937H      | 8936H      | 8935H<br>P8935 <sup>213</sup> | 8934H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> | <b>4-0</b> |
|--|-------------------|---------------|----------|------------|------------|------------|
| <b>V-7</b> 26mm<br><b>V-7</b><br>1/2 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 8978H    | 8977H      | 8976H      | 8975H      |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>6-0</b>                               | <b>7-0</b>                               |
|---|-------------------|---------------|--|--|
| <b>CC</b> 9.3mm<br><b>CC</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    |  | 8705H<br>8705H                           |
|   | Blue Monofilament | 4-18" ~ 45cm  |  | M8705 <sup>1</sup><br>M8705 <sup>1</sup> |
|   | Blue Monofilament | 24" ~ 60cm    | 8807H<br>8807H                           | 8704H<br>8704H01                         |
|   | Blue Monofilament | 4-24" ~ 60cm  | M8807 <sup>1</sup><br>M8807 <sup>1</sup> | M8704 <sup>1</sup><br>M8704 <sup>1</sup> |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>5-0</b> | <b>6-0</b> |
|---------------|-----------------|---------------|------------|------------|
|---------------|-----------------|---------------|------------|------------|

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>5-0</b>                               | <b>6-0</b>                               |
|--|-------------------|---------------|--|--|
| <b>CC-1</b> 13mm<br><b>CC-1</b><br>3/8 Circle<br> | Blue Monofilament | 18" ~ 45cm    |  | 8719H<br>8719H                           |
|  | Blue Monofilament | 4-18" ~ 45cm  |  | M8719 <sup>1</sup><br>M8719 <sup>1</sup> |
|  | Blue Monofilament | 30" ~ 75cm    |  | 8707H<br>8707H01                         |
|  | Blue Monofilament | 4-30" ~ 75cm  |  | M8707 <sup>1</sup><br>M8707 <sup>1</sup> |
|  | Blue Monofilament | 36" ~ 90cm    | 8721H                                    |  |
|  | Blue Monofilament | 4-36" ~ 90cm  | M8721 <sup>1</sup><br>M8721 <sup>1</sup> |  |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>3-0</b> | <b>4-0</b> |
|--|-------------------|---------------|------------|------------|
| <b>V-4</b> 17mm<br><b>V-4</b><br>3/8 Circle<br> | Blue Monofilament | 36" ~ 90cm    | 8926H      | 8925H      |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|--|-------------------|---------------|-------------|
| <b>CIF-4</b> 13mm<br><b>CIF-4</b><br>1/4 Circle<br> | Blue Monofilament | 8" ~ 20cm     | 788G        |

## **ULTIMA - Spatula**

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|---|-------------------|---------------|-------------|
| <b>CS160-6</b> 5.5mm<br><b>CS160-6</b><br>1/2 Circle<br> | Blue Monofilament | 12" ~ 30cm    | 7761G       |

## Sutures : Non Absorbable / PROLENE Suture

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|--|-------------------|---------------|-------------|
| <b>CS160-6</b> 5.5mm<br><b>CTC-6</b> 13mm<br>1/2 Circle<br> | Blue Monofilament | 4" ~ 10cm     | 9090G       |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|---|-------------------|---------------|-------------|
| <b>CS160-8</b> 5.5mm<br><b>CS160-8</b><br>1/2 Circle<br> | Blue Monofilament | 6" ~ 15cm     | 1755G       |

| <b>Needle</b>  | <b>Material</b>   | <b>Length</b> | <b>10-0</b>  |
|--|-------------------|---------------|--|
| <b>CTC-6L</b> 13mm<br><b>CTC-6L</b><br>1/4 Circle<br> | Blue Monofilament | 4" ~ 10cm     | 9091G<br> |

| <b>Needle</b>   | <b>Material</b>   | <b>Length</b> | <b>10-0</b> |
|---|-------------------|---------------|-------------|
| <b>CS140-6</b> 6.5mm<br><b>CS140-6</b><br>3/8 Circle<br> | Blue Monofilament | 12" ~ 30cm    | 1710G       |

<sup>1</sup> One dozen packets per box

<sup>3</sup> Three dozen packets per box

<sup>5</sup> Single strand delivery

<sup>9</sup> 1 strip per box

<sup>21</sup> 1/4 x 1/8 x 1/16 (7.0mm x 3.0mm x 1.5 mm) Firm TFE Polymer Pledgets

 **Multipass® NEEDLES**

## **Sutures : Non Absorbable / PROLENE Suture**

 CR/ CONTROL RELEASE

 HS HEMO-SEAL

 VB VISI-BLACK

# Surgical Stainless Steel Suture

## Conventional Cutting - Sternum

| Needle   | Material                 | Length       | none    | none                          | 4     | 5     | 6     | 7     |
|--|--------------------------|--------------|---------|-------------------------------|-------|-------|-------|-------|
| <b>CCS 48mm</b><br>1/2 Circle<br> |                          | 4-18" ~ 45cm | M652G33 | M653G33<br>M654G33<br>M655G33 |       |       |       |       |
|  |                          | 2-18" ~ 45cm |         | M657G33                       |       |       |       |       |
|  | Monofilament<br>B & S 22 | 4-18" ~ 45cm |         |                               | M652G |       |       |       |
|  | Monofilament<br>B & S 20 | 4-18" ~ 45cm |         |                               |       | M653G |       |       |
|  | Monofilament<br>B & S 20 | 2-18" ~ 45cm |         |                               |       | M657G |       |       |
|  | Monofilament<br>B & S 19 | 4-18" ~ 45cm |         |                               |       |       | M654G |       |
|  | Monofilament<br>B & S 18 | 4-18" ~ 45cm |         |                               |       |       |       | M655G |

| Needle   | Material                 | Length       | none    | 5     |
|--|--------------------------|--------------|---------|-------|
| <b>CCS-1 60mm</b><br>1/2 Circle<br> |                          | 2-30" ~ 75cm | M436G33 |       |
|  | Monofilament<br>B & S 20 | 2-30" ~ 75cm |         | M436G |

## Reverse Cutting

| Needle  | Material                 | Length     | none   | 3-0  |
|---|--------------------------|------------|--------|------|
| <b>CP-1 36mm</b><br>1/2 Circle<br> |                          | 20" ~ 50cm | 497G33 |      |
|   | Monofilament<br>B & S 30 | 20" ~ 50cm |        | 497G |

## Sutures : Non Absorbable / Surgical Stainless Steel Suture

---

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>3-0</b> |
|---|--------------------------|---------------|------------|
| <b>FS-1</b> 24mm<br>3/8 Circle<br> | Monofilament<br>B & S 30 | 18" ~ 45cm    | 631G       |

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>none</b> | <b>3-0</b> | <b>4-0</b> |
|---|--------------------------|---------------|-------------|------------|------------|
| <b>FS-2</b> 19mm<br>3/8 Circle<br> |                          | 18" ~ 45cm    | 603G33      |            |            |
|   | Monofilament<br>B & S 30 | 18" ~ 45cm    |             | 607G       |            |
|   | Monofilament<br>B & S 32 | 18" ~ 45cm    |             |            | 603G       |

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>none</b> | <b>3-0</b> |
|--|--------------------------|---------------|-------------|------------|
| <b>FSL</b> 30mm<br>3/8 Circle<br> |                          | 18" ~ 45cm    | 606G33      |            |
|  | Monofilament<br>B & S 30 | 18" ~ 45cm    |             | 606G       |

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>none</b> | <b>2-0</b> |
|---|--------------------------|---------------|-------------|------------|
| <b>FSLX</b> 40mm<br>3/8 Circle<br> |                          | 18" ~ 45cm    | 618G33      |            |
|   | Monofilament<br>B & S 28 | 18" ~ 45cm    |             | 618G       |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2</b> |
|---------------|-----------------|---------------|----------|
|               |                 |               |          |

## Sutures : Non Absorbable / Surgical Stainless Steel Suture

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>2</b> |
|--|--------------------------|---------------|----------|
| <b>LR</b> 75mm<br><b>LR</b><br>3/8 Circle<br> | Monofilament<br>B & S 24 | 30" ~ 75cm    | 467G     |

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>none</b> | <b>2</b> |
|---|--------------------------|---------------|-------------|----------|
| <b>LS-1</b> 45mm<br>3/8 Circle<br> |                          | 4-20" ~ 50cm  | M400G33     |          |
|   | Monofilament<br>B & S 24 | 4-20" ~ 50cm  |             | M400G    |

## Straight Cutting Needles

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>5-0</b> |
|---|--------------------------|---------------|------------|
| <b>KS</b> 60mm<br>Straight<br> | Monofilament<br>B & S 35 | 8-16" ~ 40cm  | M645G      |

## SUTUPAK Pre-Cut Sutures

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>0</b> | <b>2-0</b> | <b>3-0</b> |
|--|--------------------------|---------------|----------|------------|------------|
| <b>SUTUPAK Pre-Cut Sutures</b><br>N/A<br> | Monofilament<br>B & S 26 | 17-18" ~ 45cm | SS26G    |            |            |
|  | Monofilament<br>B & S 28 | 17-18" ~ 45cm |          | SS28G      |            |
|  | Monofilament<br>B & S 30 | 17-18" ~ 45cm |          |            | SS30G      |

## Sutures : Non Absorbable / Surgical Stainless Steel Suture

### Taper Point

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>2-0</b> |
|---|--------------------------|---------------|------------|
| <b>CT-1</b> 36mm<br>1/2 Circle<br> | Monofilament<br>B & S 28 | 8-18" ~ 45cm  | M682G      |

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>2-0</b> |
|--|--------------------------|---------------|------------|
| <b>CTX</b> 48mm<br>1/2 Circle<br> | Monofilament<br>B & S 28 | 8-24" ~ 60cm  | M414G      |

### TAPERCUT - Straight Cutting

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>4-0</b>         |
|---|--------------------------|---------------|--------------------|
| <b>TS_TS</b> 48mm<br>TS Straight<br><br>Polypropylene<br> | Monofilament<br>B & S 32 | 2-14" ~ 35cm  | 540G <sup>41</sup> |

| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>4-0</b>         |
|--|--------------------------|---------------|--------------------|
| <b>V-26</b> 26mm<br>3/8 Circle<br><br>Polypropylene | Monofilament<br>B & S 32 | 2-14" ~ 35cm  | 539G <sup>41</sup> |

## **TAPERCUT**

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>1</b> |
|---|--------------------------|---------------|----------|
| <b>V-37</b> 40mm<br>1/2 Circle<br> | Monofilament<br>B & S 25 | 4-18" ~ 45cm  | M660G    |


| <b>Needle</b>  | <b>Material</b>          | <b>Length</b> | <b>4</b> | <b>5</b> | <b>6</b> |
|--|--------------------------|---------------|----------|----------|----------|
| <b>V-40</b> 48mm<br>1/2 Circle<br> | Monofilament<br>B & S 22 | 4-18" ~ 45cm  | M651G    |          |          |
|  | Monofilament<br>B & S 20 | 4-18" ~ 45cm  |          | M650G    |          |
|  | Monofilament<br>B & S 20 | 2-18" ~ 45cm  |          | M656G    |          |
|  | Monofilament<br>B & S 19 | 4-18" ~ 45cm  |          |          | M649G    |

| <b>Needle</b>   | <b>Material</b>          | <b>Length</b> | <b>5</b> |
|---|--------------------------|---------------|----------|
| <b>V-56</b> 60mm<br>1/2 Circle<br> | Monofilament<br>B & S 20 | 4-30" ~ 75cm  | M435G    |

## **Temporary Pacing Wire**

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / Surgical Stainless Steel Suture

| <b>Needle</b>   | <b>Material</b>                  | <b>Length</b> | <b>2-0</b>         |
|---|----------------------------------|---------------|--------------------|
| <b>BB</b> 17mm<br><b>SKS</b> 60mm<br><br>3/8 Circle  | Dk. Blue Braided Stainless Steel | 24" ~ 60cm    | TPW20 <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>                  | <b>Length</b> | <b>2-0</b>         |
|--|----------------------------------|---------------|--------------------|
| <b>RB-1</b> 17mm<br><b>SKS</b> 60mm<br><br>1/2 Circle Straight Cutting  | Dk. Blue Braided Stainless Steel | 24" ~ 60cm    | TPW10 <sup>1</sup> |

| <b>Needle</b>  | <b>Material</b>                             | <b>Length</b> | <b>0</b>           | <b>2-0</b> |
|--|---|---------------|--------------------|------------|
| <b>SH</b> 26mm<br><b>SKS-3</b> 90mm<br><br>1/2 Circle Straight Cutting  | Dk. Blue Braided Stainless Steel            | 24" ~ 60cm    | TPW50 <sup>1</sup> | TPW90      |
|  | Lt. Blue & Dk. Blue Braided Stainless Steel | 2-24" ~ 60cm  | TPW52              | TPW92      |
|  | Dk. Blue Braided Stainless Steel            | 24" ~ 60cm    | TPW30<br>TPW70     |            |
|  | Lt. Blue & Dk. Blue Braided Stainless Steel | 2-24" ~ 60cm  | TPW32<br>TPW72     |            |

| <b>Needle</b> | <b>Material</b> | <b>Length</b> | <b>2-0</b> |
|---------------|-----------------|---------------|------------|
|---------------|-----------------|---------------|------------|

## Sutures : Non Absorbable / Surgical Stainless Steel Suture

| <b>Needle</b>  | <b>Material</b>                             | <b>Length</b> | <b>2-0</b>           |
|--|---|---------------|----------------------|
| <b>BB-1</b> 22mm<br><b>SKS-3</b> 90mm<br>3/8 Circle<br> | Dk. Blue Braided Stainless Steel            | 24" ~ 60cm    | TPW40 <sup>168</sup> |
|  | Lt. Blue & Dk. Blue Braided Stainless Steel | 2-24" ~ 60cm  | TPW42                |

| <b>Needle</b>  | <b>Material</b>                             | <b>Length</b> | <b>0</b>             | <b>2-0</b> |
|--|---|---------------|----------------------|------------|
| <b>TE</b> 31mm<br><b>SKS-3</b> 90mm<br>3/8 Circle<br> | Lt. Blue & Dk. Blue Braided Stainless Steel | 2-24" ~ 60cm  | TPW62 <sup>111</sup> |            |
|  | Dk. Blue Braided Stainless Steel            | 24" ~ 60cm    |                      | TPW80      |

<sup>1</sup> One dozen packets per box

<sup>11</sup> Multistrand 2" ~ 60 cm

<sup>41</sup> Pull out wire set (includes one polypropylene button)

<sup>68</sup> 1.5" long scoremarks; others .875"

| <b>Product Code</b> | <b>Material</b>   | <b>Size</b>  | <b>Packaging</b>      |
|---------------------|---|--|-----------------------|
| CLR602US            | DERMABOND PRINEO Skin Closure System 2-Octyl Cyanoacrylate with mesh applicator | 1~60cm Mesh Dispenser /<br>1~3.8ml Adhesive Applicator | 2 applicators per box |

# **DERMABOND® Adhesive**

| <b>Product Code</b> | <b>Material</b>   | <b>Size</b> | <b>Packaging</b>       |
|---------------------|---|-------------|------------------------|
| DHV12               | DERMABOND Topical Skin Adhesive<br>2-Octyl Cyanoacrylate - High Viscosity                                 | 0.5ml       | 12 Applicators Per Box |
| DPP6                | DERMABOND Propen Topical Skin Adhesive 2-Octyl Cyanoacrylate - High Viscosity with Optional Precision Tip | 0.5ml       | 6 Applicators Per Box  |
| DPPXL6              | DERMABOND Propen XL Topical Skin Adhesive 2-Octyl Cyanoacrylate - High Viscosity                          | 0.75ml      | 6 Applicators Per Box  |

| <b>Product Code</b> | <b>Material</b>   | <b>Size</b> | <b>Packaging</b>       |
|---------------------|---|-------------|------------------------|
| DNX12*              | DERMABOND ADVANCED Topical Skin Adhesive 2-Octyl Cyanoacrylate - Unique formulation with an easy-to-use applicator for precise control on wide and fine line applications | 0.7ml       | 12 Applicators Per Box |
| DNX6*               | DERMABOND ADVANCED Topical Skin Adhesive 2-Octyl Cyanoacrylate - Unique formulation with an easy-to-use applicator for precise control on wide and fine line applications | 0.7ml       | 6 Applicators Per Box  |

# BLAKE® Drains

## BLAKE Cardio Connectors

| Description | Product Code | Material                   |
|-------------|--------------|----------------------------|
|             | BCC1         | BLAKE Cardio Connector 1:1 |
|             | BCC2         | BLAKE Cardio Connector 2:1 |
|             | BCC3         | BLAKE Cardio Connector 3:1 |

## BLAKE Silicone Drains - Flat

| Description   | Product Code | Material                             |
|---|--------------|--------------------------------------|
|  | 2210         | 7mm Flat, 3/4 Fluted                 |
|   | 2211         | 7mm Flat, Full Fluted                |
|   | 2212         | 7mm Flat, Full Fluted, 3/16" Trocar  |
|   | 2213         | 10mm Flat, 3/4 Fluted                |
|   | 2214         | 10mm Flat, Full Fluted               |
|   | 2215         | 10mm Flat, Full Fluted, 3/16" Trocar |
|   | 2216         | 7mm Flat, 3/4 Fluted, 3/16" Trocar   |
|   | 2217         | 10mm Flat, 3/4 Fluted, 3/16" Trocar  |

## BLAKE Silicone Drains - Hubless

| Description | Product Code | Material |
|-------------|--------------|----------|
|             |              |          |

| Description   | Product Code | Material                                    |
|---|--------------|---|
|  | 2226         | Hubless, 10 FR Round                        |
|   | 2227         | Hubless, 10 FR Round, 1/8" Trocar           |
|   | 2228         | Hubless, 15 FR Round                        |
|   | 2229         | Hubless, 15 FR Round, 3/16" Trocar          |
|   | 2230         | Hubless, 19 FR Round                        |
|   | 2231         | Hubless, 19 FR Round, 1/4" Trocar           |
|   | 2232         | Hubless, 19 FR Round, 1/4" Bendable Trocar  |
|   | 2233         | Hubless, 15 FR Round, 3/16" Bendable Trocar |
|   | 2234         | Hubless, 24 FR Round                        |

# J-VAC Drains

## J-VAC Reservoirs - Sterile

| Description   | Product Code | Material                               |
|---|--------------|--|
|  | 2161         | 150 ml Closed Wound Drainage Reservoir |
|   | 2162         | 450 ml Closed Wound Drainage Reservoir |
|   | 2163         | 300 ml Closed Wound Drainage Reservoir |
|   | 2177         | 150 ml Closed Wound Drainage Reservoir |
|   | 2178         | 450 ml Closed Wound Drainage Reservoir |
|   | 2179         | 300 ml Closed Wound Drainage Reservoir |

## J-VAC Reservoirs - Sterile Bulbs

| Description   | Product Code | Material                      |
|---|--------------|-------------------------------|
|  | 2160         | 100 cc Bulb Suction Reservoir |

## J-VAC Silicone Drains - Sterile


| Description   | Product Code | Material                            |
|---|--------------|-------------------------------------|
|  | 2193         | J-VAC, Silicone Drain, 19 FR T Tube |
|   | 2209         | J-VAC Drain Adapters                |
|   | 2298         | J-VAC Drain Adapters                |

## ETHICON SECURESTRAP 5mm Absorbable Strap Fixation Device

| Product Code | Material  | Size | Packaging         |
|--------------|---|------|-------------------|
| STRAP12      | 5mm absorbable Fixation Device                              |      | 6 devices per box |
| STRAP25      | 5mm single-use mesh fixation device w/ 25 absorbable straps |      | 6 devices per box |

# MERSILENE Mesh

## Flat Mesh

| Product Code   | Material       | Size                        | Packaging        |
|--|----------------|-----------------------------|------------------|
| RML<br>   | MERSILENE Mesh | 12" x12"<br>(30cm x 30cm)   | 3 Pieces Per Box |
| RMS<br> | MERSILENE Mesh | 2.5" x 4.5"<br>(6cm x 11cm) | 3 Pieces Per Box |

## Tissue Separating Mesh

| <b>Product Code</b> | <b>Material</b>         | <b>Size</b>          | <b>Packaging</b> |
|---------------------|-------------------------|----------------------|------------------|
| PHY0715R            | Flexible Composite Mesh | 7 X 15 cm Rectangle  | 1 Piece Per Box  |
| PHY1015V            | Flexible Composite Mesh | 10 X 15 cm Oval      | 1 Piece Per Box  |
| PHY1515Q            | Flexible Composite Mesh | 15 X 15 cm Square    | 1 Piece Per Box  |
| PHY1520R            | Flexible Composite Mesh | 15 X 20 cm Rectangle | 1 Piece Per Box  |
| PHY1520V            | Flexible Composite Mesh | 15 X 20 cm Oval      | 1 Piece Per Box  |
| PHY2025V            | Flexible Composite Mesh | 20 X 25 cm Oval      | 1 Piece Per Box  |
| PHY2030R            | Flexible Composite Mesh | 20 X 30 cm Rectangle | 1 Piece Per Box  |
| PHY2535V            | Flexible Composite Mesh | 25 X 35 cm Oval      | 1 Piece Per Box  |
| PHY3035R            | Flexible Composite Mesh | 30 X 35 cm Rectangle | 1 Piece Per Box  |
| PHY3050R            | Flexible Composite Mesh | 30 X 50 cm Rectangle | 1 Piece Per Box  |

# PROCEED Mesh

## Tissue Separating Mesh

| Product Code | Material                                       | Size                           | Packaging       |
|--------------|--|--------------------------------|-----------------|
| PCDB1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 2" x 4"<br>(5cm x 10cm)        | 1 Piece Per Box |
| PCDD1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 4" x 8"<br>(10cm x 20cm)       | 1 Piece Per Box |
| PCDG1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 6" x 8"<br>(15cm x 20cm)       | 1 Piece Per Box |
| PCDH1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 8" x 10"<br>(20cm x 25cm)      | 1 Piece Per Box |
| PCDJ1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 8" x 12"<br>(20cm x 30cm)      | 1 Piece Per Box |
| PCDL1        | PROCEED Surgical Mesh (Tissue-Separating Mesh) | 12" x 12"<br>(30.5cm x 30.5cm) | 1 Piece Per Box |
| PCDM1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 6" x 6"<br>(15cm x 15cm)       | 1 Piece Per Box |
| PCDN1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 4" x 6"<br>(10cm x 15cm)       | 1 Piece Per Box |

## Hernia Repair: PROCEED Mesh

| Product Code | Material                                       | Size                         | Packaging       |
|--------------|--|------------------------------|-----------------|
| PCDR1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 3" x 6"<br>(7.5cm x 15cm)    | 1 Piece Per Box |
| PCDT1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 10" x 13"<br>(26cm x 34cm)   | 1 Piece Per Box |
| PCDW1        | PROCEED Surgical Mesh (Tissue Separating Mesh) | 10" x 14"<br>(25cm x 35.5cm) | 1 Piece Per Box |

# PROLENE 3D Patch Mesh

## Mesh Devices

| Product Code   | Material                                       | Size          | Packaging        |
|--|--|---------------|------------------|
| 3DPL<br>    | PROLENE 3D Patch Large with Extended Overlay   | 1 31/32"(5cm) | 3 Pieces Per Box |
| 3DPM<br>    | PROLENE 3D Patch Medium with Extended Overlay  | 1 3/8"(3.5cm) | 3 Pieces Per Box |
| P3DPL<br>  | PROLENE 3D Patch Large with Preshaped Overlay  | 1 31/32"(5cm) | 3 Pieces Per Box |
| P3DPM<br> | PROLENE 3D Patch Medium with Preshaped Overlay | 1 3/8"(3.5cm) | 3 Pieces Per Box |

# PROLENE Hernia System

## Mesh Devices



| Product Code | Material                       | Size      | Packaging         |
|--------------|--------------------------------|-----------|-------------------|
| PHSE<br>     | PROLENE Hernia System Extended | 4" (10cm) | 3 Devices Per Box |
| PHSE6<br>    | PROLENE Hernia System Extended | 4" (10cm) | 6 Devices Per Box |
| PHSL<br>     | PROLENE Hernia System Large    | 4" (10cm) | 3 Devices Per Box |
| PHSL6<br>    | PROLENE Hernia System Large    | 4" (10cm) | 6 Devices Per Box |
| PHSM<br>     | PROLENE Hernia System Medium   | 3" (7cm)  | 3 Devices Per Box |
| PHSM6<br>    | PROLENE Hernia System Medium   | 3" (7cm)  | 6 Devices Per Box |

## PROLENE Mesh & PROLENE Soft Mesh

### Flat Mesh

| Product Code | Material     | Size                     | Packaging        |
|--------------|--------------|--------------------------|------------------|
| PMH          | PROLENE Mesh | 6" x 6"<br>(15cm x 15cm) | 6 Sheets Per Box |
| PMII         | PROLENE Mesh | 3" x 6"<br>(8cm x 15cm)  | 6 Sheets Per Box |

## Hernia Repair: PROLENE Mesh & PROLENE Soft Mesh

| <b>Product Code</b>   | <b>Material</b>                               | <b>Size</b>                   | <b>Packaging</b> |
|---|---|-------------------------------|------------------|
| PML<br>    | PROLENE Mesh                                  | 12" x 12"<br>(30cm x 30cm)    | 3 Sheets Per Box |
| PMLK<br> | PROLENE Mesh<br>Pre-Shaped<br>Keyhole (Large) | 2.4" x 5.4"<br>(6cm x 13.7cm) | 6 Sheets Per Box |
| PMSK<br> | PROLENE Mesh<br>Pre-Shaped<br>Keyhole (Small) | 1.8" x 4.0"<br>(4.6cm x 10cm) | 6 Sheets Per Box |
| PMXL  | PROLENE Mesh                                  | 2" x 12"<br>(5cm x 30cm)      | 6 Sheets Per Box |
| PMXS<br> | PROLENE Mesh                                  | 1" x 4"<br>(2.54cm x 10cm)    | 6 Sheets Per Box |

## Hernia Repair: PROLENE Mesh & PROLENE Soft Mesh

| <b>Product Code</b>  | <b>Material</b>                 | <b>Size</b>                  | <b>Packaging</b> |
|--|---------------------------------|------------------------------|------------------|
| SPMH   | PROLENE Soft Polypropylene Mesh | 6" x 6"<br>(15cm x 15cm)     | 6 Sheets Per Box |
| SPMII  | PROLENE Soft Polypropylene Mesh | 3" x 6"<br>(7.6cm x 15cm)    | 6 Sheets Per Box |
| SPMLI  | PROLENE Soft Polypropylene Mesh | 10" x 10"<br>(25cm x 25cm)   | 3 Sheets Per Box |
| SPMS   | PROLENE Soft Polypropylene Mesh | 2" x 4"<br>(5cm x 10cm)      | 6 Sheets Per Box |
| SPMXS<br> | PROLENE Soft Polypropylene Mesh | 1" x 4"<br>(2.54cm x 10cm)   | 6 Sheets Per Box |
| SPMXXL   | PROLENE Soft Polypropylene Mesh | 14" x 12"<br>(35.6cm x 30cm) | 3 Sheets Per Box |

## PVP

### Tissue Separating Mesh

| Product Code | Material                     | Size          | Packaging        |
|--------------|------------------------------|---------------|------------------|
| PVPM         | PROCEED Ventral Patch Medium | 6.4cm x 6.4cm | 2 Pieces Per Box |
| PVPS         | PROCEED Ventral Patch Small  | 4.3cm x 4.3cm | 2 Pieces Per Box |

# ULTRAPRO Hernia System

## Mesh Devices

| Product Code | Material                      | Size                       | Packaging        |
|--------------|-------------------------------|----------------------------|------------------|
| UHSL         | ULTRAPRO Hernia System Large  | 4" (10cm)                  | 3 Pieces Per Box |
| UHSL6        | ULTRAPRO Hernia System Large  | 4" (10cm)                  | 6 Pieces Per Box |
| UHSM         | ULTRAPRO Hernia System Medium | 3" (7.5cm)                 | 3 Pieces Per Box |
| UHSM6        | ULTRAPRO Hernia System Medium | 3" (7.5cm)                 | 6 Pieces Per Box |
| UHSOV        | ULTRAPRO Hernia System Oval   | 4" x 4.7"<br>(10cm x 12cm) | 3 Pieces Per Box |

# **ULTRAPRO Mesh**

## **Flat Mesh**

| <b>Product Code</b> | <b>Material</b>                          | <b>Size</b>                 | <b>Packaging</b> |
|---------------------|--|-----------------------------|------------------|
| UML1                | ULTRAPRO<br>Partially<br>Absorbable Mesh | 12" x 12"<br>(30cm x 30cm)  | 1 Piece Per Box  |
| UMM3                | ULTRAPRO<br>Partially<br>Absorbable Mesh | 6" x 6"<br>(15cm x 15cm)    | 3 Pieces Per Box |
| UMR3                | ULTRAPRO<br>Partially<br>Absorbable Mesh | 3" x 6"<br>(7.6cm x 15cm)   | 3 Pieces Per Box |
| UMS3                | ULTRAPRO<br>Partially<br>Absorbable Mesh | 2.4" x 4.3"<br>(6cm x 11cm) | 3 Pieces Per Box |



# ULTRAPRO Plug

## Mesh Devices

| Product Code | Material             | Size                  | Packaging         |
|--------------|----------------------|-----------------------|-------------------|
| UPPL2        | ULTRAPRO Plug Large  | Anchor 5cm<br>Rim 5cm | 2 Devices Per Box |
| UPPL6        | ULTRAPRO Plug Large  | Anchor 5cm<br>Rim 5cm | 6 Devices Per Box |
| UPPM2        | ULTRAPRO Plug Medium | Anchor 4cm<br>Rim 5cm | 2 Devices Per Box |
| UPPM6        | ULTRAPRO Plug Medium | Anchor 4cm<br>Rim 5cm | 6 Devices Per Box |
| UPPS2        | ULTRAPRO Plug Small  | Anchor 3cm<br>Rim 5cm | 2 Devices Per Box |
| UPPS6        | ULTRAPRO Plug Small  | Anchor 3cm<br>Rim 5cm | 6 Devices Per Box |

# VICRYL Woven Mesh

## Flat Mesh

| Product Code  | Material               | Size                        | Packaging        |
|---|------------------------|-----------------------------|------------------|
| VKML<br>   | VICRYL Mesh Knitted    | 12" x 12"<br>(30cm x 30cm)  | 3 Pieces Per Box |
| VKMM<br> | VICRYL Mesh Knitted    | 6" x 6"<br>(15cm x 15cm)    | 3 Pieces Per Box |
| VM101K  | VICRYL Mesh Bag Renal  | 6" x 3"<br>(15cm x 8cm)     | 1 Piece Per Box  |
| VM102K  | VICRYL Mesh Bag Renal  | 7" x 4"<br>(18cm x 10cm)    | 1 Piece Per Box  |
| VM106K  | VICRYL Mesh Bag Spleen | 10" x 12"<br>(26cm x 30cm)  | 1 Piece Per Box  |
| VM210K  | VICRYL Mesh Bag Liver  | 17" x 12"<br>(44cm x 30cm)  | 1 Piece Per Box  |
| VM220K  | VICRYL Mesh Bag Liver  | 20" x 8"<br>(50cm x 20.5cm) | 1 Piece Per Box  |

## Hernia Repair: VICRYL Woven Mesh

| Product Code  | Material          | Size                       | Packaging        |
|---|-------------------|----------------------------|------------------|
| VWML<br> | VICRYL Mesh Woven | 12" x 12"<br>(30cm x 30cm) | 3 Pieces Per Box |
| VWMM<br> | VICRYL Mesh Woven | 6" x 6"<br>(15cm x 15cm)   | 3 Pieces Per Box |

# Nonsterile Products



## ETHI-PACK Pre-cut Sutures - Surgical 316L

| Product Code       | Material                  | Size                                 | Packaging                    |
|--------------------|---------------------------|--------------------------------------|------------------------------|
| DS18 <sup>71</sup> | B&S 18 Steel Monofilament | 18" ~45cm Long x7 (.0403") Diameter  | Box (12 Tubes - 144 Strands) |
| DS20 <sup>71</sup> | B&S 20 Steel Monofilament | 18"~45cm Long x5 (.0320") Diameter   | Box (12 Tubes - 144 Strands) |
| DS22 <sup>71</sup> | B&S 22 Steel Monofilament | 18"~45cm Long x4 (.0254") Diameter   | Box (12 Tubes - 144 Strands) |
| DS23 <sup>71</sup> | B&S 23 Steel Monofilament | 18"~45cm Long x3 (.0226") Diameter   | Box (12 Tubes - 144 Strands) |
| DS24 <sup>71</sup> | B&S 24 Steel Monofilament | 18"~45cm Long x2 (.0201") Diameter   | Box (12 Tubes - 144 Strands) |
| DS25 <sup>71</sup> | B&S 25 Steel Monofilament | 18"~45cm Long x1 (.0179") Diameter   | Box (12 Tubes - 144 Strands) |
| DS26 <sup>70</sup> | B&S 26 Steel Monofilament | 18"~45cm Long x0 (.0159") Diameter   | Box (12 Tubes - 600 Strands) |
| DS28 <sup>70</sup> | B&S 28 Steel Monofilament | 18"~45cm Long x2-0 (.0126") Diameter | Box (12 Tubes - 600 Strands) |
| DS30 <sup>70</sup> | B&S 30 Steel Monofilament | 18"~45cm Long x3-0 (.0100") Diameter | Box (12 Tubes - 600 Strands) |
| DS32 <sup>70</sup> | B&S 32 Steel Monofilament | 18"~45cm Long x4-0 (.0080") Diameter | Box (12 Tubes - 600 Strands) |
| DS35 <sup>70</sup> | B&S 35 Steel Monofilament | 18"~45cm Long x5-0 (.0056") Diameter | Box (12 Tubes - 600 Strands) |



<sup>70</sup> Each tube contains fifty 18" longstrands. Each box contains 12 tubes.

<sup>71</sup> Each tube contains twelve 18" longstrands. Each box contains 12 tubes.

# Storage Supplies

| Product Code | Description   |
|--------------|---|
| CART         | <p>CART, dimensions: 38"L x 18"W x38"H. No Suture racks included. Brackets to attach racks included. Some assembly required.</p>                       |
| LS3          | <p>3-wide LAZY SUSAN, LS3, dimensions: 15 3/8"L x 11 1/2"W x 20 3/4"H. Includes 6 modular racks. Can be expanded with additional modular racks.</p>  |
| LS4          | <p>4-wide LAZY SUSAN, LS4, dimensions: 20 1/2"L x 11 1/2"W x 20 3/4"H. Includes 8 modular racks. Can be expanded with additional modular racks.</p>  |
| MR20         | <p>DISPENSER RACK, MR20, dimensions: 25 3/4L x 5 3/4W x 27H. Includes 20 modules, 10 back plates. (each with 4 fasteners), 20 module clips, 30 divider panels. Assembly instruction sheet enclosed.</p>                                 |
| MR4          | <p>DISPENSER RACK, MR4, dimensions: 10 1/4L x 5 3/4W x 13 1/2H. Includes 4 modules, 3 back plates. (each with 4 fasteners), 4 module clips, 8 divider panels. Assembly inst. sheet encl. Shipped as comp. 4-module unit.</p>            |

## Miscellaneous: Storage Supplies

| Product Code  | Description  |
|---|--|
| PCR1<br><br>   | PORTABLE CARRYING RACK, PCR1, dimensions: 15 3/8" L x 6 1/8" W x 17 5/8" H. Includes 6 modular racks.  |
| PCR2<br><br>  | 2-SIDED PORTABLE CARRYING RACK, PCR2, dimensions: 15 3/8" L x 12 1/4" W x 16 3/4" H. Includes 12 modules racks.  |
| SPR1<br><br> | I.V. POLE WITH RACK, SPR1, dimensions: 18 1/4" L x 16 1/2" W x 56 1/2" H. Includes 18 modules racks, 15 module clips, 30 divider panels. Some assembly required. |

# Surgical Accessories

## Adjustable Retention Suture Bridge

| Product Code | Material   | Size                              | Packaging                         |
|--------------|--|-----------------------------------|-----------------------------------|
| RSB5         | Disposable Retention Suture Bridge for patients requiring retention sutures. Permits postoperative wound management. | 1/2"W x 4.375"L<br>(1.3cm x 11cm) | 3 Per Packet<br>6 Packets Per Box |

## Polypropylene Buttons

| Product Code  | Material                                 | Size | Packaging                          |
|---|--|------|------------------------------------|
| 520G<br> | For tendon repair in orthopedic surgery. | N/A  | 3 Per Packet<br>12 Packets Per Box |

## Retention Suture Bolster

| Product Code | Material                        | Size  | Packaging                          |
|--------------|---------------------------------|---|------------------------------------|
| 450G         | Surgical latex tubing, sterile. | 2.5" x 3/16" x 1/32"<br>(6cm x .48cm x .08cm) | 6 Per Packet<br>12 Packets Per Box |

## Umbilical Tape

| Product Code | Material              | Size                         | Packaging                    |
|--------------|-----------------------|------------------------------|------------------------------|
| U10T         | Cotton Umbilical Tape | 1/8" x 18"<br>(.32cm x 45cm) | 2 Strands<br>2 Dozen Packets |
| U11T         | Cotton Umbilical Tape | 1/8" x 30"<br>(.32cm x 75cm) | 2 Strands<br>2 Dozen Packets |
| U12T         | Cotton Umbilical Tape | 1/8" x 36"<br>(.32cm x 90cm) | 3 Strands<br>2 Dozen Packets |

## Miscellaneous: Surgical Accessories

| Product Code | Material                     | Size                         | Packaging                    |
|--------------|------------------------------|------------------------------|------------------------------|
| U15G         | Radiopaque<br>Umbilical Tape | 1/8" x 18"<br>(.32cm x 45cm) | 2 Strands<br>1 Dozen Packets |
| U16G         | Radiopaque<br>Umbilical Tape | 1/8" x 24"<br>(.32cm x 60cm) | 2 Strands<br>1 Dozen Packets |

# **TFE Polymer Pledgets**

| <b>Product Code</b> | <b>Material</b>             | <b>Size</b>                                     | <b>Packaging</b>                  |
|---------------------|-----------------------------|---|-----------------------------------|
| PCP20               | TFE Polymer Pledgets (Firm) | 1/4" x 1/8" 1/16"<br>(7mm x 3mm x 1.5mm)        | Box (36 Envelopes - 216 Pledgets) |
| PCP30               | TFE Polymer Pledgets (Firm) | 3/8" x 3/16" x 1/16"<br>(9.5mm x 4.8mm x 1.5mm) | Box (36 Envelopes - 216 Pledgets) |
| PCP40               | TFE Polymer Pledgets (Soft) | 3/8" x 3/16" x 1/16"<br>(9.5mm x 4.8mm x 1.5mm) | Box (36 Envelopes - 216 Pledgets) |
| PCP50               | TFE Polymer Pledgets (Soft) | 1/4" x 1/8" x 1/16"<br>(7mm x 3mm x 1.5mm)      | Box (36 Envelopes - 216 Pledgets) |