

## High Suction Tips

---

**Intended Use**

High Volume Evacuator (HVE) with 45° angle tip is intended to minimize aerosols, splashes and help maintaining a clean, clear operative field for healthcare providers during dental procedures. It provides control of fluids in the patient's mouth by quickly removing massive volumes of fluid. Fluids can be picked up without touching the surface of the tissue. The suction tip is reversible to allow the dental specialist to reach all areas of the mouth.

**Item Codes**

Item Code	Description
505	White non-vented
510	White vented one side

**Packaging**

50 units per bag  
20 bags per case  
1000 units per case

**Size**

505 = 140 mm (5.50")  
510 = 147 mm (5.77")

**Material**

Plastic

**Usage**

For single-use only

**Sterility**

Non-Sterile

**Country of Origin**

USA