



Medicom®  
SafeMask™

Feature guide  
to Medicom  
earloop face masks




1




**Non-woven spunbond outer layer material** for maximum protection against fluids resistance that will not lint, tear or shred

2



**Highest quality filter media** for best ASTM ratings (Level 1 through 3, filtration value and breathability)


3



**Inner Layer** that will not lint, tear or shred


- Non-woven spunbond inner layer made of white dye-free fluid resistant material.
- Cellulose inner layer is designed with sensitive skin in mind. This soft, dye-free, lightweight and high breathable inner layer is gentle and non-irritating. It remains cool on the skin for ultimate ease and comfort. (Available only on SafeMask® Sof Skin® line)

4




Soft, strong, sonically sealed ample-length **earloops** so there is no pulling. They fit comfortably and attached to the outside of the mask to eliminate irritation

5



Adjustable, extra-long, **enclosed aluminum nose piece** forms a strong seal for maximum protection and comfort


6



**Aluminum chin piece** creates an optimal personalized and comfortable fit around the face. It significantly reduces gapping at the bottom and sides of the mask and provides good face coverage


- Available only on the SafeMask® TailorMade line

7




**Shingle pleat design** prevents fluid pooling and cross contamination. When mask is worn, the pleats should always be facing down

8



**Embossed Medicom logo.** Logo should always appear at the bottom left hand side

9



**Soft bond ultrasonics** ensures strongest construction, eliminating holes and defects