

Medicom®
SafeMask®

NEW
ASTM
LEVEL **2**

Master Series

Feel the color



ASTM
LEVEL **1** Master
Series

● ● ● ●

● ● ● ●

ASTM
LEVEL **2** Master
Series Plus

● ● ● ●

●

#feelcolor



SafeMask® Master Series

ASTM
LEVEL 1
AND
ALL NEW
ASTM
LEVEL 2

Play With Color Not Protection

- A beautiful line of masks inspired by Augusta's striking scenery
- Designed in a variety of unique and vibrant colors so that your mask is as personalized as the dental care you provide
- These ASTM Level 1 and 2 masks offer the comfort, fit, and protection you've come to expect from SafeMask® with the added benefit of vibrant joyful colors to brighten your daily routine

NEW COLOR



2054 | Southern Bellflower
Radiant Orchid



2050 | Augusta Sky
Royal Blue



2051 | Azalea Festival
Bright Fuchsia



2052 | Lush Lawn
Vibrant Green



2056 | Azalea Festival
Bright Fuchsia

NEW ASTM LEVEL 2

Performance factors

	ASTM LEVEL 1	ASTM LEVEL 2
Bacterial filtration efficiency (BFE)	≥ 95%	≥ 98%
Sub-micron particulate efficiency at 0.1 µm (PFE)	≥ 95%	≥ 98%
Resistant to penetration by synthetic blood	@80 mm Hg	@120 mm Hg
Differential Pressure (ΔP)	< 4.0 mm H ₂ O/cm ²	< 5.0 mm H ₂ O/cm ²
Flame Spread	Class 1	Class 1

Suggested use

Packaging configuration

Daily favourite for everyday procedures
50 masks per box, 10 boxes per case



Made in
USA

Medicom® Trust in Quality

Like you, we never compromise. Our commitment is to help to make the world safer and healthier by offering consistent, reliable protection. We use carefully selected materials and engineer our products to offer protection you can count on.

As the maker of Safe brands, we take pride in the protection we offer, so you can take pride in yours.

Medicom® Pride in Protection.

www.medicom.com Rev 00

