

MEDICOM

Safe+Mask®
Sof Skin®

ASTM STANDARDS

Characteristics

Low Barrier

Bacterial filtration efficiency (BFE)(%)	≥ 95
Differential Pressure ΔP (mm H ₂ O/cm ²)	< 4.0
Sub-micron particulate efficiency at 0.1 μm (PFE)(%)	≥ 95
Resistance to penetration by synthetic blood, minimum pressure in mm Hg for pass result	80
Flame Spread	Class 1

SO SOFT
YOU FORGET
YOU'RE WEARING IT

WEAR WHAT IS COMFORTABLE

Safe+Mask® Sof Skin® is designed with sensitive skin in mind. Its soft, lightweight and highly breathable inner layer is gentle, non-irritating and remains cool on the skin for ultimate ease and comfort without forgetting ASTM approved proper protection. Wrap yourself in softness for all day comfort with Safe+Mask® Sof Skin®.



FEATURES

Cellulose-blend inner layer

ASTM Standard

Great fit

Unique shingle design

BENEFITS

White, non-irritating, extra soft and dye-free.

Leaves skin feeling clean and dry.

It will not lint, tear or shred.

Has excellent breathability.

Meets ASTM 2100 standard for everyday, low barrier protection.

Adjustable nose piece to help customize the fit.

Earloops ties are soft, stretchy, latex-free and sonically sealed.

Ease of removal without self-contamination.

Prevents fluid pooling and cross-contamination.

2080
White - Low Barrier



2081
Blue - Low Barrier



2082
Pink - Low Barrier



2083
Lavender - Low Barrier



Packaging: 50 masks per box, 10 boxes per case

MASKS

GAUZE

STERILIZATION

VARNISH

PREVENTIVES

DISINFECTANTS

REV 00

Medicom® is dedicated to manufacturing and distributing trusted high quality infection control and patient care products worldwide. Medicom® leading brand names include Safe-Touch®, Safe+Mask®, Safe-Gauze®, Duraflor® Halo, Denti-Care and Dry-Back®.

medicom.com

for free samples, visit www.safemask.ca

PROUDLY MADE IN
USA
Medicom Quality You Can Trust



MEDICOM®