

### Safe+Mask® Sof Skin® Procedure Earloop Face Mask Low Barrier Protection (Level 1)

#### Item Code:

2080 (low barrier – white), 2081 (low barrier – blue), 2082 (low barrier – pink) and 2083 (low barrier – lavender)

- Designed with sensitive skin in mind.
- Soft, lightweight, gentle, non-irritating and remains cool on the skin for ultimate ease and comfort

#### Material Description

Made in	USA
Outer material	Spunbond polypropylene
Filter media	Meltblown polypropylene
Inner material	Cellulose Polyester Blend
Number of pleats	3 shingle pleats
Ear loops	Electronically sealed knitted polyester
Nose piece	Aluminum
Embossed	Medicom LB (for Low Barrier)



ASTM 2100 <sup>1,2</sup>	LOW BARRIER
Bacterial Filtration Efficiency (BFE)	≥ 95% BFE @ 3.0 µm large bacteria
Differential Pressure (ΔP)	ΔP < 4.0 mm of H <sub>2</sub> O/cm <sup>2</sup>
Particulate Filtration Efficiency (PFE)	≥ 95% @ 0.1 µm large latex particles
Fluid Resistance	Resistant @ 80mmHg
Textile Flammability	Class I flame spread

#### Other tests<sup>3</sup>

Dermal irritation	Negligible
Dermal Sensitization	Negligible

**Storage** Store under cool, clean, and dry conditions.  
Avoid excessive heat (over 40°C or 104°F)

**Usage** For single-use only

**Packaging** 50 per box / 10 boxes per case / 500 masks per case

**NOTE 1 & 2** All tests were conducted in accordance with ASTM 2100 and the results have an AQL of 4.0%.

**NOTE 3** The values below are results from tests conducted as per ISO 10993 biocompatibility requirements