



# SafeBasics™

Procedure Earloop Face Masks



Your Everyday Necessities, Our #1 Priority.

Health Care professionals will not compromise on quality when it comes to protection. There are many options available, products with added benefits, with varying levels of protection and with many features. However, there are times when a more basic approach to disposables and every day consumables is all that is required.

For Health Care professionals who seek great quality in their essentials from a brand name they trust, Medicom extends its line of SafeBasics masks to Level 2 mask to provide the protection you expect, at a value that can't be beat.

Medicom SafeBasics Masks is a line of ASTM Level 1 & Level 2 masks focused on delivering protection at the best price and offered in a range of everyday colours.

## The perfect combination of price and protection



# SafeBasics™

Procedure Earloop Face Masks



**2130** SafeBasics Masks  
**White** - ASTM 1



**2131** SafeBasics Masks  
**Blue** - ASTM 1



**2132** SafeBasics Masks  
**Pink** - ASTM 1



**2133** SafeBasics Masks  
**Yellow** - ASTM 1



**2134** SafeBasics Masks  
**Green** - ASTM 1

## Performance factors

ASTM  
LEVEL **1**

Bacterial filtration efficiency (BFE)	≥ 95%
Sub-micron particulate efficiency at 0.1 μm (PFE)	≥ 95%
Resistant to penetration by synthetic blood	@80 mm Hg
Differential Pressure (ΔP)	< 4.0 mm H <sub>2</sub> O/cm <sup>2</sup>
Flame Spread	Class 1



**2140** SafeBasics Masks  
**White** - ASTM 2



**2141** SafeBasics Masks  
**Blue** - ASTM 2



**2142** SafeBasics Masks  
**Pink** - ASTM 2

# ALL NEW LEVEL 2

## Performance factors

ASTM  
LEVEL **2**

Bacterial filtration efficiency (BFE)	≥ 98%
Sub-micron particulate efficiency at 0.1 μm (PFE)	≥ 98%
Resistant to penetration by synthetic blood	@120 mm Hg
Differential Pressure (ΔP)	< 5.0 mm H <sub>2</sub> O/cm <sup>2</sup>
Flame Spread	Class 1



## Medicom® Trust in Quality

Like you, we never compromise. Our commitment is to help to make the world safer and healthier by offering consistent, reliable protection. We use carefully selected materials and engineer our products to offer protection you can count on.

As the maker of Safe brands, we take pride in the protection we offer, so you can take pride in yours.



Contact Your Dealer Or  
Visit [www.Medicom.com](http://www.Medicom.com)  
For Samples.